Pitalito, 4 julio de 2025

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MEDICAL ODONTOLOGICAL

OIL FREE MODELS

You have purchased another product with Schulz quality.

A company Schulz S.A., with a quality system certified by **ISO 9001** and environmental management system **ISO 14001**

The design, inspection and manufacturing of your Schulz tank is based on ASME Section VIII, Div. 1 - Code and according to NR-13 MTb of Brazil.

It is a mix of technology and easy use.

IMPORTANT

This Instruction Manual provides important installation, use, maintenance and safety instructions and must be at hand whenever needed. Before operating the equipment, or during maintenance, read this manual and make sure you understand all of its instructions in order to avoid personal injuries or damage to your compressor.

WARNING

When using electric product, basic safety precautions should always be followed to reduce risk of fire, electric shock and personal injuries.

1) INTRODUCTION

WE RECOMMEND READING THIS MANUAL PRIOR TO USING THIS PRODUCT.

It will help you to optimize product performance, guarantee safe use, and guide you in preventive maintenance. Black numbers in the text, show the Main Parts and Their Functions, Chapter 6.

If there is any problem that cannot be solved by the information provided in this manual, please contact the nearest Schulz Authorized Technical/Distributor.

In order to extend the useful life of this reciprocating piston compressor, the use of original parts Schulz are highly recommended.

It is the responsibility of the user to install, inspect, maintain, operate and provide specific documentation for the compressor and air/oil pressure vessel as required by applicable federal, state and local codes.

The pressure vessel quality certificate must be presented to authorities, in case of inspection. Keep it in a safe place along with the instruction manual.

2) COMPRESSOR INSPECTION


- Inspect for signs of obvious freight damage. Report any damage to the delivering freight carrier immediately.
- Be sure that all damaged parts are replaced and any mechanical problems are corrected prior to the operation of the air compressor.
- The compressor serial number is located on the pump. Please write the serial number in the space provided in the end this manual for future reference.

3) USES

Schulz air compressors must only be used for air compression and up to the maximum compression indicated in their Identification Tag/Information Sticker.

4) TECHNICAL DATA

CABINET (OPTIONAL)
Weight = 66.2 lbs (30 kg)

MS 3/30
MSV 6/30

CSA 6,5/30


5) SAFETY INSTRUCTIONS


When improperly used, an air compressor can cause physical and material damage. Please, follow the recommendations below to avoid them :

1. This equipment:
 - has hot and electric parts and moving pieces;
 - the models MS 3, MSV 6 and CSA 6,5 are supplied with plug 2P + GR (Ground). For your safety, do not cut the “GR” pin or use an adapter “2P”;
 - make sure the compressor is correctly grounded;
 - may start or stop automatically, when connected to the electric supply, depending on the tank's pressure or on the functions of the electric protective elements;
 - may cause mechanical or electric interference in other equipment which may be placed nearby;
 - should not operate in places where access to non-authorized people, children or animals is permitted;
 - requires an authorized person to supervise its use and maintenance, and equipment protection must be worn;
 - must be installed and operated in a well ventilated area and with protection against humidity or water.
2. Never exceed the maximum pressure indicated in the compressor's identification tag **13**.
3. Never change the adjustment of the safety valve or pressure switch, since they have already been regulated factory. If an adjustment of the pressure switch is required, look for the nearest Schulz Authorized Technical/Distributor.
4. Never perform welding jobs or repairs in the tank because they may affect its resistance or disguise more serious problems. If there is any leakage, crack or corrosion, immediately stop the equipment use and look for the nearest Schulz Authorized Technical/Distributor.
5. Never run the pressure vessel above the maximum pressure indicated in its identification tag **14**.


MSV 12/175
MSV 18/250

CSV 15/220
CSV 20/220

CSV 15/250
CSV 20/250

- NOTE:**
- 1- The tank filling up time varies from ± 10%, according to installation.
 - 2- MUNSELL 10Y 8.5/2.
 - 3- Electric motor available for working in 115 or 230V (MS 3 and MSV 6).
 - 4- Tank with internal painting (in powder) anticorrosive.
 - 5- Technical data of tank see the pressure vessel quality certificate.
 - 6- The cabinet reduce compressor noise level to less than 61 dB(A) and varies ± 3 dB(A) for models MS 3/30 and MSV 6/30. Other models MS 3/30 (77), MSV 6/30(78), MSV 6/50 (78), MSV 12/175 (84), CSV 15/220 (82), CSV 15/250 (82), MSV 18/250 (82), CSV 20/220 (84,4) and CSV 20/250 (84,4) without cabinet.
 - 7- 25,4mm = 1" (inch).

MODEL	DISPLACEMENT				MAX. PRESSURE		TANK			rpm		ELECTRIC MOTOR				DIMENSIONS (mm)						WEIGHT (kg)	NOISE dB(A) (see note 6)	PAINING (see note 2)
	50Hz		60Hz		psig	bar	Volume gal (l)	Filling up time		50Hz	60Hz	hp	kW	poles	VOLTAGE (V)	A	B	C	D	E	Ø F			
	cfm	l/min	cfm	l/min				50Hz	60Hz															
MS 3/30	2.5	70	3.0	84	120	8,3	8/29	7'10"	6'10"	1,420	1,720	1/2	0.37	4	110/220 Single-phase	500	460	645	-	-	-	34	61/77	PUMP (WHITE)
MSV 6/30	5.0	141	6.0	170			8/29	3'50"	3'10"	1,445	1,730	1	0.75			500	460	630	-	-	-	42	61/78	
MSV 6/50	5.0	141	6.0	170			12/45,7	5'25"	4'35"	1,445	1,730	1	0.75			400	400	620	-	-	-	54	78	
CSA 6,5/30	5.2	148	6.5	185			8/29	4'40"	3'30"	1,440	1,670	1.5	1.1			500	460	630	-	-	-	22	70	
MSV 12/175	10.0	282	12.0	340			47/178	11'50"	10'15"	1,445	1,730	2x1	2x0.75			520	1,080	880	550	460	18	109	84	
MSV 18/250	15.0	425	18.0	509			69/261	9'55"	8'15"	1,445	1,730	3x1	3x0.75			550	1,550	850	800	460	16	160	82	TANK (WHITE IN POWDER)
CSV 15/220	15.0	425	15.0	425			59/224	9'30"	9'30"	620	620	3	2.2	2	(110/220 Single-phase) (220/380/440 Three-phase)	500	820	1,850	-	-	16	167	85	
CSV 15/250	15.0	425	15.0	425			69/261	9'40"	9'40"	620	620	3	2.2			540	1,700	1,032	800	460	16	169	85	
CSV 20/220	20.0	566	20.0	566			59/224	4'	4'	1,050	1,050	5	3.7			500	820	1,850	-	-	16	168	89	
CSV 20/250	20.0	566	20.0	566			69/261	4'50"	4'50"	1,050	1,050	5	3.7			540	1,700	1,032	800	460	16	170	89	


TABLE 1 - TECHNICAL DATA

5) SAFETY INSTRUCTIONS

6. In case of an inspection, the user of this product must keep the pressure vessel book provided by the manufacturer available, attached to the other safety documents required by NR-13 of Labor Ministry of Brazil, for as long as the vase is used, until it is disposed of. The final user must follow the resolutions set forth in NR-13 as to the installation, maintenance, and operation of the pressure vase (compressed air receiver). The useful life of a pressure vase depends on several factors that do contribute to determine it. This aspect must be monitored and established by a certified professional, according to NR-13. Remark: The hydrostatic test carried out during the product manufacturing does not replace the initial inspection, which must be carried out at the location where the product is installed, duly followed up by a certified professional, according to NR-13 MTb. Schulz S.A., manufacturer of the product, declares that the local legislation and inspection rule upon the above provided information, and that all prudent, preventive, and wise safety procedures must prevail.

-  Air compression is a heat generating process. The pump unit **1** and tube **11** are subject to high temperatures. Therefore, one must be careful when handling them in order to avoid burns.
-  The compressed air would contain damages contaminants to the human, animal, ambient or nourishing health, among others. The compressed air must be treated with adjusted filters, as requisite of its application and use. Get more information about it at Schulz Authorized Technical/Distributor or the factory (cta.export@schulz.com.br).
-  Do not carry out any maintenance jobs with the compressor on. Do not remove accessories attached to the compressor when it's pressurized. Do not clean or touch the electric components without unplugging the compressor from the electric network. Not following the orientations may cause physical and material damage.
- The refrigeration air intake on the case and belt cover **20** and cabinet **15** must always be clean to prevent the flywheel and fan from pulling in impurities.
- Never use solvents to clean your air compressor.
- The use of solvents or paint accumulation may cause explosion risks and irreversible damage to the product in general. Cleaning products and the environment may not contain solvents.
- Always use an electric extension cord according to specifications (see Table 2) which must not be amended. If these recommendations are not followed, the user may be injured and the compressor electric components may be damaged, violating the warranty. When necessary, use a longer air hose.
- In case of any irregularities with the equipment, stop operating it immediately and look for the nearest Schulz Authorized Technical/Distributor.

6) MAIN PARTS AND THEIR FUNCTIONS

- Compressor pump unit** - sucks and compresses the atmospheric air.
- Air receiver (Tank)** - stores the compressed air.
- Motor** - makes the compressor unit run (thermally protected) (MS 3, MSV 6, CSA 6,5, MSV 12 and MSV 18).
- Electric motor** - makes the compressor unit run with to help of the pulley and belt (CSV 15 and CSV 20).
- Pressure switch** - controls the compressor's operation to avoid that the maximum working pressure is surpassed. See Installation Chapter 7, item 4 - Electrical Lead.
- ASME Safety valve** - exhausts the tank in case the pressure is above the maximum allowed (PMTA).
- Check valve** - keeps the compressed air inside the tank avoiding its return when the compressor unit stops.
- Drain petcock** - removes the condensed water inside the tank.
- Pressure gauge** - shows the pressure inside the tank in psig or bar.
- Pressure regulator** - used to adjust working pressure and control the compressed air release (see note).
- Discharge tube** - carries and cools down the compressed air.
- Intake air filter** - filters atmospheric air impurities that are sucked by the compressor, see chapter 7, item 1.
- Identification Tag/Information** - sticker shows the technical data of the pump.
- Air compressor's identification Tag** - shows the technical data in the compressor's tank.
- Cabinet** - used to lessen the compressor's noise (optional) for model MS 3/30 and MSV 6/30 only.

6) MAIN PARTS AND THEIR FUNCTIONS

- Power cable** - used to connect the compressor to electric supply (see note).
- ON/OFF button** - starts the compressor and the ventilation system of the cabinet.
- Terminal box** - used for wiring and electrical connections (MSV 12/175 and MSV 18/250). Note: not applicable for Argentinian market.
- Thermal protector** - protects the motor from overload (MS 3, MSV 6, MSV 12 and MSV 18).
- Case/belt cover** - allows better cooling of the compressor unit and protects the moving parts.
- Start key** - ON/OFF starts the compressor for models CSV 15 and CSV 20 and model MSV 18 for Argentinian market.
- Anti-vibration** - it isolates the vibration produced by the compressor to the floor.
- Hose** - transports the compressed air supplied by the compressor MSV 6/50 for the tank.
- Crankcase filter** - it helps in the breath of the air in the interior of the crankcase.

Note:

- Item **15** is optional and adaptable only for models MS 3/30 and MSV 6/30.
- For a specific compressed air quality, special filters must be installed to remove particles, condensed water and smell. Further information can be obtained at the nearest Schulz Authorized Technical/Distributor or cta.export@schulz.com.br.
- The item **16** is supplied for the models MS 3/30, CSA 6,5/30, MSV 6/30 and MSV 6/50 with plug 2P + GR (ground) see Figure 6. For the models MSV 12/175 and MSV 18/250 without the plug, its must be connected to electrical source (electrical starter and fusible) or according with the electrical laws for each country where the compressor will be installed. For the models CSV 15/220, CSV 15/250, CSV 20/220 and CSV 20/250 is responsibility of the user the installation of the compressor to electrical source, including the electrical starter.
- For the products CSV 15, MSV 18 and CSV 20 is responsibility of the user/owner installs an air control valve (register) in the connection of exit of the tank (1/2").

7) INSTALLATION

- Location:** Install the compressor in a covered, well ventilated area, free of dust, gases, toxic gases, humidity or any other kind of pollution and for model MSV 6/50 - Split, see Figure 20. The maximum ambient temperature recommended during work is 104°F (40°C).

IMPORTANT Avoid place such as: store, store-room, basements, garage and bathroom.

- Position:** Mount at least 32 inches (800 mm) from any wall or obstacle and with minimum two (2) meters of height, to allow good ventilation during operation and to make maintenance easier. When cabinet **15** is used, the minimum distance should be 4 inches (100 mm).

2.1 Compressor installation: When connected to the compressed air network, the discharge connection (1/4") and 1/2" (CSV 15, MSV 18/250 and CSV 20) must be joined by means of a hose or expansible gaskets so that the charges, thermal expansion, piping weight, mechanical and thermal shocks or obstructions are not transmitted to the air tank. Physical damage to the compressed air tank may occur if the above warnings are not followed. Instruction assembly rectangular filter (see Figure 5).

Note: Model MSV 6/50 - Split is supply of factory with a hose 5/16" x 787" (8mm x 20 m). Other distance, use a hose of 1/2".

- Foundation/floor:** It must be placed on a concrete base, which must be able to bear the equipment weight plus the liquid (water) when of the hydrostatic test. The set must be level. The compressor models is equipped with rubber foot (assembled by owner) see Figures 2, 3 and 4 and Figure 1 for model MSV 12/175, place de compressor on anti-vibration (owner's responsibility), absorbs which, in turn, must be replaced on a concrete basis. The set must be level and the absorber must be adequate to support the equipment weight and vibration. **Note:** Use a little water to assembly the absorber for models MS 3/30, CSA 6,5/30, MSV 6/30 and MSV 6/50, Figure 3.

4. Electrical lead:

- Consult a specialized technician to evaluate the electric supply and to select the adequate protection and input devices.
- Follow the technical standards of the country on Low Voltage Electric Installations.
- Figure 7, shows a electrical diagram for electric connections with pressure switch and fuse.
- The electric input wires must be in accordance to the motor power, voltage and electric energy source distance. See table 2 for guidance.
- Before plugging in the equipment, check if the voltage shown in the input wire tag is the same as the local voltage.


FIGURE 1


FIGURE 2

FIGURE 3

FIGURE 4


FIGURE 5

MOTOR POWER (hp)	ELECTRICAL VOLTAGE (V)	MOTOR CURRENT (A)	CONDUCTOR (mm²)	MAX. DIST.in (m)	F1 FUSE MAX. (A)	Diagram (Figure)	
				VOLTAGE DROP (2%)			
SINGLE-PHASE MOTOR							
1/2	110	7,6	1,5	7	10	7	
	220	3,8	1,5	31	6		
1	110	15	2,5	7	25		
	220	7,5	1,5	17	10		
1,5	110	8,4	2,5	13	10		
	220	4,4	1,5	27	6		
2 x 1	110	30	6,0	8	32	10	
	220	15	2,5	14	25	7	
3 x 1	220	22,5	4,0	15	25	10	
	3	110	32	10,0	13		50
		220	16	4,0	21		25
THREE-PHASE MOTOR							
3	220	8,1	1,5	25	16	9	
	380	4,7	1,5	75	6		
	440	4,1	1,5	101	6		
5	220	14,4	4,0	38	25		
	380	8,3	1,5	42	16		
	440	6,3	1,5	65	10		

TABLE 2 - GUIDING DATA ABOUT FUSES AND COPPER CONDUCTORS

4.7

Grounding Instructions:

This product must be grounded. In the event of an electric fault, grounding reduces the risk of electric shock by providing an escape wire for the electric current.

DANGER

The incorrect installation of the grounding wire may result in an electric shock. If it is necessary to replace or repair the cable or the connector, do not lead or join the grounding wire to the neutral wire or other. The green wire, with or without yellow stripes, is only to the grounding function. In case of doubts regarding the grounding information or whether the product is properly grounded, make sure you contact a qualified electrician to verify the connections.

110 V		220 V	
L1	3 1	L1	1
	2 7		2 3
L2	4	L2	7
			4

LEGEND WIRE OF COMPRESSOR	
Nº	COLOR
1	BLUE
2	WHITE
3	ORANGE
4	YELLOW
7/P2	BROWN
-	GREEN (WIRE GROUND)

TABLE 3 - CHANGE OF TENSION (V) AND LEGEND WIRE MS 3, MSV 6, MSV 12 AND MSV 18

7) INSTALLATION


FIGURE 6 - FOR MODELS MSV3, MSV 6 AND CSA 6,5

F1	Fuse type "D" or "NH" (see Table 2)
F2	Command fuse 2A
P	Pressure switch
M	Electric motor
FT1	Overload relay
KSFF	Sequence and phase fault relay
K1	Contactor
SH1/CH	ON/OFF switch
XS	Connector


FIGURE 7 - ELECTRIC DIAGRAM (MOTOR WITH THERMAL PROTECTOR)


FIGURE 8 - ELECTRIC DIAGRAM IN 220V - MSV 18/250



FIGURE 9 - THREE-PHASE MOTOR - DIRECT START

* Obs.: The kit K1 + FT1 can be two poles depends on product version.

FIGURE 10 - ELECTRIC DIAGRAM SINGLE-PHASE

4.8 Box connections compressor MSV 12/175


FIGURE 11 - ELECTRICAL DIAGRAM WITH BOX CONNECTIONS (CHANGE OF VOLTAGE)

4.9 For your safety, an on/off switch must be installed (before the fuse box) to switch off electric power during maintenance or equipment adjustments (must be according to current of electric motor, Table 2).

ATTENTION

For maintenance in warranty, the electrical installations or to change the voltage of the machine, it must be realized by a technical specialized electrical that must follow the diagrams of this manual. The expenses are responsibility of the owner/client, if you have any doubt, a Schulz Authorized Technical/Distributor must be consulted.

5. Cabinet installation procedure:

- Before connecting the cabinet to the electric supply, check if the indicated voltage in the input wire tag is the same as the voltage indicated in the compressor.
- Lift the cabinet.
- Place it over the compressor.
- Lower it toward the compressor until it touches the floor. Check that the rear lid of the electric motor perfectly fits the left side hole inside the cabinet.
- Connect the compressor wire 16 or the plug inside the cabinet.
- To energize the set (cabinet + compressor), connect input cabinet wire to the electric supply.

Note:

- Models MSV 12/175 and MSV 18/250 are factory supplied with electric wire without the plug.
- Select the plug or the electric connection device according to the electrical parameters of the product, see Table 2.
- The energy supply network must not present voltage changes over $\pm 10\%$.
- The voltage drop caused by the motor start must not be over 10%.
- The compressor acoustic cabinet, with vertical receiver volume of 30 liters, is factory supplied with 115 or 230 (V) electric connection. If you need to change the voltage, check the diagram placed inside the cabinet, which will guide you to change the voltage.
- The thermal protector 19 is located inside the electric motor for models MS 3, MSV 6, CSA 6,5, MSV 12 and MSV 18. If it goes off, it will re-start the compressor after the temperature drops. The installation of protection fuses is necessary, see Table 2.
- Used start key and fuse kind delayed-action.

- INSTALLATION EXPENSES ARE OF PURCHASER'S RESPONSIBILITY.


8) START PROCEDURE

FIRST START PROCEDURES

After the arrangements regarding location, the pressure vessel adequacy to local Technical Norms and Legislations, the installation of the electric network and the compressed air network are made (customer's responsibility), do the following:

- Completely open the pressure regulator 10 or air control valve.
- Turn be start key (assembly by customer for compressor without key) and check if the rotation direction is anti-clockwise, viewed from the fan.
- Let the compressor run for about 10 (ten) minutes to allow a homogeneous lubrication of the moving parts.
- Totally close the pressure regulator or air control valve, so that the compressor can fill up the tank. The compressor will automatically turn off (through the pressure switch, Figure 15) when the pressure gauge shows maximum pressure around 120 psig (8.3 bar).
- Open the pressure regulator or air control valve to release the compressed air from the tank, to decrease its internal pressure. The compressor will automatically re-start (through pressure switch) when the pressure gauge shows pressure around 80 psig (5.5 bar).
- Important:** Check the air tank's filling up time according to Table 1.
- Close the pressure regulator and unplug the equipment.
- Verify if safety valve 6 (Figure 12) is working.
- Check if the internal fan inside the acoustic cabinet continues turned on after the compressor shuts off at 8.3 bar. It must be on, in order to maintain the internal temperature equal to the environment temperature.

Note: After the job finished, the compressor's tank pressure is between 5.5 and 8.3 bar. Switch off the drive button on the cabinet 17. It will help to save electric power.

- Open the drain petcock 8 to exhaust the condensed water inside the tank and close it right after (Figures 13 and 14).
- Open the regulator or air control valve to completely exhaust the tank and close it right after.
- Your Schulz compressor is ready to operate. Connect it to the air distribution network and start the electric motor. When the compressor reaches maximum pressure, open the pressure regulator so that the air can flow into the distribution network. Check if there is any leakage along the tubing, by using water and soap lather. In case there is any, eliminate it.

ATTENTION: A compressor with good dimensions will have roughly six (6) starts per hour 70% ON and 30% OFF/unload.

Note:

- The initial start procedure must be repeated whenever maintenance occurs or when the compressor is moved to another place.
 - The pressure regulator 10 must be used as follows: Pull knob turning right or left to set desired pressure, push knob down. The desired work pressure will be shown in the regulator sight.
 - For your comfort, Schulz manufacturer of the product compressor provides an electronic drain, model PS 16, which is adaptable to the tank exhaust connection (drain) and that can be purchased at Schulz Authorized Technical/Distributor.
- Important:** When installing it, the tank must be depressurized (empty).
Attention: See item 9 of Chapter 5 - Safety Instruction.

9) PREVENTIVE MAINTENANCE

ATTENTION In order to guarantee the perfect compressor operation and to extend its lifetime, follow the recommendations below:

1. Daily

- Drain the condensed water from the tank, through petcock 8 (Figures 13 and 14).

ATTENTION

RISK OF PERSONAL INJURIES
The maximum pressure of the tank when drained must be 10 psig (0.7 bar).

- Check abnormal noise or vibration in the compressor. If the problem persists after the following corrective actions had been taken, get in touch with the nearest Schulz Authorized Technical/Distributor.

2. Weekly

- Clean the external parts of the compressor with neutral detergent or with a wet cloth.
- Verify the operation of the safety valve 6 (Figure 12).
- Inspect the intake filter 12, (Figure 5 and 16 for models CSV 15 and CSV 20) or when supply in special condition (Figure 17) when the filter element clogged is necessary change.
- Check if the belt, located inside the belt protector 20, (CSV 15 and CSV 20), is making abnormal noise or movements. If so, see item (4C) procedure below.

3. Montly

- Check the operation of the pressure switch (Figure 15) (see items 4 and 5 the Chapter 8 - Start Procedure).


FIGURE 15

FIGURE 16


FIGURE 17


FIGURE 18

4. Quarterly

- Change the intake filter element every 300 (three hundred) hours or quarterly (whichever occurs first).
- Tighten pump unit bolts (see Table 4) with a torquemeter and nuts with manual tools. For model CSA 6.5/30 check if compressor makes abnormal noise in the inner part of the case.
- Verify the tension of the belt (Figure 18) and the alignment (Figure 19) for models CSV 15 and CSV 20. The tension at the middle point yields 3/8" through 1/2" on manual force being or weight of 4,0 kg (1,8 lbs) applied at the mid-point. In order the check the deflection, remove the belt cover. After the inspection reassemble the protector, fixing the screws as per Table 4.
- Inspect the intake crankcase filter 24, when the filter clogged clean it using a maximum pressure of 2,0 bar (30 psig) or change when damaged.

5. Every six months

- Lubricate connecting rod(s) upper bearing, with Barrierta L55/2 (Klüber) for models MS 3, MSV 6, MSV 12 and MSV 18.

6. Every 1,000 hours or 9 months (whichever occurs first)

- Lubricate connecting rod(s) upper bearing, with Barrierta L55/2 (Klüber) for models CSV 15 and CSV 20.
- Inspect and clean the valves (located between the cylinder's upper part and its head).
- Replace connecting rod(s) lower bearing for models MS 3, MSV 6, CSA 6.5, MSV 12 and MSV 18.
- Replace rings on the piston [if the filling up time is 20% higher than what is specified in Table 1]
- Inspect and clean the check valve *7.

* Attention: See item 9, Chapter 5 - Safety Instruction.

7. Yearly

- Have the pressure switch, pressure gauge and safety valve, according to local technical norms. This operation must be carried in a device not attached to the air tank.


FIGURE 19

9) PREVENTIVE MAINTENANCE

NOTE
1. Items 5A and 6A/B/C and D must be carried out by the nearest Schulz Authorized Technical/Distributor .

IMPORTANT

The above instructions are based on normal operation conditions. If the compressor is installed in a polluted area, increase inspection intervals.

PLACE	FAN		CYLINDER COVER			CRANCKCASE COVER		CYLINDER TO CRANCKCASE		FLANGE	FLYWHEEL	CONNECTING ROD	CYLINDER COVER CENTER	CRANCKSHAFT COUNTERWEIGHT
MODEL	CSA 6.5	MS 3 MSV 6 MSV 12 MSV 18	CSA 6.5	MS 3 MSV 6 MSV 12 MSV 18	CSV 15 CSV 20	MS 3 MSV 6 MSV 12 MSV 18	MSV 6 CSV 15 CSV 20	MS 3 MSV 6 MSV 12 MSV 18	CSV 15 CSV 20	CSV 15 CSV 20	CSV 15 CSV 20	CSA 6.5	CSV 15 CSV 20	CSA 6.5
N x m	6.0*	8.0*	6.0	17.0	25.0	6.0	17.0	30.0	25.0	25.0	4.0	12.0	8.0	

* When refasten bolt to use in the thread “Loctite 243-Blue” .

TABLE 4 - BOLTS TORQUE SPECIFICATION

IMPORTANT

Switch off the air compressor, remove the water and exhaust the air from the tank at the end of every working day.

10) CORRECTIVE MAINTENANCE

! To assure product SAFETY and RELIABILITY, repairs, maintenance and adjustment should be performed by Schulz Authorized Technical/Distributor always using identical replacement original parts Schulz.

11) TROUBLESHOOTING TIPS

Many times, what initially looks like a defect can be solved without the help of Schulz Authorized Technical/Distributor. If the problem persists after the following corrective actions were taken, get in touch with the nearest Schulz Authorized Technical/Distributor.

PROBLEM	PROBABLE CAUSE	WHAT TO DO
Motor doesn't start or doesn't re-start. Obs.: Don't insist on starting up the motor again. Do it only after you discovered and eliminated the cause of the problem.	Voltage drop or electric supply is out.	Check the installation and/or wait for the electric supply stabilization.
	Damaged motor.	Send it to nearest Schulz Authorized Technical/Distributor.
	Thermal protector went off (MS 3, MSV 6, CSA 6.5, MSV 12 and MSV 18).	Motor will re-start as soon as temperature decreases. In case it happens again, send compressor to a Schulz Authorized Technical/Distributor.
	Tank is full of air.	Motor will re-start as soon as pressure in the tank decreases (up to 4 bar (60 psig)).
	Start key with FT1 activate.	Reset the button.
Reduced air production (pressure gauge keeps on showing a lower pressure than that of the job).	Check valve does not seal because of impurities.	Send the compressor to the nearest Schulz Authorized Technical/Distributor.
	The demand is above the compressor's capacity.	Evaluate compressor's capacity.
	Filter element clogged.	Replace it.
The compressor unit overheats.	Air leakage in the compressor.	Re-fasten bolts and/or fittings.
	High ambient temperature (max. 40°C - 104°F).	Improve installation conditions.
	Incorrect rotation direction. (Correct rotation direction isanti-clockwise, viewed from the fan and flywheel).	Take compressor to the nearest Schulz Authorized Technical/Distributor.
Premature wear of the compressor unit internal components.	Demand above compressor's capacity.	Evaluate compressor's capacity.
	Filter element clogged.	Replace it.
	Operating in a non-adequate environment.	Improve local conditions and take compressor to the nearest Schulz Authorized Technical/Distributor.
Abnormal noise or vibrations.	Loose fastening elements.	Find and re-fasten them.
	Loose fan and flywheel.	Re-fasten it.
	Compressor installed incorrectly.	Verify item 3 - Chapter 7 - Installation.
	Damaged connecting rod(s) bearing and rings.	See Chapter 9 - Preventive Maintenance, items 5 and 6.
	Loose or unaligned belt (CSV 15 and CSV 20).	Adjust it (them) according to Figures 18 and 19.
Motor overheats.	Deficient or inadequate electric supply.	Have a specialized electrician check it.
Tank's pressure goes up too fast or many frequent starts (the normal number of starts is around 6 an hour).	Water excess in the tank.	Drain the tank as indicate in Figures 13 and 14. For your comfort install an electronic drain.
High temperature of the compressor unit.	Too much dust on the pump.	Clean the pump externally.
	High ambient temperature (max. 40°C - 104°F)	Improve installation conditions.
	Filter element clogged.	Replace it.
Compressor operates uninterruptedly.	The system's demand is above compressor's capacity.	Evaluate compressor's capacity.
Insufficient pressure for the job.	Regulator not adjusted.	Adjust it, see note the item 2 on Chapter 8 - Start Procedure.
Compressor does not start	Pressure switch with mechanism actuated or damaged.	Depressurize the tank to 60 psig (4 bar). If the compressor does not restart, call a Schulz Authorized Technical/Distributor.

IMPORTANT

Further information can be obtained at the nearest Schulz Authorized Technical/Distributor or (cta.export@schulz.com.br).

12) MSV 6/50 - SPLIT INSTALLATION


FIGURE 20

13) ENVIRONMENTAL GUIDANCE AND RECOMMENDATIONS

1. Disposal of Liquid Effluents

The presence of liquid effluents or non-treated condensation from tank in rivers, lakes or in other water receiving bodies may adversely affect the aquatic life and the water quality as well. The condensation withdrawn from the tank, daily, according to the Preventive Maintenance Chapter 9, must be kept in a container and/or in an appropriate collecting network for further treatment. Schulz S.A., recommends that the liquid effluent produced inside the receiver of the compressor should be adequately treated through processes that aim at protecting the environment and the healthy quality of life of the population, complying with the country current regulation requirements. Among the treatment methods available, one may choose the physical-chemical, chemical, and biological ones. The treatment may be carried out by the company itself or by outsourcing.

2. Disposal of Solid Waste (parts in general and product packages)

The generation of solid waste is an important aspect that must be considered by the users when using and maintaining their piece of equipment. The impacts to the environment may cause meaningful changes in the quality of the soil, in surface and underground water, and in the population health, due to the inadequate disposal of the discarded residues (on streets, water springs, landfills, etc). Schulz S.A., recommends that the waste arising from the product, from its generation, handling, transportation, and treatment to its final disposal should be handled carefully. Appropriate handling should consider the following steps: quantification, qualification, classification, reduction at source, pick-ups and selective pick-ups, recycling, storage, transport, treatment and final destination. The disposal of solid waste must be carried out according to the regulation requirements of the current legislation of the country the product has been exported to.

14) LIMITED WARRANTY

Schulz S.A., assures this product's first buyer/user a warranty against any manufacturing fault that develops within 1 (one) year (including the legal warranty - first 90 (ninety) days) from date of purchase invoice issuing.

GENERAL WARRANTY CONDITIONS

- A.** The warranty servicing will be made upon purchase invoice presentation. If the client/owner of this product wishes to be attended to at home, an Schulz Authorized Technical/Distributor may charge a visit free or in accordance to local norm and legislation.
- B.** Any warranty service should be realized exclusively by Schulz Authorized Technical/Distributor.
- C.** Parts that are naturally worn out by regular use and that are affected by the installation and the way the product is used are excluded from this warranty, such as: air filter element, gaskets, valves, rings, cylinders, pistons, connecting rods, crankshaft, pressure switch, belt, star key, roller bearings, regulator valve and drain petcock. Schulz S.A., will only take responsibility on the servicing expenses of the above mentioned parts if an Schulz Authorized Technical/Distributor detects manufacturing defects.
- D.** This warranty will not apply to cleaning and installation services, lubricate bearing, neither to any damage caused to the product external parts due to incorrect use, neglect, modifications, use of inadequate accessories, inadequate use, falls, drillings, use not according to the owner manual, electrical connections in the wrong voltage, voltage conversion of the electric motor inadequate for the compressor or electrical supply nets subject to excessive voltage drops or voltage overloads.
- E.** Electric motor and electric key manufacturers assume the warranty of their products if they fail for manufacturing defects. Defects caused by bad installation are not covered by this warranty.
- F.** Representatives or dealers are not authorized to receive a product to send it to a Schulz Authorized Technical/Distributor, take it from there to return it to the client or provide information on Schulz S.A. behalf about the servicing.
Schulz S.A. or any Schulz Authorized Technical/Distributor will not be responsible for eventual damages or delays resulting from the non-accomplishment of this clause.
- G.** Only original Schulz repair for maintenance parts are used.
- H.** Will be excluded of warranty any repair or indemnification of damages occurred during the product transport.

WARRANTY TERMINATION

This warranty will be considered void in the following situations:

- A.** At the end of its validity.
- B.** When the product is delivered for servicing to non-authorized shops; signs of violation to its original features can be noticed; when it is assembled out of the manufacturer's standards.

WARNINGS

- A.** This compressor technical data were customized: displacement, tension (V) and operation pressure were specified in the purchase order at Schulz distributor.
- B.** Lubrication of your reciprocating piston compressor oil free is fundamental. Maintenance of proper is paramount and has been discussed numerous times in this manual.
- C.** The deadline of a servicing will be shown in the service report, on the date a product is delivered at a Schulz Authorized Technical/Distributor.
- D.** It's the purchaser responsibility to pay for service expenses, which are proved unfounded.
- E.** Dealers, representatives or Schulz Authorized Technical/Distributor are not authorized to change this Term or take responsibilities on Schulz S.A. behalf.
- F.** Figures and photos are only for orientation.
- G.** The technician will provide service on a come first served basis.

Note: Schulz S.A., reserves the right to make technical and design changes without prior notice.

Write here the serial number of product:

Rua Dona Francisca, 6 901
Fone (55) (47) 3451 6000 Fax (55) (47) 3451 6055
CEP 89219600 JOINVILLE SC BRASIL
schulz@schulz.com.br
www.schulz.com.br

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© 3M 2021. All rights reserved. 3M, Littmann, the L Littmann logo, Classic III, Cardiology IV, Master Classic II and Master Cardiology are trademarks of 3M. Used under license in Canada. Eko trademarks, including CORE, are marks and/or registered marks of Eko Devices Inc. All other trademarks are property of their respective owners. OMG154398