KARL STORZ — ENDOSKOPE

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## Controls, displays, connectors, and their uses

1 Power switch, lights up when in operation

2 Light outlet point (not exchangeable)
3 Button to increase brightness *
4 Button to reduce brightness *
5 Current light intensity in percent
6 Standby button **
7 SCB indicator display ***
8 to e - see next page (symbols)
r Potential equalization connector
t Fuse holder
z Power cord socket
u Two SCB ports (SCB®, see next page, (symbols))

* Keeping the buttons depressed engages continuous adjustment of the light intensity (autorepeat function).

** If you press the Standby button 6 once, the brightness is reset to the lowest value. The integrated LED illuminates. When you press this button again, the brightness is reset to the last valid value. The integrated LED extinguishes.

***If there are several power LED 175 light sources available in the same SCB system at the same time (see next page), the system assigns each power LED 175 light source its own number (1-9), provided that doing so is supported by the infrastructure. IMAGE1 S™ can only support one light source. Up to 9 light sources can be supported in the OR1™ system. This number is displayed at position 7 and makes it clear to the user on which unit he is making the changes when using the remote control via the SCB.

Example: The user changes the light intensity of LED light source no. 2 from '25 %' to '50 %' in the camera head's 'SCB menu'. This changes the light intensity of the power LED 175 for which a '2' is displayed at position 7.

If only one power LED 175 light source is connected to the same SCB system, no number is displayed. Similarly, no number appears if there is no SCB system present.

If the light source is being remote controlled by means of a CCU, an 'A' is shown.

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- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part type CF

**According to IEC 60601-1-2:**
Observe the information on electromagnetic compatibility in chapter 10 (p. 33-49).

**Directive compliance** (for **20** 1614 20-1)
**According to the Medical Device Directive (MDD):**
This medical device belongs to Class I
This medical device bears the CE mark in accordance with MDD 93/42/EEC.

**Conformidad con la norma** (para **20** 1614 20-1)
**Según CEI 60601-1, CEI 60601-1-6, ANSI/AAMI ES60601-1, CEI 60601-2-18, UL 60601-1, CAN/CSA-C22.2 No. 60601-1:**

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## Electromagnetic Compatibility (EMC) Information

The cold light fountain power LED 175 20161420-1 is suitable for use in professional healthcare facility environment.

Professional healthcare facility includes physician offices, dental offices, limited care facilities, freestanding surgical centers, freestanding birth centers, multiple treatment facilities, hospitals (emergency rooms, patient rooms, intensive care, surgical rooms, outside the RF shielded room of an ME system for magnetic resonance imaging).

***WARNING:*** *Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.*

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| Accessories and cables used for EMC Compliance | | | | |
|---|---|---|---|---|
| **Type** | **Shielded** | **Length [m]** | **Ferrite** | **Used for** |
| PA (Potential equalization) | No | > 3 | No | Potential equalization |
| Power cord | No | 3 | No | Power supply |

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