SIEMENS

OTXWG17E0501A (426)

Coagulation Factor VIII, IX, XI and XII Deficient Plasmas

FACTOR VIII DEFICIENT, FACTOR IX DEFICIENT, FACTOR XI DEFICIENT, FACTOR XII DEFICIENT

See shaded sections - updated information versus edition from August 2003

Intended Use

In vitro diagnostic reagents for the determination of the activity of coagulation factors VIII, IX, XI and XII in human plasma by coagulometric methods.

Summary and Explanation

The determination of coagulation factors VIII, IX, XI and XII in plasma is indicated in the following cases:

- for clarifying the cause of a prolonged partial thromboplastin time,
- for diagnosing congenital[1] or acquired factor deficiency states,
- for distinguishing between dysproteinemias and protein synthesis disorders (in conjunction with immunochemical methods).

Furthermore, the determination of the activity of coagulation factors VIII and IX is employed for monitoring substitution therapy with Factor VIII or Factor IX concentrates in hemophilia A or B[2]. The determination of Factor IX is also of importance in the diagnosis of consumption coagulopathy, hepatic cirrhosis and to achieve more precise monitoring of oral anticoagulant therapy[3]. Detection of an activation of the contact phase of blood coagulation can be achieved by measuring the activity of Factors XI and XII[4].

Principle of the Method

A plasma deficient in any of the factors comprising the intrinsic pathway will result in a prolonged partial thromboplastin time (APTT). Coagulation factor deficient plasma can be used to confirm a factor deficiency, in general, and to identify and quantify coagulation factor deficiency in patient plasma. A mixture of the respective factor deficient plasma and the patient plasma is tested in the APTT assay, and the result is interpreted using a reference curve obtained with dilutions of Standard Human Plasma or a normal plasma pool mixed with the deficient plasma. A patient plasma deficient in a specific factor will not be able to compensate for the absence of the factor in the corresponding coagulation factor deficient plasma and therefore result in a prolonged APTT.

Reagents

Materials provided

Coagulation Factor VIII Deficient Plasma, REF *OTXW,*
8 x ⟶ 1 mL FACTOR VIII DEFICIENT, Coagulation Factor VIII Deficient Plasma or
Coagulation Factor IX Deficient Plasma, REF, *OTXX,*
8 x ⟶ 1 mL FACTOR IX DEFICIENT, Coagulation Factor IX Deficient Plasma or
Coagulation Factor XI Deficient Plasma, REF, *OSDF,*
3 x ⟶ 1 mL FACTOR XI DEFICIENT, Coagulation Factor XI Deficient Plasma or
Coagulation Factor XII Deficient Plasma, REF, *OSDG,*
3 x ⟶ 1 mL FACTOR XII DEFICIENT, Coagulation Factor XII Deficient Plasma

Composition

Coagulation Factor Deficient Plasmas are lyophilized human plasmas with a residual Factor VIII, IX, XI or XII activity of ≤ 1 %. The deficient plasmas are manufactured by immunoadsorption from normal plasma and are free from antigen of the respective factor. Fibrinogen is present in a quantity of at least 1 g/L, and the remaining coagulation factors are present in an activity greater than 40 % of the norm. The plasmas contain mannitol (20 g/L) as a stabilizer.

Warnings and Precautions

1. For *in vitro* diagnostic use only.
2. Each donor or donor unit was tested and found to be negative for human immunodeficiency virus (HIV) 1 and 2, hepatitis B virus (HBV) and hepatitis C virus (HCV) using either tests found to be in conformance with the In Vitro Diagnostic Directive in the EU or FDA approved tests. Because no known test can offer complete assurance of the absence of infectious agents, all human derived products should be handled with appropriate caution.

Preparations of the Reagents

Deficient Plasmas: Dissolve the contents of a vial with 1 mL of distilled or deionized water. Before use, let stand at least 15 minutes at +15 to +25 °C, then shake carefully to mix (without foam formation). Mix carefully once more before using.

APTT reagents: Use according to the respective Instructions for Use.

Calcium chloride solution 0.025 mol/L: warm to +37 °C (not required for automated coagulation systems with heated reagent probes).

Storage and Stability

If stored unopened at +2 to +8 °C, Coagulation Factor Deficient Plasmas can be used up to the expiry date indicated on the label.

Stability after reconstitution:

at +15 to +25 °C	8 hours
at -20 °C	4 weeks

Coagulation Factor Deficient Plasmas can be frozen and thawed once after reconstitution without a loss of coagulation activity. The plasma must be well sealed and frozen as quickly as possible. Thawing should be accomplished at +37 °C within 10 minutes. Thawed plasma should be used within 2 hours when held at +15 to +25 °C.

Information about on-board stability is specified in the Reference Guides (Application Sheets) for the different coagulation analyzers.

Materials required, but not provided

Pathromtin* SL, REF OQGS or
Dade® Actin® Activated Cephaloplastin Reagent, REF B4218-1, -2 or
Dade® Actin® FS, REF B4218-20, -100 or
Dade® Actin® FSL, REF B4219
Calcium chloride solution 0.025 mol/L, REF ORHO
Imidazole Buffer Solution, REF OQAA or
Dade® Owren's Veronal Buffer, REF B4234 or
Dade® CA System Buffer, REF B4265 or
Physiological Saline Solution
Standard Human Plasma, REF ORKL
Control Plasma N, REF ORKE
Control Plasma P, REF OUPZ

Equipment

Coagulation Factor Deficient Plasmas can be used manually or on automated coagulation analyzers. Siemens Healthcare Diagnostics provides Reference Guides (Application Sheets) for several coagulation analyzers. The Reference Guides (Application Sheets) contain analyzer/assay specific handling and performance information which may differ from that provided in these Instructions for Use. In this case, the information contained in the Reference Guides (Application Sheets) supersedes the information in these Instructions for Use. Please also consult the instruction manual of the instrument manufacturer!

Specimens

To obtain the plasma, carefully mix 1 part sodium citrate solution (0.11 mol/L) with 9 parts venous blood, avoiding the formation of foam. Centrifuge the blood specimen at 1500 x g for no less than 15 minutes at room temperature[5].

Stability of the samples:

+15 to +25 °C	3 hours
-20 °C	4 weeks, for Factor VIII determination, 2 weeks


Plasma stored at -20 °C is to be thawed in a water bath within 10 minutes at +37 °C, mixed gently and then tested immediately. If testing cannot be performed immediately, the specimen may be held for a maximum of two hours at +4 °C until tested[5].

Procedure

Manual Testing (For Automated Testing see „Equipment")

Establishment of the reference curve:

Use either Standard Human Plasma or fresh citrated pooled plasma from no less than 10 healthy donors. Using Imidazole Buffer Solution, prepare dilutions as shown in the following scheme and determine the coagulation times as described under "Work flow". The reference curve must be re-generated if there is a change in the instrument or in the lot of reagent used, or if there is any change in the experimental conditions.

Additional dilutions can be prepared if necessary. Plot the measured coagulation times (ordinate) on log graph paper against the corresponding percentage factor activities (abscissa). The coagulation times are dependent on the measurement principle and thus also on the coagulation analyzer used. Therefore, each laboratory must establish its own reference curve.

Work flow:

Dilute the plasma sample 1:5 in Imidazole Buffer Solution (bring to room temperature (+15 to +25 °C)) immediately before the measurement.

Pipette into a test tube prewarmed to +37 °C:

Deficient Plasma	100 µL
Sample dilution	100 µL
APTT reagent	100 µL
Incubate at +37 °C	2 minutes
Calcium chloride solution (prewarmed to +37 °C)	100 µL

On adding the calcium chloride solution, start the stopwatch or timer on the coagulometer and determine the coagulation time.

Results

Read the coagulation factor content from the reference curve in % of the norm. If the given nominal value of the Standard Human Plasma is not 100 % of the norm, but e.g. only 95 %, multiply the result read from the curve by 0.95. In the case of coagulation times which correspond to a coagulation factor content of more than 100 % of the norm, further determination will be required using higher dilutions of specimen (e.g. 1:10). The percentage of norm value read from the reference curve for such a higher dilution must be multiplied by a correction factor corresponding to the dilution; e.g. for a dilution of 1:10 by a correction factor of 2.

Internal Quality Control

Normal range:	Control Plasma N
Pathological range:	Control Plasma P

Two levels of quality control material (normal and pathologic range) have to be measured at start of the test run, with each calibration, upon reagent vial changes and at least every eight hours on each day of testing.The controls should be processed just like the samples. Each laboratory should establish its own quality control range, either by means of the target values and ranges provided by the manufacturer of the controls or by means of its own confidence ranges established in the laboratory. If the control values lie outside the confidence range determined beforehand, then the reagent, calibration curve and coagulation analyzer should be checked. Patient results must not be released until the cause of the deviation is identified and corrected.

Limitations of the Procedure

Therapeutic doses of hirudin or other direct thrombin inhibitors lead to an erroneously lower factor activity[6,7].

Specific inhibitors against plasmatic coagulation factors may also modify the real factor activity[8].

Partial activation of the coagulation factors due to incorrect sample handling can lead to falsely elevated single factor results. In single factor determination, lupus anticoagulant can affect the apparent factor activity. Non-parallelism upon dilution may occur when a lupus anticoagulant is present in the test sample[9].

Siemens has validated use of these reagents on various analyzers to optimize product performance and meet product specifications. User defined modifications are not supported by Siemens as they may affect performance of the system and assay results. It is the responsibility of the user to validate modifications to these instructions or use of the reagents on analyzers other than those included in Siemens Application Sheets or these Instructions for Use.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Expected Values[10]

Factor VIII:	70 – 150 % of the norm
Factor IX:	70 – 120 % of the norm
Factor XI:	70 – 120 % of the norm
Factor XII:	70 – 150 % of the norm

Reference intervals vary from laboratory to laboratory depending on the population served and the technique, method, equipment and reagent lot used. Therefore, each laboratory must establish its own reference intervals or verify them whenever one or more of the aforementioned variables are changed.

Specific Performance Characteristics

Measuring range

The measuring range for the determination of coagulation factors VIII, IX, XI and XII extends from 1 to approximately 100 % of the norm, and can be increased to approximately 200 % by using a higher dilution of sample (see above).

Method comparison
In a study, the determination of factors VIII, IX, XI and XII on the BFA* System was compared to another coagulation analyzer. The correlations resulted in the following:

	n	r	slope	y-axis intercept (% of the norm)
Factor VIII	49	0.92	0.74	+22.6
Factor IX	50	0.96	0.88	+7.0
Factor XI	50	0.95	0.95	+5.2
Factor XII	50	0.95	1.0	-0.9

Reproducibility
Reproducibility studies for the determination of coagulation factors VIII, IX, XI and XII were run on the BFA* System with samples both in the normal and pathological range. Samples were run in 8-fold determinations over 5 days with stored calibration curves. The reproducibility studies resulted in coefficients of variation of 2 to 14 % for all factors.

Bibliography

1. Mannucci PM, Duga S, Peyvandi F. Recessively inherited coagulation disorders. Blood 2004; 104: 1243-52.
2. Levine PH. Clinical manifestations and therapy of hemophilias A and B. In: Colman RW, et al., eds. Hemostasis and Thrombosis. Philadelphia: Lippincott, 1987: 97–111.
3. Sucker C, Wenzel F, Zotz RB, et al. Monitoring oral anticoagulation may require determination of single coagulation factor activities in patients with antiphospholipid syndrome. J Rheumatol 2006; 33: 1881-2.
4. Schmaier AH, Silverberg M, Kaplan AP, Colman RW. Contact activation and its abnormalities. In: Colman RW, et al., eds. Hemostasis and Thrombosis. Philadelphia: Lippincott, 1987: 18–38.
5. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays; Approved Guideline - Fourth Edition. CLSI document H21-A4 [ISBN 1-56238-521-6]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2003.
6. Boogen C, Niederau C, Reinauer H. Assessment of the influence of r-hirudin on coagulometric factor assays. Ann Hematol 1999; 78 (Suppl): 73 (Abstract).
7. Walenga JM, Drenth A, Mayuga M, et al. Effects of Argatroban alone and combined with oral anticoagulation on coagulation parameters. Blood 2002; 100: Abstract 4006.
8. Feinstein DI. Immune coagulation disorders. In: Colman RW, et al., eds. Hemostasis and Thrombosis. Philadelphia: Lippincott, Williams & Wilkins, 2001: 1003-20.
9. Brandt JT, Triplett DA, Rock WA, et al. Effect of lupus anticoagulants on the activated partial thromboplastin time. Arch Pathol Lab Med 1991; 115: 109–14.
10. Fickenscher K. Analysis of individual coagulation factors. In: Thomas L, ed. Clinical Laboratory Diagnostics. Frankfurt: TH-Books Verlagsgesellschaft, 1998: 607-9.

Dade and Actin are trademarks of Siemens Healthcare Diagnostics.

* Pathromtin and BFA are trademarks of Siemens Healthcare Diagnostics.

Siemens Healthcare Diagnostics Products GmbH
Emil-von-Behring-Str. 76
35041 Marburg/Germany
www.siemens.com/diagnostics

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All rights reserved.
© 2008 Siemens Healthcare Diagnostics Products GmbH.