SIEMENS

Dade® Fibrinogen Calibration Graph Paper

Intended Use

For use with Dade® Fibrinogen Determination Reagents.

Use of Fibrinogen Calibration Graph

For *in vitro* diagnostic use.

The placement of the dilution lines is determined by the assay value (see vial label) for the lot number of Fibrinogen Standard used to establish the calibration curve.

The 1:4 dilution is equal to the assay value multiplied by 2.5. The 1:6.67 dilution is equal to the assay value multiplied by 1.5. The 1:10 dilution is equal to the assay value. The 1:16 dilution is equal to the 1:4 dilution divided by 4. The 1:32 dilution is equal to the 1:16 dilution divided by 2.

The resulting values are then located on the horizontal axis and a vertical line is drawn on the graph. The averaged clotting times for each dilution are then plotted using the five vertical lines calculated for the 1:4, 1:6.67, 1:10, 1:16 and 1:32 dilutions. Connecting the five points usually approximates a straight line.

For example:

If the assay for the Fibrinogen Standard is 240 mg/dL:

1. The 1:4 dilution	=	240 x 2.5	=	600 mg/dL
2. The 1:6.67 dilution	=	240 x 1.5	=	360 mg/dL
3. The 1:10 dilution	=	240 x 1	=	240 mg/dL
4. The 1:16 dilution	=	240 x 0.625	=	150 mg/dL
5. The 1:32 dilution	=	240 x 0.3125	=	75 mg/dL

The vertical lines would be drawn at 600 mg/dL, 360 mg/dL, 240 mg/dL, 150 mg/dL and 75 mg/dL.

Dade is a trademark of Siemens Healthcare Diagnostics.

Siemens Healthcare Diagnostics Products GmbH
Emil-von-Behring-Str. 76
35041 Marburg/Germany
www.siemens.com/diagnostics

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
Factor VIII Chromogenic Assay

FACTOR VIII CHROMOGENIC

See shaded sections = updated information versus edition from May 2008

Intended Use

For photometric determination of Factor VIII (antihemophilic factor) activity in human plasma.

Summary and Principle of the Method

Hemophilia is a sex-linked hemorrhagic disease caused by circulating Factor VIII (F. VIII) or Factor IX (F.IX) deficiency. Hemophilia A is the term used for a deficiency of F. VIII, and is five to seven times more common than Hemophilia B, which is due to a deficiency in Factor IX[1-4].

Both diseases are transmitted as X-linked recessive traits and occur almost exclusively in males who receive the defective gene from the X-chromosome of a carrier mother. In Hemophilia A patients, the degree of F. VIII deficiency dictates the severity of the bleeding disorder. Approximately 10 % to 40 % of normal F. VIII activity is required for normal hemostasis; below this range, a tendency towards bleeding is apparent.

Patients are generally classified by their F. VIII activity into three categories: mild, 25 % to 5 % of normal; moderate, 5 % to 1 % of normal; and severe, less than 1 % of normal. In this chromogenic assay, the Factor VIII in the sample is activated by thrombin[5].

Activated Factor VIII (F. VIIIa) then accelerates the conversion of Factor X (F.X) into Factor Xa (F.Xa) in the presence of activated Factor IX (IXa), phospholipids (PL) and calcium ions. The F.Xa activity is assessed by hydrolosis of a p-nitroanilide substrate specific to F.Xa. The initial rate of release of p-nitroaniline (pNA) measured at 405 nm is proportional to the F.Xa activity, thus to the F. VIII activity in the sample. The chemistry of the Factor VIII Chromogenic Assay is illustrated by the following equations:

$$\text{F. VIII} + \text{Thrombin} \longrightarrow \text{F. VIIIa}$$

$$\text{F. X} \xrightarrow[\text{F. IXa. PL, } Ca^{2+}]{\text{F. VIIIa}} \text{F. Xa}$$

$$CH_3OCO\text{-D-CHG-Gly-Arg-pNA} \xrightarrow{\text{F. Xa}} CH_3OCO\text{-D-CHG-Gly-Arg-OH} + \text{pNA (yellow)}$$

Reagents

Note: Factor VIII Chromogenic Assay can be used manually or on automated coagulation analyzers. Siemens Healthcare Diagnostics provides Reference Guides (Application Sheets) for several coagulation analyzers. The Reference Guides (Application Sheets) contain analyzer/assay specific handling and performance information which may differ from that provided in these instructions for Use. In such a case, the information contained in the Reference Guides (Application Sheets) supersedes the information in these Instructions for Use. In addition, please also consult the instruction manual of the instrument manufacturer.

Materials provided

Factor VIII Chromogenic, REF *B4238 with*

2 x → 2 mL REAGENT FX, Factor X Reagent
2 x → 2 mL REAGENT FIX, Factor IXa Reagent
2 x → 1 mL SUBSTRATE, Substrate Reagent
2 x 10 mL SUBSTRATE BUFFER, Stopping Buffer

Composition

Factor X Reagent: lyophilized preparation containing approximately 2 nmol of bovine F.X, tris-(hydroxymethyl)-aminomethan (Tris) Buffer at pH 8 and stabilizers.

Factor IXa Reagent: lyophilized preparation containing approximately 0.6 nmol bovine F.IXa, approximately 0.6 nmol of bovine thrombin, approximately 0.06 mmol of calcium chloride, approximately 0.12 µmol of phospholipids, Tris buffer at pH 8, and stabilizers.

Substrate Reagent: lyophilized preparation containing approximately 3.4 µmol of CH_3OCO-D-CHG-Gly-Arg-pNA.AcOH, a F.Xa substrate, Nα-(2-Naphthylsulfonylglycyl)-D,L amidino-phenylalanine piperidide (α-NAPAP)[6], a thrombin inhibitor, and stabilizers.

Stopping Buffer: solution containing Tris, ethylendiamintetraacetic acid, sodium chloride and 0.02 % sodium azide.

Warnings and Precautions

1. For *in-vitro* diagnostic use only.
2. Contains sodium azide (< 0.1 %) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations.

Preparations of the Reagents

Note: For analyzer-specific information on reagent preparation please refer to the respective Reference Guides (Application Sheet) provided by Siemens.

REAGENT FX and REAGENT FIX:
Dissolve the contents of a vial with 2.0 mL of distilled or deionized water.

SUBSTRATE:
Dissolve the contents of a vial with 1.0 mL of distilled or deionized water.

SUBSTRATE BUFFER:
For use, mix together 1.0 mL of the SUBSTRATE with 10.0 mL of the SUBSTRATE BUFFER in the brown SUBSTRATE container.

Mix the reagents carefully once more before using.

Storage and Stability

The reagents may be used up to the expiry date indicated on the label if stored at +2 to +8 °C.

Stability after reconstitution (REAGENT FX, and REAGENT FIX, SUBSTRATE/SUBSTRATE BUFFER mixture):

at +37 °C	2 hours
at +15 to +25 °C	8 hours
at +2 to +8 °C	3 days

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Signs of deterioration: No evidence of vacuum in vial upon opening; difficulty in reconstituting reagents; extreme turbidity or presence of particulate matter in SUBSTRATE BUFFER.

Specimens

Mix 9 parts of freshly collected patient blood with 1 part of 0.11 mol/L or 0.13 mol/L sodium citrate solution. Centrifuge the blood specimen at 1500 x g for no less than 15 minutes at room temperature. Immediately move plasma into a plastic tube and keep refrigerated until ready to test.

Patient plasma should be tested within 1 hour of blood collection, and can be stored up to 30 days at -20 °C. Due to the instability of Factor VIII, the plasma may be thawed only once and tested within 1 hour.

Materials required, but not provided

Sodium citrate 0.11 mol/L or 0.13 mol/L for blood collection.
Fresh plasma pool (FPP) or Standard Human Plasma (REF ORKL) for calibration.
Coagulation Factor VIII Deficient Plasma (human), REF OTXW.
Distilled or deionized water.
Spectropholometer with the following capabilities:
— Wavelength: 405 nm
— Pathlength: 1 cm
— Measuring temperatures: +25 °C, +30 °C or +37 °C.
Sodium chloride solution, 0.9 %.
For endpoint determination, 20 % acetic acid or 1 mol/L citric acid. (Mix 200 mL glacial acetic acid with 800 mL distilled or deionized water or dissolve 210 g citric acid monohydrate in distilled or deionized water and dilute to 1 liter.)
Stopwatch.
Pipettes for 0.1 mL, 0.5 mL and 3.0 mL.
Plastic tubes and 1 cm semi-micro cuvettes.
Control Plasma N, REF ORKE
Control Plasma P, REF OUPZ
Imidazole Buffer Solution, REF OQAA *or*
Dade® Owren's Veronal Buffer, REF B4234 *or*
Dade® CA System Buffer, REF B4265
Coagulation analyzer

Procedure

Manual:

Endpoint/Kinetic Assay Procedure

The assay can be performed at +25 °C, +30 °C or +37 °C. Set the wavelength of the photometer at 405 nm and zero against 0.9 % sodium chloride solution.

Sample/Calibration Plasma Dilution

The assay is performed on a 1 : 31 dilution of sample or calibration plasma. This dilution is prepared as follows: Mix 0.1 mL of sample or calibration plasma with 3.0 mL of 0.9 % sodium chloride solution in a plastic tube.

Preparation of the Calibration Curve:

The assay can be performed using a single calibration plasma dilution, but a calibration curve should be calculated for greater accuracy. For this, use a calibration plasma with known level of F. VIII (e.g. 100 %) and dilute it with 0.9 % sodium chloride solution as follows:

Dilution of calibration plasma		1 : 21	1 : 31	1 : 62	1 : 124	zero standard
0.9 % NaCl	mL	2.0	3.0	1.0	1.0	1.0
Calibrating plasma	mL	0.1	0.1			-
Mix and transfer	mL		1.0	1.0		
F. VIII*	%	148	100	50	25	0

* When a commercial assayed reference plasma is used, the percent of factor assumed present in the 1 : 31 dilution will be the assay value published. The percent of factor assumed present for the other dilutions must be adjusted accordingly. The F. VIII concentration of the 1 : 21 dilution is calculated as follows: Concentration of F. VIII in 1 : 31 dilution x 1.48 = Concentration of F. VIII in 1 : 21 dilution, e.g. 100 % x 1.48 = 148 %.

Procedure Outline

All reagents as well as the plastic tubes should be preincubated at the selected temperature. The reproducibility is optimal if the measuring temperature does not deviate by more than 0.2 °C from the selected temperature during the reaction.

Step	Kinetic Assay	Endpoint Assay
Diluted sample/calibration plasma 1 : 31**	0.1 mL	0.1 mL
Factor X Reagent	0.1 mL	0.1 mL
Factor IXa Reagent	0.1 mL	0.1 mL
Mix well and incubate	90 seconds at +37 °C***	
Add substrate Reagent and Stopping Buffer	0.5 mL	0.5 mL
Mix and measure the absorbance increase immediately (ΔA/min).	60 sec 405 nm	60 sec
Add 20 % acetic acid or 1 mol/L citric acid	-	0.1 mL
Mix and measure the absorbance against 0.9 % NaCl	-	405 nm

** See chapter titled "Notes on the Assay Procedure"
*** 90 seconds at +37 °C. 120 seconds at +30 °C or 180 seconds at +25 °C.

For the endpoint assay a plasma bank must be prepared for the samples and each calibration plasma dilution.

20 % acetic acid	0.1 mL
Add Factor X Reagent, Mix.	0.1 mL
Add Factor IXa Reagent. Mix.	0.1 mL
Add diluted Sample, Mix.	0.1 mL
Add Substrate Reagent and Stopping Buffer, Mix.	0.5 mL
Mix and measure the absorbance against 0.9 % NaCl.	405 nm

Internal Quality Control

Normal range: Control Plasma N
Pathological range: Control Plasma P

Two controls should be measured with each calibration and at least every 3 to 4 hours during each testing day (one in the normal range and one in the pathological range). The controls should be processed just like the samples. Each laboratory should determine its own quality control range, either by means of the target values and ranges provided by the manufacturer of the controls or by means of its own ranges established in the laboratory. This range is usually based on ± 2.0 to ± 2.5 standard deviations (SD) from the mean control value. If the measured control value lies outside the confidence level previously established, then the coagulation analyzer, the reagents and the calibration should be examined. Do not release patient results until the cause of deviation has been identified and corrected.

Calculation of the Results

a) Using a single calibration factor (F)

Since the calibration curve is linear between 0 - 100 %, the F. VIII activity in the sample can be calculated by using the known value of a calibration plasma (FPP, Standard Human Plasma).

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Von Willebrand Reagent
VWF REAGENT

OUBDG23C0503 (17

For the determination of ristocetin cofactor activity

See shaded sections - updated information versus edition from November 2008

Intended Use and Application

In-vitro diagnostic reagent for determination of the ristocetin cofactor activity of von Willebrand factor in human plasma using the platelet agglutination method.

Summary and Explanation

Diagnosis of von Willebrand's disease and of diseases with changes in the ristocetin cofactor such as lesions and inflammations of the vascular endothelium. Pre-operative screening of hemorrhagic diatheses: in patients with values below 40 % of the norm bleeding complications are to be expected during dental extractions and surgical interventions.

Principle of the Method

Stabilized platelets are agglutinated in the presence of von Willebrand factor and the antibiotic ristocetin A.

Reagents

Materials provided

VWF REAGENT, REF OUBD23

5 x ⟶ 2 mL REAGENT, von Willebrand Reagent

1 glass plate

1 Table of Analytical Values

Composition

The von Willebrand Reagent contains stabilized platelets, ristocetin and EDTA in lyophilized form.

Preservative: sodium azide (≤ 0.5 g/L).

Warnings and Precautions

1. For *in-vitro* diagnostic use only.
2. Contains sodium azide (< 1 g/L) as a preservative. Sodium azide can react with copper or lead pipes in drain lines to form explosive compounds. Dispose of properly in accordance with local regulations.
3. Each donor or donor unit was tested and found to be negative for human immunodeficiency virus (HIV) 1 and 2, hepatitis B virus (HBV) and hepatitis C virus (HCV) using either tests found to be in conformance with the In Vitro Diagnostic Directive in the EU or FDA approved tests. Because no known test can offer complete assurance of the absence of infectious agents, all human derived products should be handled with appropriate caution.

Preparation of the Reagent

Resuspend the reagent with the amount of distilled or deionized water indicated on the label at +15 to +25 °C by mixing (e.g. with an automatic mixer like the Heidolph mixer „Reax control" set in the range of 75 to 100 %, 2 x 5 seconds). The reagent is then immediately ready for use. Be sure to mix the reagent intensively and to an equal extent for each new reconstitution.

Mix briefly once more before using.

Storage and Stability

Store the Reagent unopened at +2 to +8 °C and use by the expiry date given on the label.

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Isotonic saline

Specimens

To separate the plasma mix 1 part sodium citrate solution 0.11 mol/L with 9 parts venous blood, taking care to avoid the formation of foam. Centrifuge the blood specimen immediately at 1500 x g for no less than 15 minutes at room temperature, remove the supernatant plasma and keep at +15 to +25 °C until required in the assay, or store deep-frozen.

Stability of the specimens:

at +15 to +25 °C	6 hours
at ≤ -20 °C	3 months

Method

Agglutination method

Dilute the plasma specimen with isotonic saline as shown in the following scheme:

Specimen or dilution	–	50 µL	500 µL	500 µL	500 µL	250 µL
Isotonic saline	250 µL	950 µL	500 µL	500 µL	250 µL	250 µL
Dilution step	blank	1 : 20	1 : 40	1 : 80	1 : 120	1 : 160
Titre	–	20	40	80	120	160

Additional dilution ratios between 1 : 80 and 1 : 120 may become necessary order to make control recovery more accurate (1 : 90, 1 : 100, and 1 : 110). Prepa additional dilutions as shown in the following scheme:

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Onto each of the 6 fields of the enclosed glass plate pipette in succession

50 µL plasma dilution or blank and

50 µL Reagent (shake vigorously before use)

and mix. Then agitate slowly for 1 minute (shaker device or by manual agitation) an then leave to stand for 1 minute. If the plates are shaken manually ensure that th test preparations are well mixed!

Assess the agglutination against a dark background with reference to the blank.

The use of an additional light source and a magnifying glass may improve the accurac of this method significantly.

Evaluation

The titre of the specimen is given by the dilution in which there is distinct agglutinatio compared to the blank. The ristocetin cofactor content in % of the norm is obtaine by multiplying the titre of the specimen with the limit of detection which is state on the enclosed Table of Analytical Values. If the content needs to be determine within narrower limits (e.g. in the decision range of 50 %) the intervals between th dilutions have to be reduced to a corresponding extent.

Internal Quality Control

Normal range: Control Plasma N

Pathological range: Control Plasma P

For each reagent filling, two controls (one in the normal range and one in th pathological range) should be measured. The controls should be treated like th samples. The results obtained must lie within the range indicated for the controls i the lot-specific Table of Assigned Values. If the measured control value lies outsid the allowable range, the reagent should be checked. Do not release patient result until the cause of deviation has been identified and corrected.

Limitations of the Procedure

Results of this test should always be interpreted in conjunction with the patient' medical history, clinical presentation and other findings.

Reference Interval

70 to 150 % of the norm[1].

Each laboratory should determine its own reference range.

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SIEMENS

Berichrom* F XIII

See shaded sections - updated information versus edition from November 2005

Intended Use

Photometric determination of Factor XIII activity in plasma samples

Monitoring of substitution therapy with Factor XIII concentrate

Summary and Explanation

The transglutaminase factor XIIIa, produced by the action of thrombin on coagulation factor XIII, produces the covalent cross-link formation in fibrin that stabilizes the fresh fibrin clot. F XIII has been discussed as an essential element in wound healing. The determination of factor XIII allows the assessment of conditions of factor XIII activity deficiency.

Such conditions may be caused by plasma levels that are too low or by functional defects. Decreased plasma levels occur in congenital F XIII deficiency or through the activation of coagulation, as in consumption coagulopathy. Levels of F XIII are typically lower following serious operations. This can lead to reduced resilience of the wound closure with reopening of the wound site.

Principle of the Method

F XIII contained in the sample is converted by the action of thrombin into F XIIIa. Fibrin formed by thrombin also accelerates this reaction. Fibrinogen is not removed prior to testing since this would cause a loss of F XIII. Instead, fibrin produced by the action of thrombin is prevented from forming clots by an aggregation inhibiting peptide and held in solution.

F XIIIa cross-links a specific peptide substrate to glycine ethyl ester, thereby releasing ammonia. The ammonia released is then determined in a parallel enzymatic reaction. The decrease in NADH is measured by monitoring the absorbance at 340 nm[1].

$$\text{fibrinogen} \xrightarrow{\text{thrombin}} \text{fibrin}_{\text{soluble}} \xrightarrow{\text{clot inhibitor}} \text{clot}$$

$$\text{F XIII} \xrightarrow[\text{fibrin}]{\text{thrombin}} \text{F XIIIa}$$

$$\text{glycine ethyl ester + peptide substrate} \xrightarrow{\text{F XIIIa}} \text{conjugate} + NH_3$$

$$NH_3 + \text{NADH} + \alpha\text{-ketoglutarate} \xrightarrow{\text{GLDH}} \text{NAD} + \text{glutamate}$$

Reagents

Materials provided

Berichrom F XIII. Test Kit,* REF *OWSU with*

3 x → 5 mL ACTIVATOR, Activator Reagent

3 x → 5 mL ACTIVATOR DILUENT, NADH Reagent

3 x → 5 mL DETECTION REAGENT, Detection Reagent

Composition

Activator Reagent: bovine thrombin (10 IU/mL), clot inhibitor (Gly-Pro-Arg-Pro-Ala-amide, 2 g/L), calcium chloride (1.2 g/L), hexadimethrine bromide (10 mg/L), bovine albumin, BICINE buffer (100 mmol/L) pH 8.3

Preservative: sodium azide (maximum 0.5 g/L)

NADH Reagent: NADH (0.5 g/L), bovine albumin

Preservative: sodium azide (maximum 0.5 g/L)

Detection Reagent: GLDH (20 IU/mL), synthetic peptide as F XIII substrate (2.4 g/L), ADP, glycine ethylester (1.4 g/L), α-ketoglutarate (2.7 g/L), bovine albumin, HEPES buffer (10 mmol/L) pH 6.5

Preservative: sodium azide (maximum 1 g/L)

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Preparation of the Reagents

ACTIVATOR DILUENT: Dissolve in 5 mL of distilled or de-ionized water.

ACTIVATOR: Reconstitute with 5 mL ACTIVATOR DILUENT.

DETECTION REAGENT: Dissolve in 5 mL of distilled or de-ionized water.

The further preparation depends on the testing method. For automatic processing with a coagulation analyzer, please refer to the information in the respective Reference Guide/Application Sheet. Manual testing is described in section **Procedure**.

Storage and Stability

The test kit may be used up to the expiry date indicated on the label if stored unopened at +2 to + 8 °C.

Stability after reconstitution:

Temperature	ACTIVATOR	DETECTION REAGENT	Reagent Mixture
+37 °C	8 hours	8 hours	4 hours
+15 to +25 °C	2 days	2 days	8 hours
+2 to +8 °C	1 week	1 week	2 days
-20 °C	6 months	6 months	2 months

Signs of expiry: Absence of vacuum when opening the vial; reagent is difficult to reconstitute; results not reproducible.

Stability of the Reagent Mixture

The stability data for the reagent mixture apply if freshly reconstituted individual reagents are used. If the reagent mixture is made from individual reagents that have already been reconstituted and then stored, it must be used immediately.

Information about on-board stability is specified in the Reference Guides (Application Sheets) for several coagulation analyzers.

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Equipment

Specimens

To obtain the plasma, carefully mix 1 part sodium citrate solution (0.11 mol/L) with 9 pa venous blood, avoiding the formation of foam. Centrifuge immediately at 1.500 x g for no le than 15 minutes at room temperature[2].

Stability of the samples:

-20 °C	2 months
+2 to +8 °C	24 hours
+15 to +25 °C	8 hours

Procedure

Manual Testing:

For the test, mix the necessary volume of reagent using equal parts of Activator Reagent a Detection Reagent.

Cuvette:	1 cm thickness
Wavelength:	340 nm
Testing temperature:	+37 °C

Prewarm the reagent mixture and the plastic cuvettes or tubes to +37 °C before testing.

Sample	100 µL
Reagent Mixture	1000 µL

Mix, then begin measuring at 340 nm.

To measure, use either a recording instrument connected to the photometer and evaluat using the speed of the change in absorbance (ΔA/min) after reaching the linear phase (afte about 5 minutes) or read the absorbance after exactly 5 and 10 minutes and calculate th change in absorbance per minute.

Results:

The speed of the change in absorbance after reaching the linear phase (after about 5 minutes) proportional to the F XIII activity.

The results can be evaluated either with a factor (F_L) or with a reference curve. To calculate tl factor, Standard Human Plasma is used as a sample in duplicate determination. To calculate tl reference curve, Standard Human Plasma is used in different dilution steps (with isotonic sali solution) in duplicate determination.

Calculation of the factor:

$$F_L = \frac{\text{Target value (\% of the norm)}_{\text{Standard Human Plasma}}}{\Delta A/\text{min}_{\text{Standard Human Plasma}}}$$

Multiply the ΔA/min of the sample by the factor obtained; this indicates F XIII in % of the norm.

$$\text{F XIII}_{\text{sample}} \text{ (\% of the norm)} = F_L \times \Delta A/\text{min}_{\text{sample}}$$

A new laboratory-internal factor F_L or reference curve must be re-calculated each time there is change in the device used or a change in the lot of Berichrom* F XIII.

Internal Quality Control

Normal range:	Control Plasma N
Pathological range:	Control Plasma P

Two levels of quality control material (normal and pathological range) have to be measured at sta of the test run, with each calibration, upon reagent vial changes and at least every 8 hours duri each day of testing. The controls should be treated like the samples. Each laboratory should d termine its own quality control range, either by means of the target values and ranges provided l the manufacturer of the controls or by means of control values determined in the laboratory. If tl measured control value lies outside of the pre-determined confidence range, then the reagents, th laboratory-internal factor F_L or the reference curve and the coagulation analyzer should be checke Do not release patient results until the cause of deviation has been identified and corrected.

Limitations of the Procedure

F XIII content may be underestimated if ammonia or ammonium levels in the samples are hig (> 0.5 mmol/L). If necessary, determine the ammonia concentration and measure for factor X again with the samples diluted in isotonic saline solution, then multiply the results obtained t the dilution factor.

Very low (< 0.8 g/L) and very high (> 6.0 g/L) fibrinogen concentrations may lead to a false lo result for F XIII activity. If the fibrinogen is over 6.0 g/L, the sample must be prediluted 1:2 or 1 with isotonic saline solution.

Unspecific reactivity of sample components can lead to overestimation of F XIII activity[3].

Siemens has validated use of these reagents on various analyzers to optimize product perfo mance and meet product specifications. User defined modifications are not supported by Siemer as they may affect performance of the system and assay results. It is the responsibility of the us to validate modifications to these instructions or use of the reagents on analyzers other than thos included in Siemens Application Sheets or these instructions for use.

Results of this test should always be interpreted in conjunction with the patient's medical histor clinical presentation and other findings.

Reference Interval

70 - 140 % of the norm[4]

Values for healthy individuals vary from laboratory to laboratory depending on the technique use Therefore, each laboratory should establish its own reference intervals based on the procedu and coagulation analyzers used.

Specific Performance Characteristics

Measuring range

The linear measuring range depends on the coagulation analyzer used. Please refer to the r spective Reference Guide (Application Sheet) or manufacturer information.

Precision and Reproducibility

The coefficient of variation within run was 5 % for Control Plasma N and 4 % for Control Plasma P. Th coefficient of variation from day to day was 6 % for Control Plasma N and 4 % for Control Plasma F

Method comparison

For other F XIII determinations (Laurell immunoelectrophoresis, clot stability test), coefficients correlation of 0.8 and 0.82 were obtained[5].

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1. Hathaway WE. Clinical aspects of antithrombin III deficiency. Semin Hematol. 1991; 28: 19-23.
2. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document H21-A5 (ISBN 1-56238-657-3). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 1908781898 USA, 2008.
3. Wendel HP, Heller W, Gallimore MJ. Aprotinin in therapeutic doses inhibits chromogenic peptide substrate assays for protein C. Thromb Res. 1994; 74: 543-8.
4. Friberger P, Egberg N, Holmer E, et al. Antithrombin assay - the use of human or bovine thrombin and the observation of a „second" heparin cofactor. Thromb Res. 1982; 25: 433-6.

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© 2010 Siemens Healthcare Diagnostics Products GmbH.

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© 2008 Siemens Healthcare Diagnostics Products GmbH.

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2014-06

Siemens Healthcare Diagnostics Products GmbH
Emil-von-Behring-Str. 76
35041 Marburg/Germany
www.siemens.com/diagnostics

SIEMENS

Dade® Actin® Activated Cephaloplastin Reagent

ACTIN

B4218 G1 E0501 (794)

See shaded sections - updated information versus edition from June 2006

Intended Use

Liquid rabbit brain cephalin with plasma activator for use in the determination of the activated partial thromboplastin time (APTT) and other coagulation procedures requiring an activated partial thromboplastin reagent.

Summary and Explanation

The activated partial thromboplastin time, a global screening procedure[1,2] used primarily to evaluate coagulation abnormalities in the intrinsic pathway, will also detect severe functional deficiencies in factors II, V, X, or fibrinogen. The APTT has also been widely advocated[3-5] as a means to monitor the effectiveness of unfractionated heparin therapy where the clotting time is prolonged in proportion to the level of heparin. In patients receiving oral anticoagulants, the circulating levels of factors II, VII, IX, and X are depressed, therefore the APTT can be expected to be prolonged. The presence of non-specific inhibitors, such as the lupus-like anticoagulant[1,7], may prolong the APTT but this effect is variable and generally recognized as being related more to the nature of the APTT reagent employed. In summary, the APTT is a clinically important screening test with wide applicability for the diagnosis of coagulant disorders and therapeutic monitoring of both hemorrhagic and thrombotic disease.

Principle of the Method

Factors of the intrinsic coagulation system are activated by incubating the plasma with the optimal amount of phospholipids and a surface activator. The addition of calcium ions triggers the coagulation process, and the clotting time is then measured.

Reagents

Materials provided

ACTIN, REF B4218-1
10 x 2 mL, ACTIN, Dade® Actin® Reagent
ACTIN, REF B4218-2
10 x 10 mL, ACTIN, Dade® Actin® Reagent

Composition

Dade® Actin® Activated Cephaloplastin Reagent: Cephalin (extract from dehydrated rabbit brain) in 1.0 x 10^{-4} M ellagic acid, buffered, stabilized and preserved. Standard activity values for rabbit brain cephalin are not available.

Warnings and Precautions

For *in-vitro* diagnostic use only.

Reagent Preparation

Dade® Actin® Reagent must be mixed gently by inversion (5 to 8 times) before use.

Storage and Stability

After use, store the closed container at +2 to +8 °C. Stored unopened at +2 to +8 °C, Dade® Actin® Reagent can be used by the labeled expiry date. Once opened, the reagent is stable for 7 days at +2 to +15 °C. Do not freeze. If the reagent is left to stand, a green deposit may form consisting of ellagic acid and lipids. Before use, mix by inverting. Avoid contamination with plasma.

Information about on-board stability is specified in the Reference Guides (Application Sheets) for the different coagulation analyzers.

Signs of expiry: Deviations from the normal laboratory value in the determination of normal plasma or controls.

Materials required but not provided

Calcium Chloride ($CaCl_2$) Solution 0.025 mol/L

Control Plasma N, or Dade® Ci-Trol® Level 1 as control for the normal range

Control Plasma P, or Dade® Ci-Trol® Level 2, or Dade® Ci-Trol® Level 3 as control for the pathological/therapeutical range

Dade® Ci-Trol® Heparin Control, low

Dade® Ci-Trol® Heparin Control, high

For blood collection, use sodium citrate (0.11 or 0.13 mol/L / 3.2 or 3.8 %) or standard commercial blood collection systems

Distilled or deionized water without preservatives

Plastic tubes

Pipettes for precise measurement of 0.1 mL

Coagulation analyzer (see "Equipment")

Equipment

Dade® Actin® Reagent can be used manually or on automated coagulation analyzers. Siemens Healthcare Diagnostics provides Reference Guides (Application Sheets) for several coagulation analyzers. The Reference Guides (Application Sheets) contain analyzer/assay specific handling and performance information which may differ from that provided in these Instructions for Use. In such a case, the information contained in the Reference Guides (Application Sheets) supersedes the information in these Instructions for Use. Please also consult the instruction manual of the instrument manufacturer!

Specimen Collection and Preparation

Mix nine parts of freshly collected patient blood with one part of 0.11 or 0.13 mol/L (3.2 or 3.8 %) sodium citrate. Evacuated tubes containing the desired anticoagulant are commercially available and may be used with caution in blood coagulation studies. For special studies, syringe technique may be preferred.

Centrifuge the blood specimen at 1500 x g for no less than 15 minutes at room temperature as soon as possible after collection. Store in an unopened tube at room temperature.

If immediate testing is to be done, the plasma may remain on the packed cells. Otherwise plasma should be separated from the cells. To separate the plasma, use a plastic transfer pipette, remove plasma to a plastic tube. Do not store on ice. Non-heparinized plasma should be tested within four (4) hours of blood collection.

Plasma containing unfractionated heparin should be centrifuged within one hour of blood collection, stored at room temperature and tested within four (4) hours.

Platelet-poor plasma may be frozen at ≤ -20 °C for up to two (2) weeks in a non-frost-free freezer. Frozen plasma should be rapidly thawed at +37 °C, gently mixed and tested immediately. Samples should not stand at +37 °C for more than five minutes. Please refer to CLSI document H21-A5[8] for detailed information on sample preparation and storage.

Procedure

Manual Testing:

Prewarm Calcium Chloride Solution at +37 °C		
Prewarm 0.1 mL Dade® Actin® Reagent for one minute at +37 °C. (Mix before use)		
Pipet into coagulation tubes as follows:		
	Test Sample	Control Sample
Dade® Actin® Reagent (prewarmed)	0.1 mL	0.1 mL
Plasma	0.1 mL	-
Control Plasma	-	0.1 mL
	Mix well. Incubate at +37 °C for 3 minutes.	
prewarmed Calcium Chloride Solution	0.1 mL	0.1 mL
	Simultaneously with addition of $CaCl_2$ start stopwatch, mix well. After 20 seconds start to observe for clot formation.	

Note: Incubation times exceeding 5 minutes may cause loss of factors V and VIII and are not recommended.

Each laboratory should determine the optimal heating-activation time for its particular assay system.

Monitoring of Unfractionated Heparin Therapy with APTT

When using the APTT for this purpose, the factors influencing the test should be kept in mind. General considerations are listed below.

A) Time of collection is important since the *in-vivo* half-life of unfractionated heparin is approximately 1.5 hours[5]. When it is administered, it has an immediate anticoagulant effect but the degree of this effect decreases rapidly with time. This is especially apparent with intermittent single intravenous injections.

B) The anticoagulant used for sample collection can alter test results.

C) Platelet factor 4, a heparin neutralizing factor in platelet alpha-granules, can be released by platelet aggregation or damage. To prevent this occurrence *in-vitro*, the specimen should be collected with a minimum of trauma. Cold temperatures are known to induce platelet aggregation and release platelet factor 4; therefore, centrifugation at room temperature is recommended for heparin studies.

D) Using APTT to monitor unfractionated heparin therapy is time-dependent. Delay in testing samples will result in prolonged APTT determinations. Therefore, it is imperative that the testing on all samples be performed as soon as possible.

E) Increased contact activation times may result in prolonged APTT in plasma containing heparin. It is imperative that the optimal heating-activation time of the Dade® Actin® Reagent-plasma mixture be rigidly standardized[9].

F) Different test systems (i.e. manual, photo-optical, etc.) will exhibit variable heparin sensitivity. Interchanging of test systems should be avoided.

G) Baseline data on the APTT of each patient before the start of therapy should be established where feasible to determine the respective patient APTT as it relates to the normal range established for the test in that laboratory.

H) Studies have shown variability in original estimates of the quality of unfractionated heparin from different sources and different manufacturers. *In-vivo* reactivity varies with the type of heparin administered, the metabolism of the individual and other co-administrated medications[5,6].

I) Because the APTT can vary with technique, method, equipment, reagent lot and heparin used, each laboratory must establish its own therapeutic ranges, or verify them whenever one or more of the aforementioned variables is changed. This can be done by simultaneously determining the APTT and the heparin concentration for samples from patients receiving heparin therapy. A dose-response curve can be calculated from the data using regression analysis, and the APTT range corresponding to a heparin concentration of 0.3 - 0.7 U/mL (by a factor Xa inhibition assay) can be derived[4-6].

Internal Quality Control

Normal range:	Dade® Ci-Trol® Level 1, or Control Plasma N
Pathological range:	Dade® Ci-Trol® Level 2, or Dade® Ci-Trol® Level 3, or Control Plasma P
Heparin monitoring:	Dade® Ci-Trol® Heparin Control, low and high

Two controls (one in the normal range and one in the pathological/therapeutical range) must be measured at the start of the test run, after each change of reagent vial, and at least once during an 8-hour shift. The control material should be prepared and processed in the same manner as the patient samples. Each laboratory should establish its own confidence intervals for the controls. This interval is generally ± 2 to ± 2.5 standard deviations (SD) from the mean control value. If the control values are outside of the confidence interval, the controls, reagents and instrument must be checked. Before reporting the patient values, it is recommended that all steps should be documented that were taken to identify and rectify the problem. New control ranges should be defined for each new lot of reagents or controls.

Results

Results of the activated partial thromboplastin time testing should be reported as the APTT in seconds. These results should be related to the normal range for APTT testing in each laboratory. It is suggested that the patient results be reported to the clinician in conjunction with the normal range. Control values for the reagent test system should never be used in place of a normal range. Furthermore, the reporting of APTT results in terms of an upper normal only may result in incorrect interpretation. Shortened APTT results may also indicate some abnormal condition in the patient's coagulation system.

Limitations of the Procedure

APTT testing encompasses the entire clotting process from contact activation to fibrin formation and is therefore more susceptible to variations than specific individual tests. The control and use of APTT is therefore subject to inherent limitations. Control of plasma sample conditions is strictly emphasized. Studies have shown that sample decomposition may occur more rapidly in stored samples that are not refrigerated. Extremely small plasma volumes (prior to testing) are to be avoided since pH changes in the plasma from physiological conditions may be encountered. Such changes may lead to the decomposition of plasma components of the blood coagulation system.

It should be noted that APTT time testing may be affected by a number of commonly administered drugs. Decrease in time of APTT determination in conjugated estrogen therapy in males and oral contraceptive administration in females has been reported[10,11]. Increase in the APTT has been seen in diphenylhydantoin, heparin, warfarin, naloxone and radiographic agent administration[12,13]. Therapeutic doses of hirudin or other direct thrombin inhibitors may prolong clotting times[14]. In addition, the choice of anticoagulant (i.e., citrate vs. oxalate) and the condition of the specimen (e.g. hemolyzed, lipemic, parenteral feeding, etc.) may affect results[4,15,16].

The latter is particularly true of optical instrumentation measurements of the APTT. Blood clotting factor deficiencies which should produce prolonged clotting times may be compensated for or made to appear normal by elevated levels of one or more different clotting factors. Similarly, the presence of active intermediates which would tend to reduce the clotting time may also mask conditions that would normally lead to prolongation of the APTT. Mild or minor deficiencies in several factors may have an additive effect on increasing the APTT. Unexpected abnormal APTT results should always be followed by additional coagulation studies to determine the source of abnormal results.

Action of heparin as an anticoagulant is related to its ability in conjunction with a plasma cofactor to interfere with several aspects of the coagulation mechanism, thus retarding the rate of fibrin formation (see section "Monitoring of Unfractionated Heparin Therapy with APTT").

Siemens has validated use of these reagents on various analyzers to optimize product performance and meet product specifications. User defined modifications are not supported by Siemens as they may affect performance of the system and assay results. It is the responsibility of the user to validate modifications to these instructions or use of the reagents on analyzers other than those included in Siemens Reference Guides (Application Sheets) or these Instructions for Use.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

Expected Values

Reference intervals vary from laboratory to laboratory depending on the population served and the technique, method, equipment and reagent lot used. Therefore, each laboratory must establish its own reference intervals or verify them whenever one or more of the aforementioned variables are changed. In a study with apparently healthy individuals using a specific lot of Dade® Actin® Reagent the following data were obtained:

	Median (seconds)	90 % Reference Interval (seconds)	
		5th Percentile	95th Percentile
111 individuals Sysmex® CA-1500	26.6	22.7	31.8
111 individuals BCS® System	26.1	22.8	31.0

Reference ranges for other populations such as pediatric groups should also be established where warranted.

Note: CLSI Document C28-A2 (cited in H47-A)[17,18] states that a parametric approach (mean ± 2 SD) can be applied. The assumption of this approach (Gaussian normal distribution) must however be checked.

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
Level / Konzentration / Niveau / Livello / Nivel / Nível

2006-09_EFIGSP

© 2008 Siemens Healthcare Diagnostics Products GmbH
All rights reserved.

SIEMENS

ORHOG37E0503 (1549)

Calcium Chloride Solution
$CaCl_2$ SOLUTION

Intended Use
Calcium Chloride Solution is used as Supplementary Reagent for various coagulation assays.

Reagent
Materials provided
Calcium Chloride Solution, REF *ORHO37*
10 x 15 mL $CaCl_2$ SOLUTION, Calc[illegible] Chloride Solution
Composition
$CaCl_2$ solution 0.025 mol[illegible]
Warnings and Precautions
For *in-vitro* diagnostic use only
Storage and Stability
Storage at 2 to 25 °C:
The expiry date is indicated on the label.
Stability once opened:
8 weeks at 2 to 25 °C
Information about on-board stability is specified in the Reference Guides (Application Sheets) for the different coagulation analyzers.
Materials required, but not provided
Coagulation Analyzer
Reagents and associated materials according to the respective Instructions for Use of the coagulation assays to be used with Calcium Chloride Solution REF ORHO37.

Procedure
- Calcium Chloride Solution is added automatically by the respective coagulation analyzer when performing the respective assay.
- Please refer to the respective Instructions for Use and Application Sheets for the assays to be used with Calcium Chloride Solution.

Siemens Healthcare Diagnostics Products GmbH
Emil-von-Behring-Str. 76
35041 Marburg/Germany
www.siemens.com/diagnostics

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© 2018 Siemens Healthcare Diagnostics Products GmbH.