Boston Scientific

EP·XT™ Dynamic Tip™ Dynamic XT™

Unidirectional Steerable Diagnostic Catheter

Directions for Use	**2**
Instrucciones de uso	**5**
Mode d'emploi	**8**
Gebrauchsanweisung	**12**
Istruzioni per l'uso	**16**
Gebruiksaanwijzing	**19**
Instruções de Utilização	**22**

90962725-01

2015-01

Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240 RevAH), DFU, MB, Steerable Catheter, Global, 90962725-01A

Black (K) ΔE ≤5.0

Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240 RevAH), DFU, MB, Steerable Catheter, Global, 90962725-01A

TABLE OF CONTENTS

WARNING **3**

DEVICE DESCRIPTION **3**

INTENDED USE/INDICATIONS FOR USE **3**

CONTRAINDICATIONS **3**

WARNINGS **3**

PRECAUTIONS **3**

HOW SUPPLIED **4**

Handling and Storage **4**

INSTRUCTIONS FOR USE **4**

WARRANTY **4**

Black (K) ΔE ≤5.0

EP•XT™ Dynamic Tip™ Dynamic XT™

Unidirectional Steerable Diagnostic Catheter

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

EP•XT, Dynamic Tip and Dynamic XT Unidirectional Steerable Diagnostic Catheters are radiopaque, flexible, insulated catheters with a polymer shaft. EP•XT Steerable Catheter has a white rotating handle mechanism which provides curve actuation of the distal tip. Dynamic Tip and Dynamic XT Steerable Catheters have a plunger mechanism, which, when moved forward or back, results in curvature of the distal tip.

INTENDED USE/INDICATIONS FOR USE

EP•XT, Dynamic Tip and Dynamic XT Unidirectional Steerable Diagnostic Catheters are intended for temporary intracardiac sensing, recording, stimulation and temporary pacing during the evaluation of cardiac arrhythmias.

CONTRAINDICATIONS

The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle patch. The retrograde transaortic approach is contraindicated in patients with aortic valve replacement.

WARNINGS

- This device should be used only by physicians thoroughly trained in the techniques of intracardiac electrophysiology studies and temporary pacing.
- The risks of using electrophysiology catheters include those risks related to heart catheterization such as thromboembolism, perforation, tamponade, and infection. The induction of an unintended arrhythmia is a known complication.
- This device is for single use only. Resterilization and reuse of this device have not been proven to be effective.
- Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.
- Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- Catheter advancement should be done under fluoroscopic guidance.

PRECAUTIONS

- The safety and effectiveness of this device as an ablation catheter have not been established. Therefore, such use is considered investigational.
- Use only sterile saline or water to wipe this catheter.
- Avoid submerging the catheter handle in any solution.
- For catheters equipped with a cable connector, use with the appropriate Boston Scientific cable.

Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240 RevAH), DFU, MB, Steerable Catheter, Global, 90962725-01A

Black (K) ΔE ≤5.0

- Excessive bending, torquing or kinking of the electrode catheter may cause damage to the catheter, including damage to the internal wires.
- After use, this product may be a potential biohazard. Handle and dispose in accordance with accepted medical practice and applicable local, state, and federal laws.

HOW SUPPLIED

EP•XT™, Dynamic Tip™ and Dynamic XT™ Unidirectional Steerable Diagnostic Catheters are supplied sterile using an ethylene oxide (EO) process and intended for single use only.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

INSTRUCTIONS FOR USE

1. Inspect the packaging carefully for overall condition. Do not use if the package is damaged.
2. Inspect the electrodes carefully for integrity, continuity of leads, and overall condition. If damaged, do not use.
3. Insert the catheter by using a standard percutaneous catheter introducer.
4. The electrode catheter should be passed from a peripheral vessel to the desired intracardiac position under fluoroscopic guidance.
5. **For EP•XT Unidirectional Steerable Diagnostic Catheters only:** The catheter tip is in a straight position when the white indicator dot is visible through the handle window. Rotate the white handle between the thumb and forefinger in a clockwise motion until the desired curve is attained. To straighten the tip, rotate in a counterclockwise motion until the white indicator dot reappears.
6. **For Dynamic Tip and Dynamic XT Unidirectional Steerable Diagnostic Catheters:** The catheter tip can be deflected by advancing or retracting the plunger mechanism on the handle. When the plunger mechanism is in the fully retracted position, the tip is approximately straight.
7. All catheter adjustments should be done under fluoroscopic guidance.
8. Upon completion of the procedure, return the catheter to the neutral position prior to removal from the patient.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240 RevAH), DFU, MB, Steerable Catheter, Global, 90962725-01A

Black (K) ΔE ≤5.0

CONTENIDO

ADVERTENCIA **6**

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CONTRAINDICACIONES **6**

ADVERTENCIAS **6**

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PRESENTACIÓN **7**

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GARANTÍA **7**

Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240 RevAH), DFU, MB, Steerable Catheter, Global, 90962725-01A

Black (K) ΔE ≤5.0

Boston Scientific (Master Brand DFUTemplate 5in x 5.25in Global, 90348240 RevAH), DFU, MB, Steerable Catheter, Global, 90962725-01A

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2015 Boston Scientific Corporation or its affiliates.