Boston Scientific

Blazer™ II
Blazer™ II HTD

Temperature Ablation Catheters

Directions for Use	**2**
Instrucciones de uso	**5**
Mode d'emploi	**9**
Gebrauchsanweisung	**12**
Istruzioni per l'uso	**16**
Gebruiksaanwijzing	**20**
Instruções de Utilização	**24**

Boston Scientific (Master Brand DFUTemplate 8.5in x 11in Global, 90106042AQ), DFU, MB, Blazer II / II HTD, Global, 50561401-01A

Black (K) ΔE ≤5.0

Blazer™ II
Blℓazer™ II HTD
Temperature Ablation Catheters

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these instructions. Failure to do so may result in patient complications. Boston Scientific relies on the physician to determine, assess and communicate to each patient all foreseeable risks of the procedure.

DEVICE DESCRIPTION

The Boston Scientific Corporation (BSC) Cardiac Ablation Systems consist of a Maestro Radiofrequency (RF) Cardiac Ablation Controller and Maestro Pod. The Maestro RF Cardiac Ablation Controller is used in conjunction with a temperature ablation catheter (henceforth, referred to as ablation catheter).

Detailed "Instructions for Use" of the BSC Cardiac Ablation Systems are provided in the appropriate Operator's Manual. Information further detailed in this Directions for Use (DFU) specifically addresses the ablation catheters. The ablation catheter is connected to a Maestro RF Cardiac Ablation Controller via a Maestro Pod. The Maestro Pod provides for additional connection to standard hospital Electrophysiology (EP) recorders. The catheter-to-pod connection may be made through the use of the Quick Connect Instrument Cable, which incorporates BSC's proprietary "Quick Connect" connectors, illustrated in Figure 1. The "Quick Connect" configuration allows for rapid connection and disconnection of the catheter.

For all ablation catheters, monopolar RF power is delivered between the catheter's distal electrode and a commercially available external Dispersive Pad. Quadripolar catheter configurations are available (three ring electrodes). Ring electrodes are utilized for either electrogram recording or for pacing.

Ablation catheters are available with temperature monitoring capabilities. Other catheter options include a variety of curve profiles, catheter lengths, distal section lengths, ring electrode spacing, torque and stiffness characteristics, and catheter shaft diameters. Similarly, tip electrodes are available in a variety of sizes (diameters, lengths) and shapes. Please refer to the Product Catalogue for a complete description of available product configurations.

User Information

An ablation catheter is to be used by physicians thoroughly trained in invasive cardiology and in the techniques of RF powered catheter mapping and ablation, in the specific approach to be used, and in a fully-equipped electrophysiology lab.

Contents

One (1) sterile ablation catheter.

INTENDED USE / INDICATIONS FOR USE

The Boston Scientific Cardiac Ablation System is indicated for interruption of accessory Atrioventricular (AV) conduction pathways associated with tachycardia, for treatment of AV nodal reentrant tachycardia and for creation of complete AV block in patients with a rapid ventricular response to an atrial arrhythmia - typically chronic, drug refractory atrial fibrillation.

CONTRAINDICATIONS

The use of this device is contraindicated in patients who have vena cava embolic protection filter devices and/or known femoral thrombus who require catheter insertion from the femoral approach.

The use of the device is contraindicated in patients with active systemic infection.

The transseptal approach is contraindicated in patients with left atrial thrombus or myxoma, or interatrial baffle or patch.

The retrograde transaortic approach is contraindicated in patients with aortic valve replacement.

WARNINGS

- Care must be taken to ensure that any equipment used in connection with the BSC catheters be type CF, be defibrillation proof, meet IEC 60601-1 electrical safety requirements, and comply with all local regulatory requirements for the specified intended use.
- Maximum Catheter Rated Voltage: 178 Vrms (251 Vpk).
- Before using, inspect for physical damage including electrical insulation on the cables and the catheter shaft. Replace damaged equipment.
- Cardiac ablation procedures should be performed only by physicians thoroughly trained in the techniques of radiofrequency powered catheter ablation and in the specific approach to be used, in a fully equipped electrophysiology laboratory.
- Catheter ablation procedures present the potential for significant x-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the x-ray beam intensity and duration of the fluoroscopic imaging. Catheter ablation should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps taken to minimize this exposure. Careful consideration must therefore be given for the use of the device in pregnant women and prepubescent children.
- Patients undergoing AV nodal modification or ablation of septal accessory pathways are at risk for inadvertent AV block. It is advisable to use lower initial power in such patients and to monitor anterior conduction closely during RF power delivery.
- Pacemakers and implantable cardioverter / defibrillators can be adversely affected by RF signals. It is important to:
 a. Have temporary external sources of pacing available during ablation.
 b. Temporarily reprogram the pacing system to minimum output or 000 mode to minimize risk of inappropriate pacing.
 c. Exercise extreme caution during ablation when in close proximity to atrial or ventricular permanent pacing leads.
 d. Perform complete pacing system analysis on all patients after ablation.
- Implanted cardioverter / defibrillators should be deactivated during delivery of RF power.
- During a transaortic approach, adequate fluoroscopic visualization is necessary to avoid placement of the ablation catheter within the coronary vasculature. Catheter placement and RF power application within the coronary artery has been associated with myocardial infarction and death.
- Patients undergoing left sided ablation procedures should be closely monitored during and post procedure for clinical manifestations of infarction, pulmonary vein injury, nerve damage, embolism and/or atrial esophageal fistula.
- Catheter entrapment within the heart or blood vessels is a possible complication of cardiac ablation procedures. The potential for catheter entrapment may be increased when the catheter is positioned in the vicinity of the chordae tendineae. The occurrence of this complication may necessitate surgical intervention and/or repair of injured tissues.
- During RF ablation, care must be taken not to deliver RF energy on or near the coronary artery even on the right side of the heart, as the resulting myocardial injury can be fatal.
- Ablation in contact with any other electrodes alters the function of the catheter and can lead to thrombus, coagulum, or char formation that may result in embolism.
- It is important to start each RF application at low power and carefully follow the power titration as specified in the instructions for use. Too rapid an increase in power during ablation may lead to arrhythmias, damage to adjacent structures, embolism and/or perforation caused by steam pop.
- Do not deliver RF energy with the catheter outside the target site. RF generators can deliver significant electrical energy and may cause patient or operator injury.
- Stimulation of cardiac tissues caused by pacing stimulus and/or RF energy may lead to inadvertent induction of arrhythmias. These arrhythmias may require defibrillation that could also result in skin burns.
- No modification of this equipment is allowed.

PRECAUTIONS

- Excessive bending or kinking of the catheter shaft may damage internal wires. Manual pre-bending of the distal curve can damage the steering mechanism and may cause patient injury.

Boston Scientific (Master Brand DFUTemplate 8.5in x 11in Global, 90106042AQ), DFU, MB, Blazer II / II HTD, Global, 50561401-01A

Black (K) ΔE ≤5.0

- The long-term risks of lesions created by RF ablation have not been established. In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown. Furthermore, the risk/benefit in asymptomatic patients has not been studied.
- The Controller is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the catheter and Dispersive Pad, particularly when operating the device.
- During power delivery, the patient should not be allowed to come in contact with grounded metal surfaces.
- Do not use the Blazer™ II/II HTD Ablation System in the proximity of Magnetic Resonance Imaging (MRI) equipment because the MRI equipment may adversely impact the function of a RF generator and the ablation system may adversely impact the image quality.
- Do not deliver RF energy with the catheter outside the target site. RF generators can deliver significant electrical energy and may cause patient or operator injury.
- The Blazer II Catheter and the Blazer II HTD Catheter are highly torqueable. Avoid overtorquing. Over-rotating the handle and catheter shaft may cause damage to the distal tip or catheter assembly. Do not rotate the handle and catheter shaft more than 1½ full rotations (540˚). If the desired catheter tip position is not achieved, adjust the catheter's curve to disengage the catheter tip from the heart wall before resuming rotation of the handle and catheter shaft.
- Peri-procedural anticoagulation therapy is recommended for patients undergoing left-sided and transeptal cardiac procedures and should be considered for selected patients undergoing right-sided procedures.
- Careful catheter manipulation must be performed in order to avoid cardiac damage, perforation or tamponade. Catheter advancement should be done under fluoroscopic guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered.
- The catheter impedance LED display of the Controller should be continuously monitored during RF power delivery. If a sudden rise in impedance is noted, power delivery should be discontinued. The catheter should be removed and the distal tip of the catheter cleaned to eliminate any coagulum.
- Adequate filtering must be used to allow continuous monitoring of the surface Electrocardiogram (ECG) during radiofrequency power applications.
- The long-term risks of lesions created by RF ablation have not been established. In particular, any long-term effects of lesions in proximity to the specialized conduction system or coronary vasculature are unknown. Furthermore, the risk / benefit in asymptomatic patients has not been studied.
- The Controller is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the catheter and Dispersive Pad, particularly when operating the device. During power delivery, the patient should not be allowed to come in contact with grounded metal surfaces.
- Apparent low power output or failure of the equipment to function correctly at normal settings may indicate faulty application of the Dispersive Pad or failure of an electrical lead. Do not increase power before checking for obvious defects or misapplication.
- The risk of igniting flammable gases or other materials is inherent in the application of RF power. Precautions must be taken to restrict flammable materials from the area where the ablation procedures are performed.
- The long-term risks of protracted fluoroscopy have not been established. Careful consideration must therefore be given for the use of the device in prepubescent children.
- Electromagnetic Interference (EMI) produced by the Controller during the delivery of RF power may adversely affect the performance of other equipment.
- When crossing the aortic valve with the ablation catheter, it is recommended that the catheter tip be deflected to resemble a "pigtail" curve to avoid damage to the valve leaflets.

ADVERSE EVENTS

The following potential risks or discomforts may be associated with cardiac ablation procedures. The frequency and severity of these adverse events can vary, and may necessitate additional medical intervention, including surgery.

- Allergic reaction
- Arrhythmias (new or exacerbation of existing arrhythmias) and complications of sedative agents / anesthesia
- Cardiac or respiratory arrest
- Cardiac valve damage
- Catheter entrapment / entanglement
- Chest pain
- Damage to vessel intima or cardiac ultrastructures
- Death
- Embolus, air embolus
- Hematoma / ecchymosis
- Hemorrhage
- Hypotension
- Infection
- Myocardial infarction
- Perforation
- Pericardial / pleural effusion
- Pericarditis / pleuritis
- Phrenic or intercostal nerve damage
- Pneumothorax
- Pseudoaneurysm
- Pulmonary edema
- Sinus or AV node injury
- Skin burns
- Stroke
- Stroke or cerebral vascular accident
- Tamponade
- Thrombosis
- Vasovagal reaction
- X-ray exposure

HOW SUPPLIED

Ablation catheters are supplied sterile using an ethylene oxide process.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Operating Environment

Ambient Temperature: 10 °C to 40 °C

Relative Humidity: 30% to 75%

Atmospheric Pressure: 70 kPa to 106 kPa

Transport Environment

Temperature: -29 °C to 60 °C

Relative Humidity: 30% to 85%

Atmospheric Pressure: Uncontrolled

Storage Environment

Ambient Temperature: 20 °C to 30 °C

Relative Humidity: Uncontrolled

Atmospheric Pressure: Uncontrolled

MATERIALS REQUIRED

Equipment required to perform cardiac ablation:

- 6F (2.0 mm), 7F (2.3 mm) or 8F (2.7 mm) hemostatic introducer sheath, sized according to the tip electrode diameter for the cardiac ablation catheter in use.
- Dispersive Pad meeting or exceeding IEC 60601-2-2 requirements for electrosurgical electrodes.

INSPECTION PRIOR TO USE

Carefully inspect the package prior to use for any breach of the sterile barrier or damage to the contents. If the sterile barrier integrity is compromised or the contents damaged, contact your Boston Scientific representative. Check the "Use By" date on the device package. Do not use the device if passed the "Use by" date.

SETTING UP THE SYSTEM

Please refer to the Operator's Manual for the appropriate BSC Cardiac Ablation Controller. The Operator's Manual describes the steps to follow to connect the system, set ablation parameters, and deliver radiofrequency power.

ATTACHING THE DISPERSIVE PAD

Plug the Dispersive Pad connector into the two-prong socket marked "INDIFFERENT ELECTRODE" on the front panel of the Maestro Pod. Make sure that the Dispersive Pad connector is firmly pressed into the socket.

DIRECTIONS FOR USE

The ablation catheter is usually inserted into a vein or artery and is then positioned into the appropriate chamber of the heart under fluoroscopic guidance. A transseptal approach may be used. Precise placement of the catheter prior to ablation is accompanied by endocardial mapping using the tip and/or ring electrodes. Once appropriate positioning has been achieved, radiofrequency power is delivered via the Controller resulting in the ablation of the targeted cardiac tissue.

Prior to insertion of the steerable ablation catheter, prepare the entry site according to standard aseptic technique practices.

1. Insert the catheter percutaneously into the appropriate artery or vein by the Seldinger technique using the hemostatic introducer sheath.
2. Once inside the vessel, the catheter tip can be deflected as necessary to facilitate advancement into the selected heart chamber.
3. The catheter curve can be straightened completely and deflected in the opposite direction against cardiac tissue, facilitating stability during ablation.
4. Connect the Maestro Pod to the Controller "Isolated Patient Connection" located on the front panel using the attached patient cable. Be sure to carefully follow the instructions in the appropriate Operator's Manual to connect the Maestro Pod.

Boston Scientific (Master Brand DFUTemplate 8.5in x 11in Global, 90106042AQ), DFU, MB, Blazer II / II HTD, Global, 50561401-01A

Black (K) ΔE ≤5.0

5. Connect the ablation catheter through the use of various adaptor cables (see Figure 1).
6. Once the ablation site has been accessed and the tip of the catheter is in contact against the endocardial surface, intracardiac electrogram signals may be obtained. Unipolar electrograms can be obtained between the distal tip electrode and any commercially available low impedance ECG electrode. Bipolar electrogram recordings can be recorded between the distal tip electrode and any ring electrode, or between any two ring electrodes.
7. Once the arrhythmogenic site has been located, the same catheter can be used therapeutically in the "Ablate" mode to deliver discrete bursts of radiofrequency power. Radiofrequency power is delivered to the tissue via the distal tip (ablation) electrode which results in thermal necrosis (ablation) of the arrhythmogenic tissue. To deliver radiofrequency power from the Cardiac Ablation Controller to the catheter, follow the instructions in the appropriate Cardiac Ablation Controller Operator's Manual.

CATHETER REMOVAL

1. Prior to removing the catheter, ensure that the distal end of the catheter is straightened completely.
2. Withdraw the catheter from the vessel.
3. Remove the introducer sheath, then follow standard practice for management of the insertion site.

REFERENCES

Available Upon Request.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose**. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments**.


Figure 1. Cable Configuration for Maestro Pod


Figure 2. Temperature Ablation Catheter


Figure 3. Distal Tip

Boston Scientific (Master Brand DFUTemplate 8.5in x 11in Global, 90106042AQ), DFU, MB, Blazer II/ II HTD, Global, 50561401-01A

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© 2017 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd