Boston Scientific

EndoVive™

Standard Balloon Replacement Kit
Straight and Right-Angle Bolsters

Directions for Use	**02**
Instrucciones de uso	**07**
Mode d'emploi	**12**
Gebrauchsanweisung	**18**
Istruzioni per l'uso	**24**
Gebruiksaanwijzing	**29**
Instruções de Utilização	**34**

ENGLISH
ESPAÑOL
FRANÇAIS
DEUTSCH
ITALIANO
NEDERLANDS
PORTUGUÊS

TABLE OF CONTENTS

WARNING **3**

DEVICE DESCRIPTION **3**

INTENDED USE **3**

INDICATIONS FOR USE **3**

CONTRAINDICATIONS **3**

WARNINGS **3**

PRECAUTIONS **4**

ADVERSE EVENTS **4**

HOW SUPPLIED **4**

HANDLING AND STORAGE **4**

OPERATIONAL INSTRUCTIONS **4**

Device Placement 4

Feeding 5

Maintenance 6

Removal 6

WARRANTY **6**

EndoVive™

Standard Balloon Replacement Kit
Straight and Right-Angle Bolsters

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

ENGLISH

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The EndoVive Standard Balloon Replacement kit includes components used to facilitate placement and use of the EndoVive Standard Balloon Replacement Tube (from herein referred to as replacement tube).

- The replacement tube is a silicone balloon-type device for feeding and/or administration of medication and decompression in conjunction with an established gastrointestinal stoma tract.
- An inflation syringe is included to facilitate balloon inflation.
- Supplied Gauze and lubricating jelly (selected markets) are intended to facilitate ease of device insertion into the stoma tract.

INTENDED USE

The replacement tube is intended to provide gastric access for Enteral feeding, medication administration, and gastrointestinal decompression.

INDICATIONS FOR USE

The replacement tube is indicated in adult and pediatric populations for use in percutaneous placement of an Enteral feeding tube for feeding and/or administration of medication in conjunction with an established GI stoma tract. Typical uses include the replacement of existing gastrostomy feeding tubes. The replacement tube may also be used for decompression.

CONTRAINDICATIONS

Placement is also contraindicated for patients with evidence of granulation tissue, infection, and/or irritation at the stoma site.

WARNINGS

- Excessive tension, tissue necrosis and tube dislodgement may cause separation of the gastric and abdominal wall.
- This product must never be used in the vasculature.

PRECAUTIONS

The replacement tube should only be used by or under the supervision of physicians trained in percutaneous gastrostomy tube placement. A thorough understanding of the technical principles, clinical applications and risks associated with percutaneous gastrostomy tube placement is necessary before using this device.

Any use of this device, other than those indicated in these instructions, is not recommended.

ADVERSE EVENTS

The following may be associated with the use of Gastrostomy Tubes: aspiration, reflux, sepsis, ascites, bleeding, granulation tissue, pressure necrosis, ulcers, severe gastroesophageal reflux or diffuse inflammatory, infectious or neoplastic disease involving the walls of the abdomen or anterior stomach, gastrointestinal obstruction, tube occlusion, tube clogging, peritonitis, perforation, malposition, leakage, kinking, tube migration, inadvertent tube removal, proximal small bowel obstruction (fistulae), irritation and infection such as redness, edema or purulent drainage.

HOW SUPPLIED

Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible. Do not reuse.

HANDLING AND STORAGE

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Device Placement

1) Remove the existing Gastrostomy tube from the stoma per the original manufacturer's instructions (noting current device diameter).
2) Choose the appropriately sized replacement tube.
3) Inspect the kit and kit contents for shipping damage, any suspected breach in sterility, foreign material, missing components, or anything that you suspect could potentially affect the replacement tube's performance.
 - Avoid sharp objects coming into contact with the replacement tube.
4) Inflate the balloon to the rated volume as indicated on the device, using sterile or distilled water and the supplied sterile syringe. See Figure 1, number 1. Do not over-inflate the balloon.


Figure 1

5) Inspect the balloon for leaks, balloon shape uniformity, occlusions, or anything that you suspect could potentially affect the replacement tube's performance.
6) Deflate the balloon using the syringe and inspect the replacement tube for occlusions, leaks, or anything that you suspect could potentially affect the replacement tube's performance. Do not attempt to use the replacement tube if the balloon cannot be deflated.
7) Apply a water-soluble lubricant to the tip of the replacement tube, deflated balloon and the shaft using sterile gauze. See Figure 1, number 2.

Note: Do not use a mineral- based oil or petroleum jelly to lubricate the replacement tube.

8) Carefully insert the balloon shaft into the stomach through the existing stoma until the balloon is completely through the stoma tract.
9) Inflate the balloon to the rated volume as indicated on the device, using sterile or distilled water and the supplied sterile syringe. Do not over-inflate the balloon.
10) Pull replacement tube outward gently, until the balloon is up against the stomach wall. See Figure 2.


Figure 2

11) Slide the external bolster down to the skin.
12) Ensure that the replacement tube can rotate freely without gastric leakage.
 - If the replacement tube cannot rotate freely, and/or there is evidence of excessive gastric leakage, reconfirm that the replacement tube chosen is the appropriate French size.
13) Ensure that the balloon has passed completely through the stoma site and is residing in the GI tract.
14) Once placement confirmed, flush with a minimum of 10 mL of sterile water.

Feeding

1) Disconnect the dust cap/plug on the replacement tube's large inline port. See Figure 3, number 1.


Figure 3

2) Flush the replacement tube's large inline port with water and inspect for leakage.
 - If there is evidence of device leakage or excessive gastric leakage, the replacement tube may need to be replaced.

3) Connect the feeding syringe or bag connector to the large inline port. See figure 3, number 1. Ensure that the feeding connector is pressed in and seated firmly into the large inline port.
4) Administer the feeding per physician's instruction.
 - Medication should be administered in liquid form. Consult with pharmacist if tablets or capsules are to be crushed and mixed with water.
 - If continuous pump feeding, flush with 10-20 mL warm water every 6 hrs through the med-port. See Figure 3, number 2.
5) When feeding is complete, flush large inline port with a minimum of 20 mL of warm water.
6) Replace the dust cap/plug on the feeding port.
7) Wash and thoroughly rinse the feeding syringe (if used).

Maintenance

1) Before each feeding, rotate the replacement tube.
 - If the replacement tube does not rotate freely, consult the physician or healthcare provider.
2) Clean the stoma site and replacement tube daily with mild soap and water, and be sure to dry both thoroughly.

 Caution: Evidence of stoma infection, irritation, or tissue granulation should be addressed medically.
3) Check the balloon volume weekly as recommended by the physician or healthcare provider.
 - Notify the physician or healthcare provider if there is any evidence of leakage.
4) Inspect the feeding set for cracks, feeding solution, or foreign material before feeding.
 - If any are present, consult the physician or healthcare provider.
5) The outside of the replacement tube should be cleaned with mild soap and water and thoroughly flushed with water after each use.

Removal

1) Rotate the balloon approximately 180 degrees.

 Note: If the balloon can not be rotated consult a physician prior to removing the replacement tube.
2) Completely deflate the balloon.
3) Gently remove the replacement tube from the Stoma site.
4) Replace as or if necessary following the chosen replacement tube's directions for use.

 Note: Spontaneous closure of the stoma may occur as early as twenty-four (24) hours after removal. Insert a new replacement tube if Enteral feeding by this route is still intended.

WARRANTY

The manufacturer warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Boston Scientific Corporation (BSC) makes no warranty, express or implied, independent of the manufacturer. Descriptions, statements or specifications appearing in promotional literature are meant solely to generally describe the product at the time of manufacture and do not constitute any express or implied warranties. Handling, storage, cleaning and sterilization of this instrument, as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond the control of the manufacturer and/or BSC, directly affect the instrument and the results obtained from its use. The manufacturer's obligation under this warranty is limited to the repair or replacement of this instrument if deemed by the manufacturer to be defective in material or workmanship, and neither the manufacturer nor BSC shall be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. The user assumes responsibility for any consequences which may arise pertaining to the use of the products. **BSC and the manufacturer neither assume, nor authorize any other person to assume for them, any other or additional liability or responsibility in connection with this instrument.**

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© 2015 Boston Scientific Corporation or its affiliates. All rights reserved.