Boston Scientific

RescueNet™

Retrieval Device

Directions for Use	**2**
Instrucciones de uso	**5**
Mode d'emploi	**8**
Gebrauchsanweisung	**11**
Istruzioni per l'uso	**14**
Gebruiksaanwijzing	**17**
Instruções de Utilização	**20**

50710957-01

2019-10

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AS), DFU, MB, RescueNet, Global, 50710957-01B

Black (K) ΔE ≤5.0

RescueNet™

Retrieval Device

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING


Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The device consists of a flexible wire cable and snare loop with net which can be extended and retracted from the flexible outer sheath using a three ring handle. The device is for use with an endoscope with a working channel that is 2.8 mm or greater.

User Information

Precaution: Read the entire Directions for Use thoroughly before using RescueNet. This device should only be used by or under the supervision of physicians thoroughly trained in the use of the device. A thorough understanding of the technical principles, clinical applications, and risks associated with endoscopy procedures is necessary before using the device.

Contents

Each packaged pouch contains the following:

- One (1) RescueNet Retrieval device

INTENDED USE/INDICATIONS FOR USE

This device is intended for use in endoscopic retrieval of foreign objects, food bolus, tissue fragments and excised tissue such as polyps.

CONTRAINDICATIONS

Contraindications include those specific to the primary endoscopic procedure to be performed in gaining access to the desired site.

WARNINGS

- Do not use this device for any purpose other than the stated intended use.
- Before using, carefully examine the unit to verify that neither the contents nor the sterile package have been damaged in shipment. Do not use if damaged. Immediately return damaged product to Boston Scientific Corporation.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AS), DFU, MB, RescueNet, Global, 50710957-01B

Black (K) ΔE ≤5.0

Precautions

- Avoid blindly passing the device past any foreign object, particularly if the entire lumen is blocked.
- Use care to avoid inadvertently grasping tissue or organs not intended for retrieval.
- Keep gentle traction on the device during retrieval so that the retrieved object does not become loosened, dislodged or aspirated into the trachea.
- This device is not recommended for the retrieval of sharp foreign objects.
- After extraction, if the retrieved material does not exit the net upon deployment, the material should be removed from the net by rinsing the net in a basin of water; do not "pick" retrieved material from the net with a finger or tool as this may damage the net.

ADVERSE EVENTS

- Bleeding
- Perforation

HOW SUPPLIED

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Read the instructions for use in its entirety prior to use.

1. Prior to advancing the device into the endoscope, uncoil the device and drape in a "U" shaped configuration. The device can be damaged if it is coiled while the handle is operated.
2. Holding the proximal end of the device in one hand and the distal sheath in the opposite hand, slowly retract the finger rings until the net is completely withdrawn into the sheath. Confirm that there are no significant deformations or abnormalities.
3. Once the foreign object, food bolus, tissue sample or excised polyp for retrieval has been identified and with the net fully retracted within the catheter, advance the catheter into the accessory channel using short, 1.0 – 1.5 cm strokes until the distal end of the sheath is endoscopically visualized.
4. Advance the distal end of the sheath slightly past the foreign object, food bolus, tissue sample or excised polyp, however, avoid blindly passing the device past any foreign object, food bolus, tissue sample or excised polyp if the entire lumen is blocked. Open the net by advancing the finger rings forward. Confirm that the net is fully open via endoscopic observation.
5. Endoscopically manipulate the net over the foreign object, food bolus, tissue sample or excised polyp.
6. Retrieve the foreign object, food bolus, tissue sample of excised polyp by slowly retracting the finger rings on the handle. Gentle traction should be applied on the handle to keep the net closed.
7. Retract the catheter retrograde toward the endoscope so that the net is endoscopically observed yet not obscuring the endoscopic view. Endoscopic observation is necessary during extubation so as not to lose sight of the foreign object, food bolus, tissue sample or excised polyp during removal.

Precaution: Extubate the endoscope and net simultaneously. Do not attempt to pull the net and captured object through the working channel.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AS), DFU, MB, RescueNet, Global, 50710957-01B

Black (K) ΔE ≤5.0

8. Once the endoscope and the net have been extubated, the retrieved foreign object, food bolus, tissue sample or excised polyp can be removed from the net by advancing the finger rings forward to open the net. Specimens should be prepared according to standard technique for histologic evaluation.

PRODUCT DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AS), DFU, MB, RescueNet, Global, 50710957-01B

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2019 Boston Scientific Corporation or its affiliates.