Right Angle Feeding Set For Use With One-Step Button and Replacement Button

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DIRECTIONS FOR USE

WARNING

Contents supplied NON STERILE. Do not use if damaged. If damage is found, please contact your Boston Scientific representative.

For single patient use only. Do not reuse, reprocess or sterilize. Reuse, reprocessing or sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DESCRIPTION:

The Microvasive® Right Angle Feeding Set is a single patient use, disposable device used to attach a pump, syringe or gravity enteral feeding device to the Microvasive® One-Step Button and Replacement Button, for direct feeding.

INDICATIONS FOR USE:

The Right Angle Feeding Set is intended to provide nutrition to a patient directly into the stomach through a gastrostomy tube. It is indicated for use on patients who are unable to consume nutrition by conventional means.

CONTRAINDICATIONS FOR USE:

None known.

PRECAUTIONS:

Do not use this product without reading and understanding the complete instructions enclosed herein.

Note: This product contains DEHP (Bis(2-ethylhexyl) phthalate). BSC has conducted a Chemical Safety Assessment of DEHP in this device based on current theoretical acceptable exposure limits for sensitive populations (infant, children and pregnant/ nursing women), and determined that the exposure to DEHP from this device is comparable to levels reported for other enteral feeding devices[1]. The risk associated with this level of exposure to DEHP has not yet been established.

POTENTIAL COMPLICATIONS:

The following complications may occur with the use of direct feeding devices:

- Fever
- Gastric Distention
- Infection
- Blockage
- Tissue Necrosis

INSTRUCTIONS FOR USE:

Prior to Use:

CAUTION: *The Microvasive® Right Angle Feeding Set is supplied in a sealed package and is intended for single patient use only. Carefully examine the unit to verify that neither the contents nor the package has been damaged in shipment.* ***DO NOT USE*** *if damaged. Immediately return damaged product to Microvasive®.*

Directions:

1. Gently push adapter into the Button. Verify that the adapter is completely inserted.

 WARNING: DO NOT USE EXCESSIVE FORCE. TOO MUCH PRESSURE MAY LOOSEN THE BUTTON, CAUSING IT TO BE PUSHED INTO THE STOMACH OR PREMATURELY PULLED OUT.

2. Attach pump, syringe or gravity enteral feeding device to adapter.
3. After feeding, use a 60 cc syringe to flush device with water.
4. Using a slight twisting motion, remove feeding set.
5. Snap Button safety plug in place to keep lumen clean.
6. Clean feeding set with soap and water. Thoroughly rinse.

HOW SUPPLIED:

The Microvasive® Right Angle Feeding Set is separately packaged as one (1) unit. The device is non-sterile.

STORAGE:

Store at controlled room temperature. Do not expose to organic solvents, ionizing radiation or ultraviolet light.

WARRANTY:

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness. Handling, storage, cleaning and sterilization of this instrument as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument.

CAUTION: *Federal law (USA) restricts this device to sale or by the order of a physician.*

[1] Takatori, S., Okamoto, Y., Kitagawa, Y., Hori, S., Izumi, S., Makino, T., and Nakazawa, H. (2008). Simulated neonatal exposure to DEHP and MEHP from PVC enteral nutrition products. Int J Pharm 352, 139-145.

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© 2010 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific International S.A.
Boston Scientific (Australia) Pty Ltd