Boston Scientific

IMAGER™ II

Angiographic Catheter

en	Instructions for Use	2
es	Instrucciones de uso	7
fr	Mode d'emploi	12
de	Gebrauchsanweisung	17
it	Istruzioni per l'uso	22
nl	Instructies voor gebruik	27
pt-EU	Instruções de Utilização	32

51457560-01

2023-03

BSC (MB IFU Template 3 x 9 Global, 92310029O), IFU, MB, IMAGER II, Global 7, 51457560-01A

MB Drawing 50573137

Black (K) ΔE ≤5.0

TABLE OF CONTENTS

REUSE WARNING **3**

DEVICE DESCRIPTION **3**

Contents **3**

Operating Principle **3**

Materials **3**

Non-Pyrogenic **3**

User Information **3**

INTENDED USE / INDICATIONS FOR USE **3**

Clinical Benefit Statement **3**

CONTRAINDICATIONS **4**

WARNINGS **4**

PRECAUTIONS **4**

ADVERSE EVENTS **4**

HOW SUPPLIED **5**

Device Details **5**

Handling and Storage **5**

OPERATIONAL INSTRUCTIONS **5**

Additional Items for Safe Use **5**

Preparation **5**

Procedure **5**

Disposal **5**

Post-Procedure **6**

PATIENT COUNSELING INFORMATION **6**

WARRANTY **6**

SYMBOL DEFINITIONS **6**

BSC (MB IFU Template 3 x 9 Global, 92310029O), IFU, MB, IMAGER II, Global 7, 51457560-01A

Black (K) ΔE ≤5.0 MB Drawing 50573137

IMAGER™ II
Angiographic Catheter

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

REUSE WARNING

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

DEVICE DESCRIPTION

Contents

Imager II Angiographic Catheter (check the label for product's characteristics prior to use).

Operating Principle

The Imager II Angiographic Catheters are sterile, single use diagnostic intravascular catheters. These catheters are available in a variety of shapes and lengths. The distal segment of these catheters is radiopaque to aid in visualization of the device under fluoroscopy. These catheters are used to create images of specific areas of interest.

The Imager II Angiographic Catheter can also be used to deliver embolic coils.

Materials

Imager II Angiographic Catheters are made of a polymer shaft braided with stainless steel, a polyamide hub, and a soft tip made of polymer - tungsten with hydrophilic coating.

Non-Pyrogenic

These devices are not pyrogenic.

User Information

Imager II Angiographic Catheters should only be used by physicians and other health care professionals experienced in angiographic procedures.

INTENDED USE / INDICATIONS FOR USE

The Imager II Angiographic Catheters are designed to provide a pathway for delivering contrast media to selected sites in the peripheral vascular system.

Additionally, the Imager II Angiographic Catheters can also be used for the controlled and selective delivery of the Boston Scientific Interlock - 35 Detachable Coils into the peripheral vasculature.

Clinical Benefit Statement

The Imager II Angiographic Catheter provides access to the area of treatment within the peripheral vasculature, a route for contrast administration and a pathway for the advancement of medical instruments. It allows the physician to maintain lesion access and can help to deliver and determine the appropriate positioning of an interventional medical device in the intravascular system.

- The Imager II Angiographic Catheters allow for the delivery of contrast which aids for a successful angiographic procedure.
- The Imager II Angiographic Catheters also provides for the successful delivery of Interlock - 35 Detachable Coils.

BSC (MB IFU Template 3 x 9 Global, 92310029O), IFU, MB, IMAGER II, Global 7, 51457560-01A

MB Drawing 50573137 Black (K) ΔE ≤5.0

CONTRAINDICATIONS

No known contraindications.

WARNINGS

- Do not store this product in an area that is subject to aggressive disinfectant exposure, such as vaporized or in-solution hydrogen peroxide. Catheter integrity can be compromised if exposed to these substances.
- The risk associated with Imager II Angiographic Catheters outside the intended use is not known.
- Do not use if the device is damaged (e.g., component detachment, breaks, kinks or deformations), this may cause adverse events.
- Due to the various sizes and shapes, and the relative stiffness of angiographic catheters, extreme care must be taken when selecting and using the Imager II Angiographic Catheters to avoid damage to the vessel walls through which this catheter passes.
- Precautions to prevent or reduce clotting should be taken when any catheter is used in the peripheral vascular system. The catheter should always be filled with either heparinized saline solution or contrast medium. Use of systemic heparinization should be considered.
- To facilitate the introduction of the catheter and avoid damage to the puncture site, the catheter should be introduced through an introducer sheath and over an appropriate guide wire.
- Do not exceed the pressure ratings printed on the device labeling. Catheter integrity may be compromised if pressure limits are exceeded.
- If resistance is encountered when removing the guide wire from the catheter, simultaneously remove the catheter and the guide wire as a single unit under fluoroscopy to prevent potential damage to the vessel and/or product damage.

PRECAUTIONS

- Use prior to the "Use By" date.
- Carefully inspect the catheter prior to use to verify that the size and shape are suitable for the specified procedure and that no damage has occurred during shipment.
- Any use of these devices, other than those indicated in the instructions is not recommended.
- The Imager II Angiographic Catheters are clinically hydratable, and the effective length may increase up to 4.4 mm after hydration. There is no effect on the outer diameter of the catheter during hydration.
- The Imager II Angiographic Catheters are designed for use with specific guide wire diameters. The recommended maximum guide wire diameter is provided on the catheter label.
- If resistance is encountered during catheter manipulation, determine the cause under fluoroscopy and take the necessary remedial action.
- If any leak is detected in the device during flush or contrast media injection, exchange the catheter.

ADVERSE EVENTS

Potential adverse events which may be associated with the use of angiographic catheters when used as indicated include but are not limited to:

- Allergic reaction (contrast, device or other)
- Cerebrovascular Accident (CVA), Transient Ischemic Attack (TIA)
- Death
- Embolism (air, plaque, thrombus, device or other)
- Hematoma
- Hemorrhage
- Infection/Sepsis

BSC (MB IFU Template 3 x 9 Global, 92310029O), IFU, MB, IMAGER II, Global 7, 51457560-01A

Black (K) ΔE ≤5.0 MB Drawing 50573137

- Ischemia
- Need for additional intervention or surgery
- Thrombosis/Thrombus
- Vasospasm
- Vessel Occlusion
- Vessel Trauma (arteriovenous fistula, dissection, perforation, pseudoaneurysm, rupture, or injury)

HOW SUPPLIED

The basic configuration of Imager II Angiographic Catheter is one (1) device held by a packaging card inside a pouch. Five (5) pouches are then placed in a shelf carton with the IFU. Each pouch and shelf carton are labeled with content identification and traceability information.

Device Details

These devices are sterilized using Ethylene Oxide (EO).

Do not use if package is damaged or unintentionally opened before use. If damage is found, notify your Boston Scientific representative.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Do not store this product in an area that is subject to aggressive disinfectant exposure, such as vaporized or in-solution hydrogen peroxide.

OPERATIONAL INSTRUCTIONS

Additional Items for Safe Use

The following supplies are not provided and need to be available and prepped prior to use of Imager II Angiographic Catheter:

- Appropriate introducer/ guide sheath.
- Contrast medium.
- Heparinized saline or other appropriate flush solution.
- Appropriate interventional/ therapeutic devices compatible with Imager II Angiographic Catheter.

Preparation

Prior to use, carefully examine the unit to verify that the catheter or sterile package has not been damaged in shipment.

Procedure

1. Obtain vascular access using standard technique and appropriate size introducer sheath and guide wire.
2. Locate the guide wire lumen hub and flush with heparinized saline solution.
3. Advance the appropriately sized guide wire and position in the vessel.
4. Advance the catheter over the guide wire through the sheath and/or puncture site and position the catheter in the vessel.
5. Remove the guide wire.
6. Engage the vessel with the catheter.
7. Perform a small test injection with contrast media to ensure the catheter is appropriately positioned.
8. Perform procedure according to standard practice.

Disposal

To minimize the risk of infection or microbial hazards after use, dispose device and packaging as follows:

After use, device may contain biohazardous substances. Devices that contain biohazardous substances are considered biohazardous waste and should be stored in a biohazard container that is labeled with the biological hazard symbol. Untreated biohazardous waste should not be disposed of in the municipal waste system. Biohazardous substances should be removed via incineration prior to disposal. Alternatively, biohazardous waste may be disposed

BSC (MB IFU Template 3 x 9 Global, 92310029O), IFU, MB, IMAGER II, Global 7, 51457560-01A

MB Drawing 50573137 Black (K) ΔE ≤5.0

of utilizing a certified facility for biohazardous waste for proper treatment in accordance with hospital, administrative, and/or local government policy.

Post-Procedure

Assess patient for hematoma and/or other signs of bleeding at the puncture site.

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

PATIENT COUNSELING INFORMATION

The physician should consider the following points while counseling patients on the use of the Imager II Angiographic Catheter in association with the interventional procedure:

Discuss the risks and benefits including review of potential adverse events listed in this document, both for the Imager II Angiographic Catheter and for other interventional treatments likely to be employed.

Discuss patient allergies, in particular the risk for patients who may be allergic to contrast.

Discuss the post procedure instructions, including any follow-up requirements, lifestyle changes, medications, and home-care or rehabilitation guidelines.

WARRANTY

For device warranty information, visit (www.bostonscientific.com/warranty).

The following are trademarks of Boston Scientific Corporation or its affiliates: Imager, Interlock.

All other trademarks are the property of their respective owners.

EU Importer: Boston Scientific International B.V., Vestastraat 6, 6468 EX Kerkrade, The Netherlands.

SYMBOL DEFINITIONS

Commonly used medical device symbols that appear on the labeling are defined at www.bostonscientific.com/SymbolsGlossary.

Additional symbols are defined at the end of this document.

BSC (MB IFU Template 3 x 9 Global, 92310029O), IFU, MB, IMAGER II, Global 7, 51457560-01A

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© 2023 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd