Boston Scientific

SAFARI 2™

PRE-SHAPED TAVR/TAVI

Guidewire

Instructions for Use	**2**
Instrucciones de uso	**8**
Mode d'emploi	**13**
Gebrauchsanweisung	**18**
Istruzioni per l'uso	**23**
Instructies voor gebruik	**28**
Instruções de Utilização	**33**

51145320-01

2020-06

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 92238512A), IFU, MB, SAFARI 2, Global, 51145320-01A

Black (K) ΔE ≤5.0

TABLE OF CONTENTS

WARNING **3**

DEVICE DESCRIPTION **3**

Contents **4**

INTENDED USE/INDICATIONS FOR USE **4**

TARGET POPULATION **4**

CONTRAINDICATIONS **4**

WARNINGS **4**

ADVERSE EVENTS **5**

HOW SUPPLIED **5**

Handling and Storage **5**

OPERATIONAL INSTRUCTIONS **5**

Preparation for Use **5**

Instructions for Use **6**

WARRANTY **6**

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 92238512A), IFU, MB, SAFARI 2, Global, 51145320-01A

Black (K) ΔE ≤5.0

SAFARI²™

PRE-SHAPED TAVR/TAVI

Guidewire

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING


Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The SAFARI² Guidewire consists of an 0.035 inch (≤0.92 mm) outside diameter by 275 cm working length guidewire comprised of a spring coil and core. The distal tip has a double-curve. The entire spring coil is coated with a LUBRIGREEN™ PTFE coating. Guidewire Length, Diameter and Curve configurations are indicated on the product label. The SAFARI² Guidewire is packaged with a J-Straightener that is used to aid in the insertion of the guidewire into the delivery catheter.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 92238512A), IFU, MB, SAFARI 2, Global, 51145320-01A

Black (K) ΔE ≤5.0

Contents

Qty	Material
1	Guidewire
1	J-Straightener

INTENDED USE/INDICATIONS FOR USE

The SAFARI²™ Guidewire is intended to facilitate the introduction and placement of interventional devices within the chambers of the heart including those used within transcatheter aortic valve procedures.

TARGET POPULATION

SAFARI² Guidewires are intended to be used in adult populations who do not meet contraindications and regardless of gender.

CONTRAINDICATIONS

This wire is not intended for use in the cerebrovasculature or coronary arteries.

WARNINGS

- The SAFARI² Guidewire should be used only by physicians trained in the introduction and placement of interventional devices including those used within transcatheter aortic valve procedures.
- Carefully read all instructions prior to use. Observe all warnings and precautions. Failure to do so may result in complications.
- Prior to use, inspect for damage. If damaged, DO NOT USE.
- Monitor wire position throughout the procedure for proper placement of curve and distal tip.
- Do not torque this guidewire.
- Use extreme caution and careful judgment in patients for whom anticoagulation is not indicated.
- Inserting the distal end of the guidewire through a Tuohy-Borst Adapter may result in damage to the tip of the guidewire.
- If contrast agents are used, use extreme caution in patients who have had a severe reaction to contrast agents and who cannot be adequately premedicated.
- When advancing or removing the guidewire, always use fluoroscopic guidance with radiographic equipment that provides high-resolution images. Never position the guidewire blindly, as this may result in misplacement, dissection or perforation.
- Exercise care in handling of the guidewire during a procedure to reduce the possibility of accidental breakage, bending, kinking, or coil separation. Resulting guidewire fractures might require additional percutaneous intervention or surgery.
- Never advance the guidewire against the resistance without first determining the reason for resistance under fluoroscopy. Excessive force against resistance may result in damage to the catheter or vessel/organ. Care should be taken when advancing a guidewire after device deployment.
- The wire should only be introduced into or withdrawn from the ventricle through a catheter already positioned in the ventricle.
- The curve of the SAFARI² Guidewire should be constrained within a catheter during insertion into or withdrawal from the body or treatment site.
- The SAFARI² Guidewire is manufactured with a double curve; attempts to modify may alter its performance. Alterations to curve may lead to complications including Perforation or Dissection, Mitral Valve Regurgitation,

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 92238512A), IFU, MB, SAFARI 2, Global, 51145320-01A

Black (K) ΔE ≤5.0

Pericardial Effusion, Cardiac Tamponade, Cardiac Arrest and Guidewire Replacement.

- Clinical data showing use of this device in pregnant/ breastfeeding women is not available.

ADVERSE EVENTS

- Access site complication
- Additional Surgical Procedure
- Air Embolism/Thromboembolism
- Allergic Reaction
- Amputation
- Aorta Complications
- Arteriovenous (AV) Fistula
- Arrhythmia
- Bleeding
- Cardiac and/or Septal Perforation
- Death
- Embolism
- Hematoma
- Hemorrhage
- Hemoglobinuria
- Hypovolemia
- Infection or Sepsis
- MACCE
- Myocardial Ischemia and/or Infarction
- Pericardial Effusion
- Pseudoaneurysm
- Renal Failure or Injury
- Stroke or other Neurologic event
- Tamponade
- Thrombus
- Valve Complications
- Vascular Complication
- Vessel Occlusion
- Vessel Perforation, Dissection, Trauma, or Damage
- Vessel Spasm
- Wire Entrapment/Entanglement
- Foreign Body/Wire Fracture

HOW SUPPLIED

This product is supplied sterile and intended for single use.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Preparation for Use

- Prior to use, carefully examine the sterile package and verify the following. Do not use if any items cannot be verified.
 - Sterile package is free of damage, unopened, and packaging is intact.
 - Labeling is complete and legible.
 - Device will be used prior to the use-by date on the label.
- SAFARI^{2TM} Guidewires should be handled carefully and inspected before use and when possible during the procedure to observe for any defect or damage that may occur. Do not use a wire that has a damaged

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 92238512A), IFU, MB, SAFARI 2, Global, 51145320-01A

Black (K) ΔE ≤5.0

coating, tip, camber (change in plane), bend or kink on the wire. Damage will hinder the performance of the SAFARI²™ Guidewire.

- Prior to use and when possible during the procedure, inspect the guidewire carefully for coil separation, bends, or kinks which may have occurred. Do not use a wire which has a damaged tip. Damage will hinder the performance of guidewire.
- If coils of a guidewire separate, DO NOT REMOVE THE CORE. Carefully remove the coils and core simultaneously.
- Inspect and prepare the catheter according to the manufacturer's instructions. This includes flushing the catheter to be used with heparinized saline solution.
- Verify compatibility of interacting devices prior to use.

Instructions for Use

1. Ensure a catheter is appropriately placed in the ventricle or other intended treatment area.
2. Carefully remove the SAFARI² Guidewire from the hoop by grasping the proximal end of the J-Straightener separating the straightener from the hoop.
3. After inspection of the SAFARI² Guidewire, advance the J-Straightener over the distal section of the SAFARI² wire to straighten the curve.
4. Insert the SAFARI² Guidewire into the hub of the catheter with the aid of the J-Straightener.
5. Carefully advance the distal tip of the SAFARI² Guidewire into the lumen of the catheter.
6. Remove the SAFARI² Guidewire J-Straightener by withdrawing it over the length of the SAFARI² Guidewire.
7. Advance the SAFARI² Guidewire through the catheter under fluoroscopic guidance.
8. Confirm SAFARI² Guidewire curve position to assure SAFARI² Guidewire placement and stability.
9. Maintain wire position while tracking or advancing a catheter or device over the wire.
10. To remove the guidewire from the treatment area:
 a. A catheter must be advanced into the treatment area over the SAFARI² Guidewire prior to withdrawing the wire.
 b. The SAFARI² Guidewire is pulled back completely into the catheter.
 c. The SAFARI² Guidewire may then be withdrawn from the catheter.
11. When finished with the procedure, dispose of guidewire by using standard methods used in a hospital environment.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes**

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 92238512A), IFU, MB, SAFARI 2, Global, 51145320-01A

Black (K) ΔE ≤5.0

no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 92238512A), IFU, MB, SAFARI 2, Global, 51145320-01A

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© 2020 Lake Region Medical.
Boston Scientific (Australia) Pty Ltd