Boston Scientific

Mach1™

Guide Catheters

Directions for Use	**2**
Instrucciones de uso	**5**
Mode d'emploi	**8**
Gebrauchsanweisung	**11**
Istruzioni per l'uso	**15**
Gebruiksaanwijzing	**18**
Instruções de Utilização	**21**

90979434-01

2019-10

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 91031819AA), DFU, MB, Mach1, Guide Cath, CV PV, Global, 90979434-01E

CV 03

Black (K) ΔE ≤5.0

Mach1™
Guide Catheters

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

Intravascular Guide Catheter

Non-Pyrogenic

INTENDED USE/INDICATIONS FOR USE

Boston Scientific Guide Catheters provide a pathway through which medical instruments, such as balloon dilatation catheters, guidewires or other therapeutic devices may be introduced. These devices are not intended for use in the cerebral vasculature.

Boston Scientific Peripheral Guide Catheters are intended for use in general intravascular and peripheral applications.

Boston Scientific Coronary Guide Catheters are intended for use in general intravascular and coronary applications.

CONTRAINDICATIONS

None Known.

WARNINGS

Guide catheters are designed for use by physicians engaged in the practice of a specialized branch of medicine. Use of these devices should be restricted to those specialists trained to perform the procedure. A thorough understanding of the technical principles, clinical applications, and risks associated with angioplasty is necessary before performing this procedure.

Due to the size and relative stiffness of guide catheters, extreme care must be taken to avoid damage to the vessel walls through which this catheter passes.

This catheter may occlude smaller vessels. Care must be taken not to completely block blood flow.

PRECAUTIONS

Manipulate the guide catheter only while under fluoroscopic observation. Care should be taken to control the position of the guide catheter tip during manipulation of the therapeutic device.

Inspect the guide catheter prior to use to verify the size, shape, and condition of the catheter as suitable for the specific procedure.

Check the labeled diameter of both the therapeutic device and guide catheter prior to use to ensure compatibility.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 91031819AA), DFU, MB, Mach1, Guide Cath, CV PV, Global, 90979434-01E

Black (K) ΔE ≤5.0

Use of a power injector is not recommended. The large internal diameter of the catheter permits rapid injection with little force being required on the syringe. Hand inject slowly whenever attempting to opacify the vessels via this catheter. The maximum recommended injection pressure is 500 psi / 3447 kPa.

Use prior to the "Use By" date.

Rotate inventory to maintain fresh stock.

To avoid damaging the tip of this catheter during removal from the package, open as follows: Grasp the hub and slowly withdraw the catheter. Rapid removal may result in damage to the catheter shaft or tip.

ADVERSE EVENTS

Potential adverse reactions which may result from the use of this device include but are not limited to:

- Access site pain
- Allergic reaction
- Death
- Hemorrhage or hematoma
- Impaired blood flow due to thrombosis, embolism, or vasospasm that could lead to myocardial infarction, transient ischemic attack (TIA), or stroke
- Infection
- Pseudoaneurysm
- Trauma that could lead to surgical intervention
- Vessel trauma including dissection or perforation

HOW SUPPLIED

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

REQUIRED MATERIALS

Appropriately sized introducer sheath; sterile saline; guidewire appropriately sized and of appropriate length; connecting tube(s); manifold; Y-adapter; and contrast medium.

DIRECTIONS FOR USE

1. Carefully inspect the guide catheter prior to use to verify that it has not been damaged during handling. Do not use if catheter is kinked or the tip is damaged.
2. Guide catheters may be introduced into the vascular system by cutdown or percutaneous techniques.
3. Flush the guide catheter with heparinized saline.
4. Insert the guide catheter through the arteriotomy or the sheath using an appropriately sized guidewire.
5. Attach a connecting tubing to a manifold and the side arm of the Y-adapter. Flush and fill the assembly with normal saline.
6. Advance the guide catheter to the vascular site for the **peripheral** indication, for the **coronary** indication advance guide catheter to the ascending aorta. Withdraw the guidewire, and connect the guide catheter to the Y-adapter.
7. Allow back bleeding through the Y-adapter and the guide catheter to remove any trapped air, and then flush with normal saline. To remove blood from the valve section of the Y-adapter, open the valve slightly while continuing to flush the device with normal saline and close the valve as soon as the blood is purged.
8. Fill the Y-adapter and guide catheter with contrast medium.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 91031819AA), DFU, MB, Mach1, Guide Cath, CV PV, Global, 90979434-01E

Black (K) ΔE ≤5.0

9. Selectively engage the catheter in the appropriate coronary ostium or other selected peripheral location. Confirm the catheter position by arteriography through the guide catheter. For **Coronary** indication, it is important to determine that coronary ostial pressure measurements are appropriate as confirmation that the guide catheter is properly positioned at the coronary orfice.
10. Insert the therapeutic device into the guide catheter through the Y-adapter valve. Perform the therapeutic procedure according to the directions provided by the manufacturer of the therapeutic device. Arteriography through the guide catheter should be utilized during the therapeutic procedure to confirm the location of the guidewire and/or therapeutic device. Care should be taken not to allow air to enter the catheter through the Y-connector or manifold.
11. Disconnect the Y-adapter from the guide catheter. Insert a guidewire into the guide catheter and remove the guide catheter over the guidewire.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 91031819AA), DFU, MB, Mach1, Guide Cath, CV PV, Global, 90979434-01E

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2019 Boston Scientific Corporation or its affiliates.