Boston Scientific

Impulse™ and Expo™

Angiographic Catheter

Directions for Use	**2**
Instrucciones de uso	**5**
Mode d'emploi	**7**
Gebrauchsanweisung	**9**
Istruzioni per l'uso	**11**
Gebruiksaanwijzing	**13**
Instruções de Utilização	**15**

50559966-01

2019-10

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 91031819AA), DFU, MB, Impulse and Expo, Global, 50559966-01B

CV 01

Black (K) ΔE ≤5.0

Impulse™ and Expo™

Angiographic Catheter

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Impulse and Expo Angiographic Catheters are single lumen catheters with a single proximal hub for attachment to a dye delivery system. The radiopaque polymer shaft is reinforced with braided stainless steel. A soft atraumatic tip protects against vessel damage. Device may have side holes. Refer to product label. Both catheters are available in three general categories of distal curves; Pigtail, Multipurpose and Selective. These preformed shapes allow access into a variety of vessels. The Impulse Angiographic Catheter is available in a variety of lengths with outer diameters of 5F (1.70 mm) or 6F (2.1 mm). The Expo Angiographic Catheter is available in a variety of lengths with outer diameters of 5F (1.70 mm), 6F (2.1 mm), or 7F (2.4 mm).

Non-pyrogenic.

INTENDED USE/INDICATIONS FOR USE

Boston Scientific Angiographic Catheters are designed to provide a pathway to be used for delivering contrast media to selected sites in the vascular system during an angiographic procedure.

CONTRAINDICATIONS

None Known.

WARNINGS

Do not exceed the pressure ratings printed on the label.

Catheter integrity may be compromised if the pressure exceeds the maximum pressure listed.

PRECAUTIONS

Angiographic catheters are designed for use by physicians engaged in the practice of a specialized branch of medicine. Use of these devices should be restricted to those specialists trained to perform the procedure. A thorough understanding of the technical principles, clinical applications and risks associated with angiography is necessary before performing this procedure.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 91031819AA), DFU, MB, Impulse and Expo, Global, 50559966-01B

Black (K) ΔE ≤5.0

Use prior to the "Use By" date.

Rotate inventory to maintain fresh stock.

To avoid damaging the tip of the catheter during removal from the package, open as follows: Grasp the hub and slowly withdraw the catheter. Rapid removal may result in damage to the catheter tip. Do not use if catheter is kinked or tip is damaged.

Carefully inspect the catheter prior to use to verify the size and shape as suitable for the specific procedure and that it has not been damaged during handling.

Angiographic catheters should only be manipulated while under fluoroscopic observation.

To facilitate its introduction and avoid damage at the puncture site, the angiographic catheter should be introduced using a guidewire.

Precautions to prevent or reduce clotting should be taken when any catheter is used. While in the vascular system, the catheter should always be filled with either flushing solution or contrast medium. Use of systemic heparinization should be considered.

ADVERSE EVENTS

Potential adverse reactions which may result from the use of this device include but are not limited to:

- Access site pain
- Allergic and neurological reactions
- Arrhythmia
- Death
- Hemorrhage or hematoma
- Impaired blood flow due to thrombosis, embolism, or vasospasm that could lead to myocardial infarction, stroke, transient ischemic attack (TIA), limb amputation, and other thrombo-embolic organ damages such as renal infarction
- Infection
- Organ failure
- Vessel trauma

In addition, the use of contrast may cause renal failure or other adverse reactions.

HOW SUPPLIED

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

REQUIRED MATERIALS

Guidewire(s) with a maximum OD of 0.038 in (0.97 mm) and of appropriate length for advancement of catheter; sterile saline; luer-lock syringe(s); appropriately sized introducer sheath; contrast medium; manifold; and hemostatic adapter.

DIRECTIONS FOR USE

1. Flush catheter with heparnized saline.
2. Advance guidewire through introducer sheath and position in or near target vessel.
3. While maintaining guidewire position, introduce the catheter over the guidewire and through the sheath.
4. Advance angiographic catheter through the sheath and position in or near target vessel.
5. Remove guidewire.
6. Connect angiographic catheter to manifold.
7. Aspirate until blood is seen.
8. Flush saline antegrade through angiographic catheter.
9. Flush contrast medium through angiographic catheter.
10. Engage target vessel.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 91031819AA), DFU, MB, Impulse and Expo, Global, 50559966-01B

Black (K) ΔE ≤5.0

11. A small test injection should be performed to insure the catheter is in a good position in the vessel.
12. Inject adequate amount of contrast media to obtain images.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 91031819AA), DFU, MB, Impulse and Expo, Global, 50559966-01B

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2019 Boston Scientific Corporation or its affiliates.