NOW !

Pulsed Biphasic Measurements

QA-45

MKIII Defibrillator-/ Transcutaneous Pacemaker Analyzer

PRODUCT HIGHLIGHTS

- Pulsed biphasic energy measurement
- Biphasic energy measurement
- Energy and cardioversion measurement
- Peak voltage and current reading
- Storage and playback of output waveform
- 12 lead ECG simulation
- ECG, Performance and Arrhythmia simulation
- Automatic defibrillator test procedures
- Pacemaker loads from 50 ohm to 2300 ohm
- Accuracy Tests; Amplitude, Rate, Pulse width
- Demand Sensitivity Tests; Amplitude, Rate, Pulse width, Square, Triangle, Haversine
- Sensed Refractory Period
- Paced Refractory Period
- 60/50 Hz Immunity Test
- Pulse Energy
- Large graphic display
- RS-232C and Centronic printer interface
- Battery operation or AC adapter

OVERVIEW

The QA-45 MKIII Defibrillator-/ Transcutaneous Pacemaker Analyzer is a precision instrument for testing defibrillators and external (transcutaneous) pacemakers. The build-in test procedures and extensive resistance loads for pacer testing, make the QA-45 a versatile test instrument and solid value.

The test result can either be written directly to a printer from the QA-45, or they may be transferred to *ansur* QA-40/45 test automation software. With *ansur*, you can make your own test protocols, store protocols and test results to disk, and provide test data to update your equipment management database.

In the interest of product improvement, Metron claims the right to alter specifications without notice.

QA-45 MKIII Specifications

DEFIBRILLATOR ANALYZER

ENERGY OUTPUT MEASUREMENT:

HIGH RANGE:

Voltage:	< 5000 V
Maximum current:	120 A
Maximum Energy:	1000 J
Accuracy:	±2 % of reading for > 100 J
	±2 J of reading for < 100 J
Trigger level:	100 volts
Playback amplitude:	1 mV/1000 V Lead I
Test pulse:	125 J ±20%

LOW RANGE:

Voltage:	< 1000 V
Maximum current:	24 A
Maximum Energy:	50 J
Accuracy:	±2 % of reading for > 20 J
	±2 J of reading for < 20 J
Trigger level:	20 volts
Playback amplitude:	1 mV/200 V Lead I
Test pulse:	5 J ±20%
Load resistance:	50 ohms ±1 %; non-inductive (< 1 µH)
Display resolution:	0.1 J
Measurement time Window:	100 ms
Absolute max. peak Voltage:	6000 V
Pulse width:	100 ms

Cardioversion: Measured time delay ±2 ms
Biphasic Measurement

OSCILLOSCOPE OUTPUT:

High Range:	1000:1 amplitude attenuation
Low Range:	200:1 amplitude attenuation

WAVEFORM STORAGE AND PLAYBACK:
Discharge can be viewed via ECG outputs and paddles

Output: 200:1 Time Base Expansion

SYNC TIME MEASUREMENTS:
Timing window: Starts 40 ms before each R wave peak
Test waveforms: All Waveform simulations available
Delay time accuracy: ±1 ms

CHARGE TIME MEASUREMENT:
From 0.1 to second to 99.9 seconds

ECG WAVE

ECG GENERAL:

Lead configuration:	12-lead simulation. RL, RA, LA, LL, V1-6
Output Impedance:	
Limb leads:	1000 ohms to RL
V Leads:	1000 ohms to RL
High level Output:	1.0 V/mV of low level (lead II)

MANUAL ECG PERFORMANCE TEST:

DC Pulse:	4.0 seconds 1.0 mV
Square wave:	2 Hz 1.0 mV p-p
Triangle wave:	2 Hz 1.0 mV
Sine waves:	0.1, 0.2, 0.5, 10, 40, 50, 60, 100Hz
Amplitudes:	0.5, 1.0, 1.5, 2.0 mV (Lead II)
Accuracy:	±5 % (Lead II 1.0 mV)

ECG PERFORMANCE TEST:

Gain/Damping:	2 Hz square wave
Frequency response:	
Low Frequency:	4 second DC pulse
Band Pass:	10 Hz sine
Monitor:	-3dB point: 40 Hz sine
Power line notch filter:	50 Hz sine
Linearity:	2 Hz triangle wave

NORMAL SINUS:

Rates:	30, 60, 80, 120, 180, 240 & 300 BPM
Accuracy:	±1 % of selection
Amplitudes:	0.5, 1.0, 1.5, 2.0 mV (Lead II)
Accuracy:	±5 % (Lead II 1.0 mV)

Automatic ECG rate Test

ARRHYTHMIA SELECTIONS:

vfib	– Ventricular Fibrillation
afib	– Atrial Fibrillation
blk II	– Second degree A-V block
RBBB	– Right Bundle Branch Block
PAC	– Premature Atrial Contraction
PVC_E	– Early PVC
PVC_STD	– PVC
PVCRonT	– R on T PVC
mfPVC	– Multifocal PVC
bigeminy	– Bigeminy
run5PVC	– Run of 5 PVCs
vtach	– Ventricular Tachycardia

SHOCK ADVISORY ALGORITHM TEST:
8 ECG signals are available to test the advisory algorithm of Automatic defibrillators:

- Asystole
- CVF: Coarse Ventricular Fibrillation
- FVF: Fine Ventricular Fibrillation
- MVT 140: Multifocal Ventricular Tachycardy @ 140 BPM
- MVT 160: Multifocal Ventricular Tachycardy @ 160 BPM
- PVT 140 : Polyfocal Ventricular Tachycardy @ 140 BPM
- PVT 160 : Polyfocal Ventricular Tachycardy @ 160 BPM
- SVTa 90 : Supra Ventricular Tachycardy @ 90 BPM

TRANSCUTANEOUS PACEMAKER ANALYZER

TEST LOAD RANGE:
50 to 2300 ohms in step of:
– 50 ohm up to 200 ohms
– 100 ohm from 200 up to 2300 ohms

Accuracy:	50 - 1300 ohm	±1 %
	1400 - 2300 ohm	±1.5 %

OSCILLOSCOPE OUTPUT:

50 - 150 ohm	10.24:1 amplitude attenuation
200 - 500 ohm	41:1 amplitude attenuation
600 - 2300 ohm	164:1 amplitude attenuation

PULSE MEASUREMENTS:

Amplitude:	4 to 300 mA (100 ohm load)
Accuracy:	±5 % or ±0.5 mA
Max. amplitude:	300 mA all loads
Rate:	30 to 800 ppm
Accuracy:	±1 % or 2 ppm
Pulse Width:	0.6 to 80 ms
Accuracy:	±1 % or ±0.3 ms

DEMAND SENSITIVITY TEST:

Waveforms: Square (SQR), Triangle (TRI) and Haversine (SSQ)

ECG OUTPUT:

Amplitude:	0 - 4 mV
Resolution:	40 µV

PACER INPUT (Load depended):

Amplitude	(50 ohm):	0 - 10 mV
Resolution	(50 ohm):	40 µV
Amplitude:	(≥500 ohm):	0 - 100 mV
Resolution:	(≥500 ohm):	1 mV

DEFIB. PADS:

Amplitude:	0 - 10 mV
Resolution:	0.1 mV

Waveform width:	10, 25, 40, 100 and 200 ms
Pacer rate:	30 to 120 ppm

IMMUNITY TEST:
50/60 Hz Interference Signal:

ECG OUTPUT:
0 - 4 mV peak in step of 0.4 mV

PACER INPUT (Load depended):
0 - 10 mV peak in step of 1 mV (50 ohm)
0 - 100 mV peak in step of 10 mV (≥500 ohm)

DEFIB. PADS:
0 - 10 mV peak in step of 1 mV

REFRACTORY PERIOD MEASUREMENT:
20 to 500 ms (both Pacing and Sensing)
Accuracy: ±2 ms

ansur QA-40/45 OFFERS THE FOLLOWING FEATURES

ansur QA-40/45 is a 'remote control' program developed to control QA-45 from a PC. The following section describes *ansur* QA-40/45, and how it works.
ansur QA-40/45 is supplied with its own user manual, which contains all information about the product.
This manual is also supplied when ordering the demonstration version of the program. *ansur* QA-40/45 is run as a Microsoft Windows application.

ansur QA-40/45 offers the following functions:

- Manual energy test
- Automatic energy test
- ECG simulation
- ECG performance test
- Manual pacer, sensitivity and refractory period test
- Automatic pacer test sequence
- Test protocols for each individual instrument

GENERAL INFORMATION

TEMPERATURE RANGE:	+15°C to +35°C operating; 0°C to +50°C storage.

DISPLAY:

Type:	LCD graphic display.
Alphanumeric format:	8 lines by 40 characters.

DATA INPUT/OUTPUTS (2):
Parallel printer port (1); bi-directional
RS-232C (1) for Computer control.

POWER:	2 x 9 Volt Alkaline Battery or 240 VAC Battery Eliminator. 115 VAC Batt. el. for US.
HOUSING:	High impact plastic case.
DIMENSIONS:	D x W x H: 28.0 cm x 24.8 cm x 9.8 cm
WEIGHT:	2.06 kg with battery.

STANDARD ACCESSORIES:
User/Service manual
Protection Cover
Internal Paddel Contact Adapter
Battery Eliminator
9 Volt Alkaline Battery

RECOMMENDED PRINTERS:
HP Desk Jet, Canon Bubble Jet or compatible

QA-45 MKIII ORDERING INFORMATION

Order no:

13040: QA-45 MKIII Defibrillator/ Transcutaneous Pacemaker Analyzer

Accessories:

17021: Battery Eliminator 240 VAC
17027: Battery Eliminator 115 VAC
13422: Carrying Case
10500: Carrying Case, ext. printer
13402 Protection cover
13403: Internal paddle contact adapter
13410: Defib. Paddle Adapter, Laerdal
13411: Defib. Paddle Adapter, Physiocontrol
13412: Defib. Paddle Adapter, Marquette
13415: Pace Adapter, HP/Agilent
13417: FR2 Paddle Contact Adapter, HP/Agilent
13428: Pace Only Test Cable, Physiocontrol
17024: Universal Banana Adapter
13600: *ansur* QA-40/45 Plug-In
13601: *ansur* QA-40/45 Plug-In, demo
13605: User manual *ansur* QA-40/45 Plug-In
13060: User/Service manual QA-45

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