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BIBLIOGRAFIA

1. Claire M. Fraser. *Complete genome sequence of Treponema Pallidum, the syphilis spirochete*. Science (1998); 281 July: 375-381.
2. Center for Disease Control. *Recommendations for diagnosing and treating syphilis in HIV-infected patients*. MMWR Morb. Mortal Wkly Rep. (1988); 37: 601
3. Aral R. Marx. *Crack, sex and STD,* Sexually Transmitted Diseases, 1991; 18:92-101
4. J.N. Wasserheit. *Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases*, Sexually Transmitted Diseases 1992; 19:61-77
5. Johnson Phillip C.*Testing for Syphilis*, Dermatologic Clinic 1994; 12 Jan: 9-17

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Abon Biopharm (Hangzhou) Co., Ltd.
#198 12th Street East,
Hangzhou Economic & Technological Development Area,
Hangzhou, 310018, P.R.China

EC REP
Wellkang Ltd
Suite B, 29 Harley Street
LONDON, W1G 9QR, U.K.

Número: 1155994401
Fecha efectiva: 2010-12-21

Abbott

Bioline™

RSV

RSV (Respiratory Syncytial Virus) Antigen Test
Test antigène du VRS(virus respiratoire syncytial)
Prueba antígeno de VSR (virus sincitial respiratorio)
Teste de antígeno de RSV(virus sincicial respiratório)
RSV-Antigen-Test (RSV = Respiratorisches Synzytial-Virus)
Test per la rilevazione dell'antigene del RSV (virus respiratorio sinciziale)

REF 40FK12

ENGLISH

About the test

[Introduction] RSV is the most common respiratory virus in infants and young children and most often infects infants under the age of two. In most infants, the virus causes symptoms resembling those of the common cold, including fever and runny nose, but can lead to more serious symptoms such as coughing, difficulty breathing, and wheezing. In infants born prematurely and/or with chronic lung disease, RSV can cause a severe or even lifethreatening disease. Prior to the introduction of Synagis, a monoclonal antibody used to prevent serious lower respiratory tract disease caused by RSV, the virus resulted in over 125,000 infantile hospitalizations each year. There was a high mortality rate among infants who were born premature. Specially selected mouse monoclonal anti-RSV antibodies are used as detector materials in Bioline™ RSV test. These enable the test to identify RSV in NPA with a high degree of accuracy.
[Test principle] Bioline™ RSV test is immobilized with mouse monoclonal anti-RSV antibody as capture material. And the specially selected antibodies are used as detector materials. These enable Bioline™ RSV test to identify RSV antigen in nasopharyngeal aspiration (NPA), with a high degree of accuracy. Bioline™ RSV test involves the extraction of RSV viral antigen by disrupting and exposing the internal viral nucleoproteins of nasopharyngeal aspirate (NPA) and the insertion of the test strip into the extraction reagent tube. Bioline™ RSV test strip has 2 pre-scored lines, "T" (Test Line) and "C" (Control Line) on the surface of the strip. Two lines in result window are not visible before applying any specimens. The Control Line is used for procedural control. Control line should always appear if the test procedure is performed properly and the test reagents of control line are working. The T line, on the result window, is the result window if there is enough RSV antigen in the specimen. If RSV antigen are not present in the specimen, there is no visible appears in "T".
[Intended use] Bioline™ RSV rapid test is a chromatographic immunoassay for the qualitative detection of Respiratory Syncytial Virus (RSV) in nasopharyngeal aspirate (NPA) samples. This kit is intended for professional use and for in vitro diagnostic use only.

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1. Bioline™ RSV test kit contains the following items to perform the assay:
 - 25 Test strips with desiccant in individual foil pouches
 - Extraction buffer (1 x 10 ml/vial)
 - 25 Disposable test tubes
 - 50 Disposable droppers
 - 1 Instructions for use
2. Active ingredients of main components
 - 1 test strip includes:
 - Gold conjugate: Mouse monoclonal anti-RSV gold-colloid (0.17±0.03 µg)
 - Test line: Mouse monoclonal anti-RSV antibody (0.40±0.08 µg)
 - Control line: Anti-mouse IgG (0.80±0.16 µg)
 - Extraction buffer includes: Tricine (0.4 M), NaCl (g.s.), Triton™ X-100 (q.s.), Sodium azide (0.02 %)

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- Protective gloves, Timer, Biohazard container

Kit storage and stability

1. The test kit should be stored between 1 °C and 30 °C. Do not freeze the kit or its components.
2. The test strip is sensitive to both heat and humidity. Perform the test immediately after removing the test strip from the foil pouch.
3. Do not use the kit or its components beyond the expiration date.
4. The shelf life of the kit is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.

Warning

1. For in vitro diagnostic use only. Do not reuse the test strip.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
3. Do not eat or smoke while handling specimens.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Do not touch inner membrane of the result window by hand. It can affect the test result.
6. Avoid splashing or aerosol formation of specimen and extraction buffer.
7. Clean up spills thoroughly using an appropriate disinfectant.
8. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials in a biohazard container as if they were infectious waste.
9. Do not mix or interchange different specimens.
10. Components must not be used after the expiration date printed on the labels. Do not mix or interchange different batches/lots of reagents.
11. Do not drink extraction buffer.
12. Bioline™ RSV extraction buffer contains a proprietary anti-microbial agent, sodium azide, which presents no hazard to the user if normal laboratory safety precautions are followed. If contact with extraction buffer to the eyes and/or skin, wash affected area with soap and water immediately. Also irritation or signs of toxicity occur, seek medical attention.

Specimen collection and handling

1. Specimen collection
 - RSV is recovered almost exclusively from the respiratory tract.
 - Use mucus trap and catheter to collect NPA specimen.
 - Nasal swab is not recommended as a specimen collection method due to great variation of RSV particles in swab. When using a nasal swab, users cannot expect the same test accuracy that can be achieved with an NPA specimen.
 - It is strongly recommended to avoid the use of sputum.
2. Specimen transport and storage
 - The NPA specimens can be tested as soon as possible after collection. If necessary, specimens may be stored at 2 - 8 °C for up to 48 hours or at -20 °C for 7 days.
 - Do not treat specimens with solutions containing formaldehyde or its derivatives.

Test procedure (refer to figure)

1. Bring the test strip and specimen to a temperature between 15-30 °C prior to testing.
2. Using a disposable dropper provided, draw the extraction buffer up to the fill line (200 - 250 µl) and add the drawn buffer into the empty test tube provided.
3. Using a different disposable dropper provided, draw the specimen up to the fill line (200 - 250 µl) and mix the drawn specimen with the test tube containing the extraction buffer to achieve a dilution ratio of 1/2 (1v:1v).
4. Mix well to homogenize the solution.
5. Let the test tube rest for 10 minutes at a temperature between 15 - 30 °C.
6. With the arrows on the test strip pointing down, insert the test strip into the test tube.
7. Do not handle or move the test strip or test tube until it is time to read the test result.
8. Read the test result at 15 minutes; do not read the test result after 15 minutes.

Test interpretation (refer to figure)

[Negative result:] Only the control line (C) appears.
[Positive result: Both the control line (C) and the test line (T) appear.

⚠ **Caution:** The presence of any test line, no matter how faint, the result is considered positive.

[Invalid result: Control line (C) fails to appear.
If the control line is not visible within the result window or only the test line is visible after performing the test, the result is considered invalid. The shortage of flow is the most likely reason for the invalid result and test may have deteriorated. It is recommended that the specimen be retested using a new test kit.

Test limitations

1. The contents of this kit are to be used for the qualitative detection of RSV from NPA.
2. Failure to follow the test procedure and instructions for interpretation of test results may adversely affect test performance and/or invalidate the test result.
3. Test results must be evaluated in conjunction with other clinical data available to the physician.
4. NPA specimens that are collected after the acute phase of infection may contain antigen titers below the reagent's sensitivity threshold. If a negative result is obtained and clinical symptoms, a culture should be started to check the specimen for the presence of RSV.

Internal quality control

Bioline™ RSV test has a test line (T) and a control line (C). Neither the test line nor the control line is visible in the result window before applying any specimens. The control line is used for procedural control only. The presence of the control line indicates that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but it does not guarantee that the specimen has been properly applied and does not rule out errors.

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TECHNICAL SUPPORT CONTACT INFORMATION

For further information, please contact your distributor or contact our support specialists:

Region	Phone	E-mail address
Europe & Middle East	+44 161 483 9032	EME.TechSupport@abbott.com
Asia Pacific	+61 7 3363 7100	APTechSupport@abbott.com
Africa	+27 10 500 9700	arcis.techsupport@abbott.com
Russia & CIS	+44 161 483 9032	arcis.techsupport@abbott.com
Latin America	+571 482 4033	LATechSupport@abbott.com
Canada	+1 800 818 8335	Canproductsupport@abbott.com
US	+1 877 441 7440	USproductsupport@abbott.com

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Matériels required but not provided

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DO NOT PRINT FROM THIS FILE - FOR REFERENCE ONLY

TRIM
FOLD
10mm clear trim to text
Gutter width is optional
Non Printing Area Keep Gutter Clear for DATA MATRIX CODE 5mm Min No Text or Graphics in RED Area.
DOCUMENT SIZE: 490 x 329 mm.

FreeStyle

Optium H β-Ketone

Blood β-Ketone Test Strips | Tiras reactivas de cuerpos cetónicos (β-HB) en sangre | Électrodes de dosage de la cétonémie | Strisce per il test della β-chetonemia | Tiras-teste para β-cetona no sangue

English

Intended Use

IVD Individually foil wrapped FreeStyle Optium H Blood β-Ketone Test Strips are for *in vitro* diagnostic use by healthcare professionals in a multi-patient setting with the FreeStyle Optium Neo H, FreeStyle Optium H, Optium, and Optium Xceed Meters (not all products are available in all countries). These systems may be used for the quantitative measurement of β-Ketone (beta-hydroxybutyrate) in fresh capillary (i.e., from the finger) and venous whole blood samples. These systems are not for use in diagnosis of diabetes mellitus, but are to be used as an aid in monitoring the effectiveness of diabetes management programmes.

What's in my test strip box?

- Test strips individually wrapped in purple foil packets
- Instructions for use
- Calibrator

Storage and Use

FreeStyle Optium H Blood β-Ketone Test Strips are for single use only. The lot number and expiry date are printed on the back of the foil packet.


- Store the test strips at temperatures between 4° and 30°C (39° and 86°F). Keep away from direct sunlight and heat. Storage outside this range may cause erroneous results.
- Do not use the test strips beyond the expiry date printed on the foil packet and outer box. If only the year and month are shown, then the expiry date is the last day of the month. For example, "EXP 2014/03" means the test strip expires on March 31, 2014.
- Use the test strip immediately after opening the foil packet.
- Do not use wet, bent, scratched or damaged test strips.
- Use the test strip only once and then discard.
- Do not use the test strip if the foil packet has a puncture or tear.
- Use FreeStyle Optium H Blood β-Ketone Test Strips only with the FreeStyle Optium Neo H, FreeStyle Optium H, Optium and Optium Xceed Meters.

Calibration

The calibration procedure programs the meter with the lot number and test strip technology. This procedure requires the calibrator supplied with this test strip package and a FreeStyle Optium Neo H, FreeStyle Optium H, Optium or Optium Xceed Meter, available separately. Please refer to the meter's User's Operator's Manual for detailed calibration procedures.

IMPORTANT: Always calibrate the meter with each new box of test strips. Failure to calibrate properly will cause incorrect results.

1. With the lot number facing toward you, insert the contact bars of the purple calibrator into the meter. The meter turns on automatically.

2. The lot number will appear in the display window with the following symbols:
 - Optium Xceed and FreeStyle Optium H displays **CODE**.
 - FreeStyle Optium Neo H displays **KET** and **Lot**.
 - Optium displays **Ketone** first and then **LOT**.
3. Check that the lot number on the meter display window matches with the number on the test strip calibrator and the first five digits on both the test strip foil packet and test strip instructions for use. Calibration is complete.

Note: Use only the calibrator supplied with the test strips. Keep the calibrator until all the test strips in the box have been used. When the box is empty, throw the calibrator away.

Sample Collection and Preparation


FreeStyle Optium H Blood β-Ketone Test Strips are designed for use with fresh whole blood samples. The minimum sample volume is 1.5 µL.

Before obtaining a fingertip blood sample, clean and dry the sample site completely. If possible hang the arm down before lancing the finger to increase the blood flow. Use the sample immediately.

Venous whole blood samples may be collected into sodium or lithium heparin tubes or EDTA tubes and used within 30 minutes. Do not use tubes containing fluoride or oxalate.

Testing Blood β-Ketone

This procedure requires a FreeStyle Optium Neo H, FreeStyle Optium H, Optium or Optium Xceed Meter, available separately. Please refer to the meter User's/Operator's Manual for detailed testing procedures.

1. Clean and dry the sample site completely.
2. Open the test strip packet by tearing at the notch.
3. Insert the contact bars into the test port.
4. Push the test strip into the test port until it stops. Check that the lot number which appears in the display window matches the first 5 digits of the lot number on the test strip foil packet.
5. Lance the finger to obtain a drop of blood.
6. Touch blood drop to white target area at tip of test strip while the Apply Sample prompt appears in the display window.
7. Hold finger in place until test begins.
8. Test will start when sample is detected.

IMPORTANT: If the test fails to start, a second drop may be applied within 30 seconds of the first blood drop. If the test fails to start again or if more than 30 seconds have passed, repeat with new strip.

9. There is a 10-second countdown before the blood β-Ketone result is displayed.
10. Dispose of the test strip correctly. The opened foil packet may be used to remove and discard the used test strip.

Refer to the User's/Operator's Manual for detailed testing procedures.

Meter Test Messages

The following messages may mean there is a blood β-Ketone result that requires immediate attention or there may be a problem with the test strip:

- **HI** means your meter has determined that your blood β-Ketone result is higher than 8.0 mmol/L.
- **Test Error 2** or **Test Error 4** (Optium Meter) or **E-3** or **E-4** (FreeStyle Optium Neo H, FreeStyle Optium H and Optium Xceed Meters) means there may be a test error.

If any of these messages show, repeat the test with a new test strip. You may also perform a control solution test to check the performance of your system. Consult the prescribing physician before making any changes to diabetes medication programme plans.

Expected Results for Non-Diabetic, Non-Pregnant Adults

The blood β-Ketone test measures β-hydroxybutyrate (β-OHB), the most important of the three ketone bodies in the blood[1]. Normally, levels of β-OHB are expected to be less than 0.6 mmol/L[1]. β-OHB levels may increase if a person fasts, exercises vigorously or has diabetes and becomes ill[1].

If the blood β-Ketone result is 0.0 mmol/L and the blood glucose result is 300 mg/dL (16.7 mmol/L) or higher, or the results are not consistent with physical symptoms, repeat both the ketone and the glucose tests with new test strips.

If the blood β-Ketone result is between 0.6 and 1.5 mmol/L and the blood glucose result is 300 mg/dL (16.7 mmol/L) or higher, this may indicate development of a problem that could require medical intervention.

If the blood β-Ketone result is higher than 1.5 mmol/L and the blood glucose result is 300 mg/dL (16.7 mmol/L) or higher, this may indicate the risk of developing diabetic ketoacidosis (DKA)[2-4].

Limitations of Procedure

Note: This test strip has not been evaluated for arterial and neonatal samples, or for alternative site testing.

- This test strip is not designed for use with serum or plasma samples.
- Use between 18° and 30°C (64° and 86°F) and 10% and 90% relative humidity for best results.
- Testing demonstrates that altitudes up to 2,195 metres (7,200 feet) above sea level do not affect results.
- Haematocrit range is 30%–60%.
- Test Results may be erroneously low if the patient is severely dehydrated, or severely hypotensive, in shock or in a hyperglycaemic-hyperosmolar state.

Quality Control

Use MediSense Glucose and Ketone Control Solutions to do a quality control check to verify the performance of the meter and test strips with each new test strip lot number, or when you question the results. In a multi-user setting, daily quality control testing is recommended to ensure the system is operating properly. Control results must be within the "Expected Results with Control Solutions" shown below on these instructions for use. Please refer to the meter User's/Operator's Manual for detailed instructions. Please contact your local Abbott Diabetes Care office or distributor to obtain control solutions.

Test Principle

When the blood sample is applied to the test strip, the β-hydroxybutyrate (β-OHB) in the blood reacts with the chemicals on the test strip, producing a small electrical current. This current is measured and then a result is displayed by the meter. The size of the current depends on the amount of β-OHB in the blood sample.

Reagents

The reagent area of each test strip contains:

β-Hydroxybutyrate Dehydrogenase (*Pseudomonas sp*)	≥ 0.03 U
NAD (Free acid form)	≥ 1.67 µg
Phenanthroline quinone	≥ 0.29 µg
Non-reactive ingredients	≥ 19.51 µg

Performance Characteristics

Performance of the test strip has been evaluated in laboratory and clinical studies.

Assay Range:	0.0–8.0 mmol/L
Sample Volume:	1.5 µL
Test Time:	10 seconds

Calibration Reference

The FreeStyle Optium H Blood β-Ketone Test Strip is calibrated to reflect plasma β-hydroxybutyrate (β-OHB) using the Randox assay kit (RB1007).

Precision

Precision testing shows that results typically vary by no more than 3.1% to 3.8%. Results were obtained in a laboratory study using venous whole blood samples (n = 20 per level). Please refer to Table 1 at the end of these instructions for use for study results.

Specificity

Systems exhibit no interference from the following substances above therapeutic levels: Acetaminophen, ascorbic acid, dopamine, ephedrine, ibuprofen, L-dopa, maltose, methyldopa, salicylate, tetracycline, tolazamide and tolbutamide.

Accuracy

Accuracy testing shows that results are comparable between trained operators and lay users. Accuracy was assessed at four hospitals by comparing whole blood β-OHB results with plasma results obtained using a reference laboratory instrument. Please refer to Table 2 at the end of these instructions for use for study results.

These studies show that blood β-Ketone test strip results compare well with laboratory method results.

Table 1 - Precision

	Low	Mid	High
Mean, mmol/L	0.34	2.36	6.32
SD, mmol/L	0.03	0.09	0.2
CV, %	—	3.8	3.1

Table 2 - Accuracy

No. of tests	288
Ketone Range, mmol/L	0.07–5.2
Slope	1.06
Intercept, mmol/L	0.07
r (corr. coef.)	0.98

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GRAV. BAR MUST BE 282 mm FROM TOP TRIM. MUST DO NOT MOVE GRAV. BAR BELOW.
shaded area for LOT, LO, MID, HI and overprint content.

LOT
Expected Results with Control Solutions/ Resultados esperados con soluciones de control/ Résultats attendus avec les solutions de contrôle/ Risultati attesi con le soluzioni di controllo/ Resultados esperados com soluções de controle
LO Low/ Bajo/ Bas/ Basso/ Baixo:
MID Mid/ Medio/ Moyen/ Medio/ Médio:
HI High/ Alto/ Haut/ Alto/ Alto:

A window of 120 mm Min clearance for OVERPRINTING is required.
This measurement begins after the colon following the lo, mid,hi text.

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Abbott Laboratories Ireland, Ltd., Abbott Diabetes Care, 4051 Kingswood Drive, Citywest Business Campus, Dublin 24, Ireland, 1800 77 66 33

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Abbott Laboratories Ltd., Abbott Diabetes Care, Abbott House, Vanwall Business Park, Vanwall Road, Maidenhead, Berkshire, SL6 4UD; UK, 0500 467 466

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References/ Referencias/ References/ Riferimenti/ Referências
1. Schade DS, Eaton RP. Metabolic and clinical significance of ketosis. Special Topics in Endocrinology and Metabolism 1982;4:1-27.
2. Wiggam MI, O'Kane MJ, Harper R, Atkinson AB, Hadden DR, Trimble ER, Bell PM. Treatment of diabetic ketoacidosis using normalization of blood 3-hydroxybutyrate concentration as the endpoint of emergency management. Diabetes Care 1997; 20:1347-52.
3. Harano Y, Kosugi K, Hyosu T, Suzuki M, Hidaka H, Kashiwagi A, Uno S, Shigeta Y. Ketone bodies as markers for Type 1 (insulin-dependent) diabetes and their value in the monitoring of diabetes control. Diabetologia 1984;26:343-8.
4. Ubukata E. Diurnal variation of blood β-ketone bodies in insulin-dependent diabetes mellitus and noninsulindependent diabetes mellitus patients: The relationship to serum c-peptide immunoreactivity and free insulin. Ann Nutr Metab 1990;34:333-42.
5. Luzi L, Barrett EJ, Groop LC, Ferrannini E, DeFronzo RA. Metabolic effects of low-dose insulin therapy on glucose metabolism in diabetic ketoacidosis. Diabetes 1988;37:1470-77.
6. Hale PJ, Crase J, Nattrass M. Metabolic effects of bicarbonate in the treatment of diabetic ketoacidosis. Br Med J 1984;289:1035-8.

Product of U.K.
This product(s) and/or its manufacture and/or use are protected by one or more of the following patents / Este producto, su fabricación y uso están protegidos por una o más de las siguientes patentes / Ce ou ces produits et ou/leur fabrication et/ou utilisation sont protégés par un ou plusieurs des brevets suivants / Questo(i) prodotto(i) e/o il relativo produttore e/o il relativo uso sono coperti da uno o più dei seguenti brevetti / Este(s) produto(s) e/ou sua fabricação e/ou uso são protegidos por uma ou mais das seguintes patentes: US5,628,890; US5,682,884; US5,727,548; US6,129,823; US6,377,894; US6,591,125; US6,600,997; US6,736,957; US6,764,581; US6,773,671; US6,939,450; US7,058,437; US7,504,019; US7,740,581; US7,905,999; US7,922,883; US7,998,337; US8,118,993; US8,182,671; US8,211,280; US8,221,612; US8,241,485; US8,241,486; US8,372,261; US8,414,485; US8,414,760; US8,425,743; US8,449,758; US8,532,731; EP1000850; EP1023455; EP1135679; EP1119637; EP1801229; CA2302448; CA2369800; CA2364415; CA2517196; CA2153670. Additional patents may be issued and/or pending. / Otras patentes pueden ser expedidas o estar en trámite. / Des brevets supplémentaires peuvent être déposés et/ou en instance. / Eventuali altri brevetti già concessi o in attesa di concessione. / Outras patentes podem ser emitidas e/ou estarem pendentes.
FreeStyle and related brand marks are trademarks of Abbott Diabetes Care Inc. in various jurisdictions.

©2014 Abbott
ART30778-005 Rev. A 06/14

Abbott Diabetes Care Ltd.
Range Road
Witney, Oxon
OX29 0YL, UK

ART30778-005_rev-A.indd 2
6/27/14 9:45 AM

DO NOT PRINT FROM THIS FILE -
FOR REFERENCE ONLY

English

IVD For external use (***in vitro*** diagnostic use)

Intended Use

MediSense Glucose and Ketone Control Solutions are designed for consumer home testing and healthcare professional use with all FreeStyle Precision, FreeStyle Precision H, Precision Xceed Pro, FreeStyle Optium, FreeStyle Optium H and all other Precision and Optium Blood Glucose and β-Ketone Monitoring Systems, except ExacTech, Optium Mini and Optium Omega. The controls are used to check the performance of your system.

When to do a control solution test:

When you question your results and want to confirm that your monitor and test strips are working properly.

For more information on how to obtain Control Solutions, contact your local Abbott Diabetes Care office or distributor.

Healthcare Professionals Please Note:

Daily control testing is recommended in a multi-user setting to ensure that the systems are operating properly.

Storage and Handling

39°F 4°C – 86°F 30°C Store the control solutions at temperatures between 39° and 86°F (4° and 30°C). They are stable for 90 days after opening.

DATE When you open a new bottle, write the date of opening on the bottle label. Always replace cap tightly.

Precautions and Warnings

- Do not use control solutions 90 days after opening or if they are expired (check the expiration date printed on each bottle and on the box they were packaged in).
- Please refer to your test strip package insert for acceptable temperature and humidity test conditions.
- If the control solutions have been refrigerated, allow them to reach the system operating temperature **before performing** control testing.
- Invert the control bottle three to four times to ensure thorough mixing before use.
- **Do not swallow or inject the solution.**
- Do not use with ExacTech, Optium Mini or Optium Omega Systems.
- Do not re-use test strips.

Procedure

Please refer to the appropriate manual for control testing instructions.

Materials Provided

Low Control Solution (white cap) Volume 4 mL
High Control Solution (blue cap) Volume 4 mL

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Monitor, Test Strips and User's or Operator's Manual.

Expected Results

Compare the control results to the "Expected Results For Control Solutions" printed on the test strip package insert sheet found in the carton of test strips or on the test strip barcode label (if applicable) being used. Check that the lot number printed on the test strip packet and insert matches. The results obtained from control testing do not reflect your personal blood glucose or ketone level.

Quality Control

If the control results are outside the expected ranges, repeat the test. Check that the instructions are followed correctly. If the results are outside the ranges again, contact your local Abbott Diabetes Care office or distributor.

Reagent Composition

LO Low Control Solution

0.04%	w/v	Glucose
0.02%	w/v	DL-Beta-hydroxybutyrate
99.94%	w/v	Other ingredients

HI High Control Solution

0.22%	w/v	Glucose
0.08%	w/v	DL-Beta-hydroxybutyrate
99.70%	w/v	Other ingredients

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QMS-DWG-069

MediSense® Glucose & Ketone Control Solutions

MediSense® Glucose & Ketone Control Solutions

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QMS-DWG-069

Manufacturer / Fabricante:
Abbott Diabetes Care Ltd.
Range Road
Witney, Oxon
OX29 0YL, UK
Product of USA
© 2010 Abbott
FreeStyle and related brand marks are trademarks of Abbott Diabetes Care, Inc. in various jurisdictions.
ART19250-101 Rev. A 12/13

Abbott

MICRO ESSENTIAL
LABORATORY
Hydrion® pH and sanitizer test kits since 1934

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TECHLAB® *C. DIFF QUIK CHEK COMPLETE®*

INTENDED USE

The TECHLAB® *C. DIFF QUIK CHEK COMPLETE®* test is a rapid membrane enzyme immunoassay for the simultaneous detection of *Clostridium difficile* glutamate dehydrogenase antigen and toxins A and B in a single reaction well. The test detects *C. difficile* antigen, glutamate dehydrogenase, as a screen for the presence of *C. difficile* and confirms the presence of toxigenic *C. difficile* by detecting toxins A and B in fecal specimens from persons suspected of having *C. difficile* disease. The test is to be used as an aid in the diagnosis of *C. difficile* disease. As with other *C. difficile* tests, results should be considered in conjunction with the patient history.

EXPLANATION

After treatment with antibiotics, many patients develop gastrointestinal problems ranging from mild diarrhea to severe pseudomembranous colitis. Many cases of the milder forms of gastrointestinal illness and most cases of pseudomembranous colitis are caused by toxigenic strains of *Clostridium difficile* (1). This organism is an opportunistic anaerobic bacterium that grows in the intestine once the normal flora has been altered by the antibiotic. Toxigenic strains of *C. difficile* carry the genes encoding the toxins while non-toxigenic strains do not carry the toxin genes. Disease onset is associated with the toxins that are produced by the toxigenic organism. The clinical symptoms associated with the disease are believed to be primarily due to toxin A, which is a tissue-damaging enterotoxin (2,3). *C. difficile* also produces a second toxin, designated toxin B. Toxin B, which has been referred to as the cytotoxin of the organism, is the toxin detected by the tissue culture assay currently used by many laboratories. Toxigenic *C. difficile* strains produce both toxins, or only toxin B (4-7). The glutamate dehydrogenase of *C. difficile* is a good antigen marker for the organism in feces because it is produced in high amounts by all strains, toxigenic or non-toxigenic (8-10). The antigen can be detected in fecal specimens by using the *C. DIFF QUIK CHEK COMPLETE®* test. A positive result in the test for the glutamate dehydrogenase of *C. difficile* confirms the presence of this organism in a fecal specimen; a negative result indicates the absence of the organism. A positive result in the test for toxins A and B confirms the presence of toxigenic *C. difficile*.

PRINCIPLE OF THE TEST

The *C. DIFF QUIK CHEK COMPLETE®* test uses antibodies specific for glutamate dehydrogenase and toxins A and B of *C. difficile*. The device contains a *Reaction Window* with three vertical lines of immobilized antibodies. The antigen test line ("Ag") contains antibodies against *C. difficile* glutamate dehydrogenase. The control line ("C") is a dotted line that contains anti-horseradish peroxidase (HRP) antibodies. The toxins A and B test line ("Tox") contains antibodies against *C. difficile* toxins A and B. The *Conjugate* consists of antibodies to glutamate dehydrogenase and antibodies to toxins A and B coupled to horseradish peroxidase. To perform the test, the sample is added to a tube containing a mixture of *Diluent* and *Conjugate*. The diluted sample-conjugate mixture is added to the *Sample Well* and the device is allowed to incubate at room temperature for 15 minutes. During the incubation, any glutamate dehydrogenase and toxins A and B in the sample bind to the antibody-peroxidase conjugates. The antigen-antibody-conjugate complexes migrate through a filter pad to a membrane where they are captured by the immobilized glutamate dehydrogenase-specific and toxins A and B-specific antibodies in the lines. The *Reaction Window* is subsequently washed with *Wash Buffer*, followed by the addition of *Substrate*. After a 10 minute incubation period, the "Ag" reaction is examined visually for the appearance of a vertical blue line on the "Ag" side of the *Reaction Window*. A blue line indicates a positive test. If the "Ag" is positive, then the "Tox" reaction should be examined visually for the appearance of a blue line on the "Tox" side of the *Reaction Window*. A blue line indicates a positive test. A positive "C" reaction, indicated by a vertical dotted blue line under the "C" portion of the *Reaction Window*, confirms that the test is working properly and the results are valid.

MATERIALS PROVIDED

MEM | DEV — ***Membrane Devices*** – each pouch contains 1 device

DIL | SPE — ***Diluent*** **(22 mL per bottle)** – Buffered protein solution with graduated dropper assembly (contains 0.05% ProClin® 300)*

WASH | REAG — ***Wash Buffer*** **(12 mL per bottle)** – Buffered solution with graduated dropper assembly (contains 0.05% ProClin® 300)*

SUBS | REAG — ***Substrate*** **(3.5 mL per bottle)** – Solution containing tetramethylbenzidine

CONJ | ENZ — ***Conjugate*** **(2.5 mL per bottle)** – Antibody specific for glutamate dehydrogenase coupled to horseradish peroxidase and antibodies specific for toxins A and B coupled to horseradish peroxidase in a buffered protein solution (contains 0.05% ProClin® 300)*

CONTROL | + — ***Positive Control*** **(2 mL)** – Antigen in a buffered protein solution

Disposable plastic transfer pipettes – graduated at 25 µL, 100 µL, 200 µL, 300 µL, 400 µL and 500 µL

IVD — For *in vitro* diagnostic use

*Signal Word: Warning
H317: May cause an allergic skin reaction
P261, P272, P280, P302, P352, P333, P313, P321, P362, P364, P501

MATERIALS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

Small test tubes (e.g., plastic microcentrifuge tubes or glass tubes)
Applicator sticks — *Timer* — *Vortex mixer*
Disposable gloves for handling fecal samples — *Pipettor and tips*

SHELF LIFE AND STORAGE

The expiration date of the kit is given on the label. Expiration dates for each component are listed on the individual labels. The kit should be stored between 2°C and 8°C.

PRECAUTIONS

1. For *in vitro* diagnostic use. For professional use only.
2. Reagents from different kits should not be mixed or interchanged. Do not use a kit past the expiration date.
3. Each component in the kit should be inspected for any signs of leakage. Upon arrival, inspect the kit to ensure that components are not frozen or warm to the touch due to improper shipping conditions.
4. Bring all components to ROOM TEMPERATURE BEFORE USE!

For Instructional Use Only

EN

5. Caps, tips and dropper assemblies are color-coded; do NOT mix or interchange!
6. Do not freeze the reagents. The kit should be stored between 2°C and 8°C.
7. The pouch containing the *Membrane Device* should be at room temperature before opening. Keep the membrane devices dry before use.
8. Use fecal specimens within 72 hours of collection to obtain optimal results. Specimens that are frozen may lose activity due to freezing and thawing. If using frozen specimens, thaw at room temperature.
9. Hold reagent bottles vertically to dispense reagents to ensure consistent drop size and correct volume.
10. Specimens and membrane devices should be handled and disposed of as potential biohazards after use.
11. Membrane devices cannot be reused.
12. The test has been optimized for sensitivity and specificity. Alterations of the specified procedure and/or test conditions may affect the sensitivity and specificity of the test. Do not deviate from the specified procedure.
13. Be attentive to the total assay time when testing more than one fecal specimen. Add *Diluent* first, and then add the *Conjugate* to each tube of *Diluent*. Then add specimen to the tube of *Diluent/Conjugate*. Thoroughly mix all of the diluted specimens, and transfer to the *Membrane Device*. The 15-minute incubation step begins after the last diluted sample-conjugate mixture has been transferred to the final *Membrane Device*.
14. If the *Substrate* reagent changes to a dark blue/violet color call technical services for replacement.
15. Fecal specimens may contain potentially infectious agents and should be handled at "Biosafety Level 2" as recommended in the CDC/NIH Manual "Biosafety in Microbiological and Biomedical Laboratories."
16. Wear disposable gloves when doing the test.
17. The *Conjugate, Diluent,* and *Wash Buffer* reagents contain 0.05% ProClin® 300 as a preservative. Although the concentration is low, ProClin® 300 is known to be harmful (skin sensitization may occur). If skin sensitization/irritation or rash occurs, get medical advice/attention. Take off contaminated clothing and wash it before reuse. Handle reagents according to existing regulations for laboratory safety and good laboratory practice. Safety Data Sheets for this product are available upon request, contact technical support.
18. Follow your national, regional, and local ordinances accordingly for waste disposal regulations.

COLLECTION, HANDLING AND STORAGE OF FECAL SPECIMENS

Acceptable Sample Types
Fresh Fecal Specimens
Frozen Fecal Specimens
Specimens in transport media (Cary Blair, C&S)

Do Not Use
Fecal specimens in Formalin-based fixative (e.g. sodium acetate formalin, 10% formalin, merthiolate formalin)
Fecal specimens in alcohol-based fixative (e.g. polyvinyl alcohol)


Sample Storage Temperature	Acceptable length of storage	Comments
2°C – 8°C	72 hours	Ideal specimens are less than 24 hours old
Frozen ≤ -10°C	Longer than 72 hours	Thaw at room temperature. Freezing and thawing multiple times may result in loss of specimen activity due to toxin degradation.

1. Standard collection and handling procedures used in-house for fecal specimens are appropriate.
2. Fecal specimens should be collected in clean, leak-proof containers.
3. Storing fecal specimens in the *Diluent* is NOT recommended.
4. Do not allow the fecal specimens to remain in the *Diluent/Conjugate* mixture for >24 hours.

SPECIMEN PREPARATION

1. Bring all reagents and the required number of devices to room temperature before use. It is recommended to remove the reagents from the foam insert to reduce the time needed to warm to room temperature.
2. Set up and label one small test tube for each specimen, and optional external controls as necessary.
3. Using the black graduated dropper assembly, add 750 µL (2nd graduation from the tip) *Diluent* to each tube for fecal specimens. **For specimens in transport media such as Cary Blair or C&S, add 650 µL (1st graduation from the tip) of *Diluent* to the tube.**

Sample Type	Volume of *Diluent*
Fresh Fecal Specimens	750 µL (2nd graduation from tip)
Frozen Fecal Specimens (frozen undiluted)	750 µL (2nd graduation from tip)
Specimens in transport media (Cary Blair, C&S)	650 µL (1st graduation from the tip)
External Controls (positive and negative)	750 µL (2nd graduation from tip)

4. Add one drop of *Conjugate* (red capped bottle) to each tube.
5. Obtain one disposable plastic transfer pipette (supplied with the kit) for each sample – the pipettes have raised graduations at 25 µL, 100 µL, 200 µL, 300 µL, 400 µL and 500 µL.

Graduated Transfer Pipette:

6. **Mix all specimens thoroughly regardless of consistency- it is essential that the specimens be evenly suspended before transferring.**
 Liquid/Semi-solid specimens – pipette 25 µL of specimen with a transfer pipette and dispense into the *Diluent/Conjugate* mixture. Use the same transfer pipette to mix the diluted specimen.
 Formed/Solid specimens – Care must be taken to add the correct amount of formed feces to the sample mixture. Mix the specimen thoroughly using a wooden applicator stick and transfer a small portion (approximately 2 mm diameter, the equivalent of 25 µL) of the specimen into the Diluent/Conjugate mixture. Emulsify the specimen using the applicator stick.

For Instructional Use Only

EN


Fecal specimens in Cary Blair or C&S transport media - pipette 100 µL (2nd graduation from tip of transfer pipette) of sample into the *Diluent/Conjugate* mixture.

7. **Optional External Control Samples:**
 External Positive Control - add one drop of *Positive Control* (gray-capped bottle) into the *Diluent/Conjugate* mixture.
 External Negative Control - add 25 µL *Diluent* into the *Diluent/Conjugate* mixture.

NOTE: Transferring too little specimen, or failure to mix and completely suspend the specimen in the Diluent/Conjugate mixture, may result in a false-negative test result. The addition of too much fecal specimen may cause invalid results due to restricted sample flow.

TEST PROCEDURE

1. Obtain one *Membrane Device* per specimen, and one device per optional external positive or negative control as necessary. The foil bags containing the devices should be brought to room temperature before opening. Use the device immediately after opening. Label each device appropriately and orient it on a flat surface so the "C. DIFF COMPLETE" print is at the bottom of the device, and the small *Sample Well* is located in the top right corner of the device.

2. Close each tube of diluted specimen and mix thoroughly. Proper mixing can be achieved by vortexing or inverting the tube. Once a patient sample or *Positive Control* has been diluted in the *Diluent/Conjugate* mixture, it may be incubated at room temperature for any period of time up to 24 hours prior to addition to the *Membrane Device*.
3. Using a new transfer pipette, transfer 500 µL of the diluted sample-conjugate mixture into the ***Sample Well*** (smaller hole in the top right corner of the device) of a *Membrane Device*, making certain to expel the liquid sample onto the wicking pad inside of the *Membrane Device*. When loading the sample into the sample well, make sure that the tip of the transfer pipette is angled towards the *Reaction Window* (larger hole in the middle of the device).
4. Incubate the device at room temperature for 15 minutes – the sample will wick through the device and a wet area will spread across the *Reaction Window*.
 NOTE FOR SAMPLES THAT FAIL TO MIGRATE:
 Occasionally, a diluted sample fails to migrate properly and the Reaction Window does not fully wet. If the Reaction Window does not appear to be completely wet within 5 minutes of adding the sample to the Sample Well, then add 100 µL (4 drops) of Diluent to the Sample Well and wait an additional 5 minutes (for a total of 20 minutes).
5. After the incubation, add 300 µL of *Wash Buffer* to the ***Reaction Window*** using the graduated white dropper assembly. Allow the *Wash Buffer* to flow through the Reaction Window membrane and be absorbed completely.
6. Add 2 drops of *Substrate* (white-capped bottle) to the ***Reaction Window***. Read and record results visually after 10 minutes.

INTERPRETATION OF RESULTS

1. Interpretation of the test is most reliable when the device is read immediately at the end of the 10 minute reaction period. Read the device at a normal working distance in a well-lit area. View with a line of vision directly over the device.
2. Observe device for the appearance of blue dots in the middle of the *Reaction Window* representing the internal positive control. The appearance of any control dot(s) represents a valid internal control. The background may appear white to light blue in color. Observe device for the appearance of blue lines on the "Ag" and "Tox" sides of the *Reaction Window* representing the test lines. The lines may appear faint to dark in intensity.
3. **Positive Antigen ("Ag") Result:** A positive antigen result may be interpreted at any time between the addition of *Substrate* and the 10-minute read time. For a positive antigen result, the blue "Ag" line and the dotted blue control line below "C" are visible (Figure 1a). The lines may appear faint to dark in intensity. An obvious partial line is interpreted as a positive result. Do not interpret membrane discoloration as a positive result. A positive result indicates the presence of *C. difficile*.
4. **Positive Antigen and Toxin ("Tox") Result:** If the antigen result is positive (i.e., a blue "Ag" line and a dotted blue control below "C" are visible), proceed to the interpretation of the toxin result. A positive toxin result may be interpreted at any time between the addition of *Substrate* and the 10-minute read time. For a positive toxin result, a blue "Tox" line is visible (Figure 1b). The line may appear faint to dark in intensity. An obvious partial line is interpreted as a positive result. Do not interpret membrane discoloration as a positive result. A positive result indicates the presence of *C. difficile* toxin.
5. **Negative Result:** A test cannot be interpreted as negative or invalid until 10 minutes following the addition of *Substrate*. A single blue dotted line is visible in the middle of the *Reaction Window,* below the "C" and no test lines are visible on the "Ag" side or the "Tox" side of the *Reaction Window,* (Figure 1c). A negative result in the antigen portion indicates *C. difficile* antigen is either absent in the specimen or is below the detection limit of the test. A negative result in the toxin portion indicates C. difficile toxin is either absent in the specimen or is below the detection limit of the test.
6. **Invalid Result:** No lines are visible in the Reaction Window (Figure 1d). The test result is invalid if a blue dotted line is not present below the "C" at the completion of the reaction period (Figures 1e, 1f, 1g).
7. **Negative Antigen ("Ag"), Positive Toxin ("Tox"):** A low percentage of specimens may test negative for antigen but positive for toxin. These samples should be considered indeterminate and retested using a fresh specimen (Figure 1h). If sample retests negative for antigen but positive for toxin, report as positive toxin result.

For Instructional Use Only

FIGURE 1: *C. DIFF QUIK CHEK COMPLETE*® INTERPRETATION OF RESULTS


Figure 1a
Positive Antigen Result


Figure 1b
Positive Antigen and Toxin Result


Figure 1c
Negative Result


Figure 1d
Invalid Result


Figure 1e
Invalid Result


Figure 1f
Invalid Result


Figure 1g
Invalid Result


Figure 1h
See #7 for Interpretation

QUALITY CONTROL

Internal: A dotted blue line must be visible in the middle of the *Reaction Window*, below the "C" on every *Membrane Device* that is tested. The appearance of the blue control dots confirms that the sample and reagents were added correctly, that the reagents were active at the time of performing the assay, and that the sample migrated properly through the *Membrane Device*. It also confirms the reactivity of the other reagents associated with the assay. A uniform background in the result area is considered an internal negative control. If the test has been performed correctly and reagents are working properly, the background will be white to give a discernible result.

External: The reactivity of the *C. DIFF QUIK CHEK COMPLETE*® kit should be verified upon receipt using the *Positive Control* and negative control (*Diluent*). The *Positive Control* is supplied with the kit (gray-capped bottle). The *Positive Control* confirms the reactivity of the other reagents associated with the assay, and is not intended to ensure precision at the analytical assay cut-off. *Diluent* is used for the negative control. Additional tests can be performed with the controls to meet the requirements of local, state and/or federal regulations and/or accrediting organizations.

LIMITATIONS

1. The *C. DIFF QUIK CHEK COMPLETE*® test is used to detect C. difficile antigen and toxin(s) in fecal specimens. The test confirms the presence of toxin in feces and this information should be taken under consideration by the physician in light of the clinical history and physical examination of the patient. The *C. DIFF QUIK CHEK COMPLETE*® test will detect levels of toxin A at ≥0.63 ng/mL, toxin B at ≥0.16 ng/mL, and glutamate dehydrogenase at ≥0.8 ng/mL.
2. Fecal specimens are extremely complex. Optimal results with the *C. DIFF QUIK CHEK COMPLETE*® test are obtained with specimens that are less than 24 hours old. Most undiluted specimens can be stored between 2°C and 8°C for 72 hours before significant degradation of the toxin is noted. If specimens are not assayed within this time period, they may be frozen and thawed. However, repeated freezing and thawing may result in loss in the immunoreactivity of antigen and toxins A and B.
3. Some specimens may give weak reactions. This may be due to a number of factors such as the presence of low levels of antigen and/or toxin, the presence of binding substances, or inactivating enzymes in the feces. The lines may appear faint to dark in intensity. These specimens should be reported as positive if any blue line, even a partial line is observed. An obvious partial blue line is interpreted as a positive result.
4. Fecal specimens preserved in 10% Formalin, merthiolate formalin, sodium acetate formalin, or polyvinyl alcohol cannot be used.
5. The *C. DIFF QUIK CHEK COMPLETE*® test is qualitative. The intensity of the color should not be interpreted quantitatively.
6. Some isolates of *C. sordellii* may react in the *C. DIFF QUIK CHEK COMPLETE*® test due to the production of immunologically related toxins (1).
7. Colonization rates of up to 50% have been reported in infants. A high rate has also been reported in cystic fibrosis patients (1,3). Results may appear positive in these groups, but should be viewed in conjunction with the potential to be a colonized carrier.
8. The only non-*C. difficile* organism to react in the toxin portion of the *C. DIFF QUIK CHEK COMPLETE*® test was *Clostridium sordellii* VPI 9048. This strain produces toxins HT and LT, which are homologous to toxins A and B, respectively.
9. No data exists on the effects of colonic washes, barium enemas, laxatives, or bowel preparations on the performance of the *C. DIFF QUIK CHEK COMPLETE*® test. All of these procedures can result in extensive dilution or the presence of additives that may affect test performance.

EXPECTED VALUES

Clostridium difficile disease is primarily a nosocomial disease of elderly patients, and the frequency of the disease is dependent on factors such as patient population, type of institution and epidemiology. The reported incidence of *C. difficile* disease in patients with antibiotic-associated diarrhea may range from 5 to 20%, and hospitals may experience rates lower or higher than this range. It is important to consider any test results in conjunction with clinical symptoms because some healthy adults and large numbers of healthy infants (up to 50%) will be positive for *C. difficile* toxin. In addition, *C. difficile* carriage rates of 22% to 32% have been reported in cystic fibrosis patients (1,3). In the studies conducted for this device, using symptomatic patients, the incidence of toxins A and B was 12% and GDH was 18%. A positive result in the antigen portion of the *C. DIFF QUIK CHEK COMPLETE*® test confirms the presence of *C. difficile* in a fecal specimen; a negative result indicates the absence of the organism. A positive result in the toxin portion of the *C. DIFF QUIK CHEK COMPLETE*® confirms the presence of *C. difficile* toxin in a fecal specimen; a negative result indicates the absence of toxin or insufficient levels of toxin for detection.

For Instructional Use Only

EN

PERFORMANCE CHARACTERISTICS

Clinical evaluation of the antigen portion of the *C. DIFF QUIK CHEK COMPLETE*® test

The antigen portion of the *C. DIFF QUIK CHEK COMPLETE*® test was compared to bacterial culture. Specimens included in the evaluation were submitted to the clinical laboratories for routine testing. The bacterial culture test was performed according to in-house procedures. The results are shown in Table 1.

Table 1. Summary of clinical performance comparing *C. DIFF QUIK CHEK COMPLETE*® test to bacterial culture

n = 1126	Bacterial Culture Positive	Bacterial Culture Negative
C. DIFF QUIK CHEK COMPLETE® Antigen Line Positive	201	62
C. DIFF QUIK CHEK COMPLETE® Antigen Line Negative	21	842

		95% Confidence Limits
Sensitivity	90.5%	85.7 – 93.9
Specificity	93.1%	91.2 – 94.7
Predictive Positive Value	76.4%	70.7 – 81.3
Predictive Negative Value	97.6%	96.2 – 98.4
Correlation	92.6%	91.8 – 93.4

Discrepant samples were evaluated using current ELISA tests for *C. difficile* glutamate dehydrogenase.
Twenty-nine of the 62 false positive samples were positive by another GDH test, and were considered true positives.
Thirteen of the 21 false negative samples were negative by another GDH test, and were considered true negatives.

The antigen portion of the *C. DIFF QUIK CHEK COMPLETE*® test was compared to the tissue culture assay for the detection of C. difficile toxin. Specimens included in the evaluation were submitted to the clinical laboratories for routine testing. The results are shown in Table 2. The antigen portion of the *C. DIFF QUIK CHEK COMPLETE*® test detected 98,7% of the tissue culture-positive samples.

Table 2. Summary of clinical performance comparing *C. DIFF QUIK CHEK COMPLETE*® test to the tissue culture assay

n = 1126	Tissue Culture Positive	Tissue Culture Negative
C. DIFF QUIK CHEK COMPLETE® Antigen Line Positive	154	109
C. DIFF QUIK CHEK COMPLETE® Antigen Line Negative	2	861

Clinical evaluation of the toxin portion of the *C. DIFF QUIK CHEK COMPLETE*® test

The toxin portion of the *C. DIFF QUIK CHEK COMPLETE*® test was compared to the tissue culture assay at two clinical laboratories and in-house at TECHLAB, Inc. Specimens included in the evaluation were submitted to the clinical laboratories for routine testing. The results are shown in Table 3.

Table 3. Summary of clinical performance comparing *C. DIFF QUIK CHEK COMPLETE*® test to tissue culture assay

n = 1126	Tissue Culture Positive	Tissue Culture Negative
C. DIFF QUIK CHEK COMPLETE® Toxin Line Positive	137	6
C. DIFF QUIK CHEK COMPLETE® Toxin Line Negative	19	964

		95% Confidence Limits
Sensitivity	87.8%	81.4 - 92.3
Specificity	99.4%	98.6 - 99.7
Predictive Positive Value	95.8%	90.7 - 98.3
Predictive Negative Value	98.1%	96.9 - 98.8
Correlation	97.8%	97.6 - 98.0

Discrepant samples were evaluated using current ELISA tests for toxins A and B.
Five of the 6 false positive samples were positive by ELISA and were considered true positives.
Twelve of the 19 false negative samples were negative by ELISA and were considered true negatives.

For Instructional Use Only

EFFECT OF FECAL SPECIMEN CONSISTENCY

Effect of fecal specimen consistency on the *C. DIFF QUIK CHEK COMPLETE*® test

The reaction of fecal specimens of varying consistencies in the antigen portion (n=978) and toxin portion (n=981) of the *C. DIFF QUIK CHEK COMPLETE*® test is shown in Tables 4 and 5. The percentages of positive reactions using either culture assay or the *C. DIFF QUIK CHEK COMPLETE*® test were similar in all three types of fecal specimens (liquid, semi-solid, and solid). All of the specimens were submitted for C. difficile testing. The basis of the submission was the clinical history of the patient and not the consistency of the specimen. In the antigen portion, the results show that the *C. DIFF QUIK CHEK COMPLETE*® test performed similarly to bacterial culture when testing samples of different consistencies. In the toxin portion, the results show the *C. DIFF QUIK CHEK COMPLETE*® test performed similarly to the tissue culture assay when testing samples of different consistencies.

Table 4. Reaction of fecal specimens of varying consistencies in the antigen portion of the *C. DIFF QUIK CHEK COMPLETE*® test

Number of specimens (n = 978)	Liquid Specimens (n = 335)	Semi-solid Specimens (n = 522)	Solid Specimens (n = 121)
Positive by bacterial culture assay	59 (17.6%)	110 (21.1%)	19 (15.7%)
C. DIFF QUIK CHEK COMPLETE® Antigen Line Positive	72 (21.5%)	128 (24.5%)	25 (20.7%)
Negative by bacterial culture assay	276 (82.4%)	412 (78.9%)	102 (84.3%)
C. DIFF QUIK CHEK COMPLETE® Antigen Line Negative	263 (78.5%)	394 (75.5%)	96 (79.3%)

Table 5. Reaction of fecal specimens of varying consistencies in the toxin portion of the *C. DIFF QUIK CHEK COMPLETE*® test

Number of specimens (n = 981)	Liquid Specimens (n = 336)	Semi-solid Specimens (n = 523)	Solid Specimens (n = 122)
Positive by tissue culture assay	43 (12.8%)	81 (15.5%)	8 (6.6%)
C. DIFF QUIK CHEK COMPLETE® Toxin Line Positive	42 (12.5%)	72 (13.8%)	7 (5.7%)
Negative by tissue culture assay	293 (87.2%)	442 (84.5%)	114 (93.4%)
C. DIFF QUIK CHEK COMPLETE® Toxin Line Negative	294 (87.5%)	451 (86.2%)	115 (94.3%)

ANALYTICAL SENSITIVITY

The cutoff for the assay was established at concentrations of 0.63 ng/mL for toxin A, 0.16 ng/mL for toxin B, and 0.8 ng/mL for glutamate dehydrogenase.

REPRODUCIBILITY

The reproducibility of the *C. DIFF QUIK CHEK COMPLETE*® test was determined using 12 fecal specimens that were coded to prevent their identification during testing. Testing was performed at 3 independent laboratories, which tested the samples for 3 days. The samples produced the expected results 100% of the time.

CROSS REACTIVITY

Fecal specimens inoculated with the following microorganisms to a final concentration of approximately 10^8 or higher organisms per mL did not react in the antigen or toxin portion of the *C. DIFF QUIK CHEK COMPLETE*® test:

Bacterium or Pathogen: *Aeromonas hydrophila, Bacillus cereus, Bacillus subtilis, Bacteroides fragilis, Campylobacter coli, Campylobacter fetus, Campylobacter jejuni, Candida albicans, Clostridium butyricum, Clostridium clostridiforme, Clostridium haemolyticum, Clostridium histolyticum, Clostridium novyi, Clostridium perfringens, Clostridium septicum, Clostridium sordellii* (nontoxigenic)*, Clostridium sporogenes, Enterobacter aerogenes, Enterobacter cloacae, Enterococcus faecalis, Escherichia coli* EIEC*, Escherichia coli, Escherichia coli* O157:H7*, Escherichia coli ETEC, Klebsiella pneumoniae, Peptostreptococcus anaerobius, Proteus vulgaris, Pseudomonas aeruginosa, Salmonella typhimurium, Serratia ligulfaciens, Shigella dysenteriae, Shigella flexneri, Shigella sonnei, Staphylococcus aureus, Staphylococcus aureus* (Cowans)*, Staphylococcus epidermidis, Vibrio cholerae, Vibrio parahaemolyticus, Yersinia enterocolitica*

The only non-*C. difficile* organism to react in the toxin portion of the *C. DIFF QUIK CHEK COMPLETE*® test was *Clostridium sordellii* VPI 9048. This strain produces toxins HT and LT, which are homologous to toxins A and B, respectively.

The following viruses of $10^{3.3}$ to $10^{8.25}$ TCID units per 0.2 mL did not react in the *C. DIFF QUIK CHEK COMPLETE*® test:

Viruses: Adenovirus types 1, 2, 3, 5, 40, 41, Human coronavirus, Coxsackievirus B2, B3, B4, B5, Echovirus 9, 11, 18, 22, 33, Enterovirus type 68, 69, 70, 71, Rotavirus.

INTERFERING SUBSTANCES

The following substances (U.S. formulation) had no effect on test results when present in feces in the concentrations indicated: mucin (3.5% w/v), human blood (40% v/v), barium sulfate (5% w/v), Imodium® (5% v/v), Kaopectate® (5% v/v), Pepto-Bismol® (5% v/v), steric/palmitic acid (40% w/v), Metronidazole (0.25% w/v), Vancomycin (0.25% w/v).

REACTION OF CLINICAL ISOLATES OBTAINED ON CYCLOSERINE-CEFOXITIN-FRUCTOSE AGAR (CCFA)

A total of 103 *C. difficile* clinical isolates, obtained by anaerobic bacterial culture on CCFA after 3 days at 37°C, were tested in the *C. DIFF QUIK CHEK COMPLETE*® test. For the analysis, individual colonies were picked and suspended in *Diluent* as recommended for fecal specimens. All 103 isolates gave a positive antigen reaction in the test.

Seventy of the 103 isolates (68%) were from fecal specimens that were positive for *C. difficile* toxin by tissue culture assay. Of these, 56 (80%) gave a positive toxin reaction when screened following anaerobic growth on CCFA for 3 days at 37°C.

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n = 1126	Bacteriekweek Positief	Bacteriekweek Negatief
C. DIFF QUIK CHEK COMPLETE® Antigenlijn positief	201	62
C. DIFF QUIK CHEK COMPLETE® Antigenlijn negatief	21	842

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n = 1126	Weefselkweek Positief	Weefselkweek Negatief
C. DIFF QUIK CHEK COMPLETE® Antigenlijn positief	154	109
C. DIFF QUIK CHEK COMPLETE® Antigenlijn negatief	2	861

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For Instructional Use Only

EFFECT VAN FECALE MONSTERCONSISTENTIE

Effect van de fecale monsterconsistentie op de *C. DIFF QUIK CHEK COMPLETE*®-test

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ENGLISH

About the test

[Introduction] Dengue viruses, transmitted by *Ae. aegypti* and *Ae. albopictus* mosquitoes, are widely distributed throughout the tropical and subtropical areas of the world. Dengue is an important arthropod-borne viral disease due to its high morbidity and mortality in humans. There are four known distinct serotypes: dengue virus 1, 2, 3 and 4. In children, infection is often sub-clinical or causes a self-limited febrile disease. However, if the patient is reinfected with a different serotype, risk of a more severe disease increases, such as dengue hemorrhagic fever or dengue shock syndrome. NS1 antigen is a highly-conserved glycoprotein that is present at high concentrations in the sera of dengue-infected patients, during the early clinical phase of the disease. NS1 antigen is found in specimens from primary or secondary dengue-infected patients between 1 and 9 days after onset of fever. IgM usually become detectable until 5 to 10 days after the onset of illness in cases of primary dengue infection and until 4 to 5 days after onset of illness in secondary infections. In primary infections, IgG appears by the 14th day and persists for life. In secondary infections, IgG levels rise for 1 - 2 days after the onset of symptoms and an IgM response is induced 20 days after infection.
[Test principle] Bioline™ Dengue IgG/IgM test is designed to simultaneously detect and differentiate IgG and IgM antibodies to dengue virus (Den 1,2,3 and 4) in human serum, plasma. This test also can detect all 4 Dengue serotypes by using a mixture of recombinant Dengue envelope proteins. Both the Mouse monoclonal anti-human IgG & IgM on the membrane and the Dengue virus envelope protein gold colloid will react specifically with IgG & IgM antibody to Dengue virus from human serum or plasma. Bioline™ Dengue IgG/IgM test device has 3 pre-coated lines, "G" (Dengue IgG Test Line), "M" (Dengue IgM Test Line) and "C" (Control Line) on the surface of the device. All three lines in result window are not visible before applying any specimen. The "Control Line" is used for procedural control. Control line should always appear if the test procedure is performed properly and the test reagents of control line are working. A purple "G" and "M" lines will be visible in the result window if there are enough IgG and/or IgM antibodies to Dengue virus in the specimen. If IgG and/or IgM antibodies to Dengue virus are not present in the specimen, there is no color appears in "G" and/or "M". When a specimen is added to the specimen well, Anti-Dengue IgGs and IgMs in the specimen will be captured by the relevant anti-human IgG and or anti-human IgM immobilized in two test lines across the test device. Mouse monoclonal anti-dengue Ab-colloidal gold conjugates will react with dengue virus in virus pad and forms a complex of antibodies-antigen. As this complex migrates along the length of the test device by capillary action, it will react with Anti-Dengue IgGs and IgMs in test lines and generate a colored line.
[Intended use] Bioline™ Dengue IgG/IgM test is a solid phase immunochromatographic assay for the rapid, qualitative and differential detection of IgG and IgM antibodies to dengue virus in human serum or plasma. This test is intended for professional use to aid in the presumptive diagnosis between primary and secondary dengue infection. This test provides only a preliminary test result. Therefore, isolation of virus, antigen detection in fixed tissues, RT-PCR and serological test like haemagglutination-inhibition test, more specific alternative diagnosis method must be used in order to obtain a confirmation of dengue virus infection.

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1. The Bioline™ Dengue IgG/IgM test kit contains the following items to perform the assay.
- 25 Test devices with desiccant in individual foil pouches
- Assay diluent (1 x 5 ml/vial)
- 25 Capillary pipettes (5 µl)
- 1 Instructions for use

2. Active ingredients of main components

1 test strip includes:	Assay diluent includes:
- Gold Conjugates: Mouse monoclonal anti-dengue Ab-gold colloid (0.467±0.014 µg), Mouse Ab-gold colloid (0.002±0.001 µg) - Test Line "G": Mouse monoclonal anti-human IgG (0.512±0.154 µg) - Test Line "M": Mouse monoclonal anti-human IgM (0.32±0.10 µg) - Control Line: Goat anti-mouse IgG (0.576±0.173 µg)	- 100 mM Tris Buffer (5 ml) - Proclin 300 (0.05 %) - Sodium azide (0.05 w/v%)

※ Classification according to Regulation (EC) No 1272/2008:

Product identifier	**Trade name**	Bioline™ Dengue IgG/IgM Assay diluent	
	Hazardous substance	Proclin 300 (5-chloro-2-methyl-4- isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1), CAS No. 55965-84-9)	
Classification	Skin sensitization Category 1	**Precautionary statements**	
Hazard Pictogram		**Prevention**	P261: Avoid breathing dust/fume/gas/mist/vapours/spray P272: Contaminated work clothing should not be allowed out of the workplace P280: Wear protective gloves/protective clothing/eye protection/ face protection
		Response	P302+P352: IF ON SKIN: Wash with plenty of soap and water P321: Specific treatment P333+P313: If skin irritation or rash occurs: Get medical advice/attention P363: Wash contaminated clothing before reuse
Signal Word	Warning		
Hazard Statement	H317: May cause an allergic skin reaction	**Disposal**	P501: Dispose of contents/container in accordance with local/regional/national/international regulation

Kit storage and stability

1. The test device is sensitive to both heat and humidity and must be used and stored properly.
- Store the test kit at between 1 °C and 30 °C. Do not freeze the kit or its components or expose them to high temperatures.
- Perform the test immediately after removing the test device from the foil pouch, to avoid exposure to moisture.

2. Do not use the test kit beyond its expiration date. The shelf life of the kit is as indicated on the outer package.
3. Do not use the test kit if the pouch is damaged or the seal is broken.

Warnings

1. For *in vitro* diagnostic use only. Do not reuse the test device and kit components.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
3. Do not eat or smoke while handling specimens.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Avoid splashing or aerosol formation of specimen and assay diluent.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials (i.e. specimen applicator, test device) in a biohazard container as if they were infectious waste.
8. Do not mix or interchange different specimens.
9. Care should be taken to avoid contamination of the end of the bottle when dropping assay diluent into a specimen well.
10. The components (test device and assay diluent) in this kit have been quality control tested as standard batch unit. Do not mix components from different lot numbers.
11. Do not drink assay diluent.
12. The assay diluent contains a proprietary anti-microbial agent, sodium azide, which presents nohazard to the user if normal laboratory safety precautions are followed. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.

Specimen collection and storage

1. **[Plasma]** Using venipuncture, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate) and then centrifuge the tube to get a plasma specimen.
2. **[Serum]** Using venipuncture, draw whole blood into the collection tube (NOT containing anticoagulant). Leave blood to coagulate for 30 minutes and then centrifuge the tube to get a serum supernatant.
- If plasma or serum specimens are not tested immediately, they must be refrigerated at 2 - 8 °C.
- For storage periods longer than 2 weeks, freezing (below -20 °C) is required. Specimens must be brought to room temperature (15 - 30 °C) prior to use.
- Plasma or serum specimens containing a precipitate may yield inconsistent test results. Use standard laboratory techniques to clarify such specimens prior to testing.

3. **Precaution**
- Repeated frozen-thawed cycle for specimen should be avoided.
- Anticoagulants including heparin, EDTA, and citrate do not affect the test result. Use of other anticoagulants have not been validated. Their use may affect the test result.
- Use a new capillary pipette or pipette tip for each specimen in order to avoid cross-contamination of specimens, which could produce erroneous results.

Test procedure (Refer to figure)

1. Bring all kit components and specimens to reach room temperature (between 15 °C and 30 °C) prior to opening and testing.
2. Remove the test device from foil pouch, and place it on a flat, dry surface. Label the test device with a patient identifier.
3. **[Using a micropipette]** Dispense 5 µl of serum or plasma specimen into the square specimen well.
 OR,
 [Using a capillary pipette] Draw serum or plasma to black line (5 µl) then dispense into the square specimen well.
4. Dispense 4 drops (about 90 - 120 µl) of assay diluent into the round assay diluent well.
5. Interpret test results in 15 - 20 minutes.
 ⚠ **Caution :** Do not read test results after 20 minutes; reading after 20 minutes can yield false results.

Test interpretation (Refer to figure)

- **Negative result:** The presence of only Control Line ("C") in the result window indicates that IgG and IgM antibodies to dengue are not present in the specimen or are present below the detectable levels. Retest in 3 - 5 days if dengue infection is suspected.
- **Positive result**
 ⚠ **Caution :** The presence of any test line, no matter how faint, the result is considered positive.

1) **IgM positive:** The presence of two colored lines ("C" and "M") in the result window, regardless of which line appears first, indicates that the speicmen is positive for IgM antibodies to dengue virus. This suggests a primary dengue infection. Test is positive even if "M" line is faint.
2) **IgG positive:** The presence of two colored lines ("C" and "G") in the result window, regardless of which line appears first, indicates that the speicmen is positive for IgG antibodies to dengue virus. This suggests a secondary or past dengue infection. Test is positive even if "G" line is faint.
3) **IgG and IgM positive:** The presence of three colored lines ("C", "M" and "G") in the result window, regardless of which line appears first, indicates that the speicmen is positive for both IgM and IgG antibodies to dengue virus. This suggests a late primary or early secondary dengue infection. Test is positive even if "M" and/or "G" line(s) are faint.

- **Invalid result:** No Control Line ("C") in the result window means the test is invalid, whether or not the Test Line ("M" and/or "G") is present. The instructions may not have been followed correctly or the test may have been compromised. When an invalid result is obtained, it is recommended that the specimen be retested using a new device.

Test limitations

1. This test detects the presence of antibodies to dengue in the specimen and should not be used as the sole criterion for the diagnosis of Dengue virus infection.
2. In early infections and some secondary infections, detectable levels of IgM antibodies may be low. Some patients may not produce detectable levels of antibody within the first 7 - 10 days after infection. If the test result is negative and clinical symptoms persist, patients should be retested 3 - 4 days after the first specimen.
3. Serological cross reactivity across Flavivirus genus (e.g. West Nile virus, St. Louis Encephalitis virus, Japanese Encephalitis virus, Yellow Fever virus, Zika virus) is common.
4. As with all diagnostic tests, results must be considered with other clinical information available to the physician.
5. If the test result is negative and clinical symptoms persist, additional follow-up testing using other clinical methods is recommended. A negative result does not preclude the possibility of early dengue virus infection.
6. The test procedure, precautions and interpretation of results for this test must be followed strictly when testing.

Internal quality control

The Bioline™ Dengue IgG/IgM test device has "G" and "M" (Test Lines) and "C" (Control Line) marked on the surface of the device. Neither the test lines nor the control line are visible in the result window before applying a specimen. The control line is used for procedural control. The presence of the Control Line shows that the diluent has been applied successfully and that the active ingredients of main components on the strip was still functional, but is not a guarantee that the specimen has been properly applied and does not represent a positive specimen control.

Expected value

Primary dengue is characterized by the presence of detectable IgM 3 - 5 days after the onset of infection. Secondary dengue is characterized by the elevation of specific IgG 1 - 2 days after the onset of infection and in the majority of cases this is generally accompanied by an elevation of IgM.

Performance characteristics

1. Bioline™ Dengue IgG/IgM test kit is designed to have 94.6 % (95 % CI : 88.8 - 97.5 %) of sensitivity and 96.5 % (95 % CI : 93.0 - 98.3 %) of specificity. A total of 312 serum specimens collected from Philippines were tested with the Bioline™ Dengue IgG/IgM in Korea and the specimens were confirmed by ELISA.

		Bioline™ Dengue IgG/IgM		Result (95 % CI)
		Positive	Negative	
Reference assay : ELISA test	Positive	106	6	Sensitivity : 94.6 % (88.8 - 97.5 %)
	Negative	7	193	Specificity : 96.5 % (93.0 - 98.3 %)

Even though our intended performance is as above, the results of individual laboratories may vary from these data because the results can be unique to the population it serves depending upon geographical, patient, dietary, environmental and other factors.
- Highest sensitivity and clinical utility is achieved by considering all results (IgG and IgM) in aggregate versus the reference methods. An individual is considered positive for dengue infection if any of the Test Lines appears. Positive results on Test Lines provide greater confidence in the diagnosis of dengue infection.
- The comparison results of sensitivity of dengue IgG and IgM specimen

Techinque Used	No. of positive Specimens	IgM positive only	IgG positive only	IgG+IgM Both positive	Total IgM positive	Total IgG positive	Negative
Commercial-ELISA	112	30	36	46	76	82	200
Bioline™ Dengue IgG/IgM	106	28	33	45	73	78	193

2. Reproducibility of the Bioline™ Dengue IgG/IgM has been demonstrated using in-house reference panels for within-run, between-run and batch-to-batch studies. All values were identical to reference panel acceptance criteria.
3. The performance of Bioline™ Dengue IgG/IgM kit was evaluated with potential relevant interfering substances. Hemolytic, lipaemic, icteric specimens and those containing rheumatoid factors were investigated. None of the specimens were found to interfere with this test kit.
4. Analytical Sensitivity : the limit of detection ; the smallest amount of the target marker that can be precisely detected ; have been equal or superior to a leading commercial dengue antibody detection rapid test.

Product Disclaimer:
While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

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1. Dengue haemorrhagic fever : Diagnosis, treatment, prevention and control. WHO New Edition 2009
2. Commercial Dengue Rapid Diagnostic Tests for Point-of-Care Application: Recent Evaluations and Future Needs? J.Biomed. and Biotech. Vol. 2012 , 12, 2012.
3. Immunochromatographic Test for the Clinical Evaluation of a Rapid Diagnosis of Dengue Virus Infection. Clin Vaccine Immunol. vol. 5 no. 3 407-409, 1998.
4. Songee L. ranch and Paul N. Levett. Evaluation of four methods for detection of immunogloulin M antibodies to dengue virus. Clin. Diagn. Lab. Immunol. Vol 6(4) p 555-557, 1999
5. Jan Groen et al. Evaluation of six immunoassays for detection of dengue-virus specific immunoglobulin M and G Antibodies. Clin. Diagn. Lab. Immunol. Vol 7(6) p 867-871, 2000
6. Shuenn-Jue L. Wu et al. Comparison of two rapid diagnostic assays for detection of immunoglobulin M antibodies to dengue virus. Clin. Diagn. Lab. Immunol. Vol 7(1) p 106-110, 2000
7. Kevin R. Porter et al. Evaluation of a commercially available immunoglobulin M capture enzyme-linked immunosorbent assay kit for diagnosing acute dengue infections. Clin. Diagn. Lab. Immunol. Vol 6(5) p 741-744, 1999
8. Ludolfs D et al. Serological differentiation of infections with dengue virus serotypes 1 to 4 by using recombinant antigens. J Clin Microbiol 40(11) p 4317-4320, 2002
9. Bangti Zhao et al. Cloning full-length dengue type 4 viral DNA sequences: analysis of genes coding for structural proteins. Viriol 155, p 77-88, 1986.

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Manufacturer
Fabricant
Fabricante

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

EN Pressing a capillary pipette lightly with finger tip, take the specimen up to fill line.

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
1

EN Dispense **5 µl** of **serum or plasma** into the square speicmen well.

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
2

EN **Immediately (within 20 seconds)**
Dispense **4 drops** of assay diluent into the round assay diluent well.

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EN ⚠ **Caution:** Do not read test results after 20 minutes. Reading test results too late can yield false results.

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EN Only Control Line ("C") in the result window.

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
INVALID / NON VALIDE / NO VÁLIDO

EN
- No Control Line ("C") in the result window.
- It is recommended that the specimen be retested using a new device.

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POSITIVE / POSITIF / POSITIVO

EN ⚠ **Caution:** The presence of any test line, no matter how faint, the result is considered positive.

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
EN
- Two Pink lines "C" and "M" in the result window.
- Test is positive even if "M" line is weak.

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
EN
- Two Pink lines "C" and "G" in the result window.
- Test is positive even if "G" line is weak.

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EN
- Three pink lines "C", "M" and "G" in the result window.
- Test is positive even if "M" and/or "G" line(s) are faint.

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ENGLISH

About the test

[Introduction] Hepatitis C virus (HCV) is recognized as a major agent of chronic hepatitis, transfusion acquired non-A, non-B hepatitis and liver disease throughout the world. HCV is an enveloped positive-sense, single-stranded RNA virus. Testing for HCV infection begins serologic testing with either a rapid or a laboratory-conducted assay for HCV antibody in blood. A reactive result indicates presumptive HCV infection. When to confirm the current HCV infection it is recommended nucleic acid testing (NAT) for the detection of HCV RNA be performed following HCV antibody reactive test result. Abbott Diagnostics Korea has constructed HCV genes for the expression of recombinant antigens in bacterium systems such as *E. coli* focused on structural and non-structural immunogenic regions of the HCV-encoded polyprotein. The major immunogenic antigens of HCV genome which have been reported as core, NS3, NS4 and NS5 regions of HCV genome, which are known to be highly immunodominant. For detection of HCV infection, these recombinant proteins were used as capture materials in this microbeads immunochromatographic test. Improving on the first generation HCV serologic antibody test using a single recombinant antigen, recombinant proteins with multiple antigens have been used to minimize non-specific cross-reactivity and to increase the sensitivity in this assay.

[Test principle] The Bioline™ HCV test contains a nitrocellulose membrane strip, which is pre-coated with recombinant HCV capture antigen (core, NS3, NS4 and NS5) in the test region (T). The protein A-colloid gold conjugate and the specimen moves along the membrane chromatographically to the test region. There the antigen-antibody protein A gold particle complex binds in a visible line with high degree of sensitivity and high specificity. This test device has letter "T" and "C" representing "Test Line" and "Control Line" on the surface of the case. Both the test line and control line in result window are not visible before applying the specimen. The control line is a procedural control. The control line should always appear if the test procedure is performed properly and the reagents in the control line are working.

[Intended use] Bioline™ HCV is an *in vitro* immunochromatographic, rapid assay designed for the qualitative detection of antibodies specific to HCV, in human serum, plasma (heparin, EDTA and sodium citrate) or whole blood. Bioline™ HCV is intended only for professional use as the initial test, as an aid to diagnosis. Reactive specimens should be reflexed for additional testing, either by nucleic acid testing (NAT) technologies for the detection of HCV RNA or HCV core antigen testing, to identify current HCV infection. This product is intended for use in a population with high HCV prevalence who have a history of HCV risk exposure behaviour including pregnant women. This test may not be suitable for diagnosis of early infection or blood donation screening. The performance of the assay has not been established for populations of infants or children.

Materials provided and active ingredients of main components

1. The Bioline™ HCV test kit contains the following items to perform the assay:

Catalog No.	Contents
02FK10	1. 30 Test devices with desiccant in individual foil pouches 2. Assay diluent (1.5 ml/vial) 3. 1 Instructions for use
02FK16	1. 25 Test devices with desiccant in individual foil pouches 2. Assay diluent (1 x 5 ml/vial) 3. 25 Capillary pipettes (10 µl), 25 Sterile lancets, 25 Alcohol swabs 4. 1 Instructions for use
02FK17	1. 25 Test devices with desiccant in individual foil pouches 2. Assay diluent (1 x 5 ml/vial) 3. 25 Capillary pipettes (10 µl), 25 Safety lancets, 25 Alcohol swabs 4. 1 Instructions for use

2. Active ingredients of main components
 - Test device includes:
 - Gold conjugates: Protein A – gold colloid (1.0±0.2 µg)
 - Test line: Recombinant HCV antigens (core, NS3, NS4, NS5) (1.5±0.3 µg)
 - Control line: Goat anti-human Immunoglobulin (2.0±0.4 µg)
 - Assay diluent includes: 50 mM Tris-HCl Buffer (5 ml), Sodium azide (0.02 %)

Materials required but not provided

- Micropipette, Protective gloves, Timer, Biohazard container

Kit storage and stability

1. The test kit should be stored at a temperature between 1°C and 30 °C. Do not freeze the kit or its components.
2. Assay diluent may be opened and resealed for each assay. Cap should be firmly sealed between each use. Assay diluent is stable until expiration date if kept at 1 - 30 °C.
3. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from foil pouch.
4. Do not use the test kit beyond its expiration date. The shelf life of the kit is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.

Warnings

1. The test devices are for *in vitro* diagnostic use only. Do not reuse the test device.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be proficient.
3. Do not use the pipette by mouth, smoke, drink, eat, apply cosmetics, or handle contact lenses in areas where specimens or kit components are being handled.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Clean up spills thoroughly using an appropriate disinfectant.
6. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials in a biohazard container as if they were infectious waste.
7. Do not mix or interchange different specimens.
8. Do not eat the desiccant from the foil pouch.
9. Avoid splashing or aerosol formation of specimen and assay diluent.
10. Do not mix or interchange components among different lots or those for other products.
11. Do not drink assay diluent.
12. Care should be taken to avoid contamination of the bottle nozzle when dropping assay diluent into the specimen well.
13. The assay diluent contains a proprietary anti-microbial agent, sodium azide, which presents no hazard to the user if normal laboratory safety precautions are followed. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.
14. The assay diluent contains sodium azide, which may react with lead or copper plumbing to form highly explosive metal azide compounds. When disposing of these reagents through plumbing fixtures, flush with a large volume of water to prevent azide build-up in the pipes.
15. Even though the color of assay diluent is changed to faint yellow, it doesn't affect the performance and stability of the kit.

Specimen collection and storage

1. Whole blood
[Collection by venipuncture]
Using venipuncture, draw whole blood into the collection tube containing anticoagulants including heparin, EDTA and sodium citrate).
- If the blood specimen is not immediately tested, it must be refrigerated at 2 - 8 °C.
- If stored at 2 - 8 °C, the blood specimen must be tested within 3 days of refrigeration.
- Do not use a blood specimen stored for more than 3 days; it can cause a nonspecific reaction.
- Bring blood specimen room temperature (15 - 30 °C) prior to use.

[Collection using a lancet]
- Clean the area to be lanced with an alcohol swab.
- Squeeze the fingertip: then prick the lateral side of the finger with a lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.
- Immerse the open end of a new capillary pipette (10 µl) in the next blood drop and release the pressure to draw blood into the capillary pipette up to the fill line.

2. Plasma or serum
(Plasma) Using venipuncture, draw the whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate) and then centrifuge blood to obtain the plasma specimen.
(Serum) Using venipuncture, draw the whole blood into the collection tube (NOT containing anticoagulants) then leave for 30 minutes to allow blood coagulation to occur. Centrifuge the tube to generate a serum specimen.
If plasma or serum specimens are not tested immediately, they must be refrigerated at 2 - 8 °C. For storage period longer than 2 weeks, freezing (below -20 °C) is required. Bring plasma or serum specimens to room temperature (15 - 30 °C) prior to use.
Plasma or serum specimens containing a precipitate may yield inconsistent test results. Such specimens must be clarified by centrifugation prior to assaying.

3. Precautions
- Repeated frozen-thawed cycle for specimen should be avoided.
- Anticoagulants including heparin, EDTA and sodium citrate do not affect the test result. Use of other anticoagulants has not been validated. Their use may affect the test result.
- Use new pipette tips for each specimen in order to avoid cross-contamination of other specimens which could cause erroneous results.

Test procedure

1. Bring all kit components and specimens to reach a temperature between 15 °C and 30 °C prior to testing.
2. Remove the test device from foil pouch and place it on a flat, dry surface. Label the test device with a patient identifier.
3. **[Using a micropipette]**
 Dispense 10 µl of plasma, serum or whole blood specimen into the specimen well marked "S".
 Or,
 [Using a capillary pipette]
 Dispense 10 µl of drawn whole blood specimen into the specimen well marked "S".
4. Dispense 4 drops of assay diluent into the specimen well "S".
 ⚠ **Caution:** If you do not hold the bottle vertically, it can lead to inaccurate results. Exactly, 4 drops should be added. Adding more than 4 drops may result in reddish color background or an invalid result.
5. As the test begins to work, you will see purple color move across the result window in the center of the test device.
6. Interpret test results 5 - 20 minutes after adding assay diluent. Do not read after 20 minutes.
 ⚠ **Caution:** If the test result is not legible after 5 minutes due to high background color, read again later but within 20 minutes of adding the diluent. Reading outside of this time frame (before 5 min or after 20 min) may result in false results.

Test interpretation

1. A colored control line will appear at "C" in the result window to show that the test is working properly. The "T" section of the test window indicates the test result.
2. **Non-reactive result:** The presence of only the control line (C) within the result window indicates a non-reactive result.
3. **Reactive result:** The presence of the test line (T) and the control line (C) within the result window, regardless of which line appears first, indicates a reactive result.
 ⚠ **Caution:** The presence of any test line, no matter how faint, the result is considered reactive.
4. **Invalid result:** If the control line (C) is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be retested using a new test device.

Test limitations

1. A non-reactive result does not preclude the possibility of infection with HCV. Other clinically available tests are required if questionable results are obtained. In addition, a definitive clinical diagnosis of infection should not be based on the result of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
2. Due to the inherent design of qualitative IVD tests, a faint or absent test line (false non-reactive) may occur in specimens containing high antibody densities; the prozone effect. In order to obtain a definitive result, all clinical and laboratory findings should be evaluated.

Internal quality control

The Bioline™ HCV test device has 2 pre-coated lines on the surface of the test: "T" (test line) and "C" (control line). Neither the test line nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the diluent has been applied successfully and that the active ingredients of the main components on the strip are functional, but it is not an assurance that the specimen has been properly applied; it is not a reactive specimen control.

Performance characteristics

1. The Bioline™ HCV test kit is designed to have 99.3 % (95 % CI: 96.1 - 99.9 %) of sensitivity and 98.1 % (95 % CI: 94.5 - 99.4 %) of specificity.
 A total of 299 serum specimens collected from Korea were tested on the Bioline™ HCV in Korea. The serum specimens were comprised of 142 positive and 157 negative clinical serum specimens confirmed by a leading commercial anti-HCV ELISA test. The results gave sensitivity of 99.3 % (141/142), a specificity of 98.1 % (154/157).

Reference assay		Bioline™ HCV		Total Results
		Reactive	Non-reactive	
Anti-HCV ELISA	Positive	141	1	142
	Negative	3	154	157
Sensitivity (95 % CI)		99.3 % (96.1 - 99.9 %)		
Specificity (95 % CI)		98.1 % (94.5 - 99.4 %)		

 Even though our intended performance is as above, the results of individual laboratories may vary from these data because the results can be unique to the population it serves depending upon geographical, patient, dietary, environmental and other factors.

2. Bioline™ HCV test kit has been evaluated in 4 different sites as below.
1) Study 1
 A total of 1,187 specimens collected from Korea were tested on the Bioline™ HCV in Korea. The serum specimens were comprised of 157 positive and 1,030 negative clinical serum specimens confirmed by a leading commercial anti-HCV ELISA test. The specimens repeatedly reactive both by ELISA and Bioline™ HCV were tested by the confirmatory assay using RT-PCR.

Reference assay		Bioline™ HCV		Total Results
		Reactive	Non-reactive	
Anti-HCV ELISA	Positive	157	0	157
	Negative	6	1,024	1,030
Sensitivity (95 % CI)		100 % (97.6 - 100 %)		
Specificity (95 % CI)		99.4 % (98.7 - 99.7 %)		

2) Study 2 : European performance evaluation according to the common technical specification (2009/886/EC)
2-1) Diagnostic sensitivity for anti-HCV detection
 - Diagnostic panel of seroconversions
 410 anti-HCV positive specimens were tested by Sanquin in the Netherlands and the German Red Cross. The diagnostic sensitivity for anti-HCV antibody detection, calculated on 410 positive specimens, was 99.3 %.

 Table1. Result obtained with the Bioline™ HCV on all Anti-HCV positive specimens

		Bioline™ HCV	
		Reactive	Non-reactive
Anti-HCV positive serum/plasma (n=213)		211	2
Anti-HCV positive whole blood/plasma couples (n=100)		100	0
Anti-HCV with known genotype (n=97)		96	1
Total (n=410)		407	3
Sensitivity (95 % CI)		99.3 % (97.9 - 99.8 %)	

 - Anti-HCV positive serum/plasma
 213 anti-HCV positive specimens confirmed with CHIRON HCV RIBA™ 3.0 SIA or INNO-LIA HCV Score without genotype information were tested at Sanquin. In two specimens the Bioline™ HCV obtained a non-reactive result.
 - Anti-HCV positive whole blood/plasma couples (Paired specimens)
 100 paired whole blood and plasma specimens anti-HCV positive with Abbott Architect anti-HCV, taken from HCV infected patients were tested at the German Red Cross. The Bioline™ HCV test was reactive on all 100 paired specimens.
 - Anti-HCV with known genotype
 94 specimens from patients with a known HCV genotype 1 to 6, based on the VERSANT HCV Genotype 2.0 Assay (LiPA) were tested by Sanquin, and 3 specimens of genotype 5 were tested by the German Red Cross.

 Table 2. Anti-HCV positive specimens with known genotype

		Bioline™ HCV	
		Reactive	Non-reactive
HCV genotype	**# Specimens**		
Genotype 1 (n=23)	1	1	0
	1a	10	0
	1b	12	1 1*
Genotype 2 (n=22)	2a/2c	13	13 0
	2b	9	9 0
Genotype 3 (n=25)	3	1	1 0
	3a	22	22 0
	3b	1	1 0
	3k	1	1 0
Genotype 4 (including non 4a) (n=20)	4a/4c/4d	4	4 0
	4c/4d	14	14 0
Genotype 5 (n=5)	5a	2	2 0
	5	3	3 0
Genotype 6 (n=2)	6	1	1 0
	6a	1	1 0
Agreement (95 % CI)		99.6 % (97.7 - 99.9 %)	

 *The specimen was weak positive of Hepatitis C confirmed by commercial anti-HCV assay.

3) Analytical performance
 Analytical specificity
 - 200/200 hospitalised patients were non-reactive on the Bioline™ HCV.
 - 200/200 specimens of pregnant women (incl. 20 specimens of multipara) were non-reactive on the Bioline™ HCV.
 - 100/100 specimens containing potentially cross reactive substances were non-reactive on the Bioline™ HCV (20 Anti-HBs positive, 20 Anti-HBc positive, 20 Anti-HIV positive, 10 Anti-HTLV I/II positive, 20 Anti-HEV positive, 10 Rheumatoid factor positive).
 The following potential cross-reacting pathogens had no effect on test results of Bioline™ HCV.

Anti-HBs	CMV	VIH	Syphilis
Borrelia burgdorferi	EBV	HTLV	Toxoplasma
Chlamydia	HBsAg	Influenza	*Trypanosoma cruzi* I / II

 - The following 9 potential interfering substances had no effect on test results of Bioline™ HCV: Pregnant women, high cholesterol (1240 mg/dl.), high bilirubin (11.4 mg/dl.), Rheumatoid factor (1228 IU/ml), lipemic, hemolyzed, autoimmune, alcoholic cirrhosis, multiparous pregnancy
 - Pharmaceutical substances
 There was no significant interference with the following 25 drugs on Bioline™ HCV. All drugs were tested at concentrations of 250 µg/mL.

Abacavir	Ergocalciferol	Naproxen	
Acetaminophen	Folic acid	Nevirapine	
Acetylsalicylic acid	Hydrochlorothiazide	Pantoprazole	
Amoxicillin	Ibuprofen	Pyrasinamide	
Aspirin	Iron Chloride	Rifampicin	
Cholecalciferol	Isoniazid	Ritonavir	
Cyclobenzaprine	L-ascorbic acid	Salicylic acid	
Darunavir	Magnesium sulfate		
Diclofenac	Metformin		

4) Sensitivity on seroconversion panels
 20 commercially available seroconversion panels (SeraCare Life Sciences and Zeptometrix) were tested with Bioline™ HCV by Sanquin.
 In total 54 specimens were tested reactive with Bioline™ HCV and 47 with the competitor tests.

Panel ID	Bioline™ HCV	Ortho HCV 3.0 or Ortho Enhanced SAVe Anti-HCV 3.0 (*)
PHV904	3/7	3/7
PHV905	3/9	4/9
PHV911	3/5	3/5
PHV913	0/4	0/4
PHV914	5/9	3/9
PHV915	2/4	2/4
PHV918	1/8	2/8
PHV919 (*)	7/7	3/7
PHV920	6/10	6/10
6213	2/12	2/12
6214	6/13	4/13
6224	2/6	0/6
6226	4/12	3/12
6228	3/12	3/12
6229	3/8	3/8
9044	2/6	2/6
9045	2/8	1/8
9047	4/10	4/10
9054	1/10	1/10 (*)
9058	2/5	1/5 (*)
Total de sangrados reactivos (*)	54	47

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2) Prozone effect
Bioline™ HCV may exhibit prozone effect (false non-reactive result) in specimens which have higher than SrCO of approx. 11.0 in Abbott Architect and SrCO approx. 3.0 in Ortho HCV version 3.0 ELISA anti HCV assays.

3) Specimen matrix
Validation on whole blood was done by testing 500 negative and 100 positive anti-HCV specimens. The performance of Bioline™ HCV on whole blood was comparable to the performance on plasma specimens. (See table 1 and table 3 presented in study 2 above.) Validation on serum and different plasma specimen types (EDTA / Heparin / Sodium citrate) was performed by German Red Cross. The results obtained on negative and positive specimens are identical in serum, EDTA plasma, heparin plasma and citrate plasma.

Table 4. Serum/plasma (EDTA / Heparin / Sodium citrate) equivalence

Specimen type	No. of Bioline™ HCV reactive/No. of positive specimens	No. of Bioline™ HCV non-reactive/ No. of negative specimens
Serum	25/25	25/25
EDTA plasma	25/25	25/25
Heparin plasma	25/25	25/25
Na-Citrate plasma	25/25	25/25

4) Complement factors interference in fresh serum specimens
In total 25 negative specimens, spiked with an anti-HCV positive specimen were tested within 24 hours after collection and retested after being stored at 4 °C for 1, 2, 3 and 4 days. No differences were observed on the results obtained on the fresh specimens and the same specimen stored for 1 to 4 days at 4 °C.

Table 5. Result obtained on spiked fresh specimens

		Bioline™ HCV	
		Reactive	Non-reactive
25 negative specimens	**day 0**	0	25
	day 0	25	0
25 negative specimens spiked with anti-HCV	**day 1**	25	0
	day 2	25	0
	day 3	25	0
	day 4	25	0

5) Reproducibility of the Bioline™ HCV has been demonstrated by within-run, between-run, and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria.

Product Disclomer:
While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

Warning:
The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect reactive or non-reactive result using this product.

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Bibliography of suggested reading / Bibliographie

1. Arsah G., Czeslaw W., Chao Lin, Stephen M. Feinstone, and Charles M. Rice : Expression and Identification of Hepatitis C Virus Polyprotein Cleavage Products. Journal of Virology, March. 1993, p.1385-1395
2. Young Gyu Cho, Min Kyung Yi, Kyung Lib jang, Chang Min Kim and Young Chul Sung : Cloning and Overexpression of the Highly Immunogenic Region of HCV Genome from Korean Patients. Mol. Cells, Vol. 3, 4-7 : 416
3. S. Osborne, E. Cecconato, S. Grira, F. Garetto, R. Calegaro, C. Rosa and F. Bonelli : Expression in *E. coli* and purification of a chimeric p22-NS3 recombinant antigen of Hepatitis C Virus(HCV). Federation of European Biochemical Societies, Volume 324, number 3, 253-257
4. A. Yoshikawa, K. Takahashi, S. Kishimoto : Serodiagnosis of hepatitis C virus infection by ELISA for antibodies against the putative core protein (p20C) expressed in Escherichia coli. Journal of Immunological Methods, 148 (1992) 143-150
5. Mitchell B.L. et al, Impact freeze-thaw cycles and storage time on plasma samples used in mass spectrometry based biomarker discovery projects. Cancer Informatics, 2005;1:98-104
6. Cuhadar S. et al, The effect of storage time and freeze-thaw cycles on the stability of serum samples. Biochemia Media, 2013;23(1):70-77

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PREPARATION / PRÉPARATION /PREPARACIÓN / PREPARAÇÃO

1

EN Open the package and look for the following:
1. Test device with desiccant in individual foil pouch
2. Assay diluent
3. Instructions for use

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Including only for Catalog No. 02FK16 and 02FK17

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PT
1. Pipeta capilar (10 µL)
2. Lanceta
3. Zaragatoa com álcool

2

EN Carefully read the instructions on how to use the Bioline™ HCV test kit.

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3

EN Look at the expiration date at the back of the foil pouch. If the expiration date has passed, use another kit. To avoid false results, ensure that the assay diluent used is from the same kit as the new test device.

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4

EN **Open the foil pouch and look for the following:**
1. Test device
2. Desiccant

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

1

EN Take 10 µl of serum, plasma or whole blood specimen using a micropipette.
Dispense 10 µl of serum, plasma or whole blood specimen into the specimen well "S".

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
2

EN **Hold assay diluent bottle vertically and dispense 4 drops of assay diluent** into the specimen well "S". **Exactly, 4 drops should be added.** Do not let bottle tip touch device in order to avoid cross-contamination.

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
3

EN Interpret test results 5 - 20 minutes after adding assay diluent. Reading outside of this time frame (before 5 min or after 20 min) may provide false results.

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
1

EN Clean the area to be lanced with an alcohol swab.

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2


EN Squeeze the fingertip then prick the lateral side of the finger with a lancet provided. Wipe away the first blood drop. Then, safely dispose of the lancet immediately after.

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
[02FK17]


3

EN Immerse the open end of a new capillary pipette (10 µl) in the next blood drop and release the pressure to draw blood into the capillary pipette up to the fill line.

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
4

EN **Dispense 10 µl of drawn whole blood specimen** in the specimen well marked "S". Lightly touch the capillary pipette to the specimen pad while dispensing.

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5

EN Hold assay diluent bottle vertically and dispense 4 drops of assay diluent into the specimen well "S". Exactly, 4 drops should be added. Do not let bottle tip touch device in order to avoid cross-contamination.

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
EN The presence of only the control line (C) within the result window indicates a non-reactive result.

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
REACTIVE / RÉACTIF / REACTIVO / REATIVO

EN The presence of the test line (T) and the control line (C) within the result window, regardless of which line appears first, indicates a reactive result.
⚠ **Caution**: The presence of any test line, no matter how faint, the result is considered reactive.

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
INVALID / NON VALIDE / NO VÁLIDO / INVÁLIDA

EN If the control line (C) is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be retested using a new test device.

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ENGLISH

About the test
[Intended use] The Bioline™ HBsAg test is an *in vitro* immunochromatographic assay designed for qualitative determination of HBsAg in human serum or plasma. Bioline™ HBsAg is intended only for professional use as an aid to diagnosis. Reactive specimens should be confirmed by a supplemental assay. This kit is for *in vitro* diagnostic use only. This product suitable for use on pregnant women, but in case of neonates, we haven't yet to perform the test. This test may not be suitable for diagnosis of early infection or blood donation screening.
[Test principle] This test device contains a membrane strip, which is pre-coated with mouse monoclonal anti-HBs antibody on test band region. The mouse monoclonal anti-HBs-colloid gold conjugate and the specimen move along the membrane chromatographically to the test region (T) and forms a visible line as the antibody-antigen-antibody gold particle complex forms.

Materials provided
The Bioline™ HBsAg test kit contains the following items to perform the assay:
1. Bioline™ HBsAg [30 Tests/kit] Cat. No. 01FK10 :
- 30 Test devices with desiccant in individual foil pouches
- 1 Instructions for use
2. Bioline™ HBsAg [100 Tests/kit] Cat. No. 01FK11 :
- 10 Multi-test devices with desiccant in individual foil pouch x 10
- 1 Instructions for use


Materials required but not provided
- Micropipette

Kit storage and stability
1. The Bioline™ HBsAg should be stored between 1°C and 30°C. Do not freeze the kit.
2. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from the foil pouch.
3. Do not use the test kit beyond its expiration date.
4. The shelf-life of the kits is as indicated on the outer package.
5. Do not use the test kit if the pouch is damaged or the seal is broken.

Specimen collection and storage
1. [Plasma] Collect the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA and sodium citrate) by venipuncture and then centrifuge blood to get plasma specimen. Anticoagulants such as heparin, EDTA, and citrate do not affect the test result.
2. [Serum] Collect the whole blood into the collection tube (NOT containing anticoagulants such as heparin, EDTA and sodium citrate) by venipuncture, leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
3. If plasma or serum specimens are not tested immediately, they should be refrigerated at 2~8°C. For storage period longer than 3 days, freezing(below -20°C) is recommended. They should be brought to 15~30°C prior to use.
4. Plasma or serum specimens containing a precipitate may yield inconsistent test results. Such specimens must be clarified prior to assaying.

Warnings
1. For *in vitro* diagnostic use only. Do not reuse the test device.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
3. Do not eat or smoke while handling specimens.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Avoid splashing or aerosol formation of specimen.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials in a biohazard container as if they were infectious waste.


Test procedure
1. Allow all kit components and specimen to room temperature (between 15°C and 30°C) prior to testing.
2. Remove the test device from the foil pouch and place it on a flat, dry surface.
3. Dispense 100µl of plasma or serum specimen into the specimen well (Figure 1).
4. As the test begins to work, you will see purple color move across the result window in the center of the test device.
5. Interpret test results at 20 minutes.


[Figure1]

SP 100 µl : Serum 100µl or Plasma 100µl

Test interpretation
1. A purple control band will appear in the left section of the results window to show that the test is working properly.
2. The right section of the results window indicates the test results. If another purple band appears in the right section of the results window, this is the test band.
- **Negative Result :** The presence of only control band (C) within the results window indicates a negative result.

- **Positive Result :** The presence of the test line (T) and the control line (C) within the results window, regardless of which band appears first, indicates a positive result.
⚠ **Caution:** The presence of any test line, no matter how faint, the result is considered positive.

- **Invalid Result :** If the control band (C) is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be re-tested using a new test kit.

Test limitations
A negative result does not preclude the possibility of infection with HBV. Other clinically available tests are required if questionable results are obtained. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

Internal quality control
The Bioline™ HBsAg test has a test line (T) and a control line (C) on the surface of the device. Neither the test line nor the control line is visible in the result window before applying a specimen. The control line is used for procedural control and shows only that the active ingredients of main components on the strip are functional, but is not a guarantee that the specimen has been properly applied and does not represent a positive specimen control.

Performance characteristics
The 276 specimens of WHO HBsAg panel that derived serum/plasma specimens from Africa, Asia, Europe and Latin America origin were tested on the Bioline™ HBsAg test kit. There were 98 HBsAg positive specimens and 178 HBsAg negative specimens. Each specimen was confirmed by commercial EIA tests, Hepanostika HBsAg Uniform II (Organon Teknika) and Monolisa Ag HBs Plus (BioRad Diagnostics).

Origin	HBsAg positive specimens	HBsAg negative specimens	Total
Africa	16	43	59
Asia	30	30	60
Europe	30	67	97
Latin America	22	38	60
Total	**98**	**178**	**276**

In this multi-site evaluation of 276 specimens, we found the relative sensitivity is 100% (98/98), the relative specificity is 100% (178/178). The results are summarized in the following tables.

Results		Bioline™ HBsAg			95% CI
		Positive	Negative	Total Results	
Reference	Positive (98)	98	0	98	96.2 - 100%
	Negative (178)	0	178	178	97.9 - 100%
Total		98	178	276	

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BIBLIOGRAFIA

1. Stallman GND. The International Committee on Systematic Bacteriology: Sub committee on the Taxonomy of Leptospira. Int J Syst Bacteriol 1987; 37:472.
2. Levett PN. Leptospirosis. Clin Microbiol Rev 2001;14:296-326
3. Faine S, ed. Guidelines for the control of leptospirasis. Geneva: World Health Organization, 1982.
4. Cumberland PC, Everard COR, Levett PN. Assessment of the efficacy of the IgM enzyme-linked immunosorbent assay (ELISA) and microscopic agglutination test (MAT) in the diagnosis of acute leptospirosis. Am J Trop Med Hyg. 1999; 61: 731–734.
5. Adler B, Murphy AM, Locarnini SA, Faine S. Detection of specific anti-leptospiral immunoglobulins M and G in human serum by solid-phase enzyme-linked immunosorbent assay. J Clin Microbiol. 1980; 11: 452–457.
6. Appassakij H, Silpapojakul K, Wansit R, et al: Evaluation of the immunofluorescent antibody test for the diagnosis of human leptospirosis. Am J Trop Med Hyg 1995; 52: 340.

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ENGLISH

About the test

[Introduction]
Malaria is a serious and sometimes fatal parasitic disease characterized by fever, chills and anemia. The disease is caused by a parasite that is transmitted from one human to another through the bite of infected *Anopheles* mosquitoes. There are four types of malaria that can infect humans: *Plasmodium falciparum, P. vivax, P. ovale* and *P. malariae*. In humans, the parasites (called sporozoites) migrate to the liver where they mature and are released in another form: merozoites. The disease occurs in more than 90 countries worldwide, and it is estimated that there are some 300 million clinical cases and 0.5 million malaria-related deaths per year. At present, malaria is diagnosed by looking for parasites in a drop of blood. Blood is placed on a microscope slide and stained so that the parasites will be visible under a microscope.

[Test principle]
Bioline™ Malaria Ag P.f/P.f/P.v test utilizes the principle of immunochromatography. Bioline™ Malaria Ag P.f/P.f/P.v test kit contains a membrane strip, which is pre-coated with mouse monoclonal antibodies specific to histidine-rich protein II (HRP2) of *P. falciparum* in test line (T1), mouse monoclonal antibodies specific to plasmodium lactate dehydrogenase (pLDH) of *P. falciparum* in test line 2(T2) and mouse monoclonal antibodies specific to pLDH of *P. vivax* in test line 3(T3) region. This test device has a letter of T1, T2, T3 and a "Test Line" and "Control Line" on the surface of the device. All the Test Lines and Control Line in result window are not visible before applying any specimens. The Control Line is used for procedural control. Control Line should always appear if the test procedure is performed properly and the test reagents of control line are working. It is indicated as follows: As the Pf (HRP2, pLDH) and Pv reactive specimen flow through the membrane assembly of the cassette after addition of the assay diluent, they form a complex with mouse monoclonal antibodies specific to Pf HRP2-gold conjugates and mouse monoclonal antibodies specific to pan pLDH -gold conjugates. This complex move further on the membrane to the test line region where it is immobilized by mouse monoclonal P.f HRPII antibodies, P.f pLDH antibodies and Pv pLDH antibodies. At this region, the membrane leading to formation of a red-purple colored band in the respective regions which confirms a reactive test result. Absence of a colored band in the test region indicates a non-reactive test result for the corresponding antigen.

[Intended Use]
Bioline™ Malaria Ag P.f/P.f/P.v test kit is an *in vitro* rapid, qualitative test for the detection of histidine-rich protein II (HRP2) antigen and plasmodium lactate dehydrogenase (pLDH) from Malaria *Plasmodium falciparum* and pLDH from *Malaria Plasmodium vivax* in human whole blood. Bioline™ Malaria Ag P.f/P.f/P.v kit is intended only for professional use and for *in vitro* diagnostic use. Positive specimens should be confirmed by a supplemental assay such as microscopic examination of blood smear.

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1. Bioline™ Malaria Ag P.f/P.f/P.v test kit contains the following items to perform the assay:
 - 25 Test devices with desiccant in individual foil pouches
 - Assay diluent (5 x 5 ml/vial)
 - 25 Disposable inverted cups (5 µl), 25 Sterile lancets, 25 Alcohol swabs
 - 1 Instructions for use
2. Active ingredients of main components
 - 1 test strip includes : Gold conjugate : Mouse monoclonal antibodies specific to *P.f* HRP2 - gold colloid (0.10±0.02 µg), Mouse monoclonal antibodies specific to pan pLDH - gold colloid (0.25±0.05 µg), Test line 1 : Mouse monoclonal antibodies specific to *P.f* HRP2 (0.5±0.1 µg), Test line 2 : Mouse monoclonal antibodies specific to *P.f* pLDH (0.5±0.1 µg), Test line 3 : Mouse monoclonal antibodies specific to *P.v* pLDH (0.5±0.1 µg), Control Line : Goat anti-mouse IgG (1.0±0.2 µg)
 - Assay diluent : Bovine serum albumin (q.s.), 200 mM Phosphate Buffer (q.s.), TRITON™ X-100 (q.s.), Sodium azide (q.s.)

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- Protective gloves, Timer, Biohazard container

Kit storage and stability

1. The test kit should be stored at a temperature between 1 °C and 40 °C. Do not freeze the kit or its components. Assay diluent can should be kept firmly sealed between each use.
 ⚠ **Caution :** When stored at refrigerator, all kit components must be brought to room temperature (15 - 40 °C) minimum 30 mins prior to the test. Do not open the pouch whilst components come to room temperature.
2. Assay diluent may be opened and resealed for each assay. It is stable until expiration date if kept at 1 - 40 °C.
3. The test device is sensitive to both heat and humidity. Perform the test immediately after removing the test device from foil pouch.
4. Do not use the test kit beyond its expiration date.
5. The shelf life of the kit is as indicated on the outer package.
6. Do not use the test kit if the pouch is damaged or the seal is broken.
7. When transporting or storing the kit, avoid exposure to high temperature (above 45 °C) for a period longer than 1 week.

Warnings

1. For *in vitro* diagnostic use only. Do not reuse the test device and kit components.
2. The instructions must be followed exactly to achieve accurate results. Any individual performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.
3. Do not eat or smoke while handling specimens.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Avoid splashing or aerosol formation of specimen and assay diluent.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials (i.e. lancet, disposable inverted cup, test device) in a biohazard container as if they were infectious waste.
8. Do not mix or interchange different specimens.
9. Do not mix reagent of different lots or those for other products.
10. Do not pipette by mouth.
11. Do not eat the desiccant in the foil pouch.
12. Do not drink assay diluent.
13. The assay diluent contains a proprietary anti-microbial agent, sodium azide, which presents no hazard to the user if normal laboratory safety procedures are followed. If contact with assay diluent to the eyes and/or skin, wash affected area with soap and water immediately. If irritation or signs of toxicity occur, seek medical attention.

Specimen collection and storage

[Collection by venipuncture]
1. Using venipuncture, draw whole blood into the collection tube (containing anticoagulants including heparin, EDTA and sodium citrate).
2. If the blood specimen is not immediately tested, it must be refrigerated at 2 - 8 °C.
3. If stored at 2 - 8 °C, the blood specimen must be tested within 3 days.
4. Do not use a blood specimen stored for more than 3 days; it can cause a nonspecific reaction.

[Collection using a lancet]
1. Clean the area to be lanced with an alcohol swab.
2. Squeeze the end of the fingertip and pierce with a sterile lancet provided.
3. Wipe away the first drop of blood with sterile gauze or cotton.
4. Using a disposable inverted cup (5 µl) provided, dip the circular end of inverted cup into the blood specimen.

[Precaution]
1. Anticoagulants including heparin, EDTA and sodium citrate do not affect the test results. Use of other anticoagulants have not been validated. Their use may affect the test results.
2. Use a new disposable inverted cup for each specimen in order to avoid cross contamination of specimens, which could produce erroneous results.
3. Discard the lancet or alcohol swab if package is pierced or damaged. The use of damaged lancet may cause any infection at the punctured skin due to cease to existing its sterility.
4. Repeated frozen-thawed cycle for specimen should be avoided.

Test procedure (Refer to figure)

1. Bring all kit components and specimen to room temperature prior to testing.
2. Remove the test device from foil pouch and place it on a flat, dry surface. Label the test device with a patient identifier.
3. Clean the area to be lanced with an alcohol swab.
4. Prick the lateral side of the finger with the sterile lancet provided. Then, safely dispose of the lancet immediately after.
5. Using a disposable inverted cup (5 µl) provided, dip the circular end of inverted cup into the blood specimen and carefully place the circular end of the inverted cup into the round specimen well.
6. Holding the assay diluent buffer bottle vertically, dispense 4 drops into the square assay diluent well.
 ⚠ **Caution :** Do not let bottle nozzle touch device in order to avoid cross-contamination. Hold bottle vertically while dispensing. If you do not hold the bottle vertically, it can lead to erroneous results.
7. Read the test result after 15 minutes.
 ⚠ **Caution :** Do not read test results after 30 minutes. ; reading after 30 minutes can yield erroneous results.

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1. Non-reactive result :
The presence of Control Line ("C") within the result window indicates a non-reactive result.
2. Reactive result :
⚠ **Caution :** The presence of any test line, no matter how faint, the result is considered reactive.
- ***P.f* Reactive result**
 "C" and "T1" (Two lines), "C" and "T2" (Two lines) or "C", "T1" and "T2" (Three lines) :
 The presence of two colored lines ("C" and "T1"), ("C" and "T2") or three colored lines ("C", "T1" and "T2") within the result window, regardless of which line appears first, indicate *P.f* reactive result.
- ***P.v* Reactive result**
 "C" and "T3" (Two lines) :
 The presence of two colored lines ("C" and "T3") within the result window, regardless of which line appears first, indicate pLDH of *P.vivax* reactive result.
- ***P.f* and *P.v* Reactive result**
 "C", "T1" and "T3"(Three lines), "C", "T2" and "T3"(Three lines) or "C", "T1", "T2" and "T3"(Four lines) : The presence of three colored lines ("C", "T1" and "T3"), ("C", "T2" and "T3") or four colored lines ("C", "T1", "T2" and "T3") within the result window, regardless of which line appears first, indicate *P.f* and *P.v* reactive result.

3. Invalid result :
If the Control Line ("C") is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be retested using a new test device.

Test limitations

1. The test procedure, precautions and interpretation of results for this test must be followed precisely.
2. The test is limited to the detection of antigen to Malaria *P. falciparum and/or P. vivax*. Although the test is very accurate in detecting HRP2 or pLDH specific to *P.f* and/or pLDH specific to *P.v*, a low incidence of false results can occur. As with all diagnostic tests, clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.
3. The tests are applicable to general population who can be taken enough volume of the claimed specimen types for the test. However, it is not verified for such a population specific to pregnant women or neonates.

Internal quality control

Bioline™ Malaria Ag P.f/P.f/P.v test device has "Test Lines (T1, T2 and T3)" and "Control Line" on the surface of the device. All the test lines and control line in result window are not visible before applying any specimens. The control line is used for procedural control. The control line used for procedural control and shows only that the diluent has been applied correctly and that the test strip and the other ingredients of main components on the strip was still functional, but is not an assurance that the specimen has been properly applied ; it is not a reactive specimen control.

Expected values

Bioline™ Malaria Ag P.f/P.f/P.v test has been compared with microscopic examination. The overall accuracy is greater or equal to 95%.

Performance characteristics

1. Evaluation interne de Bioline™ Malaria Ag P.f/P.f/P.v test
 A. Sensitivity of Bioline™ Malaria Ag P.f/P.f/P.v test for the detection of *P.f* HRP2

No. of parasites/ µl of blood	Microscopy (No. of reactive)	Bioline™ Malaria Ag P.f/ P.f/P.v (No. of reactive)	Sensitivity (95% CI)
≤50	15	14	93.3% (70.2 - 98.8%)
51-100	30	30	100% (88.7 - 100%)
101-500	83	83	100% (95.6 - 100%)
501-1,000	87	87	100% (95.8 - 100%)
1,001-5,000	62	62	100% (94.2 - 100%)
>5,000	35	35	100% (90.1 - 100%)
Total	**312**	**311**	**99.7% (98.2 - 99.9%)**

 B. Sensitivity of Bioline™ Malaria Ag P.f/P.f/P.v test for the detection of pLDH specific to *P.f*

No. of parasites/ µl of blood	Microscopy (No. of reactive)	Bioline™ Malaria Ag P.f/P.f/P.v (No. of reactive)	Sensitivity (95% CI)
≤50	6	3	50.0% (11.8 - 81.2%)
51-100	10	9	90.0% (59.6 - 98.2%)
101-500	25	24	96.0% (80.5 - 99.3%)
501-1,000	18	18	100% (82.4 - 100%)
1,001-5,000	29	29	100% (88.3 - 100%)
>5,000	23	23	100% (85.7 - 100%)
Total	**111**	**106**	**95.5% (89.9 - 98.1%)**

 C. Sensitivity of Bioline™ Malaria Ag P.f/P.f/P.v test for the detection of pLDH specific to *P.v*

No. of parasites/ µl of blood	Microscopy (No. of reactive)	Bioline™ Malaria Ag P.f/P.f/P.v (No. of reactive)	Sensitivity (95% CI)
≤50	15	8	53.3% (30.1 - 75.2%)
51-100	30	29	96.7% (83.3 - 99.4%)
101-500	83	83	100% (95.6 - 100%)
501-1,000	87	87	100% (95.8 - 100%)
1,001-5,000	62	62	100% (94.2 - 100%)
>5,000	35	35	100% (90.1 - 100%)
Total	**312**	**304**	**97.4% (95 - 98.7%)**

 D. Total Sensitivity and Specificity of Bioline™ Malaria Ag P.f/P.f/P.v test

Reference method (Microscopic examination)	Result of Bioline™ Malaria Ag P.f/P.f/P.v		Total	Result (95% CI)
	Reactive	Non-reactive		
P.f HRP2 Reactive	311	1	312	**Sensitivity** 99.7% (98.2 - 99.9%)
P.f pLDH Reactive	304	8	312	**Sensitivity** 97.4% (95 - 98.7%)
P.v pLDH Reactive	106	5	111	**Sensitivity** 95.5% (89.9 - 98.1%)
Non-reactive	2	298	300	**Specificity** 99.3% (97.6 - 99.8%)

2. Reproducibility of Bioline™ Malaria Ag P.f/P.f/P.v has been demonstrated using in-house reference panels for within-run, between-run and batch-to-batch studies. All values were identical to reference panel acceptance criteria.
3. Following 26 potential cross-reacting microorganism have no impact on the test result.

Type	Cross reaction substance	
Antigen	Dengue Virus Type 1	Dengue Virus Type 2
	Dengue Virus Type 3	Dengue Virus Type 4
	HBsAg	Yellow Fever Virus
	Chikungunya	Rubella
	West Nile Virus(WNV)	
Recombinant Antigen	*Plasmodium ovale*	*Plasmodium malariae*
Antibody	Dengue	Anti-CMV
	HIV-1	Anti-EBV
	HIV-2	Parvovirus B19
	HCV	VZV IgG/IgM
	HAMA	Syphilis
	Chagas	
Serum	*Trypanosoma brucei gambiense*	*Trypanosoma brucei rhodesiense*
Whole blood	*Plasmodium falciparum*	*Plasmodium vivax*

4. Following 36 potentially interfering substances have no impact on the test result.

Type	Potential interfering substance		
Endogenous substance	Albumin	Acetylsalicylic acid	Bilirubin (conjugate)
	Cholesterol	Hemoglobin	Triglycerides
	Nicotine		
Exogenous Substance	Acetaminophen	Amoxicillin	Cephalexin
	Doxycycline	Erythromycine	Ibuprofen
	Caffeine	Heparin	Amodiaquine
	Desethylamodiaquine	Artemether	Lumefantrine
	Artesunate	Dihydroartemesinin	Atovaquone
	Proguanil	Chloroquine	Clindamycin
	Dihydroartemisinin	Piperaquine	Doxycycline
	Mefloquine	Primaquine	Quinine
	Sulfadoxine	Pyrimethamine	
Autoimmune /Clinical specimen	Systemic lupus erythematosus	Sjogren's syndrome	Rheumatoid arthritis

5. Assay for the detection of Malaria according to the density of parasite.
 * HRP2 expressed *Plasmodium falciparum* : Sensitivity with > 50 parasite/ul of blood is 100%.
 * pLDH expressed *Plasmodium falciparum* : *Sensitivity with > 50 parasite/ul of blood is greater than 99.7%.*
 * pLDH expressed *Plasmodium vivax* : Sensitivity with > 50 parasite/ul of blood is greater than 98.2%.

Product Disclaimer :
While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

Warning :
The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect reactive or non-reactive diagnosis using this product.

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PREPARATION / PRÉPARATION / PREPARACIÓN / PREPARAÇÃO

1

EN Open the package and look for the following:
1. Test device with desiccant in individual foil pouch
2. Assay diluent
3. Instructions for use
4. Disposable inverted cup (5 µl)
5. Sterile lancet
6. Alcohol swab

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3

EN Next, look at the expiration date on the back of the foil pouch. If the expiration date has passed, use another lot. To avoid false results, ensure that the assay diluent used is from the same kit as the new test device.

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4

EN Open the foil pouch and look for the following:
1. Result window
2. Specimen well
3. Assay diluent well

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
PT Limpe a área a ser lancetada com swab de álcool.

2

EN Prick the lateral side of the finger with the sterile lancet provided. Then, safely dispose of the lancet immediately after.

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
3


EN Wipe away the first drop of blood with sterile gauze or cotton.

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

4

EN Using a disposable inverted cup (5 µl) provided, dip the circular end of a inverted cup into the blood specimen.

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
5

EN Dispense 5 µl of drawn blood into round specimen well.

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EN Use Inverted cup : Let the circular end of the inverted cup touch the pad, then press down lightly.

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
6

EN Hold assay diluent bottle vertically and dispense 4 drops of assay diluent into the assay diluent well.
Exactly, 4 drops should be added. Do not let bottle tip touch device in order to avoid cross-contamination.

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7

EN Read the test result after 15 minutes.
⚠ **Caution:** Do not read test results after 30 minutes. ; reading after 30 minutes can yield erroneous results

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REACTIVE / RÉACTIF / REACTIVO / REATIVO

EN ⚠ **Caution :** The presence of any test line, no matter how faint, the result is considered reactive.

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
EN
1. **"C", "T1"** (Two lines) : HRP2 reactive
2. **"C", "T2"** (Two lines) : pLDH reactive
3. **"C", "T1"** and **"T2"** (Three lines) : pLDH and HRP2 reactive

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EN **"C"** and **"T3"** (Two lines) : pLDH of *P.v* reactive

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
EN
1. **"C", "T1"** and **"T3"** (Three lines)
2. **"C", "T2"** and **"T3"** (Three lines)
3. **"C", "T1", "T2"** and **"T3"** (Four lines)

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INVALID / NON VALIDE / NO VALIDO / INVÁLIDO

EN If the Control Line "C" is not visible within the result window after performing the test, the result is considered invalid. Instructions may not have been followed correctly or the test may have deteriorated beyond the expiration date. It is recommended that the specimen be retested using a new test device.

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Recommendations of the method:

1.- Use 75x13 tube or similar
2.- According to the procedure of the laboratory, the blood should mix with EDTA or 3Na Citrate (1:4).
3.- A well-mixed blood to the anticoagulant after drawing and before filling the pipette must be ensure. There should be a minimum of 12 complete inversions.
4.- The volume recommended must be between 2 and 3 ml in total, in order to ensure the correct filling of the pipette (other volumes are possible according to type of tube)
5.- The test must be set up within 4 hours after the drawing, if it would be later the specimen should be in the refrigerator at 4ºC up to 12 hours and the specimen should be return to room temperature and mixing again before using it.
6.- The room must be between 18º - 25ºC.
7.- Once filling the pipette (ensuring the pipette is bubble-free and filling exactly to "0" mark), places in a vertical position (90º) in an area free of vibrations, drafts and direct sunlight.
8.- At 60 minutes read the distance (a second reading can be done at 120 minutes, but it is depend on each laboratory procedure).

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
PRODUCT INFORMATION

Dimensions in mm

Volume	0.5 – 20 µl
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	bulk
Quality	Produced under cleanroom conditions acc. to ISO 14 644-1 class 8, CE-marked acc. to EC-Directive 98/79 EC (IVD).


Subject to technical modification without notice. Errors excepted.

Status as of: October 2021

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany |T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de

0.5 – 20 µl Pipette tips
732104 | 732204 | 732244

PRODUCT INFORMATION

Dimensions in mm

Volume	0,5 – 20 µl
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, CE-marked acc. to EC-Directive 98/79 EC (IVD).

CERTIFIED BIO-CERT® QUALITY

* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml


Subject to technical modification without notice. Errors excepted.

Status as of: October 2021

BRAND GMBH + CO KG | Otto-Schott-Str. 25 | 97877 Wertheim | Germany |T +49 9342 808 0 | F +49 9342 808 98000 | info@brand.de | www.brand.de

0.5 – 20 µl Pipette tips
732124 | 732224 | 732264

PRODUCT INFORMATION

Dimensions in mm

Volume	0,5 – 20 µl
Material	PP; autoclavable at 121 °C (2 bar) acc. to DIN EN 285
Color	transparent
Packaging	racked
Quality	Cleanroom class 8 according to ISO 14 644-1, Free of DNA*, RNase and DNase, Free of pyrogens, according to LAL test **, Sterile according to ISO 11 137, Free of ATP ,CE-marked acc. to EC-Directive 98/79 EC (IVD).

LIQUID HANDLING BIO-CERT® STERILE

* human and bacterial DNA, ** according to Limulus Amebocyte Lysate (LAL) test, detection limit 0.01 EU/ml

Subject to technical modification without notice. Errors excepted.

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Revision date: 2023-09-11

Number: 40003745
V01

English: Package Insert Page 1
A rapid, one step test for the qualitative detection of Human Occult Blood in feces.

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Package Insert

INTENDED USE

The FOB One Step Fecal Occult Blood Test Device (Feces) is a rapid chromatographic immunoassay for the qualitative detection of Human Occult Blood in feces.
A rapid, one step test for the qualitative detection of Human Occult Blood in feces.
For professional *in vitro* diagnostic use only.

SUMMARY

Many diseases can cause hidden blood in the feces. This is also known as Fecal Occult Blood (FOB), Human Occult Blood, or Human Hemoglobin. In the early stages, gastrointestinal problems such as colon cancer, ulcers, polyps, colitis, diverticulitis, and fissures may not show any visible symptoms, only occult blood. Traditional guaiac-based methods lack sensitivity and specificity, and also have diet restrictions prior to testing.[1,2]

The FOB One Step Fecal Occult Blood Test Device (Feces) is a rapid test to qualitatively detect low levels of Fecal Occult Blood. The test uses a double antibody sandwich assay to selectively detect Fecal Occult Blood at 50 ng/mL or higher, or 6 µg/g feces. In addition, unlike guaiac assays, the accuracy of the test is not affected by the diet of the patients.

PRINCIPLE

The FOB One Step Fecal Occult Blood Test Device (Feces) is a qualitative, lateral flow immunoassay for the detection of Human Occult Blood in feces. The membrane is pre-coated with anti-hemoglobin antibody on the test line region of the test. During testing, the specimen reacts with the particle coated with anti-hemoglobin antibody. The mixture migrates upward on the membrane chromatographically by capillary action to react with anti-hemoglobin antibody on the membrane and generate a colored line. The presence of this colored line in the test line region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region. If the control line does not appear, the test result is not valid.

REAGENTS

The test contains anti-hemoglobin antibody particles and anti-hemoglobin antibody coated on the membrane.

PRECAUTIONS

- Each test device is for single use only.
- Standard personal hygiene measures should be taken in the case of ingestion or direct eye contact with the buffer.
- Avoid touching finger directly to the specimen pad or result window of the test as this can cause incorrect results.
- For professional *in vitro* diagnostic use only. Do not use after expiration date.
- The test must remain in the sealed pouch until use.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Do not use test if pouch is damaged.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the testing and follow standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested.
- The used test should be discarded according to local regulations.
- Humidity and temperature can adversely affect results.
- Do not mix buffer from different lots.
- Need to read the results in the required time.
- Proper storage condition is critical to product performance.
- Once the foil pouch is opened, use the product as soon as possible.
- Test the product in the prescribed environmental condition.
- Use the correct specimen type.
- Do not recycle the test device.
- Maybe the bystander is contaminated with the biological material and sample in the testing process, so stay away from the bystander.

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date.

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- Specimens should not be collected during or within three days of a menstrual period, or if the patient suffers from bleeding hemorrhoids or blood in the urine.
- Alcohol, aspirin and other medications taken in excess may cause gastrointestinal irritation resulting in occult bleeding. Such substances should be discontinued at least 48 hours prior to testing.

- No dietary restrictions are necessary before using the FOB One Step Fecal Occult Blood Test Device (Feces).

MATERIALS

Materials Provided

- Test devices
- Specimen collection tubes with extraction buffer (0.05M NaCl)
- Package insert

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
- Specimen collection container
- Timer

DIRECTIONS FOR USE

Allow test, specimen collection tube, specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

1. To collect fecal specimens:
 Collect feces in a clean, dry specimen collection container. Best results will be obtained if the assay is performed within 2 hours after collection. Specimen collected may be stored for 1 days at 2-8°C if not tested within 2 hours.
 ***NOTE:** Considering the unexpected substances and/or microorganisms exist in feces may affect the stability of hemoglobin, it's best to test the feces as soon as possible.
2. To process fecal specimens:
 Unscrew the whole cap of the specimen collection tube, then randomly **stab the specimen collection stick into the fecal specimen in at least 3 different sites.** Do not scoop the fecal specimen.
 Screw on and tighten the cap onto the specimen collection tube, then **shake the specimen collection tube vigorously** to mix the specimen and the extraction buffer. Specimens prepared in the specimen collection tube may be stored for 3 days at room temperature (15-30°C) if not tested within 1 hour after preparation.
3. Remove the test device from the sealed pouch and use it as soon as possible. Do not touch the membrane of the strip.
4. Hold the specimen collection tube upright and then unscrew and open the upper cap. Invert the specimen collection tube and **transfer 2 full drops of the extracted specimen** (approx. 90 µL) to the specimen well of the test device, then start the timer. Avoid trapping air bubbles in the specimen well. See illustration below.
5. Wait for the colored line(s) to appear. **Read results at 5 minutes.** Do not interpret the result after 10 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

POSITIVE:* Two distinct colored lines appear. One line should be in the control line region (C) and another line should be in the test line region (T).
***NOTE:** The intensity of color in the test line region (T) will vary depending on the concentration of Fecal Occult Blood present in the specimen. Therefore, any shade in the test region indicates positive result.

NEGATIVE: One colored line appears in the control line region (C). No apparent colored line appears in the test line region (T).

INVALID: No line appears in the control line region (C). If this occurs, read the directions again and repeat the test with a new test. If the result is still invalid, stop using the test kit immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control.
Control standards are not supplied with this kit; however it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The FOB One Step Fecal Occult Blood Test Device (Feces) is for *in vitro* diagnostic use only.
2. The FOB One Step Fecal Occult Blood Test Device (Feces) will only indicate the presence of Fecal Occult Blood, the presence of blood in feces does not necessarily indicate colorectal bleeding.
3. Due to the high performance characteristic the FOB One Step Fecal Occult Blood Test Device (Feces) can obtain positive results when the specimen is below 50ng/mL in extraction buffer or 6µg/g in feces.
4. As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
5. Other clinically available tests are required if questionable results are obtained.

EXPECTED VALUES

The FOB One Step Fecal Occult Blood Test Device (Feces) has been compared with another leading commercial rapid test. The correlation between these two systems is 97.95% (95%CI:96.92% - 98.71%).

PERFORMANCE CHARACTERISTICS

Accuracy

The FOB One Step Fecal Occult Blood Test Device (Feces) has been compared with another leading commercial rapid test using clinical specimens.

Method		Other Rapid Test		Total Results
FOB One Step Test Device	Results	Positive	Negative	
	Positive	206	8	214
	Negative	14	847	861
Total Results		220	855	1,075

Relative Sensitivity: 93.6% (89.6% - 96.5%)* Relative Specificity: 99.1% (98.2% - 99.6%)*
Relative Accuracy: 97.95% (96.92% - 98.71%)* * 95% Confidence Intervals

Sensitivity

The FOB One Step Fecal Occult Blood Test Device (Feces) can detect levels of Fecal Occult Blood as low as 50 ng/mL or 6 µg/g feces.

Specificity

Cross Reactivity

The FOB One Step Fecal Occult Blood Test Device (Feces) is specific to human hemoglobin. Specimens containing the following substances were diluted in the extraction buffer to a concentration of 1.0 mg/mL, and tested on both positive and negative controls with no effect on test results: Bovine hemoglobin, Chicken hemoglobin, Pork hemoglobin, Goat hemoglobin, Horse hemoglobin, Rabbit hemoglobin and Turkey hemoglobin.

Interfering Substance

The FOB One Step Fecal Occult Blood Test Device (Feces) has been tested for specimens containing up to 20mg/dl Ascoribic acid, 60mg/dl Oxalic acid, 100mg/dl Bilirubin, 60mg/dl Uric acid, 20mg/dl Aspirin, 2,000mg/dl Urea, 2,000mg/dl Glucose and 40mg/dl Caffeine and all the results were found negative.

Precision

Intra-Assay

Within-run precision has been determined by using 10 replicates of three specimens: a negative, a low titer positive, and a high titer positive. The negative, low titer positive, high titer positive values were correctly identified >99% of the time.

Inter-Assay

Between-run precision has been determined by 10 independent assays on the same three specimens: The negative, low titer positive, high titer positive values. Three different lots of the FOB One Step Fecal Occult Blood Test Device (Feces) have been tested using negative, low titer positive and high titer positive specimens. The specimens were correctly identified >99% of the time.

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BIBLIOGRAPHY/ BIBLIOGRAFÍA / BIBLIOGRAPHIE / BIBLIOGRAFIA

1. Simon J.B. *Occult Blood Screening for Colorectal Carcinoma: A Critical Review*, Gastroenterology, Vol. 1985; 88: 820.
2. Blebea J. and Ncpherson RA. *False-Positive Guaiac Testing With Iodine*, Arch Pathol Lab Med, 1985;109:437-40.

Index of Symbols/ Índice de Símbolos/ Liste des Symboles / Índice de símbolos

Consult instructions for use Consulte las instrucciones de uso Consulter la notice d'utilisation Consulte as instruções de uso	Batch code Número de lote Code de lot Código do lote
Catalogue number Número de catálogo Référence catalogue Número de catálogo	Manufacturer Fabricante Fabricant Fabricante
Use-by date Fecha de caducidad Date d'expiration Usar até	Do not reuse No reutilizar Usage unique Não reutilizar
Contains sufficient for <n> tests Contenido suficiente para realizar <n> pruebas Contenu suffisant pour <n> tests Contém o suficiente para <n> testes	Store between 2-30°C Almacenar entre 2-30°C Conserver entre 2-30°C Armazene entre 2 °C e 30 °C
In vitro diagnostic medical device Dispositivo médico para diagnóstico *In vitro* Dispositif médical de diagnostic *in vitro* Dispositivo médico de diagnóstico *in vitro*	
Authorized representative in the European Community Representante autorizado en la Comunidad Europea Représentant autorisé dans la communauté européenne Representante autorizado na comunidade europeia	

Abon Biopharm (Hangzhou) Co., Ltd.
#198 12th Street East,
Hangzhou Economic & Technological
Development Area,
Hangzhou, 310018, P.R.China
www.globalpointofcare.abbott

EC REP
MT Promedt Consulting GmbH
Ernst-Heckel-Straße 7
66386 St. Ingbert
Germany

© 2023 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

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Advice Line

For further information, please contact your distributor, or call to one of the following Abbott Product Support Care Centers:

Region	Phone	E-Mail Address
Europe & Middle East	+ (44) 161 483 9032	EMEproductsupport@abbott.com
Asia Pacific	+ (61) 7 3363 7711	APproductsupport@abbott.com
Africa, Russia & CIS	+ (972) 8 9429 683	ARCISproductsupport@abbott.com
Latin America	+ (57) 2 661 8797	LAproductsupport@abbott.com

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





EN

January 2020
241851/R9

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are deviations from the instructions in this package insert.

NAME AND INTENDED USE

Determine™ HIV Early Detect is an *in vitro,* visually read, qualitative immunoassay for the detection of antibodies (Ab) to HIV-1 and HIV-2 and the detection of non-immunocomplexed (free) HIV-1 p24 antigen (Ag) in human capillary and venous whole blood, plasma or serum. Therefore, any reactivity on 1) the Ab bar alone 2) the Ag bar alone or 3) both the Ab and Ag bars simultaneously, is considered a reactive result suggestive of infection with HIV. The test is for professional use only and intended to be used as an aid in the diagnosis of HIV-1/HIV-2 infection in adults and children over the age of 18 months.

SUMMARY AND EXPLANATION OF THE TEST

AIDS (Acquired Immunodeficiency Syndrome) is characterized by changes in the population of T-cell lymphocytes. In an infected individual, the virus causes depletion of CD4 helper T-cells, which leaves the person susceptible to opportunistic infections and some malignancies. The virus that causes AIDS exists as two related types known as HIV-1 and HIV-2. The multiplication of the HIV in the infected cells leads to cell rupture and thus the release of HIV virus particles, which are first detected in the form of HIV RNA and next in the form of HIV antigen.[1,2] This is followed by production of specific antibodies to either HIV-1 or HIV-2.[3,4,5] HIV antigen may become undetectable at this time because of the formation of antibody-antigen complexes.

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

Determine HIV Early Detect is an immunochromatographic test for the qualitative detection of free HIV-1 p24 antigen and antibodies to HIV-1 and HIV-2.
Specimen is added to the sample pad. The specimen mixes with biotinylated anti-p24 antibodies and selenium colloid-conjugates coated with recombinant HIV-1, HIV-2 and HIV-1 group O antigens, synthetic HIV-2 peptide and anti p24 mouse monoclonal antibody. This mixture continues to migrate through the solid phase to the immobilized recombinant HIV-1/HIV-1 group O antigens and synthetic HIV-1/HIV-2 peptides at the Antibody (Ab) window, immobilized avidin at the Antigen (Ag) window.
If antibodies to HIV-1 and/or HIV-2 are present in the specimen, the antibodies bind to the selenium colloid-conjugates coated with recombinant HIV-1, HIV-2 and HIV-1 group O antigens and synthetic HIV-2 peptide and to the immobilized recombinant HIV-1/HIV-1 group O antigens and synthetic HIV-1/HIV-2 peptides, forming one red bar at the Ab window site. If antibodies to HIV-1 and HIV-2 are absent, the selenium colloid-conjugates flow past the Ab window and no red bar is formed at the Ab window site.
If free HIV-1 p24 antigen is present in the specimen, the antigen binds to the biotinylated anti-p24 antibodies and the selenium colloid-conjugate coated with anti p24 mouse monoclonal antibody. This complex binds to the immobilized avidin forming a red bar at the Ag window site. If HIV-1 p24 antigen is not present, both the biotinylated anti-p24 antibodies and selenium colloid-conjugate flow past the Ag window and no red bar is formed at the Ag window site.
To ensure assay validity, a procedural control bar is incorporated in the assay device at the Control window.

CONTENTS

Determine HIV Early Detect 20 Test (7D2842) or 100 Test (7D2843)

- **Determine** HIV Early Detect Test Card, 2 or 10 cards (containing 10 tests/card) coated with HIV-1/2 recombinant antigen, synthetic peptides, anti-p24 antibodies and avidin.

Accessories (required but not provided)

For testing Whole Blood samples

- Chase Buffer (7D2243) 1 Bottle (2.5 mL) prepared in phosphate buffer. Preservatives: Antimicrobial Agents.

For testing Whole Blood samples (fingerstick assay)

- EDTA Capillary Tubes (7D2222)

Determine HIV Early Detect SET (7D2843SET) (100 Test for testing whole blood samples)

- **Determine** HIV Early Detect (7D2843)
- Chase Buffer (7D2243)
- EDTA Capillary Tubes (7D2222)
- Blood Lancet (7D2232 or 7D2233)

Materials Required But Not Provided

- Disposable gloves, timing device
- Micropipette capable of delivering 50 µL (other than fingerstick)
- Alcohol swab, gauze pad, lancet (for fingerstick)

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use.
For professional use only.
Safety data sheet available for professional user on request.

CAUTION: When handling specimens and reagents, use appropriate biosafety practices[7,8]. These precautions include, but are not limited to the following:

- Wear gloves.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens or reagents using suitable disinfectant, such as 0.5% sodium hypochlorite, or other suitable disinfectant.[7,8]
- Decontaminate and dispose of all specimens, used test strips, and other potentially contaminated materials in accordance with local regulations.[7,8]

STORAGE

Store **Determine** HIV Early Detect Test Cards and Chase Buffer at 2-30°C until expiration date.

- When handled and stored as directed, kit components are stable until the expiration date. Do not use kit components beyond expiration date.
- Immediately reseal all unused tests in the foil pouch containing the desiccant by pressing seal from end to end to close.
- Do not use devices that have become wet or if the packaging has become damaged.

SPECIMEN COLLECTION

Serum, Plasma, and Whole Blood Collection by Venipuncture.

Use EDTA collection tubes for whole blood and plasma specimens.

- Collect human serum, plasma and whole blood by venipuncture aseptically.
- To obtain serum, separate from the clot. To obtain plasma, separate from the packed cells. Separate specimens as soon as possible to avoid any hemolysis.


Whole Blood Collection by Fingerstick[6] (See Fig.1)

Use EDTA Capillary Tubes (7D2222).

CAUTION: Glass capillaries may be damaged during transportation or when in use. Handle with care in order to avoid injury when removing from the package as well as during use and during disposal.
Before collecting a fingerstick specimen, place a capillary tube on a clean dry surface.

1. Choose the fingertip of the middle, ring, or index finger (whichever is the least callused). Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.
2. Clean fingertip with alcohol; allow to air dry.
3. Position the hand palm-side up. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.
4. Wipe away the first drop of blood with a sterile gauze pad.
5. Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times. Touch the tip of the capillary tube to the drop of blood. Avoid air bubbles.
 Fill the tube with blood between the 2 lines marked on the capillary tube (50 µL).

Fig.1

SPECIMEN STORAGE

- Store serum and plasma specimens at 2-8°C and run the test within 7 days of collection. If testing is delayed more than 7 days, freeze the specimen at -20°C or colder.
- Avoid repeated freeze/thaw cycles.
- If serum or plasma specimens show particulate matter or turbidity, centrifuge at 10,000g for 5 minutes at room temperature before sampling. Carefully remove the 50 µL test sample from the supernatant. If a lipid layer is formed on the surface of the liquid, ensure to take the sample from the clear liquid below that layer.
- For whole blood collected by venipuncture, store at 2-8°C. Do not freeze whole blood specimens. Run the test within 7 days of collection.
- For whole blood collected by fingerstick, test immediately.

TEST PROCEDURE

Remove the desired number of test strips from the 10-test card by bending and tearing at the perforation.
The test is to be performed at 15-30°C.

NOTE:

- To preserve the lot number which appears on the left side of the test card, remove individual test strips starting from the right side of the test card. The lot number and expiry date are not printed on the individual test strips.
- After removing the protective foil cover from each test strip, start the assay immediately. At under 80% or less relative humidity, the test can be used within 30 minutes after removing the protective foil cover.
- For serum or plasma, ensure thorough mixing of sample prior to use. For whole blood sample, mix well by gentle inversion of the tube immediately before testing.
- Running the test in high temperature/low humidity may affect the appearance of the Ag or Ab bar. If the test strip is partially dried and difficult to read at 20 to 40 minutes, the test should be repeated using a new test strip and result read at 20 minutes. When the test strip is partially dried, it appears as mixed white spot and grayish area.
- If serum or plasma sample does not flow or shows abnormal flow, such as stopping in the middle of the window, centrifuge the specimen and repeat the test with a new test strip.

CAUTION: Abnormal flow may occur with a whole blood (fingerstick) sample, if the capillary tube is not placed in the middle of the sample pad at an upright (vertical) position. If this occurs, collect a new sample and repeat the test using a new strip.

1. Remove the protective foil cover from each test.
2. For serum or plasma samples:
 a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 b. Wait a minimum of 20 minutes from addition of the sample (up to 40 minutes maximum) and read result.
3. For whole blood (venipuncture) samples:
 a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
 b. Wait one minute to allow the sample to be absorbed, then apply one drop of Chase Buffer to the sample pad, holding the bottle vertically.
 c. Wait a minimum of 20 minutes from addition of the sample (up to 40 minutes maximum) and read result.
4. For whole blood (fingerstick) samples using an EDTA Capillary Tube:
 a. Place the capillary tube containing the blood sample onto the middle of the sample pad (marked by the arrow symbol) at an upright (vertical) position.
 b. Wait until all the blood is transferred from the capillary tube to the sample pad. Then immediately apply one drop of Chase Buffer to the sample pad, holding the bottle vertically.
 Caution: do not lift the capillary tube from the sample pad before all the blood has been transferred – a bubble may form which will prevent the complete transfer of sample and invalidate the test. It may take more than one minute for full transfer of the sample.
 c. Wait a minimum of 20 minutes from addition of the sample (up to 40 minutes maximum) and read result.

QUALITY CONTROL


To ensure assay validity, a procedural control is incorporated in the device and is labeled "C". If the control bar does not appear by assay completion, the test result is invalid. Repeat the test using a new test strip.

INTERPRETATION OF RESULTS (See pictures)

NOTES:


- Any reactivity on 1) the Ab bar alone 2) the Ag bar alone or 3) both, the Ab and Ag bars simultaneously, is considered a reactive result suggestive of infection with HIV.
- Interpret any visible red bar (even very faint) in the control window as a valid result.
- The test result is reactive even if the AG or AB bars appear lighter or darker than the control bar. Any visible red bar, no matter how faint, is interpreted as reactive.
- A test which gives very high background should be considered invalid.
- If an invalid result occurs repeatedly, or for technical assistance, contact your local distributor or call Technical Support.

ANTIBODY REACTIVE
(Two Bars - Control and AB Bars)
Red bars appear in both the control window (labeled "C") and in the Ab window (labeled "AB") of the strip. Interpret any visible red bar, no matter how faint, in the Ab window as reactive.


Antibody Reactive

ANTIGEN (p24) REACTIVE (Two Bars - Control and AG Bars)
Red bars appear in both the control window (labeled "C") and in the Ag window (labeled "AG") of the strip. Interpret any visible red bar, no matter how faint, in the Ag window as reactive. The presence of only an antigen response suggests that the infection is at an early stage. Follow up testing may be suggested in order to track the expected future detection of antibodies.


Antigen Reactive

ANTIBODY REACTIVE AND ANTIGEN (p24) REACTIVE (Three Bars - Control, AB and AG Bars)
Red bars appear in the control window (labeled "C"), the Ab window (labeled "AB") and the Ag window (labeled "AG") of the strip. Interpret any visible red bar in the Ab and Ag windows as reactive.


Antigen and Antibody Reactive

NON-REACTIVE (One Bar - Control Bar)
One red bar appears in the control window of the strip (labeled "C"), and no red bar appears in the Ab and Ag windows of the strip (labeled "AG" and "AB").

Non-Reactive

INVALID (No Control Bar)
If there is no red bar in the control window of the strip, and even if a red bar appears in the Ab or Ag window of the strip, the result is invalid. Repeat the test using a new test strip.

Invalid

LIMITATIONS OF THE PROCEDURE

- **Determine** HIV Early Detect is designed to detect antibodies (Ab) to HIV-1 and/or HIV-2 and non-immunocomplexed (free) HIV-1 p24 antigen (Ag), in human serum, plasma and capillary and venous whole blood specimens. Other body fluids or pooled specimens may not give accurate results and should not be used.
- The intensity of the AB and AG bars does not necessarily correlate to the titer of antibody and antigen in the specimen, respectively.
- A reactive result for antibodies to HIV-1/2 combined with a non-reactive result for HIV-1 p24 antigen does not preclude the possibility of acute HIV infection.
- Reactivity for both antibody and antigen lines or for antigen reactivity only may be confirmed by 4th generation EIA or NAT where the IVD is indicated as an aid for diagnosis. Antibody reactivity only can also be confirmed by using another HIV Ab RDT. If non-reactive in subsequent testing, retest at least 14 days after the original test is recommended. Follow the local guideline for confirmatory methods if applicable.
- No test provides absolute assurance that a specimen does not contain low levels of HIV-1 p24 antigen and/or antibodies to HIV-1 and HIV-2 such as those present at a very early stage or late stage of HIV infection.
- A non-reactive result for both antibodies to HIV-1/2 and HIV-1 p24 antigen does not preclude the possibility of exposure to or infection with HIV-1 or HIV-2 viruses.
- The absence of AG bar may occur when all p24 antigen is bound by antibodies. When high levels of antibodies against the p24 antigen are present in the blood after seroconversion, the antibodies tend to bind to the antigens, forming immunocomplexes. **Determine** HIV Early Detect detects only non-immunocomplexed (free) antigens; it does not detect immunocomplexed (bound) antigens.
- Some known HIV-infected persons taking antiretroviral medication have been shown to produce false negative results when tested by rapid diagnostic tests.[10,11,12]
- Where clinical presentation is inconsistent with the non-reactive **Determine** HIV Early Detect test result, the individual should be re-tested for antibodies to HIV with **Determine** HIV Early Detect at least 14 days after the original test, or with a NAT assay that is indicated for diagnosis.
- Whole blood or plasma specimens containing anticoagulants other than EDTA have not been validated for use with the **Determine** HIV Early Detect and may give incorrect results.
- Infants born to HIV-infected mothers may carry maternal antibodies and will test antibody positive until eighteen months of age, which may not necessarily indicate the true infection status of the new born. The use of HIV- 1 p24 antigen testing to exclude infection in neonates (up to around eighteen months) is not recommended by CDC, because of poor sensitivity, especially in the presence of HIV antibody. Definitive diagnosis of HIV infection in early infancy requires other assays, including HIV nucleic acid test or viral culture.[9]
- Biotin treatment higher than 20mg per day may lead to decreased Ag bar intensity. Biotin concentrations up to 200 ng/mL in serum or plasma did not impact the sensitivity. There is no impact to Ab bar by biotin.

PERFORMANCE CHARACTERISTICS
SENSITIVITY
Sensitivity was evaluated by testing confirmed HIV antibody positive specimens, commercial seroconversion panels and HIV-1 subtype virus panels.

1. HIV Antibody Positive Specimens:

Table I
HIV Antibody Positive Specimens

Types	Number of Specimens Tested	Reactive by **Determine** HIV Early Detect	Sensitivity
HIV-1	422	422	100.0%
HIV-1 non B subtypes *	56	56	100.0%
HIV-1 group O	4	4	100.0%
HIV-2	100	100	100.0%
HIV**	100	100	100.0%
Total	**682**	**682**	**100.0%**

* Subtypes: A, C, D, F, G, H, J, K, and CRF01-AE, CRF02-AG, CRF03-AB, CRF05-DF, CRF09-A/U, CRF11-cpx
** Specimens collected in Europe

A total of 682 confirmed HIV positive specimens were tested (Table I). The antibody diagnostic sensitivity of **Determine** HIV Early Detect on this population of specimens was calculated to be 100%.

2. HIV Seroconversion Panels: (See Table II on the last page.)
The sensitivity of the **Determine** HIV Early Detect was evaluated using 31 sets of seroconversion panels; each including early seroconversion panel members (98 early samples in total). The results were compared with the results of commercially available CE marked HIV 4th generation rapid immunochromatographic test kit, CE marked HIV 3rd generation rapid immunochromatographic test kit, HIV antibody EIA kit, HIV antigen EIA kit. Additionally, the results of 20 seroconversion panels were compared with the data of CE marked HIV 4th generation kit (Chemiluminescent microparticle immunoassay (CMIA)).
On 28 of 31 panels, **Determine** HIV Early Detect detected HIV infection earlier than a 3rd generation rapid immunochromatographic test kit. On 11 of 31 panels, **Determine** HIV Early Detect detected HIV infection earlier than a rapid 4th generation immunochromatographic test kit.
Detection of six of 20 panels by **Determine** HIV Early Detect was delayed by 1 bleed date when compared to a 4th generation kit (CMIA).
Detection of eight of 31 panels by **Determine** HIV Early Detect was delayed by 1 bleed date when compared to HIV Antigen Kit (EIA).

3. HIV-1 p24 Subtype Antigens:
50 HIV virus panels were prepared by diluting 50 different cultured HIV isolates including HIV-1 group M subtypes: A(4), B(8), C(7), D(5), F(5), G(2), H(1), and CRF01-AE(10), CRF02-AG(2), HIV-1 group O(4) and HIV-2(2) with HIV negative human serum. The **Determine** HIV Early Detect detected all HIV-1 panels and did not detect HIV-2 panels.

4. Analytical Sensitivity of HIV-1 p24 antigen
The analytical sensitivity of the **Determine** HIV Early Detect was evaluated by testing WHO International Standard HIV-1 p24 Antigen (NIBSC code 90/636). The results demonstrated that the test could detect a concentration of 2 IU/mL HIV-1 p24 antigen.

SPECIFICITY
A total of 2469 confirmed negative serum or plasma specimens were tested with the **Determine** HIV Early Detect and the specificity was **Determine**d for the antibody (Ab) bar and for the antigen (Ag) bar (Table III). The specificity was 99.96% for the antibody (Ab) bar and 99.76% for the antigen (Ag) bar.

Table III
Specificity of Determine HIV Early Detect

Population	Number of Tested Specimens	**Determine** HIV Early Detect					
		Ab bar		Ag bar		Combined Ab and Ag bar	
		Reactive	Non-reactive	Reactive	Non-reactive	Reactive	Non-reactive
Seronegative specimens*	1749	0	1749	2	1747	2	1747
Hospitalized specimens	218	1	217	2	216	3	215
Pregnant women	206	0	206	0	206	0	206
potentially cross-reacting specimens**	296	0	296	2	294	2	294
Total	2469	1	2468	6	2463	7	2462
Specificity (%)		99.96%		99.76%		99.72%	

* Including specimens collected in Europe (300), USA (1299) and Africa (150)
** Disease states other than HIV, and potentially interfering substances: Rheumatoid factor, antinuclear antibody, systemic lupus erythematosus, high cholesterol, high total protein, high IgM, human anti mouse IgG, other infections (HBV, HCV, HTLV I or II, CMV, Toxo IgG, Syphilis, HSV 1/2 and EBV), Flu vaccinated patients, and dialysis patient specimens

SAMPLE TYPE

Table IV
Antibody sensitivity in matched whole blood (venipuncture and fingerstick), serum and plasma specimens

No. of matched specimens tested	Type of Specimens and No. of Reactive by **Determine** HIV Early Detect (Ab detection)				Correlation between matrices
	Serum	Plasma	Whole Blood venipuncture	Whole Blood fingerstick*	
25	25	25	25	NT	100.0%
25	25	25	25	25	100.0%

NT: not tested *EDTA capillary tubes

Multiple (matched) specimens: Seropositive specimens from a total of 50 individuals from Africa were tested with the **Determine** HIV Early Detect (Table IV). Multiple (matched) specimens, 50 serum specimens, 50 plasma specimens, 50 whole blood (venipuncture) specimens and 25 whole blood (fingerstick) specimens were obtained from the donors in Africa. The results obtained from all specimen matrices showed 100% correlation, demonstrating that **Determine** HIV Early Detect gives identical results for these types of specimen matrices.

Table V
Specificity in matched whole blood (venipuncture and fingerstick), serum and plasma specimens

No. of matched specimens tested	Type of Specimens and No. of Non-reactive by **Determine** HIV Early Detect				Correlation between matrices
	Serum	Plasma	Whole Blood venipuncture	Whole Blood fingerstick*	
25	25	25	25	NT	100.0%
25	25	25	25	25	100.0%

NT: not tested *EDTA capillary tubes

Multiple (matched) specimens: Seronegative specimens from a total of 50 individuals from Africa were tested with the **Determine** HIV Early Detect (Table V). Multiple (matched) specimens, 50 serum specimens, 50 plasma specimens, 50 whole blood (venipuncture) specimens and 25 whole blood (fingerstick) specimens were obtained from the donors in Africa. The results obtained from all specimen matrices showed 100% correlation, demonstrating that **Determine** HIV Early Detect gives identical results for these types of specimen matrices.

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© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

BIBLIOGRAPHY/BIBLIOGRAFÍA/BIBLIOGRAPHIE/ BIBLIOGRAFIA

1. Pilcher C, Eron JJ, Galvin S, *et al*. (2004) Acute HIV revisited: new opportunities for treatment and prevention. *The Journal of Clinical Investigations* 113(7): 937-945.
2. Respess RA, Rayfield MA and Dondero TJ (2001) Laboratory testing and rapid HIV assays: applications for HIV surveillance in hard-to-reach populations. *AIDS* 15 Supplement 3: S49-S59.
3. Louie B, Wong E, Klausner JD, *et al.* (2008) Assessment of rapid tests for detection of Human Immunodeficiency Virus-specific antibodies in recently infected individuals. *Journal of Clinical Microbiology* 46(4): 1494-1497.
4. Lyons MS, Lindsell CJ, Hawkins DA, et al. (2008) Contributions to early HIV diagnosis among patients linked to care vary by testing venue. *BMC Public Health* 8:220
5. Gürtler LG, Hauser PH, Eberle J, et al. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. *Journal of Virology*. 1994; 68(3):1581-1585.
6. Clinical and Laboratory Standards Institute. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition GP42-A6 Vol.28 No.25 September 2008
7. Clinical and Laboratory Standards Institute. Clinical Laboratory Waste Management; Approved Guideline –Third Edition. GP05-A3 Vol.31 No.3 January 2011
8. EPA Guide for Infections Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1 ~ 5-5, R1-R3, A1-A24.
9. CDC: 2008 Case Definition; Human Immunodeficiency Virus Infection.
10. Delaney KP, Branson BM, Uniyal A, et al. Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. *Clinical Infectious Diseases*. 2011; 52(2): 257-263.
11. O'Connell RJ, Merritt TM, Malia JA, et al. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
12. O'Connell RJ, Agan BK, Anderson SA, et al. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.

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ONE-STEP Influenza Virus A, B and A(H1N1) Pandemic Rapid Test

SD BIO LINE Influenza Ag A/B/A(H1N1) Pandemic (09FK53, 19FK32)

Test rapide en une étape de détection du virus de la Grippe de Type A, B et de la souche pandémique A(H1N1).
Prueba Rápida de un paso para Influenza Tipo A,B y A(H1N1) Pandémica
Teste rápido de vírus Influenza tipo A, B e A(H1N1) Pandêmica – passo único
Einstufen-Schnelltest zum Nachweis des Influenzavirus Typ A und B sowie des pandemischen Virus A(H1N1)
Test rapid per la determinazione del virus dell'influenza di tipo A, B e (H1N1) pandemico

ENGLISH

About the test

[Introduction] In April 2009, a new strain of influenza A virus, subtype H1N1 was first detected in the United States. Influenza A (H1N1) or "swine flu" continued to spread globally, and on June 11, 2009, the WHO declared the outbreak a pandemic.
[Intended use] The **SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic** test kit is a chromatographic immunoassay for the differential and qualitative detection of influenza virus A, B and A (H1N1) antigens directly from nasal/throat/nasopharyngeal swab or nasal/nasopharyngeal aspirates/lavages/nasal wash specimens. This assay is intended for professional use and for *in vitro* diagnostic use only.
[Test principle] The **SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic** test strips are immobilized with mouse monoclonal anti-A(H1N1) hemagglutinin, anti-Influenza A and anti-Influenza B antibodies. The specially selected antibodies are used as detector materials and enable the identification of influenza virus A, B and A (H1N1) antigens directly and with a high degree of accuracy.

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The **SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic** test kit contains the following items to perform the assay:

1. SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic (10 Tests/Kit) Cat. No. 19FK31:
 10 test strips with desiccant in individual foil pouches
 10 Disposable Extraction Tubes
 Assay diluent (1 x 5ml/vial)
 10 Disposable Droppers
 10 Sterilized swabs for sample collection
 3 Control swabs: 1 Influenza A positive control swab, 1 Influenza B positive control swab, and 1 Negative control swab
 1 Instructions for use
2. SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic (25 Tests/Kit) Cat. No. 19FK32:
 25 test strips with desiccant in individual foil pouches
 25 Disposable test tubes
 Assay diluent (1 x 15ml/vial)
 25 Disposable Droppers
 25 Sterilized swabs for sample collection
 3 Control swabs: 1 Influenza A positive control swab, 1 Influenza B positive control swab, and 1 Influenza negative control swab
 1 Instructions for use

Active ingredients of main components

1. test strip includes:
 Gold conjugates: Anti-Chicken IgY-conjugated gold colloid (0.027±0.0054μg), Mouse monoclonal anti-A (H1N1) Pandemic hemagglutinin-conjugated gold colloid (0.101±0.012μg), Mouse monoclonal anti-Influenza A virus conjugated-gold colloid (0.081±0.016μg), Mouse monoclonal anti-Influenza B virus conjugated gold colloid (0.081±0.016μg)
 Test Line 'A (H1N1) Pandemic': Mouse monoclonal anti-A (H1N1) Pandemic hemagglutinin (1.068±0.214μg)
 Test Line 'A': Mouse monoclonal anti-Influenza A virus (1.068±0.214μg)
 Test Line 'B': Mouse monoclonal anti-Influenza B virus (0.212±0.142μg)
 Control Line: Mouse monoclonal anti-Chicken IgY (0.554±0.107μg)
2. Assay diluent includes: Tricine (0.4M), NaCl (q.s.), Triton X-100 (q.s.), sodium azide (0.02 w/v%)
3. Control swabs:
 Influenza A positive control swab: Recombinant haemagglutinin of influenza virus A (H1N1) (0.1ng/ml), inactivated influenza virus A (H3N2) (0.1ng/ml)
 Influenza B positive control swab: Inactivated influenza virus B (0.1ng/ml)
 Influenza negative control swab: Formalin-inactivated *Streptococcus pyogenes* (5 x 10⁷ Org)

Kit storage and stability

1. The test kit should be stored at a temperature between 1 °C and 30 °C. Do not freeze the kit or its components.
2. The test strip is sensitive to heat and humidity.
3. Check the humidity indicator on the desiccant for color change and throw the pouch if the color indicates saturation (yellow > green).
4. Perform the test immediately after removing the test strip from the foil pouch.
5. Do not use the test kit beyond its expiration date.
6. The shelf life of the kit is as indicated on the outer package.
7. Do not use the test kit if the pouch is damaged or the seal is broken.

Specimen collection and handling

[Sample collection]

I. All swab samples
- Nasopharyngeal swab specimens
 To collect a nasopharyngeal swab specimen, insert the swab into nostril parallel to the palate and leave in place for a few seconds to absorb secretions.
 ***** **Note:** It has been reported that a nasopharyngeal swab containing cellular material is preferable to nasal swabs and throat swabs in specimen quality. It is especially important to collect a sufficient quantity of virus from epithelial cells from the nasopharyngeal (mucous) membrane.
- Nasal swab specimens
 To collect a nasal swab specimen, insert the sterile swab into the nostril that presents the most secretion under visual inspection. Using sterile rotation, push the swab until resistance is met at the level of the turbinate (less than 1 inch into the nostril). Rotate the swab a few times against the nasal wall.
- Throat swab specimens
 To collect a throat swab specimen, vigorously rub a rayon throat swab on both tonsillar surfaces and on the posterior pharynx. Remove the swab from the mouth and insert the swab, tip down, into a paper wrapper.

II. Nasal/nasopharyngeal aspirate specimens
- Aspirate or wash volumes between 2ml and 2.5ml are recommended.
- Excessive wash or aspirate volumes may result in decreased test sensitivity. Transfer the specimens into a clean, dry specimen container.
- Process all specimens as described in the "Test procedure" section.

[Specimen transport and storage]
1. Specimen should be tested as soon as possible after collection.
2. Transport fresh specimens in a suitable liquid transport system maintained on ice or refrigerated at a temperature between 2 °C and 8 °C up to 3 days after collection.
3. If transport media is not used, specimens may be refrigerated (2 - 8 °C), or stored at a temperature between 1 °C and 30 °C, in a clean, dry, closed container for up to 8 hours prior to testing.
4. If the specimen should be tested later, swab specimens should be placed into 1 - 2 ml of transport media or saline by direct extraction.
5. **Transport media :** Use of the following transport media has been tested and found to be compatible with the SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic test:
 Saline EMEM+1% BSA
 PBS EMEM+0.5% BSA
 PBS+0.5% BSA Tryptose soy broth+0.5% BSA
 PBS+0.5% gelatin Tryptose soy broth +0.5% gelatin

Test procedure (Refer to figure)

[Specimen extraction]

I. All swab samples not stored in transport media should be extracted by the following method:
1. Allow the test strip and extracted sample to reach a temperature between 15 °C and 30 °C prior to testing.
2. Hold the disposable dropper vertically, draw up assay diluent in the bottle to the Fill Line as shown in the following test procedure figure (approximately 300μl). Transfer 300μl of assay diluents into the test tube.
3. Insert the patient swab sample into the test tube. Swirl the swab head against the inside of the test tube at least five times while squeezing the tube. (Refer to figure)
 Remove the swab while squeezing the tube and then dispose of the used swab in accordance with your facility's biohazard waste disposal protocol.
 Nasal/nasopharyngeal aspirate specimens and swab samples to be placed in transport media by direct extraction:
 Directly pipet 100μl of extracted specimen or aspirates and 100μl of assay diluent into the test tube. Mix specimen.

[Reaction with test strip]
1. Remove the test strip from the foil pouch prior to use.
2. Place the test strip into the test tube with the arrows on the test strip pointing down. Do not handle or move the test strip until the test is complete and ready for reading.
 Read result at 10 - 15 minutes. Some positive results may appear rapidly.
 Caution: Do not read the results after 15 minutes. Reading too late can yield false results.

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1. Negative result: Only the control line (C) appears.
2. Positive result
 Positive for influenza virus A (H1N1) Pandemic (Two or three lines appear): One purple control line (C) and the (H1N1)Pandemic or two purple colored lines (A)(H1N1)Pandemic and A) and the control line (C) appear.
 Positive for influenza virus A (Two lines appear): One purple colored line (A) and the control line (C) appears.
 Positive for influenza virus B (Two lines appear): One purple colored line (B) and the control line (C) appears.
 Positive for influenza virus A and B (co-infection suspected) (Three or four lines appear): Two purple colored lines (A and A or B and A)(H1N1) Pandemic or three purple colored lines (A, B and A)(H1N1)Pandemic) and the control line (C) appear.
 Invalid result: The control line (C) fails to appear.
 If the control line is not visible within the result window after performing the test, the result is considered invalid. It is recommended that the specimen be retested using a new test kit.

Test limitations

1. The SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic test kit is to be used for the qualitative detection of influenza virus A, B and A(H1N1) antigens from nasal/throat/nasopharyngeal swab or nasal/nasopharyngeal aspirate specimens.
2. Failure to follow the instructions for test procedure and interpretation of test results may adversely affect test performance and/or invalidate the test.
3. Test results must be evaluated in conjunction with other clinical data available to the physician.
4. A negative test result may occur if the level of antigen in a sample is below the detection limit of the test or if a sample is collected improperly. Negative test results do not preclude the possibility of infection with other non-influenza viral infections.
5. The influenza A(H1N1) virus exhibits a high frequency of mutation. Specimens from patients infected with a viral mutation can exhibit a negative test result on the "A(H1N1)Pandemic" line, while still presenting a positive result on the "A" line.

Quality control

1. Internal quality control: The test strip has test lines and a control line on the surface of the test strip. Neither the test lines nor the control line is visible in the result window before applying a sample. The control line is used for procedural control and should always appear if the test procedure is performed properly and the test reagents of the control line are working.
 External quality control:
 Control procedures: Positive and negative control swabs should be tested according to the test procedure outlined in the [Specimen extraction] and [Reaction with test strips] sections of these instructions.
 Specification:
 - The influenza A positive control swab should be interpreted as test line "A" and "A (H1N1) Pandemic" positive. The influenza B positive control swab should be interpreted as test line "B" positive. The influenza negative control swab should be interpreted as negative.
 - **Note:** If the test result of the control swab is not as expected, repeat test with a new test strip and extracted control solution. If any problems arise that cannot be resolved by this, please contact your distributor or SD Technical Service team (XKproductsupport@alere.com).

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Warning:
The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential loss, liability, claim, cost or damages arising from or related to an incorrect positive or negative diagnosis using this product.

Performance characteristics

1. Study site 1
A total of 759 samples collected from Korea were tested on the SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic test.

[Influenza A(H1N1)]

Reference Assay	SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic		
	Positive	Negative	Total
Commercial Positive	241	72	313
RT-PCR Negative	1	445	446
Sensitivity (95% CI)	77.0% (72.8%-81.3%)		
Specificity (95% CI)	99.8% (98.7-100%)		

[Influenza B]

Reference Assay	SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic		
	Positive	Negative	Total
Commercial Positive	194	119	313
RT-PCR Negative	1	445	446
Sensitivity (95% CI)	62% (56.5-67.2%)		
Specificity (95% CI)	99.8% (98.7-100%)		

2. Study site 2
A total of 124 samples collected from Korea were tested on the SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic in Korea.

[Influenza A(H1N1)]

Test of reference	SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic		
	Positive	Negative	Total
Commercial Positive	2	0	2
RT-PCR Negative	0	39	39
Sensitivity (95% CI)	100% (34.2-100%)		
Specificity (95% CI)	100% (91-100%)		

[Influenza A]

Test of reference	SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic		
	Positive	Negative	Total
Commercial Positive	4	0	4
RT-PCR Negative	0	39	39
Sensitivity (95% CI)	100% (51-100%)		
Specificity (95% CI)	100% (91-100%)		

[Influenza B]

Test of reference	SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic		
	Positive	Negative	Total
Commercial Positive	67	12	79
RT-PCR Negative	0	39	39
Sensitivity (95% CI)	84.8% (75.3-91.1%)		
Specificity (95% CI)	100% (91-100%)		

3. Cross-reactivity: Cross-reactivity of the SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic test was evaluated with various concentrations of Bacterial 17 species, viral 15 species, other sub-type influenza virus 26 species (A: 22, B: 4). These studies demonstrated no cross-reactivity on SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic test.
 Results of analytical sensitivity of the SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic test are shown below:

Virus type	Strain
Influenza A virus(H1N1)	A/WS/33 (ATCC VR-1520)
Influenza A virus(H3N2)	A/Hong Kong/8/68 (ATCC VR-544)
Influenza B virus	B/Lee/40 (ATCC VR-101)
Influenza A(Pandemic)	A/08/01/2009
Influenza A(Pandemic)	A/California/7/2009 (H1N1)rv(NIBSC, NYMC-X179A)

Analytical sensitivity was established using a total of 5 human epidemic strains of influenza viruses (4 influenza A and 1 influenza B).
5. Reproducibility of the SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic test has been demonstrated by within-run, between-run and batch-to-batch studies using in-house reference panels. All values were identical to reference panel acceptance criteria.

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
1

EN Open the package and look for the following:
1. Test strip with desiccant in individual foil pouch
2. Disposable test tube
3. Assay diluent
4. Disposable dropper
5. Sterilized swab for Sample Collection
6. 3 Control swabs
7. Instructions for use

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2

EN First, carefully read the instructions for using the SD BIOLINE Influenza Ag A/B/A(H1N1) Pandemic test kit.

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3

EN Next, look at the expiration date on the back of the foil pouch. If the expiration date has passed, use another kit.

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



EN **Caution:** Do not read the results after 15 minutes. Reading too late can yield false results.

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




Positive / Positif / Postivo

A (H1N1) Pandemic type / A (H1N1) Pandemic tipo

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


INTENDED USE

The Panbio™ Dengue IgM Capture ELISA is for the qualitative detection of IgM antibodies to dengue antigen in serum, as an aid in the clinical laboratory diagnosis of patients with clinical symptoms consistent with dengue fever. The Panbio™ Dengue IgM Capture ELISA should be used in conjunction with other dengue serology.

INTRODUCTION

Dengue, a flavivirus, is found in large areas of the tropics and sub-tropics. Transmission is by mosquito, principally *Aedes aegypti* and *Aedes albopictus.* Dengue virus infection causes a spectrum of clinical manifestations ranging from unapparent to fatal haemorrhagic disease. Classic dengue or breakbone fever is characterised by the sudden onset of fever, intense headache, myalgia, arthralgia and rash. A diphasic febrile course is common as is insomnia and anorexia with loss of taste or bitter taste. Dengue haemorrhagic fever and dengue shock syndrome are severe complications often associated with a second serotype infection.

Detection of IgM antibodies to dengue virus by ELISA is a valuable procedure, particularly in second and subsequent infections where the occurrence of complications is high. Serum IgM antibodies can be detected from dengue patients as early as three to five days after the onset of fever and generally persist for 30 - 90 days, although detectable levels may be present eight months post-infection.

PRINCIPLE

Serum antibodies of the IgM class, when present, combine with anti-human IgM antibodies attached to the polystyrene surface of the microwell test strips. A concentrated pool of Dengue 1-4 Antigens is diluted to the correct working volume with Antigen Diluent. The antigens are produced using an insect cell expression system and immunopurified utilising a specific monoclonal antibody. An equal volume of the HRP Conjugated Monoclonal Antibody (MAb) is added to the diluted antigen, which allows the formation of antigen-MAb complexes. Residual serum is removed from the assay plate by washing, and complexed antigen-MAb is added to the assay plate. After incubation, the microwells are washed and a colourless substrate system, tetramethylbenzidine / hydrogen peroxide (TMB Chromogen) is added. The substrate is hydrolysed by the enzyme and the chromogen changes to a blue colour. After stopping the reaction with acid, the TMB becomes yellow. Colour development is indicative of the presence of anti-dengue IgM antibodies in the test sample.

MATERIALS PROVIDED

Note: 01PE21 = 01PE20 x 5

1. Antibody Coated Microwells - (12 x 8 wells) Microwells are coated with anti-human IgM antibodies. Ready for use. Unused microwells should be resealed immediately and stored in the presence of the desiccant. Stable at 2-8°C until expiry.
2. Dengue 1-4 Antigen (Recombinant) - One **Clear**-capped vial, 150 µL (**Blue**) concentrated dengue viral antigens 1, 2, 3 and 4 with preservative (Proclin™ 300 and gentamycin sulphate). Unused diluted antigen must be discarded. Concentrated antigen is stable at 2-8°C until expiry. ⚠
3. Wash Buffer (20x) - One bottle, 60 mL of 20x concentrate of phosphate buffered saline with Tween 20 and preservative (Proclin™ 300). Crystallisation may occur at low temperatures. To correct, incubate at 37°C until clear. Mix well. Dilute one part Wash Buffer with 19 parts of distilled water. Diluted buffer may be stored for one week at 2-25°C. ⚠

4. Sample Diluent - Two bottles, 50 mL (**Pink**). Ready for use. Tris buffered saline with preservatives (Proclin™ 300) and additives. Stable at 2-8ºC until expiry.
5. Antigen Diluent – One bottle, 50 mL (**Clear**). Ready for use. Phosphate Buffer containing preservatives (Proclin™ 300 and gentamycin sulphate). Stable at 2-8ºC until expiry.
6. HRP Conjugated Monoclonal Antibody - One bottle, 7 mL (**Yellow**). Ready for use. Horseradish peroxidase conjugated 6B6C-1 monoclonal antibody with preservative (Proclin™ 300) and conjugate diluent. Stable at 2-8ºC until expiry.
7. TMB Chromogen (TMB) - One bottle, 15 mL. Ready for use. Working TMB. Stable at 2-8ºC until expiry.
8. Positive Control - One **Black**-capped vial, 200 µL inactivated dengue positive human serum (contains sodium azide and gentamycin sulphate). Stable at 2-8ºC until expiry.
9. Calibrator - One **Orange**-capped vial, 400 µL inactivated dengue positive human serum (contains sodium azide and gentamycin sulphate). Stable at 2-8ºC until expiry.
10. Negative Control - One **White**-capped vial, 200 µL inactivated dengue negative human serum (contains sodium azide and gentamycin sulphate). Stable at 2-8ºC until expiry.
11. Stop Solution - One **Red**-capped bottle, 15 mL. Ready for use. 1M Phosphoric acid. Stable at 2-25ºC until expiry.

Proclin™ 300 is a registered trademark of Rohm and Haas Company.

* Classification according to Regulation (EC) No 1272/2008:

Product identifier	**Trade name**	Dengue 1-4 Antigen, Wash buffer, Sample diluent, Antigen Diluent , HRP conjugate
	Hazardous substance	Proclin 300 (5-chloro-2-methyl-4- isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H -isothiazol-3-one [EC no. 220-239-6] (3:1), CAS No. 55965-84-9)
Classification		Skin corrosion/irritation Category2 Serious eye damage/irritation Category2 Skin sensitization Category1
Hazard Pictogram		
Signal Word		Warning
Hazard Statement		H315: Causes skin irritation H317: May cause an allergic skin reaction H319: Causes serious eye irritation
Precautionary statements		
Prevention		P261: Avoid breathing dust/fume/gas/mist/vapours/spray P264: Wash hands thoroughly after handling P272: Contaminated work clothing should not be allowed out of the workplace P280: Wear protective gloves/protective clothing/eye protection/face protection

ENGLISH

Response	P302+P352: IF ON SKIN: Wash with plenty of soap and water P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.Continue rinsing P321: Specific treatment P332+P313: If skin irritation occurs: Get medical advice/attention P333+P313: If skin irritation or rash occurs: Get medical advice/attention P337+P313: If eye irritation persists: Get medical advice/attention P362: Take off contaminated clothing and wash before reuse P363: Wash contaminated clothing before reuse
Disposal	P501: Dispose of contents/container in accordance with local/regional/national/international regulation

ADDITIONAL MATERIALS REQUIRED BUT NOT PROVIDED

1. Accurate adjustable micropipettors with disposable pipette tips (5-1000 µL capacity)
2. Deionised water
3. Microplate washing system
4. Microplate reader with 450 nm filter
5. Timer
6. Graduated cylinder
7. Flask
8. Test tubes or microtitre plate for serum dilutions
9. Glass or plastic tubes or vials for diluting antigen

WARNINGS

1. For *in vitro* diagnostic use.
2. All human source material used in the preparation of controls has been tested for antibody to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen and found to be negative. However no test method can offer complete assurance and all human controls and antigen should be handled as potentially infectious material. The Centers for Disease Control and Prevention and the National Institutes of Health recommend that potentially infectious agents be handled at the Biosafety Level 2[1].
3. This test should be performed on serum only. The use of whole blood, plasma or other specimen matrix has not been established.
4. Icteric or lipaemic sera, or sera exhibiting haemolysis or microbial growth should not be used.
5. Do not heat-inactivate sera.
6. All reagents must be equilibrated to room temperature (20-25ºC) before commencing the assay. The assay will be affected by temperature changes. Do not remove microwells from closed bag until they have reached room temperature (20-25ºC).
7. Dispense reagents directly from bottles using clean pipette tips. Transferring reagents may result in contamination.
8. Unused microwells should be resealed immediately and stored in the presence of desiccant. Failure to do this may cause erroneous results.
9. Substrate System:
 (a) As TMB is susceptible to contamination from metal ions, do not allow the substrate system to come into contact

with metal surfaces.

(b) Avoid prolonged exposure to direct light.

(c) Some detergents may interfere with the performance of the TMB.

(d) The TMB may have a faint blue colour. This will not affect the activity of the substrate or the results of the assay.

10. Some kit components contain sodium azide, which may react with lead or copper plumbing to form highly explosive metal azide compounds. When disposing of these reagents through plumbing fixtures, flush with a large volume of water to prevent azide build-up in drains.
11. Sodium azide inhibits conjugate activity. Clean pipette tips must be used for the conjugate addition so that sodium azide is not carried over from other reagents.

SPECIMEN COLLECTION AND PREPARATION

Blood obtained by venipuncture should be allowed to clot at room temperature (20-25°C) and then centrifuged according to the Clinical and Laboratory Standards Institute (CLSI – Approved Standard - Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture, H3).

The serum should be separated as soon as possible and refrigerated (2-8°C) or stored frozen (≤-20°C) if not tested within two days. Self-defrosting freezers are not recommended for storage. The use of icteric sera or sera exhibiting haemolysis, lipaemia or microbial growth is not recommended. The CLSI provides recommendations for storing blood specimens, (Approved Standard - Procedures for the Handling and Processing of Blood Specimens, H18).

TEST PROCEDURE

Note: Ensure all reagents are equilibrated to room temperature (20-25°C) before commencing assay. Performing the assay outside the time and temperature ranges provided may produce invalid results. Assays not falling within the established time and temperature ranges must be repeated.

Serum Predilution

1. Remove the required number of microwells from the foil sachet and insert into strip holder. Five microwells are required for Negative Control (N), Positive Control(P) and Calibrator (CAL) in triplicate. Ensure the remaining unused microwells are sealed tightly in the foil sachet.
2. Using suitable test tubes or a microtitre plate, dilute the Negative Control, Positive Control, Calibrator and patient samples:

 (a) To 10 µL serum add 990 µL of Sample Diluent. Mix well.

 Alternatively,

 (b) To 10 µL serum add 90 µL of Sample Diluent. Take 20 µL of the diluted serum and add 180 µL Sample Diluent. Mix well.


ELISA PROCEDURE

See attached figure for summary of method.

(A) Antigen

1. Determine the required number of wells for your assay. Dilute the antigen 1/250 using the Antigen Diluent. It is recommended, as a minimum, to dilute 10 µL of antigen into 2.5 mL of Antigen Diluent. This is sufficient for up to five strips (40 wells). A volume of 0.5 mL of diluted antigen is required per strip. Once the antigen is added to the Antigen Diluent, the solution becomes a pale blue colour. Ensure the remaining unused concentrated antigen remains at 2-8ºC.
2. Remove the required volume of diluted antigen and mix with an equal volume of MAb Tracer in a clean glass or polypropylene vial. Gently mix the antigen-MAb Tracer solution and leave at room temperature (20-25ºC) until required. Discard the unused diluted antigen.

(B) Assay Plate

1. Within 10 minutes after mixing the MAb Tracer and diluted antigen, pipette 100 µL diluted patient sample and Controls into their respective microwells of the assay plate.
2. Cover the plate and incubate for 1 hour at 37 ±1ºC.
3. Wash six (6) times with diluted Wash Buffer. (Refer to washing procedure).
4. Mix the antigen-MAb tracer solution before transfer. Pipette 100 µL of antigen-MAb complexes from the antigen vial to the appropriate wells of the assay plate.
5. Cover plate and incubate for 1 hour at 37 ±1ºC.
6. Wash six (6) times with diluted Wash Buffer (refer to washing procedure).
7. Pipette 100 µL TMB into each well.
8. Incubate for 10 minutes at room temperature (20-25ºC), timing from the first addition. A blue colour will develop.
9. Pipette 100 µL of Stop Solution into all wells in the same sequence and timing as the TMB addition. Mix well. The blue colour will change to yellow.
10. Within 30 minutes read the absorbance of each well at a wavelength of 450 nm with a reference filter of 600-650 nm.

Note: If a dual wavelength spectrophotometer is available, set the reference filter between 600-650 nm. Reading the microwells at 450 nm without a reference filter may result in higher absorbance values due to background.

WASHING PROCEDURE

Efficient washing to remove uncomplexed sample or components is a critical requirement of the ELISA procedure.

A. Automated Plate Washer

(1) Completely aspirate all wells.
(2) Fill all wells to rim (350 µL) during wash cycle.
(3) On completion of six (6) washes, invert plate and tap firmly on absorbent paper towel to ensure all Wash Buffer is removed.
(4) Automated plate washers must be well maintained to ensure efficient washing. Manufacturer's cleaning instructions should be followed at all times.

B. Manual Washing

(1) Discard contents of plate in appropriate waste container.
(2) Fill wells with Wash Buffer using a suitable squeeze bottle. Avoid bubbling of Wash Buffer as this may reduce wash efficiency. Discard Wash Buffer from wells immediately.

(3) Refill wells with Wash Buffer and discard immediately.
(4) Repeat step (3) another four times. This will make a total of six (6) washes with Wash Buffer.
(5) After the final wash, discard contents of wells and tap the plate on absorbent paper towel to ensure all Wash Buffer is removed.

INTERNAL QUALITY CONTROL

Each kit contains Calibrator, Positive and Negative Controls. Acceptable values for these are found on the accompanying specification sheet. The Negative and Positive Controls are intended to monitor for substantial reagent failure. The Positive Control will not ensure precision at the assay cut-off. The test is invalid and should be repeated if the absorbance readings of either the Controls or the Calibrator do not meet the specifications. If the test is invalid, patient results cannot be reported.

Quality Control (QC) requirements must be performed in conformance with local, state, and/or federal regulations or accreditation requirements and your laboratory's standard QC procedures.

It is recommended that the user refer to CLSI C24-A and 42 CFR 493.1256 for guidance on appropriate QC practices.

CALCULATIONS

IMPORTANT NOTE: The calibration factor is batch specific and is detailed in the specification sheet. Obtain the calibration factor value before commencing calculations.

1. Calculate the average absorbance of the triplicates of the Calibrator and multiply by the calibration factor. This is the Cut-off Value.
2. An index value can be calculated by dividing the sample absorbance by the Cut-off Value (calculated in step (1) above).

Alternatively,

3. Panbio™ Units can be calculated by multiplying the index value (calculated in step (2) above) by 10.

$$\textbf{Index Value} = \frac{\textbf{Sample Absorbance}}{\textbf{Cut-off Value}}$$

Example: Sample A Absorbance = 0.949
Sample B Absorbance = 0.070

Mean absorbance of Calibrator = 0.802
Calibration Factor = 0.62
Cut-off Value = 0.802 x 0.62 = 0.497

Sample A (0.949/0.497) = 1.91 Index value
Sample B (0.070/0.497) = 0.14 Index value

Panbio™ Units = Index Value X 10

Sample A 1.91 X 10 = 19.1 Panbio™ Units
Sample B 0.14 X 10 = 1.4 Panbio™ Units

TEST INTERPRETATION

The cut-off has been determined using endemic populations from South East Asia / South America and a local population from Queensland, Australia, of 208 characterised negative (208/409), 91 positive (91/409) and 110 disease controls samples (110/409). The cut-off was determined by two-graph receiver operating

characteristic analysis (TG-ROC)[2,3]. A cut-off ratio of 1.0 was selected based on the optimal F Value for sensitivity and specificity.
Diagnosis of Dengue Infection: The Dengue IgM Capture ELISA determines the level of IgM antibodies to dengue in a patient's serum. A positive result (> 11 Panbio™ Units) is indicative of either an active primary or secondary dengue infection. If differentiation between primary and secondary infection is required, the Dengue Duo (07PE10) ELISA should be used.

INDEX	Panbio™ UNITS	RESULT
<0.9	<9	Negative
0.9 – 1.1	9 – 11	Equivocal
>1.1	>11	Positive

RESULT	INTERPRETATION
Negative	No detectable IgM antibody. The result does not rule out dengue infection. An additional sample should be tested in 7-14 days if early infection is suspected. Other dengue assays should be performed to rule out acute infection.
Equivocal	Equivocal samples should be repeated. Samples that remain equivocal after repeat testing should be repeated by an alternative method or another sample should be collected.
Positive	Presence of detectable IgM antibody. Other dengue serology assays should be performed to confirm dengue infection.

The following is a recommended method for reporting the results obtained: "The following results were obtained with the Panbio™ Dengue IgM Capture ELISA. Values obtained with different methods may not be used interchangeably. The magnitude of the measured result, above the cut-off, is not indicative of the total amount of antibody present." The result should be reported as positive, negative or equivocal, and not as a numerical value.

TEST LIMITATIONS

1. The clinical diagnosis must be interpreted with clinical signs and symptoms of the patient. The results from this kit are not by themselves diagnostic and should be considered in association with other clinical data and patient symptoms.

2. Population seroepidemiology may vary over time in different geographical regions. Consequently, the cut-off may require adjustment based on local studies.
3. Screening of the general population should not be performed. The positive predictive value depends on the likelihood of the virus being present. Testing should only be performed on patients with clinical symptoms or when exposure is suspected.
4. Serological cross-reactivity across the flavivirus group is common (i.e. between dengue 1, 2, 3 & 4, Murray Valley encephalitis, Japanese encephalitis, Yellow fever and West Nile viruses). These diseases must be excluded before confirmation of diagnosis.
5. Heterophilic antibodies are a well-recognised cause of interference in immunoassays[4]. These antibodies to animal IgG may cross-react with reagent antibodies and generate a false positive signal. This must be excluded before confirmation of diagnosis.
6. The performance characteristics have not been established for visual result determination.
7. This assay employs insect expressed proteins. The cross-reactivity or interference of human anti-insect antibodies is unknown with the assay's results.
8. All sera demonstrating a positive result by the Panbio™ Dengue IgM Capture ELISA test should be referred to a reference laboratory for confirmation of positivity and epidemiological recording.
9. The Panbio™ Dengue Duo ELISA (07PE10) and Panbio™ Dengue IgG Capture-ELISA (01PE10) are most convenient for IgG determination. It uses the same capture method as the IgM ELISA, so both IgM and IgG are determined using a common method and common serum dilution.

EXPECTED VALUES

Primary dengue infection is characterised by the presence of significant or rising levels of IgM 3-5 days after the onset of infection, which can persist for 3-5 months. Secondary infection is characterised by elevation of specific IgG 1-2 days after the onset of infection and in the majority of cases (>70%) is accompanied by elevation of IgM. In early infections and some secondary infections detectable levels of IgM antibodies may be low. Some patients may not produce detectable levels of antibody within the first seven to ten days after infection. Where symptoms persist, we recommend that patients be re-tested seven days after the first specimen.

PERFORMANCE CHARACTERISTICS

255 characterised sera were tested on the Panbio™ Dengue IgM Capture ELISA as part of an in-house study. The sera included 83 endemic seronegative samples, 57 samples from patients with primary dengue infection, and 115 samples from patients with secondary dengue infection. The Panbio™ Dengue IgM Capture ELISA results were compared to the dengue status of the sera to determine sensitivity, specificity, and agreement of the assay relative to the dengue serological status. The data is summarised in Table 1.

Table 1
Dengue IgM Serological Sensitivity and Specificity of Panbio™ ELISA versus Dengue Status

Panbio™ ELISA

Dengue Status	Positive	Equivocal[a]	Negative	Total
Seronegative IgM (-) by ELISA	0	0	83	83
Primary Infection IgM (+) by HAI	41	1	1	43
Primary Infection IgM (+) by ELISA	13	0	1	14
Secondary Infection IgG (+) by HAI	41	11	34	86
Secondary Infection IgG (+) by ELISA	23	0	6	29
Total	**118**	**12**	**125**	**255**

			95% CI*
Serological Sensitivity (Primary)	= 54/57	= 94.7%	85.4 – 98.9%
Serological Sensitivity (Secondary)	= 64/115	= 55.7%	46.6 – 64.7%
Serological Specificity (Negative)	= 83/83	= 100.0%	95.7 – 100.0%
Serological Agreement (Excluding Secondary infection sera)	= 137/140	= 97.9%	93.7 – 99.6%

[a] Retesting of equivocal samples was not conducted, as the samples were unavailable.

* Confidence Interval

REPRODUCIBILITY

The reproducibility of the Panbio™ Dengue IgM Capture ELISA kit was determined by testing 7 sera 3 times each on three Panbio™ kit batch numbers on three different days. Within-run, between day, between batch and total precision were estimated by analysis of variance (ANOVA Type II) and are presented in Table 2.

Table 2
Panbio™ Dengue IgM Capture ELISA
Precision Measures (Using Index Value*)

			Within		Between Day		Between Batch		Total	
Sample	**n**	***Mean**	***SD**	**CV**	***SD**	**CV**	***SD**	**CV**	***SD**	**CV**
Positive	27	2.87	0.12	4.3%	0.05	1.7%	0.62	21.6%	0.53	18.6%
Cut-off	27	1.00	0.06	5.7%	0.00	0.0%	0.00	0.0%	0.05	5.3%
Negative	27	0.36	0.06	16.7%	0.00	0.0%	0.16	45.2%	0.15	41.0%
#1	27	6.19	0.31	5.1%	0.12	1.9%	0.41	6.7%	0.48	7.7%
#2	27	5.91	0.21	3.6%	0.11	1.9%	0.52	8.8%	0.49	8.3%
#3	27	1.25	0.05	4.1%	0.00	0.0%	0.10	7.9%	0.10	7.7%
#4	27	1.33	0.07	5.4%	0.03	1.9%	0.07	5.4%	0.10	7.2%
#5	27	1.30	0.08	6.1%	0.00	0.0%	0.12	9.3%	0.13	9.8%
#6	27	0.71	0.04	6.0%	0.00	0.0%	0.01	1.7%	0.04	6.0%
#7	27	0.81	0.06	7.7%	0.00	0.0%	0.00	0.0%	0.06	7.3%

All values are calculated from Index values (Cut-off using OD)
SD = Standard Deviation; CV = Coefficient of Variation

Note: Standard Deviation results have been rounded to two decimal places for tabulation purposes.

* Index value is calculated by dividing the sample absorbance by the cut-off

ENGLISH

CROSS-REACTIVITY

A panel of 115 specimens from patients with confirmed diseases other than dengue fever was tested to establish the analytical specificity of the Panbio™ Dengue IgM Capture ELISA. The specimens were from patients with diseases that have the potential for cross-reactivity. Each of the specimens included in the study was characterised with respect to disease diagnosis prior to analysis with the Panbio™ Dengue IgM Capture ELISA. Minimal cross-reactivity was observed for malaria, West Nile virus and rheumatoid factor specimens tested from the disease panel. Refer to Table 3 for a summary of the results.

Table 3
Panbio™ Dengue IgM Capture ELISA – Cross-reactivity Analysis


Disease Type	Total Specimens	Positive Result
Epstein-Barr virus	10	(0/10)
Malaria	10	(1/10)
Influenza A	7	(0/7)
Influenza B	3	(0/3)
Anti-nuclear antibody	30	(0/30)
Rheumatoid factor	10	(3/10)
Hepatitis A	9	(0/9)
Leptospira	7	(0/7)
Salmonella typhi	9	(0/9)
Scrub typhus	10	(0/10)
West Nile virus	10	(2/10)
Total	**115**	**(6/115)**

Product. Disclaimer:
While every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the manufacturer and distributor and test results may accordingly be affected by environmental factors and/or user error. The subject of the diagnosis should consult a doctor for further confirmation of the test result.

Warning:
The manufacturers and distributors of this product shall not be liable for any direct, indirect, or consequential losses, liability, claims, costs or damages arising from or related to an incorrect positive or negative diagnosis using this product.

Panbio™ Dengue IgM Capture ELISA 01PE20/01PE21

ASSAY PLATE
anti-human IgM

ANTIGEN-VIAL
Stabilised dengue-antigens


1 Add 100 µL of diluted samples and Controls to assay plate.

1 Add 10 µL of Antigen in 2.5 mL of Antigen Diluent and mix. Unused concentrated antigen should be stored at 2-8 °C.

2 Remove required volume of diluted antigen and mix with an equal volume of MAb HRP Conjugate in a separate glass or polypropylene vial. DISCARD UNUSED DILUTED ANTIGEN.

3 Cover plate and incubate 1 hour at 37 ± 1 °C.


3 Incubate 1 hour at 20 - 25 °C.

4 Wash the assay plate x 6. After gentle rotation to mix the antigen-MAb conjugate solution, transfer 100 µL per well to the assay plate.

5 Cover plate and incubate 1 hour at 37 ± 1 °C.

6 Wash the assay plate x 6. After the final wash, add 100 µL TMB per well and incubate at 20 - 25 °C for 10 minutes. Stop the reaction with 100 µL Stop Solution and read at 450 nm (Reference 600 - 650 nm).

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TROUBLESHOOTING GUIDE

Problem	Possible cause	Investigation / Actions
High absorbances	1. Cross-contamination from other specimens.	▶ Repeat assay taking care when washing and pipetting.
	2. Insufficient or inefficient washing or reading.	▶ Check washer efficiency
	3. Wavelength of filter not correct.	▶ Check that the wavelength is 450 nm. If a dual wavelength spectrophotometer is available, set the reference filter between 600-650 nm.
	4. High assay background.	▶ Repeat assay and include a well that contains only Sample Diluent or Sample Absorbent (i.e. a blank well).
	5. Contaminated TMB.	▶ Check that TMB is colourless or faint blue.
	6. Incubation time too long or incubation temperature too high.	▶ Check incubation time and temperature. ▶ Check incubator is at the correct temperature.
	7. Incorrect dilution of serum.	▶ Repeat assay, ensuring correct serum dilution is used.

Problem	Possible cause	Investigation / Actions
Low absorbances	1. Incubation time too short or incubation temperature too low.	▶ Ensure time and temperature of assay incubation are correct. ▶ Check incubator is set at the correct temperature.
	2. Incorrect dilution or pipetting of sera.	▶ Repeat assay ensuring correct dilutions and volumes are used. ▶ Ensure controls are sufficiently mixed.
	3. Incorrect filter wavelength.	▶ Check the wavelength is set at 450 nm. If a dual wavelength spectrophotometer is available, set the reference filter between 600-650 nm.
	4. Contaminated Conjugate solution.	▶ Dispense Conjugate directly from the bottle using a clean pipette tip; avoid transferring Conjugate to another container if possible. ▶ Do not return unused Conjugate to bottle. ▶ Ensure all pipettes and probes used to dispense the Conjugate are clean and free from serum, detergent and bleach.
	5. Kit has expired.	▶ Check expiration date of kit and do not use if expired.
	6. Air blank reading high.	▶ Investigate causes of high background absorbance.

Problem	Possible cause	Investigation / Actions
Low absorbances	7. Incorrect storage of kit.	► Ensure kit is stored at 2-8°C, plate is sealed in foil pouch and desiccant sachet is blue / purple.
	8. Kit reagents not equilibrated at room temperature.	► Allow sufficient time for reagents to equilibrate to room temperature prior to assay.
	9. Incorrect reagents used.	► Check the reagents used match those listed on the specification sheet.
	10. Over-washing of plate (e.g. inclusion of a long soak step).	► Repeat assay using recommended wash procedure.
	11. Under-washing of plate following serum incubation step (i.e. insufficient washing during wash step).	► Repeat assay using recommended wash procedure.
Poor duplicates	1. Poor mixing of samples.	► Mix reagents gently and equilibrate to room temperature.
	2. Poor pipette precision.	► Calibration may need to be checked. ► Check pipetting technique – change pipette tip for each sample and ensure excess liquid is wiped from the outside of the tip.

Problem	Possible cause	Investigation / Actions
Poor duplicates	3. Addition of reagents at inconsistent timing intervals; reagent addition takes too long, air bubbles when adding reagents.	► Use consistent timing when adding reagents. ► Ensure all dilutions are made before commencing addition to plate. ► Improve pipetting technique and speed.
	4. Inefficient washing - Wash Buffer left in wells, inconsistent washing, inadequate washing.	► Tap out Wash Buffer after washing. ► Check wells are sufficiently and uniformly filled & aspirated when washing.
	5. Reader not calibrated or warmed up prior to plate reading.	► Check reader precision. ► Check reader manual to ascertain warm-up time of instrument.
	6. Optical pathway not clean.	► Gently wipe bottom of plate. ► Check reader light source and detector are clean.
	7. Spillage of liquid from wells.	► Repeat assay, taking care not to knock the plate or splash liquid.
	8. Serum samples exhibit microbial growth, haemolysis or lipaemia.	► It is not recommended to use serum samples exhibiting microbial growth, haemolysis or lipaemia.
	9. Uneven well volumes due to evaporation.	► Cover plate with a lid or plate sealer (not provided).

Problem	Possible cause	Investigation / Actions
All wells yellow	1. Contaminated TMB.	► Check TMB is colourless or faint blue.
	2. Contaminated reagents (e.g. Conjugate, Wash Buffer).	► Check reagents for turbidity.
	3. Incorrect dilution of serum.	► Repeat assay, ensuring correct serum dilution is used.
	4. Incorrect storage of kit.	► Ensure kit is stored at 2-8°C, plate is sealed in foil pouch and desiccant sachet is blue / purple.
	5. Inefficient washing - Wash Buffer left in wells, inconsistent washing, inadequate washing.	► Tap out Wash Buffer after washing. ► Check wells are sufficiently and uniformly filled and aspirated when washing.
	6. If Conjugate reconstitution is required - Conjugate reconstituted incorrectly.	► Repeat assay ensuring Conjugate is reconstituted according to assay method.

Problem	Possible cause	Investigation / Actions
All wells negative	1. Test not performed correctly – correct reagents not added or not added in the correct sequence.	► Check procedure and check for unused reagents. ► Ensure that Stop Solution was not added before Conjugate or TMB. ► Ensure that serum was diluted in correct Sample Diluent; e.g. do not use Sample Absorbent for an IgG ELISA.
	2. Contaminated Conjugate solution.	► Dispense Conjugate directly from the bottle using a clean pipette tip; avoid transferring Conjugate to another container if possible. ► Do not return unused Conjugate to bottle. ► Ensure all pipettes and probes used to dispense the Conjugate are clean and free from serum, detergent and bleach.
	3. Over-washing of plate (e.g. inclusion of a long soak step).	► Repeat assay using recommended wash procedure.
	4. Incorrect storage of kit.	► Ensure kit is stored at 2-8°C, plate is sealed in foil pouch and desiccant sachet is blue / purple.
	5. Wash Buffer made up with Stop Solution instead of Wash Buffer Concentrate.	► Ensure Wash Buffer is made up correctly.

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1. U.S. Department of Health and Human Services: Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health. (1999). p. 8-16. *In* (ed.) Richmond JY, McKinney RW, Guidelines: Biosafety in Microbiological and Biomedical Laboratories. 4th Edition. U.S. Government Printing Office, Washington, D.C.
2. Xu, H., Lohr, J. and Greiner, M. (1997). The selection of ELISA cutoff points for testing antibody to Newcastle disease by two-graph receiver operating characteristic (TG-ROC) analysis. *J. Immunol. Meth.* **208(1)**:61-64.
3. Greiner, M., Sohr, D. and Gobel, P. (1995). A modified ROC analysis for selection of cutoff values and the definition of intermediate results in serodiagnostic tests. *J. Immunol. Meth.* **185**:123-132.
4. Seth, J. (1991). Standardisation and Quality Assurance. In: Principle and Practice of Immunoassay, Price, C.P. and Newman, D.J. (Eds). MacMillian, London

Authorized Representative
MT Promedt Consulting GmbH
Altenhofstrasse 80, 66386 St. Ingbert, Germany
Phone : +49 6894 581020, Fax : +49 6894 581021

Abbott Diagnostics Korea Inc.
65, Borahagal-ro, Giheung-gu, Yongin-si, Gyeonggi-do,
17099, Republic of Korea
Tel : 82-31-899-2800 | Fax : 82-31-899-2840
abbott.com/poct

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1. Mairuhu ATA, Wagenaar J, Brandjes DPM and van Gorp ECM. (2004). Dengue: an arthropod–borne disease of global importance. *Eur J Clin Microbiol Infect Dis.* **23**:425–433.
2. Kuhn RJ, Zhang W, Rossmann MG, Pletneve SV, Corver J, Lenches E, Jones CT, Muhkopadhyay S, Chipman PR, Strauss EG, Baker TS and Strauss JH. (2002). Structure of dengue virus: implications for flavivirus organization, maturation, and fusion. *Cell* **108**:717–725.
3. Innis BL, Nisalak A, Nimmannitya S, Kusalerdchariya S, Chongswasdi V, Suntayakorn S, Puttisri P and Hoke CH. (1989). An enzyme linked immunosorbent assay to characterize dengue infections where dengue and Japanese encephalitis co–circulate. *Am J Trop Med Hyg.* **40**:418–427.
4. Ong A, Sandar M, Chen MI and Yee Sin L. (2007). Fatal dengue hemorrhagic fever in adults during a dengue epidemic in Singapore. *Int J Infect Dis.* **11 (3)**:263–267.
5. Librarty DH, Young PR, Pickering D, Endy TP, Kalayanarooj S, Green S, Vaugh DW, Nisalak A, Ennis FA and Rothman AL. (2002). High circulating levels of the dengue virus non–structural protein NS1 EARLY in dengue illness correlate with the development of dengue haemorrhagic fever. *J Infect Dis.* **186**:1165–1168.
6. Alcon S, Talarmin A, Debruyne M, Falconar A, Deubel V,and Flammand M (2002). Enzyme–linked immunosorbent assay specific to dengue virus type 1 nonstructural protein NS1 reveals circulation of the antigen in the blood during the acute phase of disease in patients experiencing primary or secondary infections. J *Clin Microbiol*. **40**:376–381.
7. Young PR, Hilditch PA, Bletchly C and Halloran W. (2000). An antigen capture enzyme–linked immunosorbent assay reveals high levels of the dengue virus protein NS1 in the sera of infected patients. *J Clin Microbiol.* **38**:1053–1057.
8. Lindegren G, Vene S, Lundkvist A and Falk KI. (2005). Optimized diagnosis of acute dengue fever in Swedish travelers by a combination of reverse transcription–PCR and Immunoglobulin M detection. *J Clin Microbiol*. **43**:2850–2855.
9. Shu P and Huang J. (2004). Current advances in dengue diagnosis. *Clin Diagn Lab Immunol.* **11**:642–650
10. Hawkes RA, Boughton CR, Naim HM, Wild J and Chapman B. (1985). Arbovirus infections of humans in New South Wales. Seroepidemiology of the flavivirus group of togaviruses. *Med J Aust.* **143**:555–561.
11. World Health Organisation. (1997). Dengue Haemorrhagic Fever: diagnosis, treatment, prevention and control, p34–36. 2nd edition. World Health Organisation Office of Publications Geneva, Switzerland.
12. U.S. Department of Health and Human Services: Public Health Service, Centers for Disease Control and Prevention, National Institutes of Health. (1999). p. 8–16. *In* (ed.) Richmond JY, McKinney RW, Guidelines: Biosafety in Microbiological and Biomedical Laboratories. 4th Edition. U.S. Government Printing Office, Washington, D.C.

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Manufactured by
STANDARD DIAGNOSTICS, INC.
65, Borahagal–ro, Giheung–gu, Yongin–si, Gyeonggi–do, Republic of Korea
Tel : 82–31–899–2958 Email : PanbioTS@alere.com
www.alere.com

Authorized Representative
MT Promedt Consulting GmbH
Altenhofstrasse 80 D–66386 St. Ingbert Germany
Phone : +49 6894 581020, Fax : +49 6894 581021

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América Latina
+57 (1) 4824033 LAproductsupport@abbott.com

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1. Williams, KM, Jackson MA, Hamilton M. Rapid Diagnostic Testing for URIs in Children: Impact on Physician Decision Making and Cost. Infect. Med. 19(3): 109-111, 2002.
2. Bonner, A.B. et al. Impact of the Rapid Diagnosis of Influenza on Physician Decision-Making and Patient Management in the Pediatric Emergency Department: Results of a Randomized, Prospective, Controlled Trial. Pediatrics. 2003 Vol. 112 No. 2.

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CE IVD

© 2019 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Software © 2019 Axxin, used under license.
All trademarks referenced are trademarks of their respective owners.
This product is licensed and sold under agreement with Biosearch Technologies, Inc.
This product is sold under license from PHRI Properties and may be used under PHRI Properties patent rights only for human *in vitro* diagnostics.

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Latinoamérica
+57 (1) 4824033 LAproductsupport@abbott.com

REFERENCIAS

1. Shulman et al., Clinical Practice Guideline for the Diagnosis and Management of Group A Streptococcal Pharyngitis: 2012 Update by the Infectious Diseases Society of America (2012).
2. American Academy of Pediatrics. Red Book: 2012 Report of the Committee on Infectious Diseases. Pickering LK, ed. 29th ed. Elk Grove Village, IL: American Academy of Pediatrics; 2012. US Department of Health and Human Services PHS/CDC/NIH. Biosafety in microbiology and biomedical laboratories, Washington DC: US Government Printing Office, 2007.
3. Manual of Clinical Microbiology, 11th Edition, Vol. 1, ASM. (2015) pg. 279.

Abbott

Abbott Diagnostics Scarborough, Inc.
10 Southgate Road
Scarborough, Maine 04074 USA
www.abbott.com/poct

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América Latina
+57 (1) 4824033 LAproductsupport@abbott.com

REFERENCIAS

1. Manual of Clinical Microbiology, 11th Edition, Vol. 1, ASM. (2015) pg. 279.
2. https://www.iata.org/en/programs/cargo/dgr

Abbott

Abbott Diagnostics Scarborough, Inc.
10 Southgate Road
Scarborough, Maine 04074 USA
www.globalpointofcare.abbott

CE IVD

© 2020 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Software © 2020 Axxin, used under license.
All trademarks referenced are trademarks of their respective owners.
This product is licensed and sold under agreement with Biosearch Technologies, Inc.
This product is sold under license from PHRI Properties and may be used under PHRI Properties patent rights only for human *in vitro* diagnostics.

IN191000es Rev.1 2020/09

Abbott
ID NOW
COVID-19

PI - ES

Size:
Flat size: 8.375″ x 10.75″
Finished: 8.375″ x 5.375″

Printed Colors

CMYK

Incoming Inspection Colors
(For Reference Only)
Colors below are not used for printing

PMS 2995 U
Primary Blue

PMS 224 U
Magenta-Pink

PMS 303 U
Dark Blue

PN: IN191000es
Rev: 1

Date of Last Revision:
1.2 2020/09/11

Abbott

REF 62FK10

CE

Panbio™

COVID-19/Flu A&B Rapid Panel

(Nasopharyngeal)

Abbott Rapid Diagnostics Jena GmbH
Orlaweg 1, D-07743 Jena, Germany
www.globalpointofcare.abbott

Date Issued: 2022.10
62FK10A-08-A1

© 2022 Abbott. All rights reserved.
All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

Scan the QR Code for more information about this product. / Escanee el código QR para obtener más información sobre este producto. / Escaneie o código QR para mais informações deste produto. / Scanați codul QR pentru mai multe informații despre acest produs. / Skenirajte QR kod za više informacija o ovom proizvodu. / Для получения дополнительной информации об этом продукте отсканируйте QR-код. / Zeskanuj kod QR, aby uzyskać więcej informacji o produkcie. / Další informace o tomto výrobku získáte naskenováním QR kódu.

About the Test

EN

Introduction

The coronavirus disease (COVID-19) is an infectious disease caused by a newly discovered coronavirus, severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2).[1] The SARS-CoV-2 is a β-coronavirus, which is an enveloped non-segmented positive-sense RNA virus.[2] It is spread by human-to-human transmission via droplets or direct contact. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. Among patients with pneumonia caused by SARS-CoV-2, fever was the most common symptom, followed by cough.[3] Influenza, commonly called "the flu", is an infectious disease caused by influenza viruses.[4] The symptoms range from mild to severe and often include fever, runny nose, sore throat, muscle pain, headache, coughing, and fatigue. These symptoms typically begin 1-4 days after exposure to the virus and last for about 2-8 days.[5]

Test Principle

Panbio™ COVID-19/Flu A&B Rapid Panel (Nasopharyngeal) contains a membrane strip, which is pre-coated with immobilized anti-SARS-CoV-2, anti-Influenza A, and anti-Influenza B antibodies on the test line and mouse monoclonal anti-chicken IgY on the control line. Four types of gold conjugates (monoclonal antibodies specific to SARS-CoV-2, Influenza A, and Influenza B gold conjugates and chicken IgY gold conjugate) move upward on the membrane chromatographically and react with anti-SARS-CoV-2, anti-Influenza A, anti-Influenza B antibodies and pre-coated mouse monoclonal anti-chicken IgY antibody respectively. For a positive result, each pathogen of SARS-CoV-2, Influenza A and Influenza B will form a test line in the result window. Neither the test line nor the control line are visible in the result window prior to applying the patient specimen. A visible control line is required to indicate a test result is valid.

Intended Use

The Panbio™ COVID-19/Flu A&B Rapid Panel (NP) is an *in vitro* diagnostic rapid test for the qualitative detection of SARS-CoV-2 and Influenza A&B antigen (Ag) in human nasopharyngeal swab specimens from individuals, who meet COVID-19 and/or Influenza clinical and/or epidemiological criteria.

Panbio™ COVID-19/Flu A&B Rapid Panel (NP) is for professional use and is intended to be used as an aid in the diagnosis of SARS-CoV-2 and Influenza infection. The product may be used in laboratory and non-laboratory environment.

The test provides preliminary test results. Negative results don't preclude SARS-CoV-2 and Influenza A&B infection and they cannot be used as the sole basis for treatment or other management decisions. Negative results should be combined with clinical observations, patient history, and epidemiological information.

Materials Provided

- 10 Test devices with desiccant in individual foil pouches
- 10 Extraction buffer tubes in a foil pouch
- 10 Nozzle caps
- 10 Sterilized nasopharyngeal swabs for specimen collection
- 1 Tube rack
- 1 Instructions for Use with Quick Reference Guide

Materials Required but not Provided

- Personal Protective Equipment per local recommendations (i.e. gown/lab coat, face mask, face shield/eye goggles and gloves)
- Timer
- Biohazard container

Active Ingredients of Main Components

- Test device
 Gold conjugate: monoclonal antibodies to SARS-CoV-2 Ag, anti-Influenza A Ag, and B Ag - gold colloid and Chicken IgY - gold colloid, Test line: Mouse monoclonal anti-SARS-CoV-2, Mouse monoclonal anti-Influenza A and Fragment antigen-binding anti-Influenza B, Control line: Mouse monoclonal anti-Chicken IgY
- Buffer: Tricine, Sodium Chloride, Tween 20, Tergitol 15-S-9, Sodium Azide (<0.1%)

Storage and Stability of Test Kit

1. The test kit should be stored at a temperature between 2 - 30 °C. Do not freeze the kit or its components.
2. When stored in a refrigerator, all kit components must be brought to room temperature (15 - 30 °C) for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.
3. Perform the test immediately after removing the test device from the foil pouch.
4. Do not use the test kit beyond its expiration date.
5. The shelf life of the kit is indicated on the outer package.
6. Do not use the test kit if the pouch is damaged or the seal is broken.
7. Direct swab specimens should be tested immediately after collection. If immediate testing is not possible, the swab specimen can be kept in an extraction buffer tube filled with extraction buffer at room temperature (15 - 30°C) for up to 2 hours prior to testing.

Warnings

1. For *in vitro* diagnostic use only. Do not reuse the test device and kit components.
2. These instructions must be strictly followed by a trained healthcare professional to achieve accurate results. All users have to read the instructions prior to performing a test.

3. Do not eat or smoke while handling specimens.
4. Wear protective gloves while handling specimens and wash hands thoroughly afterwards.
5. Avoid splashing or aerosol formation of specimen and buffer.
6. Clean up spills thoroughly using an appropriate disinfectant.
7. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials (i.e. swab, extraction buffer tube, test device) in a biohazard container as if they were infectious waste and dispose according to applicable local regulations.
8. Do not mix or interchange different specimens.
9. Do not mix reagent of those for other product.
10. Do not store the test kit in direct sunlight.
11. To avoid contamination, do not touch the head of provided swab when opening the swab pouch.
12. The sterilized swabs should be used only for nasopharyngeal specimen collection.
13. To avoid cross-contamination, do not reuse the swabs for specimen collection.
14. Do not dilute the collected swab with any solution except for the provided buffer.
15. The buffer contains <0.1% sodium azide as a preservative which may be toxic if ingested. When disposed of through a sink, flush with a large volume of water.[6]

Test Procedure (Refer to Figures)

Nasopharyngeal Swab Specimens

Note: Healthcare professionals should comply with personal safety guidelines including the use of personal protective equipment.

Test Preparation

1. When stored in a refrigerator, all kit components must be brought to room temperature (15 - 30 °C) for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.
2. Remove the test device from the foil pouch. Place on a flat, horizontal and clean surface.
3. Keep extraction buffer tube upright and remove foil. Place the extraction buffer tube in the tube rack.

Nasopharyngeal Specimen Collection & Extraction

1. Tilt the patient's head back slightly about 45°-70° to straighten the passage from the front of the nose.
2. Insert the swab through the nostril parallel to the palate.
 ⚠ **Caution:** Use dedicated nasopharyngeal swab for specimen collection, provided in the kit.
3. Swab should reach depth equal to distance from nostrils to outer opening of the ear.

ENGLISH

⚠ **Caution:** If resistance is encountered during insertion of the swab, remove it and attempt insertion in the opposite nostril.

4. Gently rub and roll the swab, 5 times. Leave the swab in place for several seconds to absorb secretions.
5. Slowly remove swab while rotating it and insert into the extraction buffer tube.
6. Swirl the swab tip in the buffer fluid inside the extraction buffer tube, pushing into the wall of the extraction buffer tube at least five times and then squeeze out the swab by squeezing the extraction buffer tube with your fingers.
7. Dispose of the used swab in accordance with your local regulations and biohazard waste disposal protocol.
8. Press the nozzle cap tightly onto the extraction buffer tube.

Reaction with Test Device

1. Dispense 4 drops of extracted specimen vertically into the specimen well on the device.
 ⚠ **Caution:** Do not handle or move the test device until the test is complete and ready for reading.
2. Start timer. Read result at 15 minutes. Do not read result before 15 minutes or after 20 minutes.
3. Close the nozzle using cap and dispose of the extraction buffer tube according to your local regulations and biohazard waste disposal protocol.
4. Dispose of the used device according to your local regulations and biohazard waste disposal protocol.

Test Limitations

1. The contents of this kit are to be used for the professional and qualitative detection of SARS-CoV-2 and Influenza A&B antigens from nasopharyngeal swab. Other specimen types may lead to incorrect results and must not be used.
2. Failure to follow the instructions for test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
3. A negative test result may occur if the specimen was collected, extracted or transported improperly. A negative test result does not eliminate the possibility of SARS-CoV-2 and Influenza A&B infection and should be confirmed by viral culture or a molecular assay.
4. Positive test results do not rule out co-infections with other pathogens.
5. Test results must be evaluated in conjunction with other clinical data available to the physician.
6. Reading the test results earlier than 15 minutes or later than 20 minutes may give incorrect results.
7. Panbio™ COVID-19/Flu A&B Rapid Panel (NP) is not intended to detect from defective (non-infectious) virus during the later stages of viral shedding that might be detected by PCR molecular tests.[7]

8. Positive results may occur in cases of infection with SARS-CoV or Albumin.
9. Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was not evaluated for the following influenza strains as referenced by WHO guideline: A/Darwin/9/2021 (H3N2), B/Austria/1359417/2021.
10. Wait 4 hours before repeating the test following an invalid result.
11. Clinical performance was evaluated with frozen samples, and test performance may be different with fresh samples.

Quality Control

Procedural Controls

The test device has a control line ("C") and test lines ("FLU B", "FLU A", "COV 19") on the surface of the test device. Neither the test line nor the control line are visible in the result window before applying a specimen. The control line is used for procedural control and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

External Controls

Good laboratory practice recommends the use of external control materials. Users should follow federal, state and local guidelines concerning the frequency of external control testing. Specific positive controls can be ordered by the following catalogue number: 62FK11C.

Performance Characteristics

1. **Clinical performance with NP Specimen of Panbio™ COVID-19/Flu A&B Rapid Panel (NP)**

 The Panbio™ COVID-19/Flu A&B Rapid Panel (NP) is designed for nasopharyngeal specimen collection from one nostril. The performance of Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was established with Nasopharyngeal swabs tested in Korea in September 2021 with 100 COVID-19 positives, 50 Influenza A positives, 50 Influenza B positives, and 200 negative specimens for COVID-19 and Influenza A&B antigen.

 Samples were collected from subjects who presented with symptoms of Influenza A&B and/or COVID-19. One specimen was collected for testing on Panbio™ COVID-19/Flu A&B Rapid Panel (NP) and another for PCR. Swabs were collected and extracted into viral transport media. Samples were frozen after collection and thawed before testing. The performance of the Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was evaluated by comparison to a CE-marked RT-PCR method.

 The study showed a sensitivity of 88.0% (95% CI: 80.0-93.6%) and specificity of 100.0% (95% CI: 98.2-100.0%) for SARS-CoV-2, a sensitivity of 92.0% (95% CI: 80.8-97.8%) and specificity of 100.0% (95% CI: 98.2-100.0%) for Influenza A, and a sensitivity of 100.0% (95% CI: 92.9-100.0%) and specificity of 100.0% (95% CI: 98.2-100.0%) for Influenza B.

Panbio™ COVID-19/Flu A&B Rapid Panel (NP) Results for NP COVID-19

SARS-CoV-2		RT-PCR			Overall Sensitivity (95% CI)	Overall Specificity (95% CI)
		Positive	Negative	Total		
Panbio™ COVID-19/ Flu A&B Rapid Panel (NP)	Positive	88	0	88	88.0% (80.0-93.6%)	100.0% (98.2-100.0%)
	Negative	12	200	212		
	Total	100	200	300		

- Positive agreement of the Panbio™ COVID-19/Flu A&B Rapid Panel (NP) is higher with specimens of Ct values ≤30 with a sensitivity of 94.4% (95%CI:87.5-98.3%) and lower with specimens of Ct values >30 with a sensitivity of 30.0% (95%CI:6.7-65.3%) for SARS-CoV-2. As suggested in references 7-9, patients with Ct values>30 are no longer contagious.[7, 8, 9]

Panbio™ COVID-19/Flu A&B Rapid Panel (NP) Results for NP Influenza A

Influenza A		RT-PCR			Overall Sensitivity (95% CI)	Overall Specificity (95% CI)
		Positive	Negative	Total		
Panbio™ COVID-19/ Flu A&B Rapid Panel (NP)	Positive	46	0	46	92.0% (80.8-97.8%)	100.0% (98.2-100.0%)
	Negative	4	200	204		
	Total	50	200	250		

Panbio™ COVID-19/Flu A&B Rapid Panel (NP) Results for NP Influenza B

Influenza B		RT-PCR			Overall Sensitivity (95% CI)	Overall Specificity (95% CI)
		Positive	Negative	Total		
Panbio™ COVID-19/ Flu A&B Rapid Panel (NP)	Positive	50	0	50	100.0% (92.9-100.0%)	100.0% (98.2-100.0%)
	Negative	0	200	200		
	Total	50	200	250		

2. Detection Limit

The detection limit of the Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was determined on the following virus strains:

- 1.54 X 10^2 $TCID_{50}$/ml of SARS-CoV-2 which was isolated from a COVID-19 confirmed patient in Korea.
- 4.07 X 10^2 $TCID_{50}$/ml of Influenza A H1N1/California/07/09 strain.

- 1.14 X 10^2 $TCID_{50}$/ml of Influenza A H3N2/Perth/16/09.
- 2.46 X 10^4 IU/ml of Influenza B/Wisconsin/1/2010.

This influenza A viral strains were obtained from ZeptoMetrix with titer information. The performance characteristics for influenza A virus can vary depending on the subtypes, region, and year of occurrence since influenza A viruses have undergone minor amino acid changes in the target epitope region.

3. Hook Effect

There is no hook effect at:

- 1.0 X $10^{5.8}$ $TCID_{50}$/ml with SARS-CoV-2 culture fluid which was isolated from a COVID-19 confirmed patient in Korea.
- 1.17 X 10^5 $TCID_{50}$/ml with Influenza A H1N1/California/07/09.
- 1.17 X 10^5 $TCID_{50}$/ml with Influenza A H3N2/Perth/16/09.
- 1.40 X 10^7 IU/ml with Influenza B/Wisconsin/1/2010.
- 5.0 X $10^{7.25}$ $TCID_{50}$/ml with Influenza A H1N1/New Jersey/8/76.
- 3.40 X 10^9 $TCID_{50}$/ml with Influenza A/Denver/1/57.
- 5.0 X $10^{6.5}$ $TCID_{50}$/ml with Influenza A/Hong Kong/8/68.

4. Cross Reactivity

Cross-reactivity of Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was evaluated by testing 76 viruses and 37 other microorganisms (5 replicates of each). The human Influenza A strains, as expected, reacted with the Influenza A test line and did not cross react with the influenza B or COVID-19 test lines. The human Influenza B strains, as expected, reacted with the Influenza B test line and did not cross react with the influenza A or COVID-19 test lines. The non-human influenza strains did not react with the Influenza A, B or COVID-19 test lines. Panbio™ COVID-19/Flu A&B Rapid Panel (NP) has cross-reactivity with Human-SARS-coronavirus Nucleoprotein at a concentration of ≥200 pg/ml due to SARS-CoV having a high homology to the SARS-CoV-2 (79.6%).

Tested Viruses		
SARS-CoV Nucleoprotein	Influenza A H1N1 (California/07/09)	Influenza A H1N1 (A/Beijing/262/1995)
Influenza A H1N1 (A/New Jersey/8/76)	Influenza A H1N1 (A/Denver/1/57)	Influenza H1N1 (A/Puerto Rico/8/34)
Influenza A H1N1 (A/ FM/1/47)	Influenza H3N2 (A/Hong Kong/8/68)	Influenza A H3N2 (A/Victoria/3/75)
Influenza A H3N2 (A/Port Chalmers/1/73)	Influenza A H3N2 (Perth/16/09)	Influenza A H3N2 (Wisconsin/67/05)
Influenza A H5N1 (A/ China/2006)	Influenza A H5N2 (A/chicken/ Iowa/ 04-20/2015)*	Influenza A H7N1 (A/turkey/ Italy/4602/99)*
Influenza A H7N2 (A/ ruddy turnstone/ New Jersey/563/2006)*	Influenza A H7N7 (A/chicken/ Netherlands/1/03)*	Influenza B (B/Flori- da/78/2015, Victoria Lineage)
Influenza B (B/ Florida/4/2006, Yamagata Lineage)	Influenza B (B/Maryland/1/59, Victoria Lineage)	Influenza B (B/GL/1739/54)
Influenza B (B /Taiwan/2/62)	Influenza B (B/Malaysia/2506/04, Victoria Lineage)	Influenza B (B/Wisconsin/1/2010, Yamagata Lineage)
Influenza A H1N1 (A/Wisconsin/588/2019)/ (A/Victoria/2570/2019)	Influenza A H3N2 (A/Cambodia/ e0826360/2020)	Influenza B (B/ Phuket/3073/2013, Yamagata Lineage)
Influenza B (B/ Washington/02/2019, Victoria Lineage)		

* non-human influenza strain

There was no cross-reaction with the viruses and microorganisms listed below.

Tested Viruses		
HCoV-HKU1	HCoV-NL63	HCoV-OC43
MERS-CoV	HCoV-229E	Human adenovirus 1
Human adenovirus 2	Human adenovirus 3	Human adenovirus 4
Human adenovirus 5	Adenovirus 6	Adenovirus 7
Adenovirus 10	Human adenovirus 18	Adenovirus 21
Norovirus	Epstein Barr Virus	Human herpesvirus (HSV) 1
Human herpesvirus (HSV) 2	Human cytomegalovirus	Mumps virus
Parainfluenza virus type 4A	Parainfluenza virus type 4B	Parainfluenza type 1
Parainfluenza type 2	Parainfluenza type 3	Human respiratory syncytial virus A
Human respiratory syncytial virus B	Measles virus	hMPV 16 type A1
Human rhinovirus type 1B	Rhinovirus type A2	Human rhinovirus 14
Rhinovirus A16	Rotavirus A	Sendai virus
Rubella virus	Coxsackievirus A9	Coxsackievirus B4
Coxsackievirus B5	Human echovirus 2	Human echovirus 6
Echovirus 7	Echovirus 11	Human enterovirus 71
Poliovirus 1	Enterovirus D68	

Tested Microorganisms		
Bordetella pertussis	*Mycobacterium avium*	*Staphylococcus epidermidis*
Candida albicans	*Mycobacterium tuberculosis*	*Staphylococcus haemolyticus*
Chlamydophila pneumoniae	*Mycoplasma pneumoniae*	*Staphylococcus saprophyticus*
Enterococcus faecalis	*Neisseria lactamica*	*Streptococcus pneumoniae*
Escherichia coli	*Neisseria gonorrhoeae*	*Streptococcus pyogenes*
Fusobacterium necrophorum	*Neisseria meningitidis*	*Streptococcus salivarius*
Haemophilus influenzae	*Neisseria sp. (Neisseria lactamica)*	*Streptococcus anginosus*
Klebsiella pneumoniae	*Proteus vulgaris*	*Pneumocytis jirovecci (PJP)*
Legionella pneumophila	*Pseudomonas aeruginosa*	*Pooled human nasal wash*
Moraxella catarrhalis	*Corynebacterium diphtheriae*	*Enterococcus faecium*
Lactobacillus plantarum	*Neisseria sicca*	*Streptococcus agalactiae*
Streptococcus bovis	*Streptococcus oralis*	*Arcanobacterium haemolyticum*
Staphylococcus aureus		

5. Interfering Substances

The following 53 potentially interfering substances (5 replicates of each) have no impact on Panbio™ COVID-19/Flu A&B Rapid Panel (NP) with the exception of Albumin at concentrations >55 g/L.

Type	Tested Substances			
Endogenous	Mucin	Human blood	Bilirubin	Antinuclear antibody (ANA)
	Blood lipid	Hemoglobin	Rheumatoid factor	Human Anti-Mouse Antibody
	Albumin	Mometasone	Lopinavir	Acetylsalicylic acid
Exogenous	Phenylephrine	Fluticasone	Ritonavir	Ibuprofen
	Oxymetazoline	Chlorpheniramine	Arbidol	Sulphur (Sulfur)
	Sodium chloride	Diphenhydramine	Levofloxacin	Luffa opperculata
	Nasal spray (NasalCrom)	Histamine dihydrochloride	Azithromycin	Benzocaine
	Nasal spray (Equate)	Histamine	Meropenem	Menthol
	Beclomethasone	Galphimia glauca	Tobramycin	Acetaminophen
	Dexamethasone	Histaminum hydrochloricum	Mupirocin	Mouthwash (Crest)
	Flunisolide	Zanamivir	Albuterol	Mouthwash (Listerine)
	Triamcinolone acetonide	Ribavirin (Rebetol)	Ephedrine	Mouthwash (Therabreath)
	Budesonide	Tamiflu (Oseltamivir)	Guaiacol glyceryl ether	Heparin
	Phenylpropanolamine			

6. Repeatability & Reproducibility

Repeatability & Reproducibility of Panbio™ COVID-19/Flu A&B Rapid Panel (NP) was established using 3 positive specimens (moderate, weak, LoD proximity concentration) per each strain (SARS-CoV-2, Influenza A H1N1, Influenza A H3N2, and Influenza B) and negative specimen. There were no differences observed within-run, between-run, between-lots, between-sites, and between-days.

PREPARATION

1 When stored in a refrigerator, all kit components must be brought to room temperature (15 - 30 °C) for a minimum of 30 minutes prior to performing the test. Do not open the pouch while components come to room temperature.

Note: Healthcare professionals should comply with personal safety guidelines including the use of personal protective equipment.

2 **Open the package and look for the following:**

1. Test device with desiccant in individual foil pouch
2. Extraction buffer tube in the foil pouch
3. Nozzle cap
4. Sterilized nasopharyngeal swab for specimen collection
5. Tube rack
6. Instructions for Use with Quick Reference Guide


3 Carefully read these instructions prior to using Panbio™ COVID-19/Flu A&B Rapid Panel (NP).

4 Look at the expiration date of the kit box. If the expiration date has passed, select another kit for use.

5 **Open the foil pouch of test device and look for the following:**


1. Result window
2. Specimen well

Then, label the device with the patient identifier.


TEST PROCEDURE


1 Keep extraction buffer tube upright and remove foil. Place the extraction buffer tube in the tube rack.

2 Tilt the patient's head back (about 45°-70°). Insert the swab through the nostril. Gently rub and roll the swab, 5 times. Leave the swab in place for several seconds. Slowly remove swab.

3 Insert the swab specimen into the extraction buffer tube. Swirl the swab tip in the buffer fluid inside the extraction buffer tube, pushing into the wall of the extraction buffer tube at least five times and then squeeze out the swab by squeezing the extraction buffer tube with your fingers.


4 Dispose of the used swab in accordance with your local regulations and biohazard waste disposal protocol.


TEST PROCEDURE

5 Press the nozzle cap tightly onto the extraction buffer tube.

6 Dispense 4 drops of extracted specimen vertically into the specimen well on the device.


⚠ **Caution:** Do not handle or move the test device until the test is complete and ready for reading.


7 Start timer. Read result at 15 minutes.

Do not read result after 20 minutes.

Note: A Control ("C") line may appear in the result window in a few minutes but a test line ("FLU B", "FLU A", "COV 19") may take as long as 15 minutes to appear.

8 Close the nozzle using cap and dispose of the extraction buffer tube according to your local regulations and biohazard waste disposal protocol.


ENGLISH


TEST INTERPRETATION

POSITIVE


- **COVID-19 positive:** The presence of two colored lines (test line “COV 19” and control line “C”) within the result window.


- **Influenza A positive:** The presence of two colored lines (test line “FLU A” and control line “C”) within the result window.

- **Influenza B positive:** The presence of two colored lines (test line “FLU B” and control line “C”) within the result window.

- **Multi positives:** The presence of three or four colored lines (test line “FLU B”, “FLU A”, “COV 19” and control line “C”) within the result window.


⚠ **Caution:** The presence of any test line (“FLU B”, “FLU A”, “COV 19”), no matter how faint, indicates a positive result.

Note: In the case of multi-positives, User should be aware of this unique situation as both can present with similar symptoms but vary in treatment.

Note that the test provides preliminary test results.

NEGATIVE

The presence of only the control line (“C”) and no test line (“FLU B”, “FLU A”, “COV 19”) within the result window indicates a negative result.


TEST INTERPRETATION

INVALID

If the control line (“C”) is not visible within the result window after performing the test, the result is considered invalid (even though the test line may be visible).

Instructions may not have been followed correctly. It is recommended to wait 4 hours, read the IFU again and collect a new specimen and re-test using a new test device.

9 Dispose of the used device according to your local regulations and biohazard waste disposal protocol.

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
1. Rothan HA, Byrareddy SN. The epidemiology and pathogenesis of coronavirus disease (COVID-19) outbreak. J Autoimmun. 2020 May;109:102433. doi: 10.1016/j.jaut.2020.102433.
2. Guo YR, Cao QD, Hong ZS, et al. The origin, transmission and clinical therapies on coronavirus disease 2019 (COVID-19) outbreak - an update on the status. Mil Med Res. 2020 Mar 13;7(1):11. doi: 10.1186/s40779-020-00240-0.
3. Lai CC, Shih TP, Ko WC, et al. Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) and coronavirus disease-2019 (COVID-19): The epidemic and the challenges. J Antimicrob Agents. 2020 Mar;55(3):105924. doi: 10.1016/j.ijantimicag.2020.105924. Epub 2020 Feb 17.
4. "Flu Symptoms & Diagnosis". Centers for Disease Control and Prevention (CDC). 10 July 2019. Retrieved 24 January 2020.
5. "Flu Symptoms & Complications". Centers for Disease Control and Prevention (CDC). 26 February 2019. Retrieved 6 July 2019.
6. Current Intelligence Bulletin 13: Explosive Azide Hazard DHHS (NIOSH) Publication Number 78-127 August 16, 1976.
7. CDC. Discontinuation of Transmission-Based Precautions and Disposition of Patients with COVID-19 in Healthcare Settings (Interim Guidance). (2020).
8. CDC. Duration of Isolation of Precautions for Adults with COVID-19. (2020).
9. Bullard, et al. Predicting Infectious Severe Acute Respiratory Syndrome Coronavirus 2 From Diagnostic Samples. Clin Infect Dis. 2020 Dec 17;71(10):2663-2666. doi: 10.1093/cid/ciaa638.

GLOSSARY OF SYMBOLS / GLOSARIO DE SÍMBOLOS / GLOSSÁRIO DE SÍMBOLOS / GLOSAR DE SIMBOLURI / REČNIK SIMBOLA / УСЛОВНЫЕ ОБОЗНАЧЕНИЯ / OBJAŚNIENIE SYMBOLI / SLOVNÍK SYMBOLŮ

	Temperature limitation Limitación de temperatura Limitação de temperatura Limite de temperatură Ograničenje temperature Ограничение температуры Ograniczenie temperatury Omezení týkající se teploty
IVD	For *in vitro* diagnostic use only Para uso diagnóstico *in vitro* únicamente Somente para uso em diagnóstico *in vitro* Destinat utilizării doar pentru diagnostic *in vitro* Samo za *in vitro* dijagnostičku upotrebu Только для диагностики *in vitro* Wyłącznie do użytku w diagnostyce *in vitro* Pouze k použití pro diagnostiku *in vitro*
	Do not reuse No reutilizar Não reutilizar A nu se reutiliza Ne koristiti ponovo Повторное использование запрещено Nie używać ponownie Nepoužívejte opakovaně
	Do not use if package is damaged No lo use si el envase está dañado Não usar se a embalagem estiver danificada A nu se utiliza dacă ambalajul este deteriorat Ne koristiti ako je pakovanje oštećeno Не использовать при повреждении упаковки Nie używać, jeśli opakowanie jest uszkodzone Nepoužívejte, je-li je obal poškozen

LOT	Lot Number Número de lote Número do lote Număr de lot Broj serije Номер партии Numer partii Číslo šarže
REF	Catalog Number Número de catálogo Número de catálogo Număr de catalog Kataloški broj Номер по каталогу Numer katalogowy Katalogové číslo
	Consult instructions for use Consulte las instrucciones de uso Consulte as instruções de uso A se consulta instrucțiunile de utilizare Pogledajte uputstvo za upotrebu Обратитесь к инструкциям по применению Sprawdzić w instrukcji użycia Prostudujte si návod k použití
	Keep dry Manténgase seco Mantenha seco A se păstra uscat Održavajte suvim Беречь от влаги Chronić przed wilgocią Uchovávejte v suchu
	Biological Risks Riesgos biológicos Riscos biológicos Riscuri biologice Biološki rizici Биологическая опасность Zagrożenie biologiczne Biologická rizika

	Use By Fecha de caducidad Usar até A se utiliza până la Upotrebiti do Срок годности Użyć przed Spotřebujte do
	Manufacturer Fabricante Fabricante Producător Proizvođač Производитель Producent Výrobce
	Date of manufacture Fecha de fabricación Data de fabricação Data fabricației Datum proizvodnje Дата изготовления Data produkcji Datum výroby
	Keep away from sunlight Manténgase alejado de la luz solar Mantenha longe da luz solar A se feri de lumina solară Držati dalje od sunčeve svetlosti Беречь от солнечных лучей Chronić przed światłem słonecznym Udržujte mimo dosah slunečního záření
	CE mark Marca CE Marca CE Marcaj CE CE oznaka Маркировка CE Znak CE Označení CE

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English

Explanation of the test

[Introduction] Dengue viruses, transmitted by the mosquito, *Aedes aegypti* and *Aedes albopictus* mosquitoes, are widely distributed throughout the tropical and subtropical areas of the world. There are four known distinct serotypes (dengue virus 1, 2, 3 and 4). In children, infection is often subclinical or causes a self-limited febrile disease. However, if the patient is infected second times with a different serotype, a more severe disease, dengue hemorrhagic fever or dengue shock syndrome, is more likely to occur. Dengue is considered to be the most important arthropod-borne viral disease due to the human morbidity and mortality it causes. NS1 is a highly-conserved glycoprotein that is present at high concentrations in the sera of dengue-infected patients during the early clinical phase of the disease. NS1 antigen is found from the first day and up to 9 days after onset of fever in sample of primary or secondary dengue infected patients. Usually IgM does not become detectable until 5 to 10 days after the onset of disease in cases of primary dengue infection and until 4 to 5 days after onset of illness in secondary infections. In primary infections, IgG appear the 14th day and persist for life. Secondary infections show that IgGs rise within 1~2 days after the onset of symptoms and induce IgM response after 20 days of infection.

[Intended Use] The SD BIOLINE Dengue Duo rapid test is an *in-vitro* immunochromatographic, one step assay designed to detect both dengue virus NS1 antigen and differential IgG / IgM antibodies to dengue virus in human serum, plasma or whole blood. This SD BIOLINE Dengue Duo rapid test contains two devices (left side ; Dengue NS1 Ag test, right side ; Dengue IgG/IgM test). The Dengue NS1 Ag rapid test in the left side is an in-vitro immunochromatographic, one step assay designed for the qualitative determination of dengue virus NS1 antigen in human serum, plasma or whole blood. It is for the diagnosis of early acute dengue infection. This test device contains a membrane strip, which is pre-coated with anti-dengue NS1 Ag capture on test band region. The anti-dengue NS1 Ag -colloid gold conjugate and serum, plasma or whole blood sample migrate along the membrane strip by capillary action and react with the reagent (?) and forms a visible line on the antibody-antigen-antibody gold particle complex forms. The Dengue IgG / IgM Rapid Test in the right side is a solid phase immunochromatographic assay for the rapid, qualitative and differential detection of IgG and IgM antibodies to dengue virus in human serum, plasma or whole blood. This test is intended for professional use to aid in the presumptive diagnosis between primary and secondary dengue infection. This SD BIOLINE Dengue Duo provides only a preliminary test result. Therefore, isolation of virus, antigen detection in fixed tissues, RT-PCR and serological test like haemagglutination-inhibition test, more specific alternative diagnosis method must be used in order to obtain a confirmation of dengue virus infection.

[Principle] The Dengue NS1 Ag test device result window has 2 pre-coated lines, "I" (NS1 Ag Test Line) and "C" (Control Line). Both the Test Line and the Control Line in result window are not visible before applying any samples. The Control Line is used for procedural control and should always appear if the test procedure is performed correctly. The Dengue NS1 Ag can identify dengue virus NS1 antigen in serum, plasma or whole blood specimens with a high degree of sensitivity and specificity. Dengue IgG / IgM rapid test is designed to simultaneously detect and differenciate IgG and IgM antibodies to dengue virus in human serum, plasma or whole blood. This test also can detect all 4 dengue serotypes by using a mixture of recombinant dengue envelope proteins: Dengue IgG / IgM test device has 3 pre-coated lines, "G" (Dengue IgG Test Line), "M" (Dengue IgM Test Line) and "C" (Control Line) on the surface of the device. All three lines in result window are not visible before applying any samples. The "Control Line" is used for procedural control. Control line should always appear if the test procedure is performed properly and the test reagents of control line are working. A purple "G" and "M" lines will be visible in the result window if there are enough IgG and/or IgM antibodies to Dengue virus in the sample. If IgG and/or IgM antibodies to Dengue virus are not present in the sample, there is no color appearance in "G" and/or "M". When a specimen is added to the sample well, anti-dengue IgG and IgMs in the specimen will react with recombinant dengue virus envelope proteins-colloidal gold conjugates and forms a complex of antibodies-antigen. As this complex migrates along the length of the test device by capillary action, it will be captured by the relevant anti-human IgG and/or anti-human IgM immobilized in two test lines across the test device and generate a colored line.

Materials provided / Active ingredients of main component

1. SD BIOLINE Dengue Duo kit contains the following items to perform the assay.
- Dengue NS1 Ag and Dengue IgG/IgM Combo device
- Assay diluent for Dengue IgG/IgM test
- Sol capillary pipette for Dengue IgG/IgM test
- Disposable dropper for Dengue NS1 Ag test
- Instruction for use
2. Active ingredients of main components
[SD BIOLINE Dengue NS1 Ag]
- 1 test strip included ; Gold Conjugates (as main component) : Mouse monoclonal anti-dengue NS1 – gold colloid (0.27±0.05μg), Test Line (as main component) : Mouse monoclonal antidengue NS1 (0.72±0.14μg), Control Line (as main component) : Goat anti-mouse IgG (0.72±0.14μg)

[SD BIOLINE Dengue IgG / IgM]
- 1 test strip included ; Gold Conjugates (as main component) : Recombinant Dengue virus envelope protein – gold colloid (1± 0.2 μg), Test Line "G" (as main component) : Mouse monoclonal anti-human IgG (5± 1μg), Test Line "M" (as main component) : Mouse monoclonal anti-human IgM (5± 1μg), Control Line (as main component) : Rabbit anti-Dengue IgG (2.5±0.5μg),
- Assay buffer included ; 100 mM Phosphate buffer (5 ml), Sodium azide (0.01% w/w)

Kit Storage and Stability

1. The SD BIOLINE Dengue Duo test should be stored at 1~30°C.
2. The test strip is sensitive to humidity and as well as to heat.
3. Perform the test immediately after removing the test strip from the foil pouch.
4. Do not use it beyond the expiration.

Precautions / Storage and kit stability

Whole blood
- Collect the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA and sodium citrate) by venipuncture.
- If blood specimens are not immediately tested, they should be refrigerated at 2~8°C.
- When stored at 2~8°C, the blood specimens should be used within 3 days.
- For storage period longer than 3 days, freezing is recommended. They should be brought to room temperature (15~30°C) prior to test.
- Using the blood specimens in the long-term keeping more than 3 days can cause non-specific reaction.

Plasma or Serum
- **[Plasma]** Collect the whole blood into the collection tube (containing anticoagulants such as heparin, EDTA and sodium citrate) by venipuncture and then centrifuge blood to get plasma specimen.
[Serum] Collect the whole blood into the collection tube (NOT containing anticoagulants such as heparin, EDTA and sodium citrate) by venipuncture, leave to settle for 30 minutes for blood coagulation and then centrifuge blood to get serum specimen of supernatant.
- If plasma or serum specimens are not tested immediately, they should be refrigerated at 2~8°C. For storage period longer than 2 weeks, freezing is recommended. They should be brought to room temperature (15~30°C) prior to use.
- Plasma or serum specimens containing a precipitate may yield inconsistent test results. Such specimens must be clarified prior to assaying.

Precaution
- Anticoagulants such as heparin, EDTA, and citrate do not affect the test result.
- As known relevant interference, hemolytic samples, rheumatoid factors-contained samples and lipaemic, icteric samples can lead to impair the test results.
- Use separately disposable capillary pipettes or pipette tips for each sample in order to avoid cross-contamination of either samples which could cause erroneous results.

Warnings

1. For *in vitro* diagnostic use only. DO NOT RE-USE test device.
2. Do not eat or smoke while handling specimens.
3. Wear protective gloves while handling specimens. Wash hands thoroughly afterwards.
4. Avoid splashing or aerosol formation.
5. Clean up spills thoroughly using an appropriate disinfectant.
6. Decontaminate and dispose of all specimens, reaction kits and potentially contaminated materials, as if they were infectious waste, in a biohazard container.
7. Do not use the test kit if the pouch is damaged or the seal is broken.
8. The instruction must be followed exactly to get accurate results.

Procedure of the test (Refer to figure of Test Procedure)

[SD BIOLINE Dengue NS1 Ag]
1. Remove the test device from the foil pouch, and place it on a flat, dry surface.
2. With a disposable dropper, add 3 drops (about 100μl) of specimen into the sample well (S).
3. As the test begins to work, you will see purple color move across the result window in the center of the test device.
4. Interpret test results at 15~20 minutes.
5. A positive result will not change once it has been established at 15 ~ 20 minutes. However, in order to prevent any incorrect results, the test result should not be interpreted after 20 minutes.

[SD BIOLINE Dengue IgG / IgM]
1. Allow all kit components and specimen to room temperature prior to testing.
2. Remove the test device from foil pouch, place it on a flat, dry surface.
3. With 10μl capillary pipette, add 10μl of specimen drawn to black line into the square sample well marked "S".
4. Add 4 drops (about 90~120μl) of assay diluent to the assay diluent well round shaped.
5. Interpret test results at 15~20 minutes.
Caution : Do not read test results after 20 minutes. Reading too late can give false results.

Interpretation of the test (Refer to figure)

[SD BIOLINE Dengue NS1 Ag]
- **Negative Result :** The presence of each one color line within the result window indicates a negative result.
- **Positive result :** The presence of two color lines ("I" band and "C" band) within the result window no matter which line appears first, indicates a positive result.
- **Invalid Result :** If the color line is not visible within the result window after performing the test, the result is considered invalid. The directions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be re-tested.

[SD BIOLINE Dengue IgG / IgM]
- **Negative :** The control line is only visible on the test device. No IgG and IgM antibodies were detected. Retest in 5-7 days if dengue infection is suspected.
- **IgM Positive :** The control line (C) and IgM line (M) are visible on the test device. This is positive for IgM antibodies to Dengue virus. This is indicative of a primary dengue infection.
- **IgG Positive :** The control line and IgG line (G) are visible on the test device. This is positive for IgG antibodies. This is indicative of secondary or past dengue infection.
- **IgG and IgM Positive :** The control line, IgM (M) and IgG line (G) are visible on the test device. This is positive for both IgM and IgG antibodies. This is indicative of a late primary or early secondary dengue infection.
- **Invalid :** The control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Repeat the test using a new test device.

Limitations of the test

1. A negative result can occur if the quantity of Dengue virus NS1 antigen present in the specimen is below the detection limits of the assay, or the antigens that are detected are not present during the stage of disease in which a sample is collected.
2. A negative test result cannot exclude a recent infection.
3. The presence of detectable Dengue virus NS1 Ag may mean positive for early Dengue Infection. As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
4. In early infections and some secondary infections, detectable levels of IgM antibodies may be low. Some patients may not produce detectable levels of antibody within the first seven to ten days after infection. Where symptoms persist, patients should be re-tested 3-4 days after the first species.
5. Serological cross-reactivity across the Flavivirus group (Dengue virus, St. Louis encephalitis, Japanese encephalitis, West Nile and yellow fever virus) is common.

Internal Quality Control

The SD BIOLINE Dengue NS1 Ag test device has "Test Line" and "Control Line" on the surface of the cassette. And the SD BIOLINE Dengue IgG / IgM test device has "G(Dengue IgG Test Line)", "M(Dengue IgM Test Line)" and "C(Control Line)". All the Test Lines and Control Lines in result window are not visible before applying any samples. The Control Line is used for procedural control.Control line will appear if the test reagents of control line are working and assay diluent has been applied properly. Control line does not guarantee that the sample has been applied properly, or that the sample was correctly stored particularly if it was a whole blood sample.

Expected value

The NS1 is expected to be detected 1 day after the onset of fever and persist up to 9 days in both primary and secondary dengue infection. But if anti-NS1 antibodies produced, the detection of NS1 is inhibited. Primary dengue is characterized by the presence of detectable IgM 3-5 days after the onset of infection. Secondary dengue is characterized by the elevation of specific IgG 1-2 days after the onset of infection and in the majority of cases this is generally accompanied by an elevation of IgM.

Product Disclaimer:
Whilst every precaution has been taken to ensure the diagnostic ability and accuracy of this product, the product is used outside of the control of the Manufacturer and Distributor and the result may accordingly be affected by environmental factors and / or user error. A person who is the subject of the diagnosis should consult a doctor for further confirmation of the result.

Warning:
The Manufacturers and Distributors of this product shall not be liable for any losses, liability, claims, costs or damages whether direct or indirect or consequential arising out of or related to an incorrect diagnosis, whether positive or negative, in the use of this product.

Date Issued : 2011. 09
11FK45-En-1

Performance Characteristics

1. Dengue NS1 Ag Sensitivity and Specificity
The specimens used in this test were confirmed by RT-PCR and/or Viral culture.

		RT-PCR and/or Viral culture		Total
		Positive	Negative	
SD Dengue NS1 Ag	Positive	104	3	107
	Negative	8	186	194
	Total	112	189	301

* Sensitivity = 104/112 x 100 = 92.8 %, Specificity = 186/189 x 100 = 98.4 %

2. Dengue IgG/IgM Sensitivity and specificity
The comparison results of SD BIOLINE Dengue IgG/IgM test with Haemagglutination inhibition (HI) test showed that SD BIOLINE Dengue IgG/IgM test had good correlation with HI test.

		Reference assay : HI test		Total
		Positive	Negative	
SD Dengue IgG/IgM	Positive	159	3	162
	Negative	1	39	40
	Total	160	42	202

* Sensitivity = 159/160 x 100 = 99.4 %, Specificity = 39/42 x 100 = 93 %

Technique used	No. of samples teste	IgM positive only	IgG positive only	IgG+IgM Both positive	Total IgM positive	Total IgG positive	Negative
MAC-ELISA	202	39	45	76	115	121	42
Commercial Rapid	202	41	88	32	73	120	41
SD Dengue IgG/IgM	202	33	21	108	141	129	40

3. Sensitivity of SD Dengue Duo related to clinical signs appriation

Days after symptom onset	No. of samples	No. of positive samples (%)					
		IgG antibody	IgM antibody	Total Antibody(*1)	NS1 Ag	Final Result(*2)	
1st	1~7	52	9 (17.3)	21 (40.4)	23 (44.2)	43 (82.7)	49 (94.3)
	8~14	30	26 (86.7)	30 (100)	30 (100)	15 (50)	30 (100)
	15~21	36	35 (97.2)	36 (100)	36 (100)	5 (13.9)	36 (100)
2nd	1~7	36	24 (66.7)	8 (22.2)	24 (66.7)	24 (66.7)	32 (88.9)
	8~14	34	34 (100)	18 (53)	34 (100)	14 (41.2)	34 (100)
	15~21	42	42 (100)	24 (57.2)	42 (100)	3 (7.2)	42 (100)

* 1 : Dengue IgG and/or IgM * 2 : Dengue NS1 and/or IgG&IgM antibody
4. Cross-reactivity test with other Flavivirus mediated and mosquitoes-borne disease

Disease	Dengue IgM Negative/Total	Dengue IgG Negative/Total	Dengue NS1 Ag Negative/Total
Japanese Encephalitis	25/25	25/25	25/25
Yellow Fever	25/25	25/25	25/25
Malaria *P. falciparum*	25/25	25/25	25/25
Malaria *P. vivax*	25/25	25/25	25/25
Total	100/100	100/100	100/100

5. Precision : Within-run and between-run precisions have been determined by the testing of ten specimens three times : 6 of negative, 3 of low positive, 3 of medium positive and 3 of strong positive. All values were collectly identified 100% of the time.
6. To evaluate the interference of Dengue Duo rapid kit with known relevant interfering specimens, the hemolytic samples, rheumatoid factorscontained samples and lipaemic, icteric samples were investigated. In these studies, those specimens did not interfere with this test kit.
7. Analytical sensitivity : the limit of detection ; the smallest amount of the target marker that can be precisely detected ; have been equal or superior to a leading commercial dengue NS1 Ag and antibody detection rapid test.

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1

EN Now, open the package and look for the following;
1. Test device individually foil pouched with a desiccant
2. Assay diluent
3. Instructions for use
4. Capillary Pipette 10µl
5. Disposable dropper

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
EN Open the foil pouch and look for the following.
1. Result window
2. Sample well
3. Assay diluent well

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2

EN FIRST, read carefully the instruction on how to use the SD BIOLINE Dengue Duo test kit.

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3


EN Next, look at the expiry date at the back of the foil pouch. Use another kit, if expiry date has passed.

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
4


EN Using disposable dropper provided, add 3 drops (about 100µl) of serum, plasma or whole blood into the sample well marked "S".

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
Or

EN Using capillary pipette provided, add 10µl of serum, plasma or whole blood into the sample well marked "S".

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EN Put 4 drops of assay diluent into the assay diluent well round shaped.

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

Dengue NS1 Ag

Negative / Negatif / Negativo

EN The presence of only one color line within the result window indicates a negative result.

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Positive / Positif / Positivo

EN The presence of two color lines ("T" band and "C" line) within the result window in left-side, no matter which line appears first, indicates a positive result.

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EN -No control (C) line in result window.
-It is recomended that the specimen be re-tested.

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EN (primary dengue infection)
-Two lines "C" and "M" in result window.
-It is positive even if "M" line is weak.

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2. IgG positive / IgG positivo


EN (Secondary or past dengue infection)
-Two lines "C" and "G" in result window.
-It is positive even if "G" line is weak.

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3. IgG and IgM positive / IgM e IgG positivo

EN (Late primary or early secondary dengue infection)
-Three lines "C", "M" and "G" in result window.

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4. Negative / Négatif / Negativo

EN (No dengue infection)
-One line "C" in result window at right

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EN - No control (C) line in result window.
- It is recommended that the specimen be re-tested.

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© 2021 TECHLAB, Inc. All rights reserved.
C. DIFF QUIK CHEK COMPLETE, the TECHLAB Logo, and TECHLAB are trademarks of TECHLAB, Inc.
ProClin is a trademark of Rohm and Haas Company.
All other trademarks referenced are trademarks of their respective owners.