Boston Scientific

CRE™ PRO Wireguided

Balloon Dilatation Catheter

Directions for Use	**2**
Instrucciones de uso	**8**
Mode d'emploi	**14**
Gebrauchsanweisung	**20**
Istruzioni per l'uso	**26**
Gebruiksaanwijzing	**32**
Instruções de Utilização	**38**

91042082-01

2015-04

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AP), DFU, MB, CRE PRO Wireguided, Global, 91042082-01A

Black (K) ΔE ≤5.0

TABLE OF CONTENTS

WARNING **3**

DEVICE DESCRIPTION **3**

INTENDED USE/INDICATIONS FOR USE **4**

CONTRAINDICATIONS **4**

PRECAUTIONS **4**

ADVERSE EVENTS **4**

HOW SUPPLIED **4**

HANDLING AND STORAGE **4**

OPERATIONAL INSTRUCTIONS **4**

- Inspection and Preparation **4**
- Catheter Insertion for Dilatation of the Alimentary Tract **5**
- Catheter Insertion for Dilatation of the Sphincter of Oddi with or without Prior Sphincterotomy **6**
- Balloon Inflation **6**
- Catheter Withdrawal **7**

WARRANTY **7**

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AP), DFU, MB, CRE PRO Wireguided, Global, 91042082-01A

Black (K) ΔE ≤5.0

CRE™ PRO Wireguided
Balloon Dilatation Catheter

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The CRE PRO Wireguided Balloon Dilatation Catheter is capable of 3 distinct and progressively larger size diameters via controlled radial expansion. Specific balloon sizes are printed on each package and hub label.

	Balloon OD ∅		Inflation Pressure	
UPN	**mm**	**F**	**ATM**	**kPa**
180 cm				
M00558600	6-7-8	18-21-24	3-6-10	304-608-1013
M00558610	8-9-10	24-27-30	3-5.5-9	304-557-912
M00558620	10-11-12	30-33-36	3-5-8	304-507-811
M00558630	12-13.5-15	36-40.5-45	3-4.5-8	304-456-811
M00558640	15-16.5-18	45-49.5-54	3-4.5-7	304-456-709
M00558650	18-19-20	54-57-60	3-4.5-6	304-456-608
240 cm				
M00558660	6-7-8	18-21-24	3-6-10	304-608-1013
M00558670	8-9-10	24-27-30	3-5.5-9	304-557-912
M00558680	10-11-12	30-33-36	3-5-8	304-507-811
M00558690	12-13.5-15	36-40.5-45	3-4.5-8	304-456-811
M00558700	15-16.5-18	45-49.5-54	3-4.5-7	304-456-709
M00558710	18-19-20	54-57-60	3-4.5-6	304-456-608

The CRE PRO Wireguided Balloon Dilatation Catheter is designed to pass through the working channel of an endoscope and accept a 0.035 in (0.89 mm) guidewire through its guidewire lumen. This catheter comes packaged with a 0.035 in (0.89 mm), floppy tip guidewire preloaded in the guidewire lumen. The guidewire is 25 cm longer than the balloon catheter with the excess length extending from the hub end of the catheter.

For all sizes two radiopaque markers are placed under the balloon segment of the catheter to provide visual reference points fluoroscopically for balloon positioning within the stricture.

A guidewire locking device is attached to the guidewire hub of the catheter. The locking device will be packaged in the “OFF” or unlocked position. The guidewire may only be advanced or removed from the catheter when the switch on the locking device is in the “OFF” position. The guidewire may be held in place within the catheter by moving the switch to the “ON” position.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AP), DFU, MB, CRE PRO Wireguided, Global, 91042082-01A

Black (K) ΔE ≤5.0

INTENDED USE/INDICATIONS FOR USE

The CRE™ PRO Wireguided Balloon Dilatation Catheters are indicated for use in adult and adolescent populations to endoscopically dilate strictures of the alimentary tract. Also indicated in adults for endoscopic dilatation of the Sphincter of Oddi with or without prior sphincterotomy. The recommended application is printed on the package label.

CONTRAINDICATIONS

None known.

PRECAUTIONS

- Balloon dilatation catheters should be used by or under the supervision of physicians thoroughly trained in endoscopic balloon dilatation. A thorough understanding of the technical principles, clinical application, and risks associated with balloon dilatation of the alimentary and biliary tracts is necessary before using these devices.
- If resistance is met during the procedure, do not advance the catheter without first determining the cause of resistance and taking remedial action.
- The CRE PRO Wireguided Balloon Dilatation Catheter is supplied sterile by method of ethylene oxide (EO). If catheter package is opened or damaged prior to use, do not use the catheter and contact Boston Scientific for replacement.
- Any use for procedures, other than those indicated in these instructions, is not recommended.

ADVERSE EVENTS

Possible adverse events that may result from alimentary and biliary tracts balloon dilatation procedure include, but may not be limited to:

- perforation
- hemorrhage
- hematoma
- sepsis/infection
- allergic reaction to contrast medium
- pain and tenderness
- inflammation

HOW SUPPLIED

Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible. Do not reuse.

HANDLING AND STORAGE

Store in a cool, dry, dark place. Rotate inventory so that catheters are used prior to the expiration date on the package label.

OPERATIONAL INSTRUCTIONS

Inspection and Preparation

Precaution: A thorough understanding of the technical principles, clinical applications, and risks associated with endoscopic balloon dilatation is necessary before using this product.

Precaution: The CRE PRO Wireguided Balloon Dilatation Catheter is packaged sterile in a pouch. Before using, inspect the pouch for any breach of the package to ensure a sterile product.

1. Open the package and remove the CRE PRO Wireguided Balloon Dilatation Catheter from the pouch.
2. Carefully extend the catheter from its coiled configuration.
3. Prior to use, carefully examine the unit to verify that the catheter and sterile package have not been damaged during shipment. Confirm that the blue guidewire tip is positioned inside the transparent catheter tip and move the Locking Device switch to the “ON” position. This will prevent

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AP), DFU, MB, CRE PRO Wireguided, Global, 91042082-01A

Black (K) ΔE ≤5.0

guidewire movement during scope introduction. Do not use if tip is frayed, shaft or core wire is kinked, or balloon material is damaged.

4. Attach the CRE™ PRO Wireguided Balloon Dilatation Catheter hub to an integrated inflation system such as the Alliance™ Inflation System or other 60 cc inflation device with gauge to monitor balloon pressure.
5. To maximize endoscope passage, apply a vacuum to the catheter before removing protective sleeve.
6. Remove protective sleeve from balloon catheter prior to use.
7. Maintain vacuum to the catheter during insertion through the scope.

Precaution: Do not preinflate, pretest balloon, or attempt to refold balloon into protective sleeve.

Catheter Insertion for Dilatation of the Alimentary Tract

The CRE PRO Wireguided Balloon Dilatation Catheter is designed to pass through the endoscope channel size specified on the package label. Removal of the rubber valve covering the working channel may facilitate the insertion of larger balloon sizes (i.e., 15 mm or larger).

Warning: Check for proper position of the balloon catheter using endoscopic visualization. Balloon inflation in an improper location may lead to patient injury.

Precaution: The CRE PRO Wireguided Balloon Dilatation Catheter should be advanced through the endoscope using short, deliberate 2-3 cm movements to prevent inadvertent damage to the device.

1. Attach the "BALLOON" hub to an integrated inflation system, such as the Alliance Inflation System, or other 60 cc inflation device with gauge to monitor balloon pressure.
2. To maximize endoscope passage, apply and maintain a vacuum to the catheter during insertion.
3. Advance the catheter into the endoscope using short, deliberate, 2-3 cm movements. Due to variations in endoscope construction, some resistance may be experienced immediately upon entering the endoscope and again 2-3 cm before exiting the distal end of working channel.
4. Once the balloon has exited the distal end of the endoscope and is within endoscopic view, the guidewire may be advanced beyond the distal end of the catheter. To use the guidewire as a catheter guide:
 - Confirm position of switch on the locking device is in "OFF" position.
 - Advance guidewire into desired position beyond catheter tip (fluoroscopy is recommended if advancing without direct visualization).
 - Length of guidewire advancement can be approximated using the yellow marking bands on the segment of guidewire extending from the locking device (marking bands are placed 5 cm apart).
 - Advance catheter over extended portion of guidewire until balloon segment is in desired position. The CRE PRO Wireguided Balloon Dilatation Catheter includes 2 radiopaque markers, located under the balloon to help proper placement across the stricture. Match the distal and proximal radiopaque markers of the balloon with the distal and proximal ends of the stricture.

 This technique may be useful in very tight or tortuous areas of the anatomy.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AP), DFU, MB, CRE PRO Wireguided, Global, 91042082-01A

Black (K) ΔE ≤5.0

5. If desired, the guidewire may be removed from the catheter prior to endoscope insertion. A standard 0.035 in (0.89 mm) guidewire may be placed through the endoscope across the strictured area, and the catheter backloaded over the guidewire.

Precaution: Endoscopy should be used to confirm proper placement of the catheter. Verify that the shaft segment of the catheter is within endoscopic view. This ensures that the balloon has exited the endoscope completely. Fluoroscopy may also be used to confirm balloon placement. The CRE™ PRO Wireguided Balloon Dilatation Catheter includes 2 radiopaque markers, located under the balloon to help proper placement across the stricture.

6. Once the balloon is in position across the stricture, inflate the balloon per the Balloon Inflation instructions.

Catheter Insertion for Dilatation of the Sphincter of Oddi with or without Prior Sphincterotomy

Precaution: The safety and effectiveness of the use of this device to dilate regions of the biliary tree other than the Sphincter of Oddi have not been established.

The balloon size should be selected based on the size of the stone and must not exceed the diameter of the distal bile duct. Inflate the balloon to the desired diameter by utilizing the balloon OD/Inflation Pressure Chart. Watch the waist of the balloon to determine complete balloon expansion. Once inflated, maintain position for 10 to 60 seconds until desired effect is reached.

If upon completion of the initial procedure an additional procedure is desired, it may be performed as long as the balloon diameter does not exceed the diameter of the distal bile duct.

1. Remove the preloaded guidewire from the catheter prior to scope insertion.
2. With the standard 0.035 in (0.89 mm) guidewire already positioned in the common bile duct, backload the CRE PRO Wireguided Balloon Dilatation Catheter over guidewire.
3. Advance the catheter through the scope until it exits the distal end of the scope.
4. Once the balloon has exited the distal end of the duodenoscope and within endoscopic view, advance the balloon catheter over the guidewire.
5. Under endoscopic or fluoroscopic visualization position the balloon across the Sphincter of Oddi. The CRE PRO Wireguided Balloon Dilatation Catheter includes 2 radiopaque markers, located under the balloon to help proper placement across the stricture. Match the distal and proximal radiopaque markers of the balloon with the distal and proximal ends of the stricture.
6. Once the balloon is in the proper position, inflate the balloon per the Balloon Inflation instructions.

Balloon Inflation

Balloon must be filled with fluid. Depending on technique, the balloon can be filled with either sterile water, sterile saline, or a contrast mixture (e.g., 50/50 saline and contrast medium).

Warning: NEVER USE AIR OR A GAS MEDIUM TO INFLATE THE BALLOON.

1. Each CRE PRO Wireguided Balloon Dilatation Catheter is capable of the 3 distinct diameters listed on the package and hub labels. Inflate balloon to the pressure corresponding to the smallest balloon diameter and maintain until desired dilatation is achieved. To achieve larger balloon diameters, increase pressure as indicated up to the maximum inflation pressure listed on the catheter hub and package label.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AP), DFU, MB, CRE PRO Wireguided, Global, 91042082-01A

Black (K) ΔE ≤5.0

Warning: To prevent balloon burst, do not exceed the inflation pressure given for the largest diameter on the catheter's hub and package label. If the balloon does rupture or a significant loss of pressure within the balloon occurs, deflate the balloon completely and carefully remove the balloon and endoscope together as a unit. Do not attempt to withdraw a ruptured balloon through the endoscope. Continue procedure with a new catheter.

2. Monitor pressure using an inflation device and pressure gauge system, such as the Alliance™ Inflation System. As dilatation takes place, the pressure reading may fluctuate. Adjust the balloon pressure as necessary to meet the pressure given for the desired diameter (a slight drop in pressure at each diameter is normal).

Catheter Withdrawal

1. The CRE™ PRO Wireguided Balloon Dilatation Catheters are designed for rapid deflation. To deflate completely, maintain endoscopic view of proximal end of balloon as a vacuum is applied using the inflation device.

Precaution: Do not pull back on the catheter until the balloon is deflated completely

Warning: THE BALLOON MUST BE THOROUGHLY DEFLATED AND ALL FLUID REMOVED PRIOR TO WITHDRAWAL (approximately 10-30 seconds depending on balloon size and inflation medium).

Precaution: For improved withdrawal, straighten the distal end of the endoscope as much as possible. Any excess bend in the working channel will increase the force needed to withdraw the catheter through the endoscope.

2. Slowly remove the catheter from the endoscope.

Precaution: If excessive resistance is felt, remove the endoscope and balloon catheter together as a complete unit to prevent damage to body tissue, the catheter, or the endoscope.

3. The catheter should be disposed of in accordance with accepted hospital guidelines.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AP), DFU, MB, CRE PRO Wireguided, Global, 91042082-01A

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All rights reserved.
Boston Scientific Limited
Boston Scientific (Australia) Pty Ltd
© 2015 Boston Scientific Corporation or its affiliates.