Boston Scientific

WallFlex™ Duodenal

Stent System with Anchor Lock Delivery System

Directions for Use	**2**
Instrucciones de uso	**8**
Mode d'emploi	**14**
Gebrauchsanweisung	**20**
Istruzioni per l'uso	**26**
Gebruiksaanwijzing	**32**
Instruções de Utilização	**38**

90676311-01

2011-05

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Duodenal, Global, 90676311-01A

WallFlex™ Duodenal

Stent System with Anchor Lock Delivery System

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The WallFlex Duodenal Stent System with Anchor Lock Delivery System consists of two components: the implantable metallic stent and the Anchor Lock Delivery System (See Figure A). The stent is composed of nitinol wire braided in a tubular mesh configuration. This design configuration results in a stent that is both flexible and self-expanding. The delivery system consists in part of coaxial tubes. The exterior tube serves to constrain the stent until retracted during deployment. The interior tube contains a central lumen that accommodates a 0.035 in (0.89 mm) guidewire. The device may be inserted through the working channel of an endoscope (minimum channel diameter 3.7 mm).

There are 3 radiopaque marker bands; The Exterior Tube Marker Band, the Deployment Limit Marker Band and the Post Deployment Marker Band. The Exterior Tube Marker Band is located adjacent to the leading end of the stent and the Deployment Limit Marker Band is located adjacent to the trailing end of the stent. The Deployment Limit Marker Band identifies the stent deployment limit, the point beyond which the stent cannot be reconstrained. The Post Deployment Marker Band helps to facilitate accurate placement of the stent as it represents the approximate final position of the trailing end of the stent when fully expanded.


Figure A

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Duodenal, Global, 90676311-01A

DIRECTIONS FOR USE

(See Figure A for Device Diagram)

The exterior tube is pulled back by immobilizing the hub handle in one hand, grasping the exterior tube handle with the other hand, and gently sliding the exterior tube handle along the stainless steel tube towards the hub handle. Retraction of the exterior tube releases the stent.

The deployment process can be reversed if repositioning is desired. The stent can be reconstrained by the exterior tube if the stent Deployment Limit Marker Band has not been passed during deployment. To reposition, first fully reconstrain the stent by gently pushing the exterior tube handle forward while holding the hub handle completely immobile. Once reconstrained, the stent can be repositioned either distally or proximally and the deployment process restarted. Reversing the deployment process can be repeated twice, allowing a total of three deployment attempts.

INDICATIONS FOR USE

The device is indicated for the palliative treatment of gastroduodenal obstructions produced by malignant neoplasms.

CONTRAINDICATIONS

Contraindications associated with the use of the WallFlex™ Duodenal Stent System with Anchor Lock Delivery System include:

- Enteral ischemia
- Suspected or impending perforation
- Intra-abdominal abscess/perforation
- Strictures that do not allow passage of a guidewire
- Patients for whom endoscopic techniques are contraindicated
- Any use other than those specifically outlined under indications for use

WARNINGS

- The device should be used with caution and only after careful consideration in patients with elevated bleeding times or coagulopathies.
- Stents cannot be repositioned after complete deployment.
- The safety and effectiveness of this device for use in benign strictures have not been established.
- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent stent migration.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.

PRECAUTIONS

- Read the entire directions for use thoroughly before using the WallFlex Duodenal Stent System with Anchor Lock Delivery System. The WallFlex Duodenal Stent System with Anchor Lock Delivery System should only be used by or under the supervision of physicians thoroughly trained in the placement of duodenal stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device.
- The system must not be resterilized.
- The packaging and the device should be inspected prior to use. Do not use the device if the product is damaged in shipping.
- The device is intended for single use only. Do not attempt to reload deployed stents onto the delivery system.
- Use of fluoroscopy is recommended. Not using fluoroscopy can result in misplacement of the stent.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Duodenal, Global, 90676311-01A

MR CONDITIONAL

Non-clinical testing has demonstrated that the WallFlex™ Duodenal Stent System with Anchor Lock Delivery System is MR Conditional. A patient with this device may undergo MRI immediately under the following conditions:

- static magnetic field of 3 Tesla or less
- maximum spatial magnetic gradient field of 720 Gauss/cm or less
- maximum whole body averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning

In non-clinical testing, the WallFlex Duodenal Stent System with Anchor Lock Delivery System produced a temperature rise of less than or equal to 0.6°C at a maximum whole body averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of MR scanning in a 3 Tesla Excite®, G3.0-052B, General Electric Medical Systems, Milwaukee, WI; active-shielded, horizontal field MR scanner. MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the WallFlex Duodenal Stent System with Anchor Lock Delivery System.

ADVERSE EVENTS

Complications associated with the use of the WallFlex Duodenal Stent System with Anchor Lock Delivery System may include:

Procedural Complications

- Bleeding
- Pain
- Stent misplacement or inadequate expansion
- Intestinal perforation
- Death

Post Stent Placement Complications

- Bleeding
- Perforation
- Pain
- Stent migration
- Stent occlusion due to tumor in-growth through stent
- Stent occlusion due to tumor over-growth around ends of stent
- Stent occlusion
- Foreign body sensation
- Bowel impaction
- Ulceration
- Fever
- Septicemia
- Death (other than that due to normal disease progression)
- Infection

EQUIPMENT REQUIRED

- Endoscope (therapeutic gastroscope or duodenoscope with a channel size of 3.7 mm or greater)
- Fluoroscope
- 0.035 in (0.89 mm) guidewire

PROCEDURE

The WallFlex Duodenal Stent System with Anchor Lock Delivery System should be placed with the combined aid of fluoroscopy and direct visualization with an endoscope.

1. Prepare and sedate the patient as per endoscopy protocol for the site.
2. Pass the endoscope to the site of the stricture. Under direct visualization and fluoroscopic guidance, pass a 0.035 in (0.89 mm) guidewire through the working channel of the endoscope. (Radiographic contrast can also be injected through a catheter to define the contour and dimensions of

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Duodenal, Global, 90676311-01A

the stricture.) Maneuver the guidewire across the site of the stricture in preparation for stent placement.

3. Identify the stricture location, and the proximal and distal stricture margins, using fluoroscopic and endoscopic visualization. Fluoroscopic markers (e.g. external metal markers taped to the skin or internal injection of radiopaque dye) may be used to indicate the stricture margins.
4. Estimate stricture length by one of several methods (eg withdrawal of catheter with radiopaque marker at its tip or withdrawal of marked guidewire).
5. Allowing for possible further tumor development and post implant stent shortening (due to continued expansion), determine the length and number of stents necessary to cross the stricture. (Note: The fully expanded stent length should be at least 4 cm longer than the length of the stricture.) Should multiple stents be required, the leading stent should be placed first followed by the trailing stent, allowing for adequate overlapping.

 If the obstruction happens to fall near a bend, it is recommended that the stent(s) total length should be such that the two free ends fall beyond the bend and open freely into the lumen.
6. Carefully remove the delivery system from its protective packaging.

Caution: Visually inspect the entire system for damage. Visually check that the exterior tube covers the leading end of the stent. Ensure that the delivery system does not have any kinked sections.

7. Pass the WallFlex™ Duodenal Stent System with Anchor Lock Delivery System over the guidewire, through the endoscope working channel and to the site of the stricture.
8. Advance the stent across the stricture until the Post Deployment Marker Band is at the outermost point of the proximal end of the stricture. (See Figure B)

Note: Accounting for maximum shortening, the Post Deployment Marker Band identifies the approximate furthest position of the proximal end of the deployed stent. Because shortening will vary based on anatomy, the Post Deployment Marker Band should only be used to approximate the final proximal end of the stent. The Exterior Tube Marker Band should be at least 4 cm beyond the distal end of the stricture (See Figure B). If the Exterior Tube Marker Band is not at least 4 cm beyond the distal end of the stricture, a longer stent may be required or a second overlapping stent may be used to cross the stricture adequately.


Figure B

9. To begin stent deployment, immobilize the hub handle in one hand and grasp the exterior tube handle with the other hand. Gently slide the exterior tube handle back along the stainless steel tube towards the hub handle. At the same time and the same rate, the endoscopist should gently withdraw the exterior tube to prevent inadvertent stent advancement.

Note: The stent can be reconstrained at any time up to when the Exterior Tube Marker Band reaches the Deployment Limit Marker band. (See Figure C).

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Duodenal, Global, 90676311-01A


Figure C

Caution: Do not push forward or pull backward on the hub handle with the stent partially deployed. The hub handle must be securely immobilized. Inadvertent movement of the hub handle may cause misalignment of the stent and possible intestinal wall damage.

Caution: The stent should deploy easily. Do not deploy the stent if unusual force is required. To remove the device, see the Caution in Step 12 below.

Caution: Do not reconstrain around tortuous anatomy as it may cause damage to the device.

10. Assess stent position and reposition if desired. To reposition, first fully reconstrain the stent by gently pushing the exterior tube handle while holding the hub handle completely immobile. When fully reconstrained, the delivery system can be moved either proximally or distally and the deployment process restarted. Repositioning can be repeated twice, allowing a total of three deployment attempts.

Caution: A stent cannot be repositioned after the Deployment Limit Marker Band has been passed.

11. To complete stent deployment immobilize the hub handle with one hand, grasp the exterior tube handle with the other hand, and gently slide the exterior tube handle along the stainless steel tube towards the hub handle until the stent is fully deployed. (See Figure D) As before, the endoscopist should gently withdraw the exterior tube at the same time and at the same rate to prevent inadvertent stent advancement.


Figure D

Caution: If the stent can not be fully deployed or fully reconstrained, the entire delivery system can be pulled into the endoscope, using the working channel of the endoscope to reconstrain the stent.

Warning: this may cause damage to the working channel of the endoscope. In addition, the stent could inadvertently be forced off of the delivery system.

12. After the stent is correctly positioned and fully deployed, the delivery system should be closed, by pushing the exterior tube handle forward, and removed. The implanted stent length should allow for adequate overlapping into the non-obstructed anatomy to compensate for further tumor progression and stent shortening. In the event the stent does not adequately cover the stricture, a second stent should be implanted providing adequate overlapping of the initially placed stent.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Duodenal, Global, 90676311-01A

Caution: Do not place stents end to end (without overlap) as this may cause kinking. Do not use in combination with stents from other manufacturers.

Caution: If the delivery system is not fully closed prior to withdrawal, there is a possibility that the tip of the delivery system will get caught in the stent.

13. Perform routine post implant procedures.

Caution: Do not dilate the stent after placement, this may result in perforation.

HOW SUPPLIED

Supplied sterile by method of ethylene oxide.

Do not use if package is opened or damaged.

Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Excite is a trademark of General Electric Company Corporation.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Duodenal, Global, 90676311-01A

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All rights reserved.
Boston Scientific International S.A.
Boston Scientific (Australia) Pty Ltd
© 2011 Boston Scientific Corporation or its affiliates.