Boston Scientific

WallFlex™ Biliary

RX FULLY COVERED

Stent System

Directions for Use	**2**
Instrucciones de uso	**11**
Mode d'emploi	**21**
Gebrauchsanweisung	**31**
Istruzioni per l'uso	**42**
Gebruiksaanwijzing	**52**
Instruções de Utilização	**62**

90598950-01

2010-07

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Biliary RX FC, Global, 90598950-01A

TABLE OF CONTENTS

INTENDED USE/INDICATIONS FOR USE **3**

CONTRAINDICATIONS **3**

WARNING **3**

WARNINGS **3**

DEVICE DESCRIPTION **4**

Figure 1. WallFlex™ Biliary RX Fully Covered Stent **4**

Figure 2. WallFlex Biliary RX Stent System **5**

MR Conditional **5**

POTENTIAL COMPLICATIONS **5**

CAUTIONS **6**

DIRECTIONS FOR USE **6**

PRE-PROCEDURE NOTES **6**

PROCEDURE **6**

1. Choosing the Appropriate Stent **6**
2. Accessing the Stricture Using a Guidewire **7**
3. Positioning the Stent **7**

Figure 3. WallFlex Biliary RX Stent System Deployment **8**

4. Deploying the Stent **8**
5. Assessing Stent Position **8**
6. Repositioning of the Stent **8**
7. Removing a Partially Deployed Stent **9**
8. Post-deployment **10**
9. Removing Endoscope **10**

HOW SUPPLIED **10**

Handling and Storage **10**

WARRANTY **10**

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Biliary RX FC, Global, 90598950-01 A

WallFlex™ Biliary

RX FULLY COVERED

Stent System

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

INTENDED USE/INDICATIONS FOR USE

The WallFlex Biliary RX Fully Covered Stent System is indicated for use in the palliative treatment of biliary strictures produced by malignant neoplasms.

CONTRAINDICATIONS

The WallFlex Biliary RX Fully Covered Stent System is contraindicated for:

- Placement in biliary strictures caused by benign tumors, as the long-term effects of the stent in the bile duct is unknown.
- Placement in strictures that cannot be dilated enough to pass the delivery system.
- Placement in a perforated duct.
- Placement in very small intrahepatic ducts.
- Those patients for whom endoscopic techniques are contraindicated.
- Any use other than those specifically outlined under indications for use.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

WARNINGS

Visually inspect the system for any signs of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this warning may result in patient injury.

NO WARRANTY IS MADE WITH REGARD TO REMOVABILITY OF THIS DEVICE BY ENDOSCOPIC MEANS OR OTHERWISE. Careful consideration must be taken when removing a stent from an intrinsic malignant tumor. Removal may result in perforation, bleeding or tissue abrasion.

The safety and effectiveness of this device for use in the vascular system has not been established.

Passing a second stent delivery system through a just deployed stent is not recommended and could cause the stent to dislodge.

Use caution when placing stent near ductal branches to avoid obstruction of duct. Placement of a fully covered biliary stent across a branch duct or major bifurcation may result in complications due to blockage of flow from the branch duct and prevent endoscopic or transhepatic access for future procedures.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Biliary RX FC, Global, 90598950-01A

A stent cannot be reconstrained after the reconstrainment limit has been exceeded.

No warranty is made with regard to removability of this device by endoscopic means or otherwise. Careful consideration must be taken when removing a stent from an intrinsic malignant tumor. Removal may result in perforation, bleeding or tissue abrasion.

The WallFlex™ Biliary RX Fully Covered Stent should not be moved or removed after completion of the initial stent placement procedure. Manipulating, repositioning or removal of the stent may result in perforation, bleeding, tissue abrasion or other patient injury.

DEVICE DESCRIPTION

The WallFlex Biliary RX Fully Covered Stent System consists of a flexible delivery system preloaded with a self-expanding biliary metal stent. The stent is made from a metallic radiopaque material that is formed into a cylindrical mesh. The WallFlex Biliary RX Fully Covered Stent is offered fully covered with Permalume™ Coating, a translucent silicone polymer, to reduce the potential for tumor ingrowth through the stent (Figure 1). The WallFlex Biliary RX Fully Covered Stent has a retrieval loop for removal during the initial stent placement procedure, to be used in the event of incorrect placement. The stent has a flare at both ends to aid in preventing migration after the stent has been placed in the bile duct. The WallFlex Biliary RX Fully Covered Stent System is an RX compatible system only. The WallFlex Biliary RX Fully Covered Stent is provided sterile using ethylene oxide and is a single use device.


Figure 1. WallFlex Biliary RX Fully Covered Stent

The delivery system is a coaxial tube design. The exterior tube is used to constrain the stent before deployment and reconstrain the stent, if stent repositioning is necessary, after partial deployment. The exterior tube has a clear section so that the constrained stent is visible. A yellow transition zone on the inner tube of the delivery system is visible between the stent and the blue outer sheath. There are four radiopaque (RO) markers to aid in the deployment of the stent while using fluoroscopy (Figure 2). There are two RO markers on the inner tube of the delivery system identifying the ends of the constrained stent (Figure 2, marker 1 and 3). Between these RO markers is an additional RO marker that indicates at what point reconstrainment is no longer possible (Figure 2, marker 2). The fourth RO marker at the leading end of the exterior tube indicates how far the stent has been deployed (Figure 2, marker 4). There is one visual marker on the interior tube between the handles to aid in the deployment of the stent (Figure 2, marker 5). The visual marker indicates the point at which reconstrainment is no longer possible. The interior tube has a single central lumen to accommodate a 0.035 in (0.89 mm) guidewire.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Biliary RX FC, Global, 90598950-01A


Figure 2. WallFlex™ Biliary RX Stent System

Caution: Read the entire Directions for Use thoroughly before using the WallFlex Biliary RX Fully Covered Stent System. The WallFlex Biliary RX Fully Covered Stent System should only be used by or under the supervision of physicians thoroughly trained in biliary prosthesis placement. A thorough understanding of the technical principles, clinical applications, and risks associated with this procedure is necessary before using this device.

MR Conditional

Non-clinical testing has demonstrated that the Covered WallFlex Biliary RX Stent System is MR Conditional. A patient with this device may undergo MRI immediately under the following conditions:

- static magnetic field of 3 Tesla or less
- maximum spatial magnetic gradient field of 720 Gauss/cm or less
- maximum whole body averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning

In non-clinical testing, the Covered WallFlex Biliary RX Stent System produced a temperature rise of less than or equal to 1.4°C at a maximum whole body averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of MR scanning in a 3 Tesla Excite™, G3.0-052B, General Electric Medical Systems, Milwaukee, WI; active-shielded, horizontal field MR scanner.

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the Covered WallFlex Biliary RX Stent System.

POTENTIAL COMPLICATIONS

The following complications have been reported in the literature for biliary prostheses or have been observed in the Boston Scientific clinical trial of this device.

These include, but are not limited to:

- Pain
- Bleeding
- Fever
- Nausea
- Vomiting
- Infection
- Inflammation
- Recurrent obstructive jaundice
- Stent occlusion
- Tumor overgrowth around ends of stent
- Tumor ingrowth through the stent
- Mucosal hyperplasia
- Cholangitis

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Biliary RX FC, Global, 90598950-01A

- Cholecystitis*
- Pancreatitis
- Bile duct ulceration
- Perforation of duodenum or bile duct
- Stent migration
- Death (other than that due to normal disease progression)
- Stent misplacement
- Perforation of the gallbladder due to the stent covering the cystic duct*

***Note:** In a small clinical trial of this device, two out of four (50%) subjects who had a stent placed across the cystic duct developed cholecystitis. One of these subjects suffered a perforated gallbladder due to the stent covering the cystic duct, requiring a drain to be placed.

CAUTIONS

Excessive force should not be used to position or deploy the stent. This may cause inadvertent damage to the device and/or endoscope.

The sterile packaging and device should be inspected prior to use. If sterility or performance of the device is suspected to be compromised, it should not be used.

DIRECTIONS FOR USE

Required Equipment:

- RX Locking device
- Endoscope with minimum 3.2 mm working channel
- 0.035 in (0.89 mm) stiff-bodied guidewire
- WallFlex™ Biliary RX Fully Covered Stent System containing a stent of the appropriate length and diameter
- Fluoroscopic capability for pre-stent placement and stent placement confirmation

PRE-PROCEDURE NOTES

Radiography of pertinent anatomy performed no more than 10 days before the procedure should be available.

Initial preparation of delivery system

- Carefully remove the delivery system from the protective packaging.
- Visually inspect the device for damage or defects.
- Visually check that the leading end of the stent is covered by the exterior tube.
- Ensure that no stent wires have perforated the exterior tube.

Note: DO NOT remove the shipping mandrel from the leading end of the device (Figure 2), this will help facilitate guidewire access.

RO markers are used to aid in positioning the stent across the stricture. During deployment, these RO markers indicate when the reconstrainment limit is reached and when the stent is fully deployed. The RO markers are fully described in the Device Description section of these directions.

PROCEDURE

1. **Choosing the Appropriate Stent**

 Fluoroscopy can be used to locate the stricture with the aid of a contrast medium.

 Measure the stricture using fluoroscopy. Determine the stent length necessary to adequately cover the lesion. The stent length should allow for further lesion development.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Biliary RX FC, Global, 90598950-01A

In the event a single stent will not adequately cover the stricture, a second stent of the same diameter should be implanted providing adequate overlapping (minimum 10 mm) of the initially placed stent to ensure a smooth transition between the stents. The second stent should be placed to ensure complete tumor coverage. It is recommended that the proximal stent be placed first followed by the distal stent.

Warning: Passing a second stent delivery system through a just deployed stent is not recommended and could cause the stent to dislodge.

2. Accessing the Stricture Using a Guidewire

Endoscopic and fluoroscopic placement of the guidewire is recommended.

Ensure that the RX Locking Device (supplied separately) is firmly seated on the endoscope against the side of the working channel port as outlined in its directions for use.

Place the 0.035 in (0.89 mm) guidewire through the endoscope, into the ampula and across the biliary stricture. Lock the guidewire in place with the locking device.

A sphincterotomy and/or predilatation of the biliary stricture may be performed prior to stent implantation at the discretion of the physician.

Insert the trailing end of the guidewire through the tip of the delivery system and advance in short strokes until it exits the guidewire access port. The guidewire access port allows the guidewire to exit the delivery system of the WallFlex™ Biliary Fully Covered Stent System approximately 27 cm from the tip of the device. The guidewire then lies along the side of the delivery system.

The shipping mandrel will be pushed out as the guidewire is inserted. Do not remove the shipping mandrel before loading the guidewire.

Advance the system over the guidewire. Unlock the wire to pass the stent into the scope.

Lock the guidewire into place and continue advancing the system through the scope until the stent is in the desired location.

Caution: Attempting to place the WallFlex Biliary RX Stent in patients with severe anatomical angulation may prevent the stent from deploying or cause damage to the device.

3. Positioning the Stent

The system should be advanced through the endoscope using short, controlled, 2-3 cm movements.

Use the leading RO marker to position the stent a minimum of 1 cm beyond the stricture.

As the stent is deployed, the stent will foreshorten. The amount of foreshortening depends on the anatomy of the stricture.

For stent placements through the papilla, keep direct endoscopic visualization of the yellow transition zone at the end of the stent. Additionally, at least two diamonds of the deployed stent should be showing outside the papilla at all times.

There are RO markers to aid in the deployment of the stent while using fluoroscopy (Figure 3).

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Biliary RX FC, Global, 90598950-01A


Figure 3A. Markers 1 and 3 identify the ends of the constrained stent.


Figure 3B. Marker 4 indicates to what extent the stent has been deployed.


Figure 3C. Marker 4 meets marker 2 when the reconstrainment limit has been reached. Marker 5 (on the handle) is additional indication.

Figure 3D. Marker 4 moves distally past marker 1 when the stent is fully deployed.

Figure 3. WallFlex™ Biliary RX Stent System Deployment

4. Deploying the Stent

Monitor the stent position fluoroscopically during the deployment process.

To begin stent deployment, immobilize the stainless steel tube in one hand, grasp the leading handle with the other hand, and gently slide the handle back along the stainless steel tube until the reconstrainment limit is reached. The reconstrainment limit is marked by a visual marker on the stainless steel tube of the deployment system (Figure 3C, marker 5). Fluoroscopically, this point is reached when marker 4 meets marker 2.

Warning: Use caution when placing stent near ductal branches to avoid obstruction of duct. Placement of a fully covered biliary stent across a branch duct or major bifurcation may result in complications due to blockage of flow from the branch duct and prevent endoscopic or transhepatic access for future procedures.

5. Assessing Stent Position

When the reconstrainment limit is reached, check the position of the stent endoscopically and fluoroscopically. If the position of the RO markers and stent are correct, then complete deployment (Figure 3D). If the position is not correct, reposition the stent system per the directions in step 6.

6. Repositioning of the Stent

To reposition the stent, first reconstrain the stent by holding the leading handle stationary and moving the stainless steel tube away from the leading handle. It may be necessary to guide the delivery system into the endoscope.

Caution: Do not push forward on the delivery system with the stent partially deployed. The stainless steel tube must be immobilized securely. Pushing on the delivery system may cause misalignment of the stent and possible duct damage. The stent should deploy easily. Do not deploy the stent if unusual force is required, since this may indicate a failed device.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Biliary RX FC, Global, 90598950-01A

Under fluoroscopy, the stent will not be completely reconstrained until the leading marker 4 is even with the exterior tube marker 3 (Figure 3A). When fully constrained, the delivery system can be moved either proximally or distally and the deployment process re-started.

The stent can be reconstrained and the position adjusted proximally or distally at any time before the reconstrainment limit is exceeded. The reconstrainment limit is marked by a visual marker on the stainless steel tube of the deployment system. Fluoroscopically, this point is reached when marker 4 meets marker 2.

As an alternative method for distal repositioning only, immobilize both the stainless steel tube and the leading handle and pull the entire delivery system back.

Caution: Do not reconstrain around tortuous anatomy as it may cause damage to the device.

Caution: Do not allow the unconstrained stent to re-enter the endoscope during the reconstrainment process.

Warning: A stent cannot be reconstrained after the reconstrainment limit has been exceeded.

Stent reconstrainment can be completed twice, allowing a total of three deployment attempts.

Warning: No warranty is made with regard to removability of this device by endoscopic means or otherwise. Careful consideration must be taken when removing a stent from an intrinsic malignant tumor. Removal may result in perforation, bleeding or tissue abrasion.

If the positioning of a WallFlex™ Biliary RX Fully Covered stent is not correct and one of the following has occurred, continue to fully deploy the stent:

A. The stent has already been deployed past the reconstrainment limit

 OR

B. The stent has already been reconstrained twice

Then in either case, using forceps, grasp the retrieval loop on the end of the stent (Figure 1). Gently pull the stent back with the scope to remove.

Note: Pre-clinical testing to evaluate stent removal during the initial stent placement procedure was assessed using the WallFlex Biliary RX Partially Covered stent in a porcine model. The WallFlex Biliary RX Partially Covered Stent is identical in design to the WallFlex Biliary RX Fully Covered Stent except for the covering length. In a limited evaluation of the WallFlex Biliary RX Partially Covered Stent in a porcine model, 6 stents were successfully removed intra-procedurally from 6 separate animals using the methods described above for removal of a fully deployed stent during the initial stent placement procedure.

Warning: The WallFlex Biliary RX Fully Covered Stent should not be moved or removed after completion of the initial stent placement procedure. Manipulating, repositioning or removal of the stent may result in perforation, bleeding, tissue abrasion or other patient injury.

7. Removing a Partially Deployed Stent

To remove a partially deployed stent, first reconstrain the entire stent and then pull back on the entire delivery system.

If reconstrainment is not possible, immobilize both the stainless steel tube and the leading handle and pull the entire delivery system back. The unconstrained stent should be removed through the body external of the endoscope.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Biliary RX FC, Global, 90598950-01A

8. **Post-deployment**

After the stent is correctly positioned and fully deployed, while monitoring under fluoroscopic guidance, hold the leading handle stationary and carefully retract the trailing handle until the tip is flush with the end of the outer sheath. Then retract the delivery system and guidewire through the endoscope.

Note: If during delivery system withdrawal, the delivery system is not free of the stent and/or the stent begins to move distally down the bile duct, immediately stop retraction of the delivery system. Advance the inner sheath of the delivery system forward by advancing the trailing handle of the delivery system while holding the leading handle (outer sheath) stationary. Carefully advance the inner sheath forward approximately 1 cm and commence retraction of the delivery system again. Repeat until the inner sheath can be retracted without interfering with the position of the deployed stent.

Using standard operative procedure, perform routine post implant radiographic procedures to demonstrate location and patency of the stent.

The stent length should allow for adequate overlapping into the non-strictured duct to compensate further tumor progression and stent foreshortening. In the event the stent does not adequately cover the stricture, a second stent of the same diameter should be implanted providing adequate overlapping of the initially placed stent. Directions on placing two overlapping stents are provided in Step 1.

9. **Removing Endoscope**

Withdraw endoscope from the patient.

HOW SUPPLIED

Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Excite is a Trademark of General Electric Company Corporation.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Biliary RX FC, Global, 90598950-01 A

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All rights reserved.
Boston Scientific International S.A.
Boston Scientific (Australia) Pty Ltd
© 2010 Boston Scientific Corporation or its affiliates.