# RIGIFLEX™ II SINGLE USE ACHALASIA BALLOON DILATOR

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Boston Scientific

**℞ ONLY CAUTION: FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.**

**Warning:** Contents supplied STERILE using an ethylene oxide process. Do not use if sterile barrier is damaged. If damage is found call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious desease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

## DESCRIPTION

The Rigiflex™ II Achalasia Balloon Dilator is used for dilatation of the cardia Lower Esophageal Sphincter (LES) in patients with Achalasia. It is comprised of a catheter and a balloon. The balloon is composed of an inelastic material and may be filled with air to exert radial force on a narrowed segment in the alimentary tract resulting from Achalasia. This balloon can be inflated with the Boston Scientific Achalasia Pneumatic Hand Pump and Monitor.

## INTENDED USE/INDICATIONS FOR USE

The Rigiflex II Achalasia Balloon Dilator is indicated for dilatation of the cardia (LES) in patients with achalasia.

## CONTRAINDICATIONS

- Acute or incompletely healed esophageal perforation
- Bleeding disorders
- Severe pulmonary disease
- Recent myocardial infarction
- Recent esophageal perforation or surgery
- Pharyngeal or cervical deformity
- Recent laparotomy
- Large thoracic aortic aneurysm

## WARNINGS

Complications may include:

- Chest pain
- Reflux
- Perforation
- Bleeding
- Infection
- Asymptotic Esophageal Tear
- Esophageal Intramural Hematoma

## CONTENTS

1 - Rigiflex™ II Achalasia Balloon Dilator

1 - 0.038 in./ 0.97 mm floppy tip guidewire

## HOW SUPPLIED

The Single Use Achalasia Balloon Dilator is supplied sterile by means of ethylene oxide (EO). If the catheter package is opened or damaged prior to use, do not use the catheter or try to resterilize. Do not use if labeling is incomplete or illegible.

## PRECAUTIONS

- Exercise care when handling the floppy tip of the guidewire to avoid damage.
- Insert the floppy tip of the guidewire first into the working channel of the endoscope.
- Do not exceed the maximum inflation pressure of the balloon.

## INSPECTION

1. Prior to opening package, inspect seal integrity and package surface integrity to ensure sterile barrier has not been compromised.
2. Open the package and remove the Single Use Achalasia Balloon Dilator and guidewire hoop from the pouch.
3. Prior to use, carefully examine the catheter, balloon, and guidewire to verify that they were not damaged in shipment. Verify that the catheter tip is not frayed and that the catheter and guidewire are free of sharp ends and kinks.

***NOTE:*** The guidewire is packaged with the floppy tip of the guidewire at the inner diameter of the hoop.

## INSTRUCTIONS FOR USE

## PREPARATION AND INSERTION OF THE GUIDEWIRE AND BALLOON

1. Attach a hospital supplied 60 ml syringe to the catheter.
2. Draw back on syringe and create a negative vacuum, drawing the balloon to its lowest profile.
3. Remove the guidewire introducer from the guidewire hoop.
4. Insert the guidewire introducer into the biopsy cap of the endoscope.
5. After full endoscopic exam, insert the floppy tip of the provided guidewire, indicated by the flag on the guidewire hoop, through the working channel of the scope and cross the LES.
6. Once the guidewire is in place, remove the endoscope, leaving the guidewire in place.
7. The balloon catheter is then backloaded over the guidewire and the balloon is centered in the LES.
8. The scope may be reinserted for endoscopic visualization.
9. The Rigiflex II Achalasia Balloon Dilator may also be used under fluoroscopic guidance. The catheter has 4 radiopaque markers: one at the distal end of the balloon, one at the proximal end of the balloon and two in the center of the balloon. These positioning markers are easily seen under fluoroscopy as is the radiopaque tip. The soft atraumatic tip at the end of the catheter acts as a guide for the balloon catheter and gently dilates the esophagus and sphincter as the catheter is advanced.

## BALLOON INFLATION

When the balloon is positioned in the esophogogastric junction, inflate, using the Achalasia Pneumatic Hand Pump and Monitor. The hand pump is attached to the balloon lumen on the catheter and is operated like the pump on a blood pressure cuff. Pump the balloon until no balloon waisting is observed. Do not exceed the maximum inflation pressure of 20 psi (1.4 ATM, 138 kPa). Depending on the patient, you may wish to deflate and reinflate the balloon several times until the criteria for sufficient dilatation are met.

## CATHETER WITHDRAWAL

1. When dilation is complete, deflate the balloon with a 60 ml syringe to its lowest profile by applying suction to the balloon lumen.

***NOTE:*** The greater the suction applied and held during withdrawal the smaller the deflated profile of the balloon. For maximum suction a 60 ml (cc) syringe is recommended.

2. Once the balloon is at its lowest profile, remove the balloon and guidewire from the patient.

## HANDLING AND STORAGE

Store in a cool, dry, dark place. Rotate inventory so that catheters are used prior to the expiration date on the package.

## WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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© 2010 Boston Scientific Corporation or its affiliates. All rights reserved.
**Boston Scientific (Australia) Pty Ltd**
**Boston Scientific International S.A.**