Boston Scientific

WallFlex™ Esophageal

FULLY COVERED

Stent System

Directions for Use	**2**
Instrucciones de uso	**11**
Mode d'emploi	**20**
Gebrauchsanweisung	**29**
Istruzioni per l'uso	**38**
Gebruiksaanwijzing	**47**
Instruções de Utilização	**56**

90599988-01

2010-06

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Esophageal FC, Global, 90599988-01A

TABLE OF CONTENTS

WARNING **3**

DEVICE DESCRIPTION **3**

Figure 1. WallFlex™ Esophageal Fully Covered Stent **4**

INTENDED USE/INDICATIONS FOR USE **4**

CONTRAINDICATIONS **4**

ADVERSE EVENTS **5**

WARNINGS AND CAUTIONS **5**

HOW SUPPLIED **6**

DIRECTIONS FOR USE **6**

Required Equipment **6**

Pre-Procedure **6**

Initial Preparation of Delivery System **6**

OPERATIONAL INSTRUCTIONS **6**

Start of Initial Stent Placement Procedure **6**

Figure 2. Delivery System and Radiopaque (RO) Markers **7**

Figure 3. Delivery system, Visual Markers and Handles **8**

POST-PROCEDURE **10**

WARRANTY **10**

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Esophageal FC, Global, 90599988-01A

WallFlex™ Esophageal

FULLY COVERED

Stent System

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied NON-STERILE. Do not use if damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or sterilize. Reuse, reprocessing or sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The WallFlex Esophageal Fully Covered Stent System consists of a flexible delivery system preloaded with a self-expanding esophageal metal stent. The stent is made from a metallic radiopaque material that is formed into a cylindrical mesh (Figure 1). It is fully covered with a translucent silicone polymer to restrict tumor ingrowth through the wire mesh and to occlude concurrent esophageal fistulas. A suture is threaded through the proximal end of the stent and is intended to aid in removal during the initial stent placement procedure, to be used in the event of incorrect placement. The stent has flares at both extremities to aid in minimizing migration after the stent has been placed in the esophagus. The delivery system is a coaxial tube design. The exterior tube is used to constrain the stent before deployment and reconstrain the stent, if desired, after partial deployment. The exterior tube has a clear section so that the constrained stent is visible. There are four radiopaque (RO) markers to aid in the deployment of the stent while using fluoroscopy (Figure 2). There are RO marker bands on the inner tube of the delivery system identifying the ends of the constrained stent (Figure 2, #1 and #3). Between these marker bands is an additional RO marker band to indicate at what point reconstrainment is no longer possible (Figure 2, #2). The fourth radiopaque marker band at the distal end of the exterior tube indicates up to what point the stent has been deployed (Figure 2, #4). The tip and interior tube are also radiopaque for use with fluoroscopy. There are two visual markers on the interior tube between the handles to aid in the deployment of the stent (Figure 3, #5 and #6). Visual marker #5 indicates the start position and marker #6 indicates the point at which reconstrainment is no longer possible. The interior tube has a single central lumen to accommodate a 0.038 in (0.97 mm) guidewire.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Esophageal FC, Global, 90599988-01A


Figure 1. WallFlex™ Esophageal Fully Covered Stent

Caution: Read the entire Directions for Use thoroughly before using the WallFlex Esophageal Fully Covered Stent System. The WallFlex Esophageal Fully Covered Stent System should only be used by or under the supervision of physicians thoroughly trained in esophageal prosthesis placement. A thorough understanding of the technical principles, clinical applications, and risks associated with this procedure is necessary before using this device.

Note: MRI Conditional Labeling

Non-clinical testing has demonstrated that the WallFlex Esophageal Stent is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient field of 720 Gauss/cm or less
- Maximum MR system reported whole-body-averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of scanning.

In non-clinical testing, the WallFlex Esophageal Stent produced a temperature rise of less than 0.6 °C at a maximum whole body averaged specific absorption rate (SAR) of 3 W/kg for 15 minutes of MR scanning in a 3 Tesla MR scanner (Excite™, General Electric).

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the WallFlex Esophageal Stent.

INTENDED USE/INDICATIONS FOR USE

The WallFlex Esophageal Fully Covered Stent System is intended for maintaining esophageal luminal patency in esophageal strictures caused by intrinsic and/or extrinsic malignant tumors, and occlusion of concurrent esophageal fistulas.

CONTRAINDICATIONS

The WallFlex Esophageal Fully Covered Stent System is contraindicated for:

- Placement in esophageal strictures caused by benign tumors, as the long-term effects of the stent in the esophagus are unknown.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Esophageal FC, Global, 90599988-01A

- Placement in strictures that cannot be dilated enough to pass the endoscope or the delivery system.
- Placement of the proximal end of stent within 2 cm of the cricopharyngeal muscle.
- Placement in an esophago-jejunostomy (following gastrectomy), as peristalsis and altered anatomy may displace stent.
- Placement in necrotic chronically bleeding tumors, if bleeding is active at the time of placement.
- Placement in polypoid lesions.
- Those patients for whom endoscopic techniques are contraindicated.
- Any use other than those specifically outlined under indications for use.
- Placement in patients who have an underlying bleeding diathesis

ADVERSE EVENTS

The potential adverse effects associated with esophageal stent placement may include:

- Bleeding
- Perforation
- Pain
- Aspiration
- Stent migration
- Tumor overgrowth around stent ends
- Foreign body sensation
- Food bolus impaction
- Reflux
- Esophagitis
- Edema
- Ulceration
- Fever
- Infection
- Sepsis
- Septicemia
- Recurrent dysphagia
- Fistula formation
- Tracheal compression/obstruction (or acute airway compression)
- Hematemesis
- Death (other than that due to normal disease progression)
- Stent fracture

Possible Post Stent Complications

- Sensitivity to the metal component of the stent
- Mediastinitis
- Aspiration
- Intestinal obstruction (secondary to stent migration)
- Granulation of tissue around stent ends

WARNINGS AND CAUTIONS

The WallFlex™ Esophageal Fully Covered Stent System should be used with caution and only after careful consideration in patients with:

- Strictures exceeding 12 cm in length
- Significant preexisting pulmonary or cardiac disease

Warning: Stent is considered to be a permanent device. Once stent placement is permanently achieved, stent removal or repositioning is not recommended.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Esophageal FC, Global, 90599988-01A

Warning: Visually inspect the system for any signs of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this warning may result in patient injury.

HOW SUPPLIED

The WallFlex™ Esophageal Fully Covered Stents are supplied in single pack configuration. Store in a cool, dry, dark place. See product label for expiration date. Do not use if labeling is incomplete or illegible. Do not use if package is opened or damaged.

DIRECTIONS FOR USE

Required Equipment

- Endoscope (if desired)
- 0.038 in (0.97 mm), stiff-bodied 260 cm guidewire with floppy tip
- WallFlex Esophageal Fully Covered Stent System containing a stent of the appropriate length and diameter
- Fluoroscopic capability for pre-stent placement and stent placement confirmation
- Rat tooth forceps

Pre-Procedure

Radiography of the esophagus performed no more than 10 days before the procedure should be available.

Prepare for the procedure as you would an upper endoscopy. Administer a mild sedative, if necessary and a topical anesthesia to the throat to be repeated during the procedure as patient comfort warrants.

Initial Preparation of Delivery System

- Carefully remove the delivery system from the protective packaging.
- Visually inspect the device for damage or defects.

OPERATIONAL INSTRUCTIONS

Start of Initial Stent Placement Procedure

1. Locate the Stricture

If using an endoscope, intubate the patient using a standard endoscope per standard technique. Access the stricture location upon direct visualization. Fluoroscopy can also be used to locate the stricture with the aid of a contrast medium.

2. Examine the Stricture (Endoscopically and/or Fluoroscopically)

A. Examine Stricture Endoscopically

Endoscopically examine both the proximal and distal segments of the stricture. Using the external ruler on the endoscope, measure the distance between the distal margin of the stricture to the patient's incisors. Withdraw the endoscope to the proximal margin of the stricture and measure the distance to the patient's incisors. The stricture length is calculated as the difference between those two distances. To minimize the potential of stent migration, dilate the stricture ONLY if passage of the endoscope or the delivery system through the stricture lumen is not possible.

Caution: In some patients, tumor encroachment will make dilation of the stricture challenging. Physicians should use judgment based on experience in dilating esophageal strictures. Perforation or bleeding of an esophageal tumor is a risk during a tumor dilation procedure.

Warning: Placement of the WallFlex Esophageal Fully Covered Stent should not be attempted in patients with esophageal strictures that cannot be dilated wide enough for passage of the endoscope or delivery system.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Esophageal FC, Global, 90599988-01A

B. Examine the Stricture Fluoroscopically

The stricture may also be examined fluoroscopically. Leaving the endoscope in place, observe both the proximal and distal margins of the tumor fluoroscopically. Mark the locations with either radiopaque markers or use anatomical landmarks such as ribs or vertebras. It is recommended to re-measure the stricture length by measuring the distance between the radiopaque markers.

3. Choose the Stent Size

The size of the stricture must be accurately calculated to ensure the ideal stent size is used. The WallFlex™ Esophageal Fully Covered Stent should bridge the tumor and/or the fistula and should extend >1 cm above and below the stricture or fistula. For stent use with a fistula, it is critical to ensure that the stent completely covers the fistula to avoid leakage and facilitate healing. If stent length choice is questionable, always use the longer stent. A second stent of the same diameter may be placed if the first stent does not cover the entire stricture length. The second stent should be placed to ensure complete tumor coverage and a smooth transition between the stents. It is recommended that the proximal stent be placed first followed by the distal stent in order to maximize luminal diameter of the interlocked stents. Care should be used when passing the delivery system through the first stent.

Warning: Passing the scope through a just deployed stent is not recommended and could cause the stent to dislodge.

Warning: Use only in combination with WallFlex Esophageal Fully Covered Stents.

4. Insert Guidewire and Place Through Stricture

Pass a guidewire through the working channel of the endoscope and then through the stricture and into the stomach. A floppy tip guidewire is recommended in order to reduce potential trauma from the tip of the wire. Endoscopic and fluoroscopic placement of the guidewire is also recommended to assure proper passage through the stricture and proper placement in the stomach. Maintain guidewire position throughout procedure.

Caution: A stiff bodied 0.038 in (0.97 mm) guidewire with a floppy tip is recommended to facilitate passage through tortuous anatomy. The Jagwire™ Guidewire M00556621 is recommended.

5. Advance the Delivery System over the Guidewire and Position Stent


Figure 2. Delivery System and Radiopaque (RO) Markers

There are four radiopaque (RO) markers to aid in the deployment of the stent while using fluoroscopy. Radiopaque marker bands attached to the inner tube indicate the length of the constrained stent on the delivery system (Figure 2, #1 and #3). A marker band on the inner catheter between the constrained stent markers indicates at what point reconstrainment is no longer possible, Figure 2, #2. Up to this point, the stent can be reconstrained and repositioned if desired, no more than two times. The fourth radiopaque marker band at the distal end of the exterior catheter indicates up to what point the stent has been deployed (Figure 2, #4). The tip and inner tube are also radiopaque for use with fluoroscopy. After placing the guidewire, remove

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Esophageal FC, Global, 90599988-01A

the endoscope from the patient leaving the guidewire in place. Reinsert the endoscope alongside the guidewire if direct visualization of stent deployment is desired. Tilt the patient's head to provide the straightest pathway for advancement of the delivery system over the guidewire. Lubrication can be applied if desired.

Using fluoroscopic guidance, position the stent with the proximal RO marker above the proximal tumor margin and the distal RO marker below the distal tumor margin, centering the tumor between the markers (Figure 2, #1 and #3). This ensures that the stent will properly bridge the tumor. If it is not necessary to cross the Lower Esophageal Sphincter (LES), the distal end of the stent should remain above the LES in order to leave the LES functional and reduce gastric reflux. The stent can cross the LES if necessary due to tumor involvement and stricture.

6. Deploy Stent


Figure 3. Delivery system, Visual Markers and Handles

Caution: Do not twist the delivery system or use a boring motion during stent deployment as this may affect stent positioning and ultimately, the stent function.

Begin stent deployment by holding the distal handle (farthest from the operator, Figure 3, #7) of the delivery system with one hand, and using the other hand, grasp the proximal handle (closest to the operator, Figure 3, #8) and hold this handle stationary. Between the handles is a yellow segment with visual markers. These visual markers are intended to aid in stent deployment. Prior to deploying the stent, you can visualize a marker at the fully constrained/ undeployed position (Figure 3, #5).

In order to deploy the stent, hold the distal handle (farthest from the operator, Figure 3, #7) with one hand and the proximal handle (handle closest to the operator, Figure 3, #8) with the other hand. To deploy the stent, slowly pull the distal handle toward the proximal handle while holding the proximal handle stationary. Monitor the stent release fluoroscopically and/or endoscopically, keeping the fluoroscopic markers on the delivery system between the identified stricture margins. If necessary, it is possible to stop deployment and adjust the stent position proximally prior to passing the reconstrainment marker (Figure 3, #6). Refer to the reconstrainment technique section. If satisfied with stent placement, proceed with full deployment.

Reconstrainment Technique

- The stent can be reconstrained at any point up to the reconstrainment markers (radiopaque Figure 2, #2 and visual Figure 3, #6).

Note: Once the visual reconstrainment marker cannot be seen during deployment the stent cannot be reconstrained.

- Reconstrainment is done by reversing the direction of deployment, by holding the proximal handle (closest to the operator) steady while pushing the distal handle (farthest from the operator) away.
- The stent has been designed to be reconstrained no more than two times.
- Prior to full deployment, if repositioning is desired, the stent can be pulled proximally by slowly pulling back on

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Esophageal FC, Global, 90599988-01A

the delivery system. The ability to pull proximally will be limited by the amount of stent deployed and the tightness of the stricture. Full reconstrainment, when possible, is always preferred and recommended over pulling the device proximally.

Note: The stent is fully constrained if visual marker (Figure 3, #5) is fully visible.

Caution: Pulling proximally when partially deployed could further deploy the stent if there is resistance on the stent.

Caution: Do not push the delivery system forward once deployment has begun. The delivery system can be pulled proximally if necessary. The ability to pull proximally will be limited by the amount of stent deployed and the tightness of the stricture.

Warning: Stent is considered to be a permanent device. Once stent placement is permanently achieved, stent removal or repositioning is not recommended.

If the positioning of a WallFlex™ Esophageal Fully Covered Stent is not correct and one of the following has occurred, continue to fully deploy the stent:

A. The stent has already been deployed past the reconstrainment limit

OR

B. The stent has already been reconstrained twice

Then in either case, using rat tooth forceps, grasp the suture on the proximal end of the stent (Figure 1). Gently pull the stent back with the scope to remove the stent during the initial stent placement procedure.

Note: In a limited evaluation of the WallFlex Esophageal Partially Covered Stent in a porcine model, 6 stents were successfully removed during the initial stent placement procedure from 6 separate animals using the methods described above.

Warning: Stent is considered to be a permanent device. Once stent placement is permanently achieved, stent removal or repositioning is not recommended.

7. Assess Deployed Stent Position and Remove the Delivery System

Following stent deployment, view the stent endoscopically to confirm stent expansion as tumor impingement may prevent the stent from achieving its maximum diameter immediately.

Carefully remove the delivery system and the guidewire.

Caution: An attempt to remove the delivery system and guidewire prior to stent expansion or when a stent is partially deployed may dislodge the stent. If excessive resistance is felt during delivery system removal, wait 3-5 minutes to allow further stent expansion then proceed with the following steps:

A. Slowly withdraw the delivery system and guidewire.

B. If removal is still not possible, use a balloon dilation catheter to dilate the stent. It should not be necessary for the balloon diameter/size to be equal to the stent diameter. Judgement should be used when selecting the balloon size. Carefully position the balloon catheter within the stent. Inflate the balloon to its recommended pressure.

C. Deflate the balloon catheter and withdraw into the endoscope. Slowly withdraw the delivery system and guidewire.

Boston Scientific (Master Brand DFUTemplate 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Esophageal FC, Global, 90599888-01A

Warning: A stent may require 24-72 hours to expand fully. Physicians should use judgement based on experience when dilating esophageal strictures. Perforation or bleeding of the esophageal tumor is a risk during any dilation of the tumor. Never use a rigid-type dilator for post stent placement dilation because the axial force may dislodge the stent.

Warning: Once stent is in desired location, passing the scope through a just deployed stent is not recommended and could cause the stent to dislodge.

8. Remove Endoscope

Withdraw endoscope from the patient.

This completes the initial stent placement procedure. Stent placement is considered permanent upon completion of initial stent placement procedure.

POST-PROCEDURE

Patients should have P-A (postoanterior) and lateral chest films taken as a permanent record of the stent position. Observe the patient for development of any complications of endoscopy, esophageal dilation and stent placement. Vital signs should be monitored and clear liquids given in an upright position during the first 24 hours post-stent placement. Patients being treated for fistula should have no fluids or foods orally until successful sealing of the fistula has been confirmed. Post 24 hours, the patient should be instructed to eat only in an upright sitting position, to chew food thoroughly, to avoid certain foods (such as meats, raw vegetables, and breads), and to drink fluids during and after meals. Patients with stents placed in the distal esophagus or across the LES should be instructed to elevate the head of the bed, and should be prescribed acid suppression therapy to minimize gastric reflux into the stent. Subsequent follow-up at 1 week and at 3 month intervals, thereafter, or for symptomatic dysphagia may be performed to verify patency and placement.

Note: Recurrence or worsening of dysphagia may occur after stent placement due to tumor ingrowth or overgrowth over time, severe hyperplasia reaction or stent migration. Repeat endoscopy may be required.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or sterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

Excite is a Trademark of General Electric Company Corporation.

Boston Scientific (Master Brand DFU Template 3in x 9in Global, 90106040AK), DFU, MB, WallFlex Esophageal FC, Global, 90599988-01A

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All rights reserved.
Boston Scientific International S.A.
Boston Scientific (Australia) Pty Ltd
© 2010 Boston Scientific Corporation or its affiliates.