Boston Scientific

Trapezoid™ RX

Wireguided Retrieval Basket

Directions for Use	2
Instrucciones de uso	7
Mode d'emploi	10
Gebrauchsanweisung	13
Istruzioni per l'uso	16
Gebruiksaanwijzing	19
Instruções de Utilização	22

90571837-01

2010-06

Boston Scientific, (Master Brand DFU Template 5in x 5.25in Global, 90348240AD), DFU, MB, Trapezoid, Global, 90571837-01A

Trapezoid™ RX

Wireguided Retrieval Basket

℞ ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Trapezoid RX Wireguided Retrieval Basket is a device used to facilitate the removal of biliary calculi. Once advanced through an endoscope and into the common bile duct, this device serves as a retrieval apparatus with crushing capability whereby complete stone fragments may be removed. Please refer to the package label for dimensions specifics of the basket and working length of the sheath.

INDICATIONS FOR USE

The Trapezoid RX Wireguided Retrieval Basket is indicated for endoscopic crushing and removal of biliary calculi. The lithotripter compatible basket should not be used for any other purpose other than its intended application.

Important: In the event the biliary calculi cannot be crushed, the Trapezoid RX Wireguided Retrieval Basket has been designed for the basket tip to disengage minimizing the potential for stone entrapment in the basket. Flushing may follow.

CONTRAINDICATIONS

The Trapezoid RX Wireguided Retrieval Basket is contraindicated whenever Endoscopic Retrograde Cholangiopancreatography (ERCP) or Endoscopic Retrograde Sphincterotomy (ERS) would be contraindicated.

Caution: Trapezoid RX Wireguided Retrieval Basket is not intended for use in the pancreas.

Boston Scientific, (Master Brand DFU Template 6in x 5.25in Global, 90348240AD), DFU, MB, Trapezoid, Global, 90571837-01A

ADVERSE EVENTS

Possible complications include, but may not be limited to:

- Pancreatitis
- Bleeding
- Stone impaction
- Cholangitis
- Septicemia
- Perforation by catheter tip

PRECAUTIONS

Do not use this product without reading and understanding the complete instructions enclosed herein.

PRIOR TO USE

Precautions:

This device should only be used by or under the supervision of physicians trained in ERCP with mechanical lithotripsy. A thorough understanding of the technical principles, clinical applications and risks associated with ERCP and with mechanical lithotripsy in necessary before using this device.

This device is supplied sterile. Carefully examine the unit to verify that neither the contents nor the sterile package has been damaged in shipment. Do not use if damaged. Immediately return damaged product to Boston Scientific Corporation.

Any use of this device, other than those indicated in these instructions, is not recommended.

This device is for use with an endoscope with a working channel that is 3.2 mm or greater.

Preparation:

Prepare the designated endoscope and the Trapezoid™ RX Wireguided Retrieval Basket.

1. Open the package and inspect the device.
2. Verify smooth open/close movement of the basket by actuating the handle.
3. The Trapezoid RX Wireguided Retrieval Basket allows for contrast injection during operation. With the basket tip protruding approximately 5 mm from the sheath, attach a contrast filled syringe to the luer lock. Ensure that the contrast flows from the sheath's distal end when the syringe plunger is pushed.
4. Visually inspect the device's external surfaces for broken, crushed, kinked, or excessively bent areas.

DIRECTIONS FOR USE

Insertion:

1. Once the endoscope is in position, fully raise the elevator of the endoscope.
2. Retract the basket tip.
3. Carefully insert the device into the endoscope instrument channel.

Note: This may be accomplished as an over the wire procedure or without a guidewire. Due to various tolerances among manufacturers, it is recommended that the Trapezoid RX Wireguided Retrieval Basket be used with a Boston Scientific 0.035 in (0.89 mm) guidewire.

Boston Scientific, (Master Brand DFU Template 5in x 5.25in Global, 90348240AD), DFU, MB,Trapezoid, Global, 90571837-01A

4. Lower the elevator when the basket tip touches the elevator.
5. Carefully advance the catheter out of the endoscope approximately 3 mm. Raise the elevator to bring the basket tip into view. Advance the catheter carefully to the duodenal papillary opening.
6. While observing the duodenal papilla endoscopically, insert the catheter tip through the papilla into the common bile duct.

Caution: If excessive resistance is encountered and insertion becomes difficult, lower the elevator, advance the catheter and then operate the elevator again. Monitor the catheter under fluoroscopy to ensure the tip is at the proper angle to minimize the force against the wall of the bile duct.

Injecting Contrast:

1. For cholangiography, attach a 10 ml (cc) syringe filled with contrast medium to the female luer on the handle. Advance the basket approximately 5 mm from the tip of the sheath and inject contrast medium according to standard practices.

Basketing the Calculus:

1. Advance the catheter beyond the calculus to be captured.
2. Slowly open the basket.
3. Entrap the calculus by pulling the basket over the stone.

Note: Small movements of the device may be necessary to capture the calculus.

Note: For wire-guided procedure, it may be necessary to withdraw the wire into the wire guide to facilitate the capture of the calculus.

Crushing the Calculus:

In the event that the physician has decided that the calculus in the basket cannot be removed through the papilla, a mechanical advantage may be applied with the Alliance™ II System to fragment the calculus.

Note: The Trapezoid™ RX Wireguided Retrieval Basket 3 cm Basket is designed for crushing calculus larger than 1.5 cm (15 mm) in diameter.

1. Attach the Trapezoid RX Wireguided Retrieval Basket to the Alliance II System by aligning the thumb ring and the handle ring with the heel plate and the body portion of the Alliance II System respectively.

Note: For easy attachment, the Alliance II System directional switch should be pointed upward to allow free movement of the handle lever.

2. Push the Trapezoid RX Wireguided Retrieval Basket handle down into the Alliance II System until secure.

Boston Scientific, (Master Brand DFU Template 6in x 5.25in Global, 90348240AD), DFU, MB,Trapezoid, Global, 90571837-01A

Note: The Alliance™ II Heel Clip is not intended for use with the Trapezoid™ RX Wireguided Retrieval Basket.

3. Turn the Alliance II System directional switch toward the front of the device. Slowly squeeze the Alliance II System handle until the calculus is fragmented.

Note: Once the calculus is fragmented, stop squeezing the handle of the Alliance II System to avoid disengaging the Trapezoid RX Wireguided Retrieval Basket tip.

Note: In the event the biliary calculi cannot be crushed, the Trapezoid RX Wireguided Retrieval Basket has been designed for the basket tip to disengage minimizing the potential for the unreleasable stone entrapment. Flushing may follow.

4. To extend the Trapezoid RX Wireguided Retrieval Basket, turn the directional switch towards the back of the Alliance II System and slowly squeeze the handle until the basket is fully extended.
5. Detach the Trapezoid RX Wireguided Retrieval Basket from the Alliance II System by returning the directional switch to the upward position and pull upwards on the Trapezoid RX Wireguided Retrieval Basket handle.

Withdrawing the Trapezoid RX Wireguided Retrieval Basket

Lower the endoscope elevator and remove the device from the endoscope with the basket in the closed position.

EMERGENCY PRECAUTIONS

In the event that the **Trapezoid RX Wireguided Retrieval Basket will not crush the calculus** and the basket tip has not disengaged, use the following guidelines:

- If the calculus **can** be removed:
 Release the calculus from the basket, retract the basket into the catheter and remove the Trapezoid device.
- If calculus can not be removed from the basket, follow clinical and/or surgical standards of practice.

 The Sohendra® Lithotripsy System may be utilized for endoscopic calculus removal. Using wire cutting pliers, cut the handle off the basket sheath, leaving the basket wire behind with the calculus in situ followed by the Sohendra Lithotripsy System or other lithotripsy.

If the *basket wires are broken:*

Slowly withdraw the basket into the catheter and withdraw the catheter from the scope.

- In the event the **damaged part of the wire cannot be withdrawn** into the catheter, be careful not to injure the duct wall with the basket wire.

Boston Scientific, (Master Brand DFU Template 5in x 5.25in Global, 90348240AD), DFU, MB,Trapezoid, Global, 90571837-01A

- In the event that **all four basket wires are broken** and the **distal tip of the basket has come off,** withdraw basket wires into catheter and remove. Flushing may follow.

Precaution: The Trapezoid™ RX Wireguided Retrieval Basket should be advanced through the endoscope using short, deliberate 2-3 cm (20-30 mm) strokes to prevent inadvertent damage to the device.

Note: The Emergency Precautions section outlines a general technique for mechanical lithotripsy. The endoscopist should follow standard surgical considerations to decide on technical details of the procedure from a professional medical viewpoint.

DEVICE DISPOSAL

After procedure, dispose of the Trapezoid RX Wireguided Retrieval Basket in accordance with hospital protocol.

HOW SUPPLIED

The Trapezoid RX Wireguided Retrieval Basket is packaged one (1) unit per sealed package. The package is supplied sterile by ethylene oxide gas. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Handling and Storage

Store in a cool, dry, dark place.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose**. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments**.

Soehendra is a trademark of Wilson-Cook Medical, Inc.

Boston Scientific, (Master Brand DFU Template 6in x 5.25in Global, 90348240AD), DFU, MB,Trapezoid, Global, 90571837-01A

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All rights reserved.
Boston Scientific International S.A.
Boston Scientific (Australia) Pty Ltd
© 2010 Boston Scientific Corporation or its affiliates.