BD CD22 (S-HCL-1)

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Referencias

1. Foon KA, Gale RP, Todd RF III. Recent advances in the immunologic classification of leukemia. *Semin Hematol.* 1986;23:257-283.
2. Mason DY, Stein H, Gerdes J, et al. Value of monoclonal anti-CD22 (p135) antibodies for the detection of normal and neoplastic B lymphoid cells. *Blood*. 1987;69:836-840.
3. Schwarting R, Stein H, Wang CY. The monoclonal antibodies alpha S-HCL 1 (alpha Leu-14) and alpha S-HCL 3 (alpha Leu-M5) allow the diagnosis of hairy cell leukemia. *Blood*. 1985;65:974-983.
4. van Dongen JJ, Lhermitte L, Böttcher S, et al. EuroFlow antibody panels for standardized n-dimensional flow cytometric immunophenotyping of normal, reactive and malignant leukocytes. *Leukemia*. 2012;26:1908-1975.
5. Loken MR, Shah VO, Hollander Z, Civin CI. Flow cytometric analysis of normal B lymphoid development. *Pathol Immunopathol Res*. 1988;7:357-370.
6. Nadler LM. B Cell/Leukemia Panel Workshop: Summary and Comments. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human B Lymphocytes*. Vol 2. New York, NY: Springer-Verlag; 1986:3-43.
7. Zola H, Swart B, Nicholson I, Voss E. *Leukocyte and Stromal Cell Molecules: The CD Markers. Hoboken*, NJ: John Wiley & Sons, Inc.; 2007.
8. Terstappen LW, Johnsen S, Segers-Nolten IM, Loken MR. Identification and characterization of plasma cells in normal human bone marrow by high-resolution flow cytometry. *Blood*. 1990;76:1739-1747.
9. Borowitz MJ, Bonsvaros A, Brynes RK, et al. Monoclonal antibody phenotyping of B-cell non-Hodgkin's lymphomas: The Southeastern Cancer Study Group experience. *Am J Pathol*. 1985;121:514-521.

10. Law CL, Sidorenko SP, Clark EA. Regulation of lymphocyte activation by the cell-surface molecule CD22. *Immunol Today*. 1994;15:442-449.

11. Doody GM, Dempsey PW, Fearon DT. Activation of B lymphocytes: integrating signals from CD19, CD22 and Fc gamma RIIb1. *Curr Opin Immunol*. 1996;8:378-382.

12. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI document M29-A4.

13. Centers for Disease Control and Prevention. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html. Accessed March 12, 2019.

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Información de contacto

Becton, Dickinson and Company
BD Biosciences
155 North McCarthy Boulevard
Milpitas, California 95035 EE. UU.

bdbiosciences.com
ResearchApplications@bd.com

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BD y el logotipo de BD son marcas comerciales de Becton, Dickinson and Company o sus filiales. Todas las demás marcas comerciales son propiedad de sus respectivos dueños. © 2024 BD. Todos los derechos reservados.