[LOGO] BD

CD11b (D12)

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Referencias

1. Bray R, Gottschalk L, Landay A, Gebel H. Differential surface marker expression in patients with $CD16^+$ lymphoproliferative disorders: In vivo model for NK differentiation. *Human Immunol*. 1987;19:105.

2. Clement LT, Grossi CE, Gartland FL. Morphologic and phenotypic features of the subpopulation of Leu-2+ cells which suppress T cell proliferative responses. *J Immunol*. 1984;131:2757.

3. Gebel HM, Kaizer H, Landay AL. Characterization of circulating suppressor T lymphocytes in bone marrow transplant recipients. *Transplantation*. 1987;43:258.

4. Ross GD, Cain JA, Lachmann PJ. Membrane complement receptor type three (CR_3) has lectin-like properties analogous to bovine conglutinin and functions as a receptor for zymosan and rabbit erythrocytes as well as a receptor for iC3b. *J Immunol*. 1985;134:3307.

5. Patarroyo M, Makgoba MW. Leucocyte adhesion to cells: molecular basis, physiological relevance, and abnormalities. *Scand J Immunol*. 1989;30:129-164.

6. Lanier LL, Phillips JH. A map of the cell surface antigens expressed on resting and activated human natural killer cells. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human Myeloid and Hematopoietic Cells*. New York, NY: Springer-Verlag; 1986;3:157-170.

7. Landay A, Gartland L, Clement LT. Characterization of a phenotypically distinct subpopulation of $Leu\text{-}2^+$ cells which suppress T cell proliferative responses. *J Immunol*. 1983;131:2757.

8. Bernstein ID, Self S. Joint report of the myeloid section of the Second International Workshop on Human Leukocyte Differentiation Antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human Myeloid and Hematopoietic Cells*. New York, NY: Springer-Verlag; 1986;3:1-25.

9. Shalekoff S, Page-Shipp L, Tiemessen CT. Effects of anticoagulants and temperature on expression of activation markers CD11b and HLA-DR on human leukocytes. *Clin Diagn Lab Immunol*. 1998;5:695-702.

10. Repo H, Jansson S-E, Leirisalo-Repo M. Anticoagulant selection influences flow cytometric determination of CD11b upregulation in vivo and ex vivo. *J Immunol*. 1995;185:65-79.

11. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI document M29-A4.

12. Centers for Disease Control and Prevention. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html. Accessed March 12, 2019.

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Información de Contacto

Becton, Dickinson and Company
BD Biosciences
155 North McCarthy Boulevard
Milpitas, California 95035, EE. UU.

bdbiosciences.com
ResearchApplications@bd.com

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BD, el logotipo de BD y de Horizon son marcas comerciales de Becton, Dickinson and Company o de sus filiales. Todas las demás marcas comerciales pertenecen a sus respectivos propietarios. © 2024 BD. Todos los derechos reservados.