BD CD3 (SK7)

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Referencias

1. Kaneoka H, Perez-Rojas G, Sasasuki T, Benike CJ, Engleman EG. Human T lymphocyte proliferation induced by a pan-T monoclonal antibody (anti-Leu-4): heterogeneity of response is a function of monocytes. *J Immunol*. 1983;131:158-164.

2. Allison JP, Lanier LL. Structure, function, and serology of the T-cell antigen receptor complex. *Ann Rev Immunol*. 1987;5:503-540.

3. Friedrich W, O'Reilly RJ, Kozinger B, Gebhard DF, Good RA, Evans RL. T lymphocyte reconstitution in recipients of bone marrow transplants with and without GVHD: imbalances of T-cell subpopulations having unique regulatory and cognitive functions. *Blood*. 1982;59:696-701.

4. Kurrle R. Cluster report: CD3. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:290-293.

5. Maecker HT, Maino VC. Analyzing T-cell responses to cytomegalovirus by cytokine flow cytometry. *Hum Immunol*. 2004;65:493-499.

6. Maino VC, Picker LJ. Identification of functional subsets by flow cytometry: intracellular detection of cytokine expression. *Cytometry*. 1998;34:207-215.

7. Campbell MJ, Scott J, Maecker HT, Park JW, Esserman LJ. Immune dysfunction and micrometastases in women with breast cancer. *Breast Cancer Res Treat*. 2005;91:163-171.

8. Tu W, Chen S, Sharp M, et al. Persistent and selective deficiency of $CD4^{+}$T cell immunity to cytomegalovirus in immunocompetent young children. *J Immunol*. 2004;172:3260-3267.

9. van Dongen JJM, Krissansen GW, Wolvers-Tettero ILM, et al. Cytoplasmic expression of the CD3 antigen as a diagnostic marker for immature T-cell malignancies. *Blood*. 1988;71:603-612.

10. Brenner MB, Groh V, Porcelli SA, et al. Structure and distribution of the human γ/δ T-cell receptor. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:1049-1053.

11. Clevers H, Alarcón B, Wileman T, Terhorst C. The T cell receptor/CD3 complex: a dynamic protein ensemble. *Annu Rev Immunol*. 1988;6:629-662.

12. Reichert T, DeBruyère M, Deneys V, et al. Lymphocyte subset reference ranges in adult Caucasians. *Clin Immunol Immunopathol*. 1991;60:190-208.

13. Lanier LL, Allison JP, Phillips JH. Correlation of cell surface antigen expression on human thymocytes by multi-color flow cytometric analysis: implications for differentiation. *J Immunol*. 1986;137:2501-2507.

14. Garson JA, Beverley PCL, Coakham HB, Harper EI. Monoclonal antibodies against human T lymphocytes label Purkinje neurones of many species. *Nature*. 1982;298:375-377.

15. Haynes BF. Summary of T-cell studies performed during the Second International Workshop and Conference on Human Leukocyte Differentiation Antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human T Lymphocytes*. Vol 1. New York, NY: Springer-Verlag; 1986:3-30.

16. Knowles RW. Immunochemical analysis of the T-cell–specific antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human T Lymphocytes*. Vol 1. New York, NY: Springer-Verlag; 1986:259-288.

17. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI document M29-A4.

18. Centers for Disease Control and Prevention. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html. Accessed March 12, 2019.

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Información de Contacto

Becton, Dickinson and Company
BD Biosciences
155 North McCarthy Boulevard
Milpitas, California 95035 EUA

bdbiosciences.com
ResearchApplications@bd.com

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BD, el logotipo de BD y Horizon son marcas registradas de Becton, Dickinson and Company o sus filiales. Todas las demás marcas registradas son propiedad de sus respectivos dueños. © 2024 BD. Todos los derechos reservados.