BD CD5 (L17F12)

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Referencias

1. Antin JH, Emerson SG, Martin P, Gadol N, Ault KA. Leu-1^{+} (CD5^{+}) B cells. A major lymphoid subpopulation in human fetal spleen: phenotypic and functional studies. *J Immunol*. 1986;136:505-510.

2. Warnke R, Miller R, Grogan T, Pederson M, Dilley J, Levy R. Immunologic phenotype in 30 patients with diffuse large-cell lymphoma. *N Eng J Med*. 1980;303:293-300.

3. Zipf TF, Fox RI, Dilley J, Levy R. Definition of the high risk acute lymphoblastic leukemia patient by immunological phenotyping with monoclonal antibodies. *Cancer Res*. 1981;41:4786-4789.

4. Royston I, Majda JA, Baird SM, Meserve BL, Griffiths JC. Human T-cell antigens defined by monoclonal antibodies: the 65,000-dalton antigen of T cells (T65) is also found on chronic lymphocytic leukemia cells bearing surface immunoglobulin. *J Immunol*. 1980;125:725-731.

5. Gadol N, Ault KA. Phenotypic and functional characterization of human Leu1 (CD5) B cells. *Immunol Rev*. 1986;93:23-24.

6. Hardy RR, Hayakawa K, Shimizu M, Yamasaki K, Kishimoto T. Rheumatoid factor secretion from human Leu-1^{+} B cells. *Science*. 1987;236:81-83.

7. Knowles RW. Immunochemical analysis of the T-cell–specific antigens. In: Reinherz EL, Haynes BF, Nadler LM, Bernstein ID, eds. *Leukocyte Typing II: Human T Lymphocytes*. Vol 1. New York, NY: Springer-Verlag; 1986:259-288.

8. Engleman EG, Warnke R, Fox RI, Dilley J, Benike CJ, Levy R. Studies of a human T lymphocyte antigen recognized by a monoclonal antibody. *Proc Natl Acad Sci USA*. 1981;78:1791-1795.

9. Ledbetter JA, Evans RL, Lipinski M, Cunningham-Rundles C, Good RA, Herzenberg LA. Evolutionary conservation of surface molecules that distinguish T lymphocyte helper/inducer and cytotoxic/suppressor subpopulations in mouse and man. *J Exp Med*. 1981;153:310-323.

10. Ledbetter JA, Frankel AE, Herzenberg LA, Herzenberg LA. Human Leu T-cell differentiation antigens: quantitative expression on normal lymphoid cells and cell lines. In: Hämmerling G, Hämmerling U, Kearney J, eds. *Monoclonal Antibodies and T-Cell Hybridomas: Perspectives and Technical Advances*. New York, NY: Elsevier/North Holland; 1981:16-22.

11. Warnke R, Levy R. Detection of T and B cell antigens hybridoma monoclonal antibodies: a biotin-avidin-horseradish peroxidase method. *J Histochem Cytochem*. 1980;28:771-776.

12. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI document M29-A4.

13. Centers for Disease Control and Prevention. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html. Accessed March 12, 2019.

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Información de Contacto

Becton, Dickinson and Company
BD Biosciences
155 North McCarthy Boulevard
Milpitas, California 95035 EUA

bdbiosciences.com
ResearchApplications@bd.com

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BD, el logotipo de BD y Horizon son marcas registradas de Becton, Dickinson and Company o sus filiales. Todas las demás marcas registradas son propiedad de sus respectivos dueños. © 2024 BD. Todos los derechos reservados.