BD CD34 (8G12)

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2. Leary AG, Strauss LC, Civin CI, Ogawa M. Disparate differentiation in hemopoietic colonies derived from human paired progenitors. *Blood*. 1985;6:327-332.

3. Terstappen LW, Huang S, Safford M, Lansdorp PM, Loken MR. Sequential generations of hematopoietic colonies derived from single nonlineage-committed $CD34^{+}CD38^{-}$ progenitor cells. *Blood*. 1991;77:1218-1227.

4. Lansdorp PM, Dougherty GJ, Humphries RK. CD34 epitopes. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:826-827.

5. Loken MR, Shah VO, Dattilio KL, Civin CI. Flow cytometric analysis of human bone marrow. II. Normal Blymphocyte development. *Blood*. 1987;70:1316-1324.

6. Peschel C, Köller U. Cluster report: CD34. In: Knapp W, Dörken B, Gilks W, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:817-818.

7. Siena S, Bregni M, Brando B, et al. Flow cytometry for clinical estimation of circulating hematopoietic progenitors for autologous transplantation in cancer patients. *Blood*. 1991;77:400-409.

8. Brocklebank AM, Sparrow RL. Enumeration of CD34+ cells in cord blood: a variation on a single-platform flow cytometric method based on the ISHAGE gating strategy. *Cytometry*. 2001;46:254-261.

9. Hurwitz CA, Loken MR, Graham ML, et al. Asynchronous antigen expression in B lineage acute lymphoblastic leukemia. *Blood*. 1988;72:299-307.

10. Ryan D, Kossover S, Mitchell S, Frantz C, Hennessy L, Cohen H. Subpopulations of common acute lymphoblastic leukemia antigen-positive lymphoid cells in normal bone marrow identified by hematopoietic differentiation antigens. *Blood*. 1986;68:417-425.

11. Terstappen LWMM, Safford M, Könemann S, et al. Flow cytometric characterization of acute myeloid leukemia. Part II. Phenotypic heterogeneity at diagnosis. *Leukemia*. 1992;6:70-80.

12. Gore SD, Kastan MB, Civin CI. Normal human bone marrow precursors that express terminal deoxynucleotidyl transferase include T-cell precursors and possible lymphoid stem cells. *Blood*. 1991;77:1681-1690.

13. Andrews RG, Singer JW, Bernstein ID. Precursors of colony-forming cells in humans can be distinguished from colony-forming cells by expression of the CD33 and CD34 antigens and light scattering properties. *J Exp Med*. 1989;169:1721-1731.

14. Loken MR, Shah VO, Dattilio KL, Civin CI. Flow cytometric analysis of human bone marrow. I. Normal erythroid development. *Blood*. 1987;69:255-263.

15. Kurtzberg J, Denning SM, Nycum LM, Singer KH, Haynes BF. Immature human thymocytes can be driven to differentiate into nonlymphoid lineages by cytokines from thymic epithelial cells. *Proc Natl Acad Sci USA*. 1989;86:7575-7579.

16. Civin CI, Trischmann TM, Fackler MJ, et al. Report on the CD34 cluster workshop. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. New York, NY: Oxford University Press; 1989:818-825.

17. Greaves MF, Titley I, Colman SM, et al. CD34 cluster workshop report. In: Schlossman SF, Boumsell L, Gilks W, et al, eds. *Leucocyte Typing V: White Cell Differentiation Antigens*. Vol 1. New York, NY: Oxford University Press; 1995:840-846.

18. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI document M29-A4.

19. Centers for Disease Control and Prevention. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html. Accessed March 12, 2019.

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Información de Contacto

Becton, Dickinson and Company
BD Biosciences
155 North McCarthy Boulevard
Milpitas, California 95035 EUA

bdbiosciences.com
ResearchApplications@bd.com

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BD, el logotipo de BD y Horizon son marcas registradas de Becton, Dickinson and Company o sus filiales. Todas las demás marcas registradas son propiedad de sus respectivos dueños. © 2024 BD. Todos los derechos reservados.