CD57 (HNK-1)

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1. Lanier LL, Le AM, Civin C, Loken MR, Phillips JH. The relationship of CD16 (Leu-11) and Leu-19 (NKH-1) antigen expression on human peripheral blood NK cells and cytotoxic T lymphocytes. *J Immunol*. 1986;136(4):4480-4486.
2. Lanier LL, Le AM, Phillips JH, Warner NL, Babcock GF. Subpopulations of human natural killer cells defined by expression of the Leu-7 (HNK-1) and Leu-11 (NK-15) antigens. *J Immunol*. 1983;131:1789-1796.
3. Manara DC, de Panfilis B, Ferrari C. Ultrastructural characterization of human large granular lymphocyte subsets defined by the expression of HNK-1 (Leu-7), Leu-11, or both HNK-1 and Leu-11 antigens. *J Histochem Cytochem*. 1985;33:1129-1133.
4. Mechtersheimer G, Staudter M, Möller P. Expression of the natural killer cell-associated antigens CD56 and CD57 in human neural and striated muscle cells and in their tumors. *Cancer Res*. 1991;51:1300-1307.
5. Kruse J, Mailhammer R, Wennecke H, et al. Neural cell-adhesion molecules and MAG share a common carbohydrate moiety recognized by monoclonal antibodies L2 and HNK-1. *Nature*. 1984;311:153-155.
6. McGarry RC, Helfand SL, Quarles RH, Roder JC. Recognition of myelin-associated glycoprotein by the monoclonal antibody HNK-1. *Nature*. 1983;306:376.
7. Bunn PA, Gazdar AF, Carney DN, Minna J. Small cell lung carcinoma and natural killer cells share an antigenic determinant, Leu-7. *Clin Res*. 1984;32:416A.
8. Lewis DE, Puck JM, Babcock GM, Rich RR. Disproportionate expansion of a minor T cell subset in patients with lymphadenopathy syndrome and acquired immunodeficiency syndrome. *J Infect Disease*. 1985;151:555.
9. Prince HE, Kreiss JK, Kasper CK, et al. Distinctive lymphocyte subpopulation abnormalities in patients with congenital coagulation disorders who exhibit lymph node enlargement. *Blood*. 1985;66:64.
10. Stites D, Casavant C, McHugh T, et al. Flow cytometric analysis of lymphocyte phenotypes in AIDS using monoclonal antibodies and simultaneous dual immunofluorescence. *Clin Immunol Immunopathol*. 1986;38:161-177.
11. Abo T, Balch CMA. A differentiation antigen of human NK and K cells identified by a monoclonal antibody (HNK-1). *J Immunol*. 1981;127:1024-1029.
12. Schubert J, Lanier LL, Schmidt RE. Cluster report: CD57. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:711-714.

13. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI document M29-A4.
14. Centers for Disease Control and Prevention. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html. Accessed March 12, 2019.

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Información de contacto

Becton, Dickinson and Company
BD Biosciences
155 North McCarthy Boulevard
Milpitas, California 95035, EE. UU.

bdbiosciences.com
ResearchApplications@bd.com

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BD, el logotipo de BD y FACScan son marcas comerciales de Becton, Dickinson and Company o sus filiales. Todas las demás marcas comerciales son propiedad de sus respectivos propietarios. © 2024 BD. Todos los derechos reservados.