CD13 (L138)

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7. Lo Coco F, Pasqualetti D, Lopez M, et al. Immunophenotyping of acute myeloid leukemia: relevance of analysing different lineage-associated markers. *Blut.* 1989;58:235-240.

8. Sakai K, Hattori T, Sagawa K, Yokoyama M, Takatsuki K. Biochemical and functional characterization of MCS-2 antigen (CD13) on myeloid leukemic cells and polymorphonuclear leukocytes. *Cancer Res.* 1987;47:5572-5576.

9. van Dongen JJ, Lhermitte L, Böttcher S, et al. EuroFlow antibody panels for standardized n-dimensional flow cytometric immunophenotyping of normal, reactive and malignant leukocytes. *Leukemia.* 2012;26:1908-1975.

10. Lee SH, Crocker PR, Westaby S, et al. Isolation and immunocytochemical characterization of human bone marrow stromal macrophages in hemopoietic clusters. *J Exp Med.* 1988;168:1193-1198.

11. Look AT, Ashmun RA, Shapiro LH, et al. Report on the CD13 (aminopeptidase N) cluster Workshop. In: Knapp W, Dörken B, Gilks WR, et al., eds. *Leucocyte Typing IV: White Cell Differentiation Antigens.* New York, NY: Oxford University Press; 1989:784-787.

12. Pombo de Oliveira MS, Matutes E, Rani S, Morilla R, Catovsky D. Early expression of MCS2 (CD13) in the cytoplasm of blast cells from acute myeloid leukaemia. *Acta Haematol.* 1988;80:61-64.

13. Terstappen LW, Hollander Z, Meiners H, Loken MR. Quantitative comparison of myeloid antigens on five lineages of mature peripheral blood cells. *J Leukoc Biol.* 1990;48:138-148.

14. Zola H, Swart B, Nicholson I, Voss E. *Leukocyte and Stromal Cell Molecules: The CD Markers.* Hoboken, NJ: John Wiley & Sons, Inc.; 2007.

15. Bradstock KF, Favaloro EJ, Kabral A, Kerr A, Hughes WG, Musgrove E. Myeloid progenitor surface antigen identified by monoclonal antibody. *Br J Haematol.* 1985;61:11-20.

16. Kirshenbaum AS, Goff JP, Semere T, Foster B, Scott LM, Metcalfe DD. Demonstration that human mast cells arise from a progenitor cell population that is CD34(+), c-kit(+), and expresses aminopeptidase N (CD13). *Blood.* 1999;94:2333-2342.

17. Mina-Osorio P, Winnicka B, O'Conor C, et al. CD13 is a novel mediator of monocytic/endothelial cell adhesion. *J Leukoc Biol.* 2008;84:448-459.

18. Yeager CL, Ashmun RA, Williams RK, et al. Human aminopeptidase N is a receptor for human coronavirus 229E. *Nature.* 1992;357:420-422.

19. Ashmun RA, Holmes KV, Shapiro LH, et al. CD13 (aminopeptidase N) cluster workshop report. In: Schlossman SF, Boumsell L, Gilks W, et al., eds. *Leucocyte Typing V: White Cell Differentiation Antigens.* New York, NY: Oxford University Press; 1995:771-775.

20. Howard MR, Thomas L, Reid MM. Variable detection of myeloid antigens in childhood acute lymphoblastic leukaemia. *J Clin Pathol.* 1994;47:1006-1009.

21. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline — Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI document M29-A4.

22. Centers for Disease Control and Prevention. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html. Accessed March 12, 2019.

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Información de contacto

Becton, Dickinson and Company
BD Biosciences
155 North McCarthy Boulevard
Milpitas, California 95035 EE. UU.

bdbiosciences.com
ResearchApplications@bd.com

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BD, el logotipo de BD y Horizon son marcas comerciales de Becton, Dickinson and Company o sus filiales. Todas las demás marcas comerciales son propiedad de sus respectivos dueños. © 2024 BD