BD Anti–TCRγ/δ-1 (11F2)

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Referencias

1. Lanier LL, Ruitenberg JJ, Bolhuis RL, Borst J, Phillips JH, Testi R. Structural and serological heterogeneity of γ/δ T cell antigen receptor expression in thymus and peripheral blood. *Eur J Immunol*. 1988;18:1985-1992.
2. Urban EM, Chapoval AI, Pauza CD. Repertoire development and the control of cytotoxic/effector function in human gammadelta T cells. *Clin Dev Immunol*. 2010;2010:732893.
3. Voogt PJ, Falkenburg JH, Fibbe WE, et al. Normal hematopoietic progenitor cells and malignant lymphohematopoietic cells show different susceptibility to direct cell-mediated MHC-non-restricted lysis by T cell receptor^{-}/CD3^{-}, T cell receptor γ/δ^{+}/CD3^{+} and T cell receptor-α/β^{+}/CD3^{+} lymphocytes. *J Immunol*. 1989;142:1774-1780.
4. Blink SE, Miller SD. The contribution of gammadelta T cells to the pathogenesis of EAE and MS. Curr Mol Med. 2009;9:15-22.
5. Cairo C, Hebbeler AM, Propp N, Bryant JL, Colizzi V, Pauza CD. Innate-like gammadelta T cell responses to mycobacterium Bacille Calmette-Guerin using the public V gamma 2 repertoire in Macaca fascicularis. *Tuberculosis (Edinb).* 2007;87:373-383.
6. Chen ZW. Immune regulation of gammadelta T cell responses in mycobacterial infections. *Clin Immunol*. 2005;116:202-207.
7. Garcia VE, Sieling PA, Gong J, et al. Single-cell cytokine analysis of gamma delta T cell responses to nonpeptide mycobacterial antigens. *J Immunol*. 1997;159:1328-1335.
8. Huang D, Chen CY, Zhang M, et al. Clonal immune responses of Mycobacterium-specific gammadelta T cells in tuberculous and non-tuberculous tissues during M. tuberculosis infection. *PLoS One.* 2012;7:e30631.
9. Bonneville M, O'Brien RL, Born WK. Gammadelta T cell effector functions: a blend of innate programming and acquired plasticity. *Nat Rev Immunol*. 2010;10:467-478.
10. Carding SR, Egan PJ. The importance of gamma delta T cells in the resolution of pathogen-induced inflammatory immune responses. *Immunol Rev*. 2000;173:98-108.
11. Borst J, van Dongen J, Bolhuis R, et al. Distinct molecular forms of human T cell receptor γ/δ detected on viable T cells by a monoclonal antibody. *J Exp Med*. 1988;167:1625-1644.
12. Lanier LL, Serafini AT, Ruitenberg JJ, et al. The gamma T-cell antigen receptor. *J Clin Immunol*. 1987;7:429-440.
13. Testi R, Lanier LL. Functional expression of CD28 on T cell antigen receptor gamma/delta-bearing T lymphocytes. *Eur J Immunol*. 1989;19:185-188.
14. Ichinohasama R, Miura I, Takahashi T, et al. Peripheral CD4+ CD8- γδ T cell lymphoma: a case report with multiparameter analyses. *Hum Pathol*. 1996;27:1370-1377.
15. Lanier LL, Federspiel N, Ruitenberg JJ, et al. The T-cell antigen receptor complex expressed on normal peripheral blood CD4, CD8 T lymphocytes: A CD3-associated disulfide-linked γ chain heterodimer. *J Exp Med.* 1987;165:1076-1094.
16. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI document M29-A4.

17. Centers for Disease Control and Prevention. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html. Accessed March 12, 2019.

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Información de contacto

Becton, Dickinson and Company
BD Biosciences
155 North McCarthy Boulevard
Milpitas, California 95035 EE. UU.

bdbiosciences.com
ResearchApplications@bd.com

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BD, el logotipo de BD y Horizon son marcas comerciales de Becton, Dickinson and Company o sus filiales. Todas las demás marcas comerciales son propiedad de sus respectivos dueños. © 2024 BD. Todos los derechos reservados.