CD42a (Beb1)

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Referencias

1. von dem Borne AEGKr, Modderman PW. Cluster report: CD42. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:1000-1002.
2. Ruggeri ZM. The platelet glycoprotein Ib-IX complex. In: Coller B, ed. *Progress in Hemostasis and Thrombosis*. Philadelphia: WB Saunders Co. 1991:10:35-68.
3. Hickey MJ, Williams SA, Roth GJ. Human platelet glycoprotein IX: An adhesive prototype of leucine-rich glycoproteins with flank-center-flank structures. *Proc Natl Acad Sci USA*. 1989;86:6773-6777.
4. Du X, Beutler L, Ruan C, Castaldi PA, Berndt MC. Glycoprotein Ib and glycoprotein IX are fully complexed in the intact platelet membrane. *Blood*. 1987;69:1524-1527.
5. Hourdillé P, Heilmann E, Combrié R, et al. Thrombin induces a rapid redistribution of glycoprotein Ib-IX complexes within the membrane systems of activated human platelets. *Blood*. 1990;76:1503-1513.
6. Modderman PW, Admiraal LG, Sonnenberg A, von dem Borne AEGKr. Glycoproteins V and Ib-IX form a noncovalent complex in the platelet membrane. *J Biol Chem*. 1992;267:364-369.
7. Michelson AD, Benoit SE, Furman MI, et al. The platelet surface expression of glycoprotein V is regulated by two independent mechanisms: Proteolysis and a reversible cytoskeletal-mediated redistribution to the surface-connected canalicular system. *Blood*. 1996;87:1396-1408.
8. Michelson A, Ellis P, Barnard M, et al. Downregulation of the platelet surface glycoprotein Ib-IX complex in whole blood stimulated by thrombin, adenosine diphosphate, or an in vivo wound. *Blood*. 1991;77:770- 779.
9. Michelson AD, Benoit SE, Kroll MH, et al. The activation-induced decrease in the platelet surface expression of the glycoprotein Ib-IX complex is reversible. *Blood*. 1994;83:3562-3573.
10. Fox JEB, Aggerbeck LP, Berndt MC. Structure of the glycoprotein Ib-IX complex from platelet membranes. *J Biol Chem*. 1988;263:4882-4890.

11. Nieuwenhuis HK. Report on functional studies. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV: White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:1002-1003.
12. Yamamoto N, Greco NJ, Barnard MR, et al. Glycoprotein Ib (GPIb)-dependent and GPIb-independent pathways of thrombin-induced platelet activation. *Blood*. 1991;77:1740-1748.
13. Silverstein RL. Platelet antigens: Section report. In: Schlossman S, Boumsell L, Gilks W, ed. *Leucocyte Typing V: White Cell Differentiation Antigens*. Oxford: Oxford University Press. 1995:1197.
14. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline — Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI document M29-A4.
15. Centers for Disease Control and Prevention. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html. Accessed March 12, 2019.

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Información de contacto

Becton, Dickinson and Company
BD Biosciences
155 North McCarthy Boulevard,
Milpitas, California 95035, EE. UU.

bdbiosciences.com
ResearchApplications@bd.com

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BD y el logotipo de BD son marcas comerciales de Becton, Dickinson and Company o sus filiales. Todas las demás marcas comerciales son propiedad de sus respectivos dueños. © 2024 BD. Reservados todos los derechos.