BD CD61 (RUU-PL 7F12)

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Referencias

1. Abrams CS, Ellison N, Budzynski AZ, Shattil SJ. Direct detection of activated platelets and platelet-derived microparticles in humans. *Blood*. 1990;75:128-138.
2. Modderman PW. Cluster report: CD61. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV, White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:1025.
3. Parmentier S, Kaplan C, Catimel B, McGregor JL. New families of adhesion molecules play a vital role in platelet functions. *Immunol Today*. 1990;11:225-227.
4. von dem Borne AEGKr, Modderman PW, Admiraal LG, Nieuwenhuis KH. Platelet antibodies, the overall results. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV, White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:950-966.
5. Modderman PW. New clusters: CD29/CDw49, CD47, CD51, CD55, and CD61. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV, White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:1017-1019.
6. Springer TA. Adhesion receptors of the immune system. *Nature*. 1990;346:425-433.
7. Hynes RO. Integrins: Versatility, modulation, and signaling in cell adhesion. *Cell*. 1992;69:11-25.
8. Fitzgerald LA, Steiner B, Rall Jr SC, Lo S-S, Phillips DR. Protein sequence of endothelial glycoprotein IIIa derived from a cDNA clone. *J Biol Chem*. 1987;262:3936-3939.
9. Fijnheer R, Modderman PW, Veldman H, et al. Detection of platelet activation with monoclonal antibodies and flow cytometry. *Transfusion*. 1990;30:20-25.
10. von dem Borne AEGKr, Modderman PW. Cluster report: CD41. In: Knapp W, Dörken B, Gilks WR, et al, eds. *Leucocyte Typing IV, White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1989:997-999.
11. Jennings LK, Ashmun RA, Wang WF, Dockter ME. Analysis of human platelet glycoproteins IIb-IIIa and Glanzmann's thrombasthenia in whole blood by flow cytometry. *Blood*. 1986;68:173-179.

12. Shattil SJ, Hoxie JA, Cunningham M, Brass L. Changes in the platelet membrane glycoprotein IIb-IIIa complex during platelet activation. *J Biol Chem*. 1985;260:11107-11114.

13. Wong DA, Springer TA. CD61 (b3) cluster report. In: Schlossman SF, Boumsell L, Gilks W, et al, eds. *Leucocyte Typing V, White Cell Differentiation Antigens*. Oxford: Oxford University Press; 1995:1664-1665.

14. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline — Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI document M29-A4.

15. Centers for Disease Control and Prevention. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html. Accessed March 12, 2019.

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Información de contacto

Becton, Dickinson and Company
BD Biosciences
155 North McCarthy Boulevard
Milpitas, California 95035 EE. UU.

bdbiosciences.com
ResearchApplications@bd.com

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BD, el logotipo de BD y FACScan son marcas comerciales de Becton, Dickinson and Company o sus filiales. Todas las demás marcas comerciales son propiedad de sus respectivos dueños. © 2024 BD. Todos los derechos reservados.