INSTRUCTION MANUAL

OPERATION & MAINTENANCE
TECHNICAL MANUAL

HIBLOW SERIES COMPRESSORS

MODELS

HP-40, HP-60, HP-80,
HP-100, HP-100LL, HP-120, HP-120LL,
HP-150, HP-200

This is a hazard alert symbol: ⚠ When you see this symbol, be aware that personal injury or property damage is possible. The hazard is explained in the text following the symbol. Read the information carefully before proceeding.

The following is an explanation of the three different types of hazards:

⚠ **DANGER:** Severe personal injury or death will occur if hazard is ignored.
⚠ **WARNING:** Severe personal injury or death can occur if hazard is ignored.
⚠ **CAUTION:** Minor injury or property damage can occur if hazard is ignored.

GENERAL INFORMATION

⚠ **DANGER:** Do not pump flammable or explosive gases or operate the unit in an atmosphere containing them.
⚠ **CAUTION:** The pump is designed for air only. Do not allow corrosive gases or particulate material to enter the pump. Water vapor, oil-based contaminants, or other liquids must be filtered out.
⚠ **CAUTION:** Ambient temperature should not exceed 40℃(104°F). For operation at higher temperatures, consult the factory.
⚠ **CAUTION:** Close supervision is necessary when any appliance is used by or near children.

This pump's function is to pump air and under no circumstances should it be used to pump any other gases. The pump must not be used for the pumping of fluids, particles, solids or any substance mixed with air, particularly combustible substances likely to cause explosions.

Never lubricate this oil-less pump, the moving parts consist of an actuating rod supported by two special synthetic rubber diaphragms which vibrate laterally, permitting long-term continuous operation. Operation is based upon the principle of electromagnetic vibration which eliminates the need for sliding parts, thereby minimizing power consumption and offering high efficiency.

INSTALLATION

⚠ **WARNING:** To avoid risk of electrocution do not use this product in an area where it could become submerged in water or other liquids. The Canadian (c-UL) pumps need to be stored in a dry location and it is not intended for outdoor use.
⚠ **CAUTION:** Do not block the flow of cooling air around the pump in any way. The life of pump may be reduced or malfunction could occur, if hazard is ignored.


GROUNDING INSTRUCTIONS

FOR ALL GROUNDED, CORD-CONNECTED PRODUCTS

This product should be grounded. In the event of an electrical short circuit, grounding reduces the risk of electric shock by providing an escape wire for the electric current. This is equipped with a cord having a grounding wire with an appropriate grounding plug.
The plug must be plugged into an outlet that is properly installed and grounded in accordance with all local codes and ordinances.

⚠ **DANGER:** Improper installation of the grounding plug can result in a risk of electric shock. If repair or replacement of the cord or plug is necessary, do not connect the grounding wire to either flat blade terminal. The wire with insulation having an outer surface that is green, with or without yellow stripes, is the grounding wire.

Check with a qualified electrician or serviceman if the grounding instructions are not completely understood, or if in doubt as to whether the product is properly grounded. Do not modify the plug provided; if it will not fit the outlet, have the proper outlet installed by a qualified electrician.

For a grounded, cord-connected product rated less than 15 amperes and intended for use on a nominal 120 volt supply circuit.


This product is for use on a nominal 120 volt circuit, and has grounding plug. Make sure that the product is connected to an outlet having the same configuration as the plug. No adapter should be used with the product.

FOR PERMANENTLY CONNECTED PRODUCT

This product should be connected to a grounded, metallic, permanent wiring system, or an equipment-grounding terminal or lead on the product(refer to wiring diagram).

WIRING

120V AC Wiring Diagram

EXTENSION CORDS

Use only a 3-wire extension cord that has a 3-blade grounding plug, and a 3-slot receptacle that will accept the plug on the product. Make sure your extension cord is in good condition. When using an extension cord, be sure to use one heavy enough to carry the current that your product will draw. A cord no longer than 100 feet, No.16 AWG extension cord should be used. An undersized cord will cause a drop in line voltage resulting in loss of power and overheating. For outdoor use, use an extension cord that is approved for outdoor use, marked W-or W-A.

MOUNTING

⚠ **CAUTION:** When using the pump to inject air into a liquid, make sure that the pump is higher than the surface level of the liquid, otherwise that liquid may run back into the pump if the power is interrupted.


The pump should be in an upright, level position and remain on a stable, rigid operating surface for proper operation.

PLUMBING

Improper connection will cause leakage of air.

HP Series : Connect the pumps exhaust port and PVC(plastic)pipe with an L-shaped rubber hose. Be sure to secure place of connection with a hose tightening clip(Refer to Fig.1).

Fig. 1

STORAGE

Store indoors out of the weather in a dry area. Make certain the electrical cord is placed so there will be no chance of damage.

⚠ **CAUTION:** To avoid risk of electrocution or damage to the pump do not store this unit where it could come into contact with water or other liquids. Store indoors. Inspect compressor and cord for damage prior to putting compressor back into service.

・After the Diaphragm mounting block have been installed, replace with the new Head (Casing) block, reattaching the L-tube, and reinstall the 4 screws holding the Head (Casing) on. (Refer to Fig.4-③or⑤)


〔3〕 Replacement of Safety Screw (Refer to Fig.6)

SP Switch (For HP-60/HP-80/HP-100/HP-100LL/HP-120/HP-120LL/HP-150/HP-200 Series)

The body of the switch is set on the upper part of the Rod. If the Diaphragm is damaged or unit is dropped, the unusual movement of the Rod will cause the upper part of the Rod to touch the L-shaped lever, contacting the spring terminal, breaking the safety screw. This turns the power off, stopping the pump automatically. This is a built in safety device to prevent any further damage to the pump. To replace SP Screw:

REPLACEMENT SAFETY SCREW PART#
HP-60/HP-80/HP-100/HP-100LL/HP-120/HP-120LL/HP-150/HP-200 ····· PASPSW0200

Fig. 6

・First remove the pump plug from the electrical outlet.

⚠ **WARNING:** All power must be disconnected or deenergized before servicing unit. If hazard is ignored personal injury or property damage is possible.

・Dispose of broken screw. Be sure all debris is removed (may be necessary to turn pump upside down) from unit.
・Fasten screw with plastic nut. The screw is designed so that the nut will turn freely when it is properly fastened, stop tightening when this happens.

〔4〕 Reinstall the Upper housing (Refer to Fig.4)

・Reinstall the 4 screws holding the Frame cover on. (For HP-100/120, Refer to Fig.4-④)
・Reinstall the Sound absorber (except for HP-40, HP-100LL, HP-120LL).
・Reinstall the Upper housing and tighten corner screws. (Refer to Fig.4-①)

OPERATION

⚠ **WARNING:** Lower housing and pump surfaces can become very hot during operation. Do not touch these surfaces until unit has been shut off and allowed to cool.

⚠ **CAUTION:** Carefully examine the air pump after installation. It should not be plugged in if there is water on parts not intended to be wet.

⚠ **CAUTION:** Do not operate if the cord or plug is damaged, or if the pump is malfunctioning, dropped, or damaged in anyway.

⚠ **CAUTION:** Ensure pump is securely mounted prior to operation.

When alternating current is applied to electromagnet as in the figures, the actuating rod moves first in the arrow direction shown in Fig. 2A and then in the arrow direction shown in Fig. 2B, by the magnetic attraction and repellent forces exerted between the electromagnet and the permanent magnets attached to the rod.

The rod vibrates at the power supply(60Hz/50Hz)frequency and changes the volume of the space enclosed between the head(casing) and the diaphragm. Thus, air intake, compression, and exhaust, can be performed as the rod vibrates.

Fig. 2A



Fig. 2B

MAINTENANCE

⚠ **WARNING:** Always disconnect power before installing or servicing. Failure to do so could result in electrical shock, personal injury or death. The motor may be thermally protected and will automatically restart when it cools, if the thermal protection switch is tripped.

⚠ **CAUTION:** The filters must be periodically cleaned or replaced. A clogged filter can cause overheating or pump failure.

Any servicing of the unit other than that recommended in the instruction manual must be performed by an Authorized Service Facility. Inspection of filters is suggested every three months.
When inspecting filter, also clean inlet port area of filter and housing. Filters can be cleaned with soap and water. Be sure to dry all parts and filter well before reinstalling. If the pump makes an abnormal noise, or the amount of discharged air is greatly decreased, turn off the power immediately.


For repairs or replacement parts, please contact the factory, or send unit to an Authorized Service Facility.

Cleaning the filter (Refer to Fig.3)

REPLACEMENT FILTERPAD PART#	
HP-40	40PA000010
HP-60/HP-80	80PA000040
HP-100/HP-100LL/HP-120/HP-120LL/HP-150/HP-200	120PA20010

- Undo the truss screw securing the cover to the pump housing, and remove the Filter pad.
- Gently dust the Filter pad then hand-wash in mild soapy water and be sure to dry well before reinstalling
- Replace the Filter pad and cover. Secure with the truss screw.

Fig. 3


Changing the Pump (Chamber) Assembly-HP Series

REPLACEMENT PUMP ASSEMBLY PART#		CONTENTS
HP-40	40PC000030	Head (Casing) block
HP-60/HP-80	80PC000041	Diaphragm mounting block
HP-100/HP-120	120PC20011	Hex nut (U-lock nut) / washer
HP-100LL/HP-120LL	120PC20021	Safety screw / plastic nut (HP-60 and over)
HP-150/HP-200	200PC20011	

〔1〕 Remove the Head (Casing) block (Refer to Fig.4)

- Remove the four corner screws(Refer to Fig.4-①), and the Upper housing. (Refer to Fig.4-②) If stuck, insert screwdriver between exhaust outlet and Upper housing, and gently pry apart.
- Remove the Sound absorber. (Refer to Fig.4-②, except for HP-40, HP-100LL, HP-120LL)
- Remove the screws from all the corners of Frame cover for subsequent removal of the Frame cover. (For HP-100/120, Refer to Fig.4-④)
- Undo the hose band and remove the L-tube from the Head (Casing) block.
- Remove screws attaching the Head (Casing) to the pump mechanism, and remove Head (Casing) block. (Refer to Fig.4-③. For HP-100/120, Refer to Fig.4-⑤)


Fig. 4

〔2〕 **Replacing Diaphragm mounting block** (Refer to Fig.5)

- After removing Head (Casing) block, remove one of the hex nut (U-lock nut) and washer in the center of the diaphragm.
- Remove one of the Diaphragm mounting blocks from the Rod, and pull out the other Diaphragm mounting block from the pump body together with the Rod without loosening its hex nut (U-lock nut) and washer. (Refer to Fig.5-①)
- Then remove the other Diaphragm mounting block from the Rod. (Refer to Fig.5-②)
- Reattach a new Diaphragm mounting block on one side of the Rod with new hex nut (U-lock nut) and washer, insert the Rod in accordance with the procedure reverse to that performed when the Rod was removed the pump.
- Reattach washer and new hex nut (U-lock nut) to Rod screw protruding through other Diaphragm mounting block and tighten.

Fig. 5

Operation of Compressor Alarm Module


General Information


⚠ **Danger:** Do not attempt to open the alarm module. If the unit is connected to electricity, opening of the alarm module can result in a risk of electric shock.

⚠ **Danger:** Service of the alarm module should only be done by a qualified electrician or serviceman.

⚠ **Danger:** If the outside of the alarm module or light lens becomes cracked or damaged, unplug or shut-off the electric power immediately and contact an electrician or serviceman.

⚠ **Caution:** Do not carry the unit by the alarm module or electrical cord. It could damage the alarm module.

The three position toggle switch is clearly labeled on the alarm module.

RUN : During normal operation the switch is in the "RUN" or center position.


MUTE : The "MUTE" position will silence the audible alarm while leaving the visual alarm operational.

TEST : The "TEST" position will activate both the audible and visual alarms to check to see if they are working.

In Case of an Alarm :

The alarm module audible alarm (buzzer) and visual alarm (light) are intended to signal a system malfunction. Contact your service provider whenever you hear or see the alarm when the switch is in the "RUN" or center position. The "MUTE" position can be used to silence the audible alarm until the service provider arrives.

TROUBLE SHOOTING

YOUR WARRANTY

Hiblow products, when properly installed and operated under normal conditions of use, are warranted by Hiblow to be free from defects in material and workmanship for a period of twelve (12) months from the date of purchase from Hiblow. In order to obtain performance under this warranty, the buyer must give written notice within thirty (30) days after discovery of the defect to Hiblow USA, Inc., 1300 Tefft Court Suite 8, Saline Michigan USA 48176. Buyer is responsible for freight charges both to and from Hiblow in all cases.

Hiblow's warranties do not extend to any goods or parts which have been subjected to misuse, lack of maintenance, neglect, damage by accident or transit damage.

This warranty excludes all other warranties or representations expressed or implied by any literature, data or person. Hiblow's maximum liability under this exclusive remedy shall never exceed the cost of the subject product and Hiblow reserves the right, at its sole discretion, to refund the purchase price in lieu of repair or replacement.

Hiblow will not be responsible or liable for indirect or consequential damages of any kind, however arising, including but not limited to those for use of any products, loss of time, inconvenience, lost profit, labor charges, or other incidental or consequential damages with respect to persons, business or property, whether as a result of breach of warranty, negligence or otherwise.

Buyer's remedy against Hiblow for goods supplied or for non-delivered goods or failure to furnish goods, whether or not based on negligence, strict liability or breach of expressed or implied warranty is limited solely, at Hiblow's option, to replacement of or cure of such nonconforming goods or non delivered goods or return of the purchase price for such goods and in no event shall exceed the price or charge for such goods. Hiblow expressly disclaims any warranty or merchantability or fitness for a particular use or purpose with respect to goods sold. There are no warranties which extend beyond the descriptions set forth in this warranty, notwithstanding any knowledge of Hiblow regarding the use or uses intended to be made of goods, proposed changes or additions to goods, or any assistance or suggestions that may have been made by Hiblow personnel.

Unauthorized extensions of warranties by the customer shall remain the customer's responsibility.

Customer is responsible for determining the suitability of Hiblow products for customer's use or resale or for incorporating them into objects or applications which customer designs, assembles, constructs or manufactures.

This warranty can be modified only by authorized Hiblow personnel by signing a specific, written description of any modifications.

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FLA-EU Series

Features

- Built-in fire protection system.
- Scalable expansion up to 15 pcs in parallel.
- Built-in WiFi for remote monitor the battery pack data.
- LiFePO4: Higher safe performance and longer cycle life.
- Multiple Protection with built-in smart BMS, breaker and fuse.
- Compatible with leading inverter brands.
- Up to 10-year long warranty.

FLA48100-EU~FLA48250-EU

Guangzhou Felicity Solar Technology Co.,Ltd.

www.felicitysolar.com sales@felicitysolar.com

No. 2, Donghua Huaye Road, Renhe Town, Baiyun District, Guangzhou, China **(Only 10 Mins from Baiyun Airport)**

Official Website Facebook

Model	FLA48100-EU	FLA48171-EU	FLA48230-EU	FLA48250-EU
Battery Type	LiFePO4			
Nominal Energy	5.12kWh	8.75kWh	11.8kWh	12.5kWh
Nominal Voltage	51.2V	51.2V	51.2V	51.2V
Operating Voltage	44.8-57.6V	44.8-57.6V	44.8-57.6V	44.8-57.6V
Max. Continuous Charge/Discharge Current[1]	100A	120A	150A	150A
Peak Charge/Discharge Current(15s)	150A	150A	200A	200A
Max. Charge/Discharge Power(15s)	7,500W	6,000W	7,500W	10,000W
Scalability	Max. 15 pcs in Parallel(76.8kWh)	Max. 15 pcs in Parallel(131.5kWh)	Max. 15 pcs in Parallel(177kWh)	Max. 15 pcs in Parallel(187.5kWh)
Depth of Discharge(DOD)	≥95%			
Communication	RS485/CAN			
Protection Degree	IP21			
Cycle Life	≥ 6000 Cycles			
Charging Temperature Range	0°C~+55 °C			
Discharging Temperature Range	-20°C~+55 °C			
Display Type	LCD+LED			
Installation	Wall-Mounted / Floor-Mounted			
Protection	Built-in smart BMS, Breaker, Fuse			
Warranty Period	10 Years			
Certification	\	\	UN38.3	UN38.3、CE-EMC
Product Weight	55KG	68KG	97KG	97KG
Package Weight	77KG	88KG	115KG	115KG
Product Dimension	775X425X209MM	735X433X188MM	783X450X274MM	783X450X274MM
Package Dimension	892X545X385MM	849X550X362MM	900X570X450MM	900X570X450MM

[1] Max. Continuous Charge/Discharge Current is Affected by Temperature and SOC.

*Product specifications may vary during continuous improvements. Please contact us for latest details.
*Refer to Felicitysolar's Warranty Policy for applicable conditions.

Product Display

Longer Lasting LiFePO4 Battery

Built-in WiFi/Bluetooth

Max. 15 Units in Parallel

Aerosol Fire Extinguishing System

Easy Fuse Replacement Design

Real-time Status via LCD Display

Guangzhou Felicity Solar Technology Co.,Ltd.
www.felicitysolar.com sales@felicitysolar.com
No. 2, Donghua Huaye Road, Renhe Town, Baiyun District, Guangzhou, China **(Only 10 Mins from Baiyun Airport)**

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
Instruction Manual

ZEISS Primostar 3

Upright Microscope for Education and Routine

ZEISS Primostar 3

Original Manual

EC REP

Carl Zeiss Microscopy GmbH
Carl-Zeiss-Promenade 10
07745 Jena
Germany
info.microscopy.de@zeiss.com
www.zeiss.com/microscopy

Carl Zeiss Suzhou Co., Ltd.
Modern Industrial Square 3-B, No.333 XingPu Road SIP
215126 Suzhou
China

Document Name: Instruction Manual ZEISS Primostar 3
Order Number: 415501-7011-111
Revision: 3
Language: en-US
Effective from: 05/2022

IVD CE

© 2022 Without the prior written consent of ZEISS, this document or any part of it must neither be translated nor reproduced or transmitted in any form or by any means - including electronic or mechanic methods, by photocopying, recording or by any information or filing system. The right to make backup-copies for archiving purposes shall remain unaffected thereby. Any violations may be prosecuted as copyright infringements.

The use of general descriptive names, registered names, trademarks, etc. in this document does not imply that such names are exempt from the relevant intellectual property laws and regulations and therefore free for general use. This shall also apply if this is not specifically referred to. Software programs shall entirely remain the property of ZEISS. No program or subsequent upgrade thereof may be disclosed to any third party, copied or reproduced in any other form without the prior written consent of ZEISS, even if these copies or reproductions are destined for internal use at the customer's only, the only exception being one single back-up copy for archiving purposes.

Table of Contents

1 About this Instruction Manual 6

1.1 Text Conventions and Link Types 6
1.2 Explanation of Warning Messages and Additional Information 6
1.3 Further Applicable Documents 7
1.4 Contact 7

2 Safety 9

2.1 Intended Use 9
2.1.1 Intended Purpose 9
2.1.2 Lifetime 9
2.1.3 EMC Information 9
2.2 General Safety Information 10
2.2.1 Requirements for Operators 10
2.2.2 Safe Operating Condition 10
2.2.3 Order and Use of Spare Parts 10
2.3 Prevention of Hazards 11
2.3.1 Mechanical Hazards 11
2.3.2 Electrical Hazards 11
2.3.3 Hazards Generated with the Operating Environment 11
2.3.4 Ergonomic Hazards 11
2.3.5 Hazards Generated by Materials and Substances 12
2.3.6 Hazards Generated by Radiation 12
2.4 Labels and Lights 12
2.4.1 Labels on the Primostar 3 13
2.5 Safety Devices and Interlocks 15
2.5.1 Safety Device of the Reflected-Light FL iLED 15

3 Product and Functional Description 16

3.1 Main Components of the Primostar 3 16
3.1.1 Components of Primostar 3 Fixed-Köhler 17
3.1.2 Components of Primostar 3 iLED (Fixed-Köhler) 18
3.1.3 Components of Primostar 3 (Fixed-Köhler) with Binocular Tube 19
3.1.4 Components of Primostar 3 Full-Köhler 20
3.2 Controls and Connections 21
3.2.1 Fixed-Köhler Stand 21
3.2.2 Full-Köhler Stand 24
3.2.3 Reflected-Light FL iLED 27
3.2.4 Binocular Tube 25°/22 w/int Cam 8 MPx 28
3.2.5 Abbe Condenser 0.9/1.25 Field 20 37
3.2.6 Eyepiece 38
3.2.7 Cable and Screw Tool Storage 38
3.2.8 Accessories for Connecting the Microscope 39
3.3 Objective Labeling 40

4 Installation 42

4.1 Safety During Installation 42
4.2 Preparing Installation 42
4.3 Installing Additional Components 43
4.3.1 Installing or Changing the Tube 43

4.3.2 Installing the Reflected-Light FL iLED onto the Stand 44
4.3.3 Assembling the Yellow Filter or the TL Cover Plate (Fixed Köhler Stand) 47
4.3.4 Inserting Color Filters (Full Köhler Stand) 47
4.3.5 Installing the Polarizer Plate (Full Köhler Stand) 48
4.3.6 Installing the Analyzer Plate 49
4.3.7 Installing the Binocular Tube 25°/22 w/int Cam 8 MPx 49
4.3.8 Installing a Camera to the Trinocular Tube 51
4.3.9 Exchanging Objectives 53
4.3.10 Installing/Removing the Condenser 53
4.3.11 Installing/Removing the Mirror 54
4.3.12 Folding Over the Eyecups 55
4.3.13 Changing the Eyecups 56
4.4 Connecting the Microscope to the Power Supply 56
4.4.1 Connecting the Microscope to the Mains Power Supply 56
4.4.2 Connecting the Microscope to a Power Bank 56
4.5 Connecting the Binocular Tube 25°/22 w/int Cam 8 MPx to External Devices 57
4.5.1 Connecting a USB Flash Drive 57
4.5.2 Connecting to the PC via USB 3.0 Port 58
4.5.3 Connecting to a Display (without PC) 58
4.5.4 Integrating the Binocular Tube 25°/22 w/int Cam 8 MPx into a Network 59
4.6 Installing Labscope Software for Using the Binocular Tube 25°/22 w/int Cam 8 MPx 61
4.6.1 Installing the Labscope Software on PC 61
4.6.2 Installing the Labscope App on iPad 61

5 Operation 62

5.1 Prerequisites for Commissioning and Operation 62
5.2 Switching On the Microscope 62
5.2.1 Switching On the Microscope with Mains Power Connection 62
5.2.2 Switching On the Microscope with Power Bank Connection 62
5.3 Adjusting the Tube 63
5.3.1 Adjusting the Position of the Eyepieces 63
5.3.2 Adjusting the Eyepiece with Eyepiece Pointer or Eyepiece Micrometer .. 64
5.4 Adjusting the Transmitted-Light Brightfield on the Full-Köhler Microscope 65
5.5 Adjusting the Transmitted-Light Brightfield on the Fixed-Köhler Microscope 68
5.6 Adjusting the Transmitted-Light Phase Contrast or the Transmitted-Light Darkfield Using a Contrast Slider 70
5.7 Adjusting the Transmitted-Light Phase Contrast or the Transmitted-Light Darkfield Using the Turret Condenser 73
5.8 Adjusting the Transmitted-Light Simple Polarization Contrast 75
5.9 Adjusting the Reflected-Light Fluorescence 75
5.10 Capturing Images and Videos Using the Integrated Smart 8 MP Color Camera 77
5.10.1 Capturing Images and Videos Using the Operator Panel 77
5.10.2 Capturing Images and Videos Using the OSD Menu 78
5.10.3 Capturing Images with Labscope 78
5.11 Switching Off the Microscope 79

6 Care and Maintenance 80

6.1 Safety During Cleaning and Maintenance 80
6.2 Maintenance Schedule 80
6.3 Maintenance Work 81
6.3.1 Cleaning an Optical Surface 81

6.3.2 Removing Water-soluble Contamination 81
6.3.3 Updating the Firmware of the Binocular Tube 25°/22 w/int Cam 8 MPx 81
6.3.4 Exchanging the Halogen Lamp of the Halogen Module 82

7 Troubleshooting 83

7.1 Troubleshooting Binocular Tube 25°/22 w/int Cam 8 MPx..... 85

8 Decommissioning and Disposal..... 87

8.1 Decommissioning..... 87
8.2 Transport and Storage 87
8.3 Disposal 87
8.4 Decontamination 88

9 Technical Data and Conformity..... 89

9.1 Performance Data and Specifications 89
9.2 Performance Data and Specifications of the Optional Components 90
9.3 Performance Data and Specifications of the Binocular Tube 25°/22 w/int Cam 8 MPx 92
9.4 Applicable Standards and Regulations..... 93

10 Accessories and System Expansions..... 94

Revision History 96

Glossary 97

Index..... 98

1 About this Instruction Manual

This Instruction Manual (further called "document") is considered to be part of the Primostar 3, herein after referred to as "microscope".

Primostar 3 microscopes include:

- Primostar 3
- Primostar 3 iLED

This document contains basic steps and safety information that must be observed during operation and maintenance. Therefore, the document must be read by the operator prior to commissioning and must always be available at the place of use of the microscope.

This document is an essential part of the microscope and, if the microscope is resold, the document must remain with the microscope or be handed over to the new owner.

1.1 Text Conventions and Link Types

Example	Explanation
	The names of controls and important information are shown in bold letters, for instance:
Click **Start**.	Software controls and GUI elements.
Press the **Standby** button.	Hardware controls and elements.
Press **Enter** on the keyboard.	Key on the keyboard.
Press **Ctrl + Alt + Del**.	Press several keys on the keyboard simultaneously.
Select **Tools > Goto Control Panel > Airlock**.	Follow a path in the software.
Enter *example.pdf* in this field.	Text to be entered by the user.
`Programming and Macros`	Anything typed in literally during programming, including, for example, macro codes, keywords, data types, method names, variables, class names, and interface names.

Tab. 1: Text conventions

Example	Explanation
See: *Text Conventions and Link Types [▸ 6]*.	Link to further information for this topic.
https://www.zeiss.com/corporate/int/home.html	Link to a website on the internet.

Tab. 2: Link types

1.2 Explanation of Warning Messages and Additional Information

DANGER, WARNING, CAUTION, and NOTICE are standard signal words used to determine the levels of hazards and risks of personal injury and property damage. Not only the safety and warning messages in the **Safety** chapter are to be considered also all safety and warning messages in other chapters. Failure to comply with these instructions and warnings can result in both personal injury and property damage and involve the loss of any claims for damages.

The following warning messages indicating dangerous situations and hazards are used in this document.

⚠ DANGER

Type and source of danger

DANGER indicates an imminently hazardous situation which, if not avoided, will result in death or serious injury.

⚠ WARNING

Type and source of danger

WARNING indicates a potentially hazardous situation which, if not avoided, may result in death or serious injury.

⚠ CAUTION

Type and source of danger

CAUTION indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury.

NOTICE

Type and source of danger

NOTICE indicates a potentially harmful situation which, if not avoided, may result in property damage.

Info

Provides additional information or explanations to help operator better understand the contents of this document.

1.3 Further Applicable Documents

Also take note of the following documents:

Brochures and Certificates For brochures, ISO certificates, CSA certificates, and EU declarations of conformity ask your ZEISS Sales & Service Partner.

System and third-party Components, Accessories Information about the individual components, enhancements, and accessories can be obtained from your ZEISS Sales & Service Partner. Also refer to the documentation of third-party manufacturers.

1.4 Contact

If you have any questions or problems, contact your local ZEISS Sales & Service Partner or one of the following addresses:

Headquarters

Phone:	+49 1803 33 63 34
Fax:	+49 3641 64 3439
Email:	info.microscopy.de@zeiss.com

Microscopy Courses, Training, and Education

For information on microscopy courses, training, and education contact us on our homepage (https://www.zeiss.com/microscopy/int/service-support/training-and-education.html#contact).

ZEISS Portal

The ZEISS Portal (https://portal.zeiss.com/) offers various services that simplify the daily work with your ZEISS systems (machines and software). It is constantly improved and extended to meet your needs and requirements better.

ZEISS Sales & Service Partner

You can find a ZEISS Sales & Service Partner in your area under https://www.zeiss.de/mikroskopie/website/forms/sales-and-service-contacts.html.

Service Germany

Phone:	+49 7364 20 3800
Fax:	+49 7364 20 3226
Email:	service.microscopy.de@zeiss.com

2 Safety

This chapter contains general requirements for safe working practices. Any person using the microscope or commissioned with installation or maintenance must read and observe these general safety instructions. Knowledge of basic safety instructions and requirements is a precondition for safe and fault-free operation. Operational safety of the supplied microscope is only ensured if it is operated according to its intended use.

If any work is associated with residual risks, this is mentioned in the relevant parts of this document in a specific note. When components must be handled with special caution, they are marked with a warning label. These warnings must always be observed.

Any serious incident that has occurred in relation to the microscope and its components shall be reported to these institutions:

- the competent authority of the Member State in which the user is established
- ZEISS
 - for users within the EU:
 Carl Zeiss Microscopy GmbH, Jena, Germany
 - for users outside the EU:
 Carl Zeiss Suzhou Co., Ltd., Suzhou, China

2.1 Intended Use

Improper use of the microscope and its components can easily lead to impairment of their function or even damage them. Damage caused by incorrect operation, negligence, or unauthorized intervention, in particular by removing, modifying, or replacing parts of the microscope or its components, cannot be held liable by the device manufacturer. Third-party devices or components that are not expressly approved by ZEISS may not be used.

2.1.1 Intended Purpose

The microscopes Primostar 3, Primostar 3 iLED are instruments for general microscopic imaging for the in vitro-examination of various biological samples including samples collected from humans or animals. This imaging provides information to further assess physiological and pathological conditions.

The microscopes are intended to be used by trained professionals only.

2.1.2 Lifetime

A microscope is an opto-electronic device. Its availability for use is significantly determined by the performed maintenance. ZEISS guarantees the ability for maintenance and repair within eight years after initial operation. This is ensured by a corresponding service and spare parts concept, thus enabling the intended purpose within this duration.

2.1.3 EMC Information

The microscope is intended to be used in a basic electromagnetic environment. The electromagnetic environment should be evaluated prior to operation of the microscope. Do not use the microscope in proximity to sources of strong electromagnetic radiation, as these can interfere with proper operation.

Use of this microscope in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets, etc.), may cause electrostatic discharges that may cause erroneous results.

Electromagnetic interference (EMI) according to CISPR 11 Group 1:

- Class B (without Binocular tube 25°/22 w/int cam 8 MPx)
- Class A (with Binocular tube 25°/22 w/int cam 8 MPx)

If in doubt, contact a ZEISS service representative.

2.2 General Safety Information

This document must be read before commissioning in order to ensure safe and uninterrupted operation. Pay particular attention to all listed safety notes. Make sure, that

- the operating personnel has read and understood this manual, associated documents and particularly all safety regulations and instructions, and applies them.
- the local and national safety and accident prevention regulations must be observed, as well as the applicable laws and regulations in your country.
- this document is always available at the place of use of the microscope.
- the microscope is always in perfect condition.
- the microscope is secured against access by unauthorized persons.
- maintenance and repair work, retrofitting, removal or replacement of components, as well as any other intervention in the microscope not described in this document, may only be carried out by the manufacturer ZEISS or persons expressly authorized by ZEISS to do so.

2.2.1 Requirements for Operators

The microscope, components, and accessories may only be operated and maintained by authorized and trained personnel. The microscope may only be used in accordance with this document. If the microscope is not used as described, the safety of the user may be impaired and/or the microscope may be damaged.

Any unauthorized intervention or use other than within the scope of the intended use shall void all rights to warranty claims. The regional regulations on health protection and accident prevention must be observed at all times and during all work on and with the microscope.

Training Authorized ZEISS personnel will provide basic training in operating the microscope. As well as information on equipment safety and maintenance work that can be conducted by the operator. The training will be documented by ZEISS and its completion is to be confirmed by the operator.

Special application training is offered for a fee. Current training dates, additional information and the registration form can be found at https://www.zeiss.com/microscopy/int/service-support/training-and-education.html.

2.2.2 Safe Operating Condition

If circumstances occur which impair safety and cause changes in operating behavior, the microscope must be shut down immediately and a ZEISS service representative should be informed.

The microscope may only be operated after correct installation by a ZEISS service representative and if the operating conditions are adhered to.

- Do not operate the microscope until you have completely read and understood the entire documentation.
- Make sure that all protective cover panels are installed and all warning labels are available and legible.
- Ensure conditions and take measures to prevent the build up of electrostatic charge on the workplace.

2.2.3 Order and Use of Spare Parts

Using spare parts that are not provided by ZEISS can be hazardous or can lead to property damage.

- Unless authorized by ZEISS, all spare parts should be installed by a ZEISS service representative.
- Contact your ZEISS service representative for information on spare parts order.
- Only genuine parts supplied by ZEISS are to be used in servicing the microscope.

2.3 Prevention of Hazards

This section summarizes potential hazards and recommended safety precautions. Failure to follow the safety instructions and instructions may result in personal injury and property damage.

2.3.1 Mechanical Hazards

Property Damage due to Transport

There is a risk of injury and property damage if the microscope is improperly handled and transported.

- Only use the handle, if applicable, for transport of the microscope. Otherwise hold the microscope with one hand and the base plate with the other hand.

2.3.2 Electrical Hazards

Voltage Hazards

Risk of electric shock in case of contact with live parts.

Always use the power adapter supplied by ZEISS. When an unsuitable power adapter is used, ZEISS can no longer guarantee the electrical safety and functionality of the microscope.

- Shut down the microscope.
- Disconnect the power supply before cleaning.
- Set up and operate the microscope so that the connectors are easily accessible.
- Position the microscope stand in a way so that you can easily unplug the power cable at any time.

Safe disconnection from the mains is ensured exclusively by removing the mains plug. The switch on the rear side of the microscope only switches into standby mode.

2.3.3 Hazards Generated with the Operating Environment

Dirt, Dust, and Moisture

Dirt, dust, and moisture can impair the microscope's functionality.

- Shut down the microscope whenever it is not used and cover it with a dust protection cover.
- Always cover unused openings/ports.
- Perform regular maintenance and cleaning according to the instructions in this manual.
- Make sure that no cleaning liquid or moisture gets inside the microscope.
- Make sure that the electrical parts never come into contact with moisture.
- Never expose the microscope to inadmissible climate conditions (high humidity and temperature).

2.3.4 Ergonomic Hazards

Prevention of Musculoskeletal Disorders

Musculoskeletal disorders (MSDs) affect the muscles, nerves, blood vessels, ligaments and tendons. Workers in many different industries and occupations can be exposed to risk factors at work, such as lifting heavy items, bending, reaching overhead, pushing and pulling heavy loads, working in awkward body postures and performing the same or similar tasks repetitively. Employers are responsible for providing a safe and healthful workplace for their workers.

2.3.5 Hazards Generated by Materials and Substances

Infection Hazards Direct contact with the eyepieces can be a potential way of passing on bacterial and viral infections.

- The risk can be lowered by using personal eyepieces or eyecups. If eyepieces need to be disinfected frequently, ZEISS recommends to use the eyepieces without eyecups.
- To avoid infections, the use of personal protective equipment (PPE), e.g. gloves, for operation, cleaning, and decontamination is highly recommended. Disposable gloves can be decontaminated with alcohol for example, if necessary, or should be changed frequently to minimize the risk of contamination.

Consumable Hazards Incorrect handling of consumables and cleaning agents can lead to property damage or skin and eye injuries. Consumables that are not approved by ZEISS can lead to property damage. Consult your ZEISS Sales & Service Partner to learn what consumables you can order and how to handle them.

Disinfectant Hazards Ensure adequate ventilation in closed rooms. In case of insufficient ventilation, wear respiratory protective equipment. Remove any harmful residue. Allow the device to dry off after disinfection, particularly after disinfection of eyepieces. Do not inhale vapors. When using disinfectants, do not eat, drink or smoke. Avoid contact with eyes and skin. Remove contaminated clothing and wash before reuse.

Eye, Skin, Respiratory Tract Irritation Exposure to chemicals and their aerosols can cause eye, skin and respiratory tract irritation. Use appropriate personal protective equipment (PPE).

2.3.6 Hazards Generated by Radiation

Optical Radiation Hazards Gas discharge lights, LED lights and other sources of white light emit strong optical radiation (e.g. UV, VIS, IR). Optical radiation may cause damage to the skin and eyes. The extent of the damage depends on the parameters such as wavelength, exposure time, mode of operation (continuous or pulsed), etc.

- Avoid exposure of eyes and skin to radiation.
- Do not introduce reflective objects into the beam path.
- Never remove covers or cover panels during operation.
- Do not disable any interlock system elements.
- Use suitable protective equipment / protective clothing if required.

Electromagnetic Radiation Hazards In a domestic environment the microscope may cause radio interference, which may be mitigated by relocating or re-orienting the equipment. The use of non-specified accessories, cables, or other auxiliary parts from the field of information technology may lead to increased electromagnetic emissions and reduced immunity to interference. Any integration into the system may result in a degradation of the EMC performance.

2.4 Labels and Lights

This chapter shows labels and, where applicable, indicator lights.

All parts that may pose specific hazards are marked with warning labels.

Always observe **all** warning labels!

- Check all warning labels for availability and legibility.
- Immediately replace damaged or illegible warning labels.

In case a label is missing please contact your ZEISS service representative for free of charge replacement.

2.4.1 Labels on the Primostar 3


Fig. 1: Labels on the microscope rear side


Fig. 2: Label behind the nosepiece


Fig. 3: Labels on the fluorescence unit




Fig. 4: Label on the halogen module of the Full-Köhler stand

Pos.	Label or light	Explanation
1		CAUTION Hot surface below Operate microscope only on a stable, solid, smooth and not tinderlike (non-flammable) surface.
2		Microscope type label
3		Microscope type label
4		Primostar 3 EU representative label Carl Zeiss Microscopy GmbH Carl-Zeiss-Promenade 10 07745 Jena, Germany
5		CAUTION LED Radiation Do not stare at operating lamp. May be harmful to the eyes.
6		Fluorescence unit type label

Pos.	Label or light	Explanation
7		Hot surface
*		EU representative and importer label * Label is placed on the microscope packaging.

2.5 Safety Devices and Interlocks

In order to prevent injuries and/or property damage, the microscope is equipped with several safety devices and interlocks. In case of defect or damage, the affected parts and the microscope must be taken out of operation immediately and must be secured against unintentional use.

To verify the safety of the microscope, contact your ZEISS service representative and keep the service logs and logbooks.

2.5.1 Safety Device of the Reflected-Light FL iLED

To prevent personal injury and/or property damage, the reflected-light FL iLED is equipped with an interlock function, which causes the built-in LED of the reflected-light FL iLED to be switched off as soon as the reflected-light FL iLED is turned with regard to the stand or detached.

Defective and damaged safety devices can lead to injuries and damage. In the event of damage or defect, the affected parts and the entire microscope must be taken out of operation immediately and secured against unintentional use.

To check the microscope safety, contact your ZEISS Sales & Service Partner and keep the service protocols for your microscope.

3 Product and Functional Description

The Primostar 3 is a transmitted-light microscope of compact design with a small footprint. The microscope provides high resolution, infinity-corrected objectives for all important transmitted-light techniques, such as brightfield, darkfield, simple polarization, and phase contrast microscopy.

The microscope is optionally available with a camera port for photo and video documentation. For special camera applications the optional binocular tube 25°/22 is equipped with an integrated smart 8 MP color camera.

Typical Applications

- examination of blood and tissue samples taken from the human body, from plants, or animals
- medical examinations in laboratories, hospitals, and doctors' offices
- academic and practical education in medicine and biology
- industrial applications, e.g. in pharmacology, food technology, and wastewater examination

Info

Additional information about the hardware configuration and optional enhancements can be obtained from your ZEISS Sales & Service Partner.

3.1 Main Components of the Primostar 3

The Primostar 3 has the following main components:


Fig. 5: Main components of the Fixed-Köhler microscope

1	*Eyepiece [▸ 38]*	2	Objective
3	Stage	4	*Condenser [▸ 37]*
5	Illumination	6	Tube
7	Stand	8	Focusing drive

3.1.1 Components of Primostar 3 Fixed-Köhler

The Fixed-Köhler microscope has the following components:


Fig. 6: Components of the Fixed-Köhler microscope

1	Nosepiece with four positions	2	Rackless stage 75x40, drive right or drive left
3	Condenser carrier	4	Luminous-field diaphragm
5	Trinocular tube 25°/20 (50:50)	6	Small hand lift
7	Specimen holder left for rackless stage	8	Stage carrier

3.1.2 Components of Primostar 3 iLED (Fixed-Köhler)

The Fixed-Köhler microscope with reflected light FL iLED has the following main components:


Fig. 7: Components of the Fixed-Köhler microscope with iLED

1	Special eyecups with light protection mounted	**2**	Nosepiece with four positions
3	Rackless stage 75x40, drive right or drive left	**4**	Condenser carrier
5	Yellow filter inserted (with filter position for adapting the color temperature in transmitted light and with position for blocking the transmitted-light path in case of reflected-light fluorescence applications)	**6**	Luminous-field diaphragm
7	Binocular tube 25°/20	**8**	Reflected-light FL iLED
9	Specimen holder left for rackless stage	**10**	Stage carrier

3.1.3 Components of Primostar 3 (Fixed-Köhler) with Binocular Tube

The Fixed-Köhler microscope with Binocular tube 25°/22 w/int cam 8 MPx has the following components:


Fig. 8: Main components of Fixed-Köhler microscope with Binocular tube 25°/22 w/int cam 8 MPx

1 Nosepiece with four positions	**2** Rackless stage 75x40, drive right or drive left
3 Condenser carrier	**4** Luminous-field diaphragm
5 *Binocular Tube 25°/22 w/int Cam 8 MPx [▸ 28]*	**6** Specimen holder left for rackless stage
7 Stage carrier	

3.1.4 Components of Primostar 3 Full-Köhler

The Full-Köhler microscope has the following components:


Fig. 9: Components of the Full-Köhler microscope

1	Nosepiece with five positions	2	Rackless stage 75x50, drive right or drive left
3	Turret condenser 0.9/1.25 field 22 or *Abbe Condenser 0.9/1.25 Field 20* *[▸ 37]*	4	Condenser carrier with centering screws for condenser adjustment
5	Luminous-field diaphragm	6	Trinocular tube 25°/22 (50:50)
7	Specimen holder for 2 sliders	8	Stage carrier
9	Transmitted-light illuminator, LED or HAL		

3.2 Controls and Connections

3.2.1 Fixed-Köhler Stand

Purpose The controls on the stand control the main functions of the microscope.

Position The following controls are located on the right side of the stand.


Fig. 10: Stand controls on the right side

1 Clamping screw for tube

2 Knurled ring for turning the nosepiece

3 Lever for adjusting the aperture diaphragm of the condenser, see also *Abbe Condenser 0.9/1.25 Field 20* *[▸ 37]*

4 Centering screw (Allen screw) for condenser on condenser carrier (left side)

5 Centering screw (Allen screw) for condenser on condenser carrier (right side)

6 Cover cap of the luminous-field diaphragm

7 Clamping screw for condenser

8 Small hand lift (only Fixed-Köhler stand)

9 Illumination-intensity indicators for transmitted light (right side)
The selected intensity is indicated in five layers by the Cyan light-emitting diodes belt.

10 Vernier and scale, displaying the Y position of the stage

11 Control knob for Y travel of rackless stage

12 Rotary knob to adjust the illumination intensity for transmitted light

13 Control knob for X travel of rackless stage

14 Coarse focusing drive (right side)

15 Fine focusing drive (right side)

Position The following controls are located on the left side of the stand.


Fig. 11: Stand controls on the left side

1	Illumination-intensity indicators for transmitted light (left side) The selected intensity is indicated in five layers by the Cyan light-emitting diodes belt.	**2**	Clamping screw for the specimen holder
3	Vernier and scale, displaying the X position of the stage	**4**	Lever of the specimen holder for fixing the specimen
5	Knurled ring for condenser height adjustment	**6**	Knurled ring for adjusting the smoothness of the coarse focusing drive
7	Coarse focusing drive (left side)	**8**	Fine focusing drive (left side)

Position The following controls and connections are located on the rear side of the stand.


Fig. 12: Stand controls and connections on the rear side

1	Carrying handle	2	Holder for power cable
3	**Main power ON/OFF** button	4	5V USB port for mobile power bank
5	Storage place for phase plate adjustment tool	6	**DC 5V** switch
7	Mains socket	8	Connection port for external DC power supply (12 V)

Info

When the **DC 5V** switch 6 is set to the **OUT** position, the maximum output of the USB port 4 is 1 A.

3.2.2 Full-Köhler Stand

Purpose The controls on the stand control the main functions of the microscope.

Position The following controls are located on the right side of the stand.


Fig. 13: Stand controls on the right side

1 Clamping screw for tube

2 Knurled ring for turning the nosepiece

3 Turret condenser, controls see *Abbe Condenser 0.9/1.25 Field 20 [▸ 37]*

4 Centering screw (knurled screw) for condenser on condenser carrier (right side)

5 Cover cap of the luminous-field diaphragm

6 Knurled ring for adjusting the luminous-field diaphragm

7 Clamping screw for condenser

8 Illumination-intensity indicators for transmitted light (right side)
The selected intensity is indicated in five layers by the Cyan light-emitting diodes belt.

9 Rotary knob for varying the brightness of the microscopic image, for switching on / off the Light manager and the ECO mode

10 Control knob for Y travel of rackless stage

11 Control knob for X travel of rackless stage

12 Fine focusing drive (right side)

13 Coarse focusing drive (right side)

Position The following controls are located on the left side of the stand.


Fig. 14: Stand controls on the left side

1 Illumination-intensity indicators for transmitted light (left side)
The selected intensity is indicated in five layers by the Cyan light-emitting diodes belt.

2 Clamping screw for the specimen holder

3 Vernier and scale, displaying the X position of the stage

Knurled ring for condenser height adjustment

4 Vernier and scale, displaying the Y position of the stage

5 Lever of the specimen holder for fixing the specimen

6 Knurled ring for condenser height adjustment

7 Knurled ring for adjusting the smoothness of the coarse focusing drive

8 Coarse focusing drive (left side)

9 Fine focusing drive (left side)

10 Centering screw (knurled screw) for condenser on condenser carrier (left side)

Position The following controls and connectors are located on the rear side of the stand.


Fig. 15: Stand controls and connections on the rear side

1	Carrying handle	**2**	Holder for power cable
3	**Main power ON/OFF** button	**4**	Storage place for phase plate adjustment tool
5	Mains socket		

Function For Full-Köhler stands, Light Manager, encoding function and sleeping mode are factory set up:

- The Light Manager enables the microscope to remember the best illumination setting from the previous use. When the microscope is set at the best illumination with objective A, light intensity will be remembered. You can switch to objective B and adjust the illumination. When you change back to objective A, light intensity will return to the last setting. The light intensity of all five objectives can be remembered.
- The Light Manager feature can be disabled by pushing the intensity knob for five seconds. You will see that the Cyan LED Belt reduces its intensity from Max to Min, indicating that the Light Manager is disabled. You can re-activate the Light Manager by the same process.
- If you do not touch the microscope for 30 minutes it will go into ECO mode. The lowest layer of the Cyan LED belt will "breathe" (blink slowly). Any movement of the nosepiece or the intensity knob will re-activate the microscope.
- The ECO mode can be disabled by triple clicking the intensity knob (for example in case of long time video shooting). You will see the Cyan LED belt illuminating once as a sign. Triple clicking the intensity knob again will reactivate ECO mode.

3.2.3 Reflected-Light FL iLED

Purpose To carry out examinations according to the reflected-light fluorescence method, the reflected-light FL iLED is required.

Position The reflected-light FL iLED is mounted on the Fixed-Köhler stand (without small handle lift).


Fig. 16: Controls of the reflected-light FL iLED

1 Pilot lamp for reflected-light FL iLED: lights up blue when switched on; brightness corresponds to intensity

2 Clamping screw for tube

3 Rotary knob for switching ON/OFF and adjusting the intensity for reflected light

4 **fluorescence/brightfield** changeover switch
Always turn the **fluorescence/bright-field** changeover switch first upward and then to the desired position. Using force to turn it downward will damage the reflected-light FL iLED.


Fig. 17: Connections and special components on the rear side of the reflected-light FL iLED

1 RJ12 port
For power supply of the binocular tube 25°/22 w/int cam 8 MPX, if mounted additionally

2 RJ12 port
For power supply of the reflected-light FL iLED

3.2.4 Binocular Tube 25°/22 w/int Cam 8 MPx

Purpose The Binocular tube 25°/22 w/int cam 8 MPx is used to visualize and transfer the microscopic image to an external medium, such as a separate monitor, PC or tablet PC, via data line or Wi-Fi connection.

Position The Binocular tube 25°/22 w/int cam 8 MPx can only be mounted on the stand without small handle lift.


Fig. 18: Controls of the Binocular tube 25°/22 w/int cam 8 MPx

1 LED indicator:
Shows the current camera working status.

2 **Reset** button:
Restarts the camera by continuously pressing the button for 5 seconds using a clip.

3 **Snap** button:
Press briefly to trigger an exposure, press and hold to start a video recording, press briefly to stop a video recording

The states of the Binocular tube 25°/22 w/int cam 8 MPx are displayed by LED indicator.

Color	State	Status
pink	continuously	Start sequence
blue	flashing	Snapping/recording in process, and saving data to USB flash drive
blue	continuously	Power supplied and camera ready
red	flashing	Firmware updating/factory resetting
pink	flashing	No USB flash drive or the drive is full


Fig. 19: Cover plate on the rear side of the Binocular tube 25°/22 w/int cam 8 MPx

1 Cover plate
For covering the USB TYPE A port and the RJ12 port

2 Cover plate lock screw
For locking the cover plate


Fig. 20: Connections and special components of the Binocular tube 25°/22 w/int cam 8 MPx on the rear side

1 USB TYPE B
For image data transfer to a PC

2 Fast Ethernet port (RJ45)
For communication and image data transfer

3 USB TYPE A
Especially for Wi-Fi adapter

4 USB TYPE A
Optional: For storing image data to memory stick or for mouse or keyboard for OSD control

5 USB TYPE A
Optional: For storing image data to memory stick or mouse or keyboard for OSD control

6 HDMI port
For image data transfer to a monitor

7 RJ12 port
For power supply from stand

8 Cover plate
For covering the USB TYPE A port and the RJ12 port

Function The visualized images of the Binocular tube 25°/22 w/int cam 8 MPx may only be used for training and research.

Info

Direct generation of diagnostic results from these images is not recommended.

3.2.4.1 On Screen Display (OSD) Menu of the Binocular Tube 25°/22 w/int Cam 8 MPx

When the Binocular tube 25°/22 w/int cam 8 MPx is powered and connected to a display device via HDMI, the OSD menu will be shown automatically on the connected display device.

3.2.4.1.1 Home Menu


*Fig. 21: **Home** menu, control elements*

1 **Home** button
Opens the **Home** menu.

2 **Snap** button
Snaps a single image.

3 **Record** button

Starts a video recording. A time counter is displayed.

When pushing the button the following controls expand:

Enables you to manually re-adjust the exposure time, e.g. when the **Auto Exposure** mode does not deliver appropriate results.

When pushing the button the following controls expand:

The recording can be stopped by clicking the red square symbol above the time display.

4 **Auto Exposure** button

Ensures consistent brightness of the image by continuously calculating the correct exposure time based on the current light intensity.

5 **Target Intensity** controls

In **Auto Exposure** mode, the light intensity can be fine-tuned using the slider or the input field.

6 **Auto White Balance** button

Ensures a consistent color temperature of the image by continuously calculating the white balance.
Note: Proper functioning of the **Auto White Balance** mode requires enough empty space in the field of view.

7 **White** button

If active, in the **Auto White Balance** mode the color temperature is calculated based on the assumption that the light source is white.


8 **Settings** button

Opens the **Settings** menu.

9 **Manual Exposure** button

Enables you to manually re-adjust the exposure time, e.g. when the Auto Exposure mode does not deliver appropriate results.

When pushing the button, the following controls expand:

The exposure time can be fine-tuned using the upper slider or input field.

The gain value can be set using the lower slider or input field.

10 **Manual White Balance** button

Enables you to manually re-adjust the color temperature, e.g. when the Auto White Balance mode does not deliver appropriate results.

When pushing the button, the following controls expand:

The color temperature can be set warmer (more red) or colder (more blue) using the slider.

The white balance can be calculated once by clicking the **Push** button.

11 **Thumbnail** area

Shows a thumbnail of the last image taken, the first frame of the last video captured, or the merged image of the last multi-channel acquisition with overlaid pseudo-colors.

The thumbnail can be closed by clicking the **Close** symbol.

The thumbnail can be maximized to the full display by clicking the **Expand** symbol.

3.2.4.1.2 Settings Menu


*Fig. 22: **Settings** menu, control elements*

1 **Image Settings** icon

Opens the **Image Settings** menu.

2 **Microscope System Settings** icon

Opens the **Microscope System Settings** menu.

3 **Operating System Settings** icon

Opens the **Operating System Settings** menu.

4 **System Information** button

Shows the status of the entire camera/ microscope system.

5 **Back** button

Opens the previous menu.

3.2.4.1.2.1 Image Settings Menu


*Fig. 23: **Image Settings** menu, control elements*

1 **Gamma** controls

Enable you to adjust the gamma value.

The gamma value can be fine-tuned using the slider or the input field.

The gamma value can be reset to default by clicking the Reset button.

Note: Default gamma value is 0.45.

2 **Optimize Denoise** button

Reduces noise.

3 **Image Orientation** controls

Sets the image orientation.

0°: original image orientation

Mirror vertical: image mirrored in vertical direction

Mirror horizontal: image mirrored in horizontal direction

180°: image rotated by 180°

4 **Show ... Scale Bar** button

Shows a scale bar on the live display.

Note: The scale bar will also appear in the captured image.

5 **Show ... Image after Snap** button

If activated, the image is shown on the full display for five seconds after being captured. The full display image can be closed or minimized to a thumbnail if required.

3.2.4.1.2.2 Microscope System Menu


*Fig. 24: **Microscope system** menu, control elements*

1 **Objectives** button (only for Full-Köhler stands)

Opens a drop-down list to define the item installed at the objective turret's current position.

2 **Shading Correction** button

Opens a menu to define the shading correction for each combination of reflector, objective and light source.

3 **Configuration Wizard** button

Opens a wizard to support you in setting up all microscope configurations listed above.

3.2.4.1.2.3 Operating System Menu

The **Operating System** menu provides options such as setting the language, defining the file-name format, and updating the system firmware.


*Fig. 25: **Operating System** menu, control elements*

1 **Language** button

Opens a menu to set English, Chinese or German as operating system language.

2 **Date & Time** button

Opens a menu to set date and time.

3 **File Options** button

Opens a menu to set the filename format and the file type for captured images, see *File Options Menu [▸ 36]*.

4 **HDMI Resolution** button

Opens a menu to set the desired HDMI resolution for live view.

Notes:

Switching from 1080p to 4K must be confirmed by pressing the **Confirm** button.

The camera restarts automatically for the new resolution to take effect.

If the connected monitor does not support 4K it automatically switches back to 1080p after 15 seconds.

5 **Wi-Fi Settings** button (for USB Wi-Fi adapter)

Opens a workflow to set up a wireless connection, see *Wi-Fi Settings Menu [▸ 36]*.

6 **Bright Pixel Correction** button

Opens the pixel correction dialog.

7 **Firmware Update** button

Starts a firmware update when the update file is available via the USB interface, see *Updating the Firmware of the Binocular Tube 25°/22 w/int Cam 8 MPx [▸ 81]*.

3.2.4.1.2.3.1 File Options Menu


*Fig. 26: **File Options** menu, control elements*

1 Filename format example

Shows the template of the filenames based on the currently selected components.

2 Filename format components

Allows you to adjust the filename template.

Static text components (e.g. for sample names etc.) can be added by typing text into the Custom fields. Empty fields will be removed from the template.

The date and time format can be set.

A counter number is appended to each filename by default.

3 **Set Filename after each Snap** button

Activates a prompt after each image capturing to manually set the filename.

4 File type selection

Allows you to set the file type (**JPEG** or **TIFF**) for captured images.

3.2.4.1.2.3.2 Wi-Fi Settings Menu

When a compatible USB Wi-Fi adapter is inserted in the Binocular tube 25°/22 w/int cam 8 MPx, the **Wi-Fi Settings** menu will be displayed.

It offers the option to connect the camera to a Wi-Fi device (e.g. iPad or PC).

3.2.5 Abbe Condenser 0.9/1.25 Field 20

Purpose The condenser focuses the light from the light source to illuminate the object as evenly as possible, and to guide it as accurately as possible into the imaging beam path.

Position The condenser is mounted on the condenser carrier behind the specimen stage.


Fig. 27: Abbe condenser and turret condenser

1	Abbe condenser	**2**	Magnification marker for aperture setting
3	Lever for adjusting the aperture diaphragm	**4**	Slot for phase contrast or darkfield slider
5	Dovetail ring mount	**6**	Front lens
7	Turret condenser	**8**	Display field for selected contrast stop
9	Turret disk with five positions for: Brightfield **H** Phase contrast **Ph1**, **Ph2**, **Ph3** Darkfield **D**	**10**	Phase ring adjustment holes

3.2.6 Eyepiece

Purpose The eyepieces (10x/20 Br. Foc. and 10x/22 Br. Foc.) serve to observe the microscopic image.

Position The eyepieces are inserted into the tube.


Fig. 28: Eyepiece

1 Foldover rubber eyecup

2 Focusing ring
for compensation of defective vision

3 Diopter scale
to facilitate finding the correct setting

Function Both eyepieces are suitable for spectacle wearers. Additionally, they contain a focusing ring for compensation of defective vision. The provided diopter scale helps to find the correct setting. When using the microscope with the Reflected-Light FL iLED for fluorescence applications, the special eyecups with light protection can be used. However, they cannot be folded over and are not suitable for spectacle wearers.

3.2.7 Cable and Screw Tool Storage

The main power cable and phase plate adjustment tools can be stored on the rear side of the microscope.
If reflected light FL iLED is used, the hex key 3 mm can be stored on its rear side.



Fig. 29: Storage places


1 Holder for power cable

2 Storage place for phase plate adjustment tool

3.2.8 Accessories for Connecting the Microscope


The following connecting cables and accessories are required for voltage supply and use of the ports:

Name	Figure	Remarks
Power cable		Connection between microscope and power outlet
RJ 12 cable		Connection between Binocular tube 25°/22 w/int cam 8 MPx and stand for power supply
HDMI cable (can be ordered in demo kit)		Connection between camera and monitor, TV, or projector
USB flash drive, Type A (not supplied)		Connection to camera for immediate image and video storage
USB hub, Type A (not supplied)		Connection between camera and multiple USB Type-A devices, such as keyboard, mouse and flash drive
Mouse (not supplied)		For control and navigation in the **OSD** menu
Keyboard (not supplied)		For typing in the **OSD** menu
USB Wi-Fi adapter (has to be ordered separately)		Connection between the camera and a network or iPad

Name	Figure	Remarks
USB 3.0 cable, Type B to Type A (not supplied)		Connection between camera and PC
Ethernet cable (can be ordered in demo kit)		Connection between camera and network or WLAN router

3.3 Objective Labeling

The functional characteristics of an objective are indicated by the respective labeling, e.g. **A-Plan 10x/0,25 ∞/-**.

Obligatory label components

1 **Magnification of the objective** (**10x**)

The magnification factor specifies the magnification of the objective. Larger magnifications allow to resolve smaller features on the sample.

The objective's magnification factor is also expressed by the objective's color ring.

The objective magnification multiplied by the eyepiece magnification gives the overall magnification.

2 **Numerical aperture** (NA, **0,25**)

The numerical aperture multiplied by 1000 represents the maximum useful magnification - above that limit no additional resolution is provided.

3 **Applicable cover glass thickness** (**-**)

- : The objective can be used without cover glass or with a cover glass with a thickness of 0.17 mm.

0: The objective should be used without cover glass.

0.17: The objective should be used with a cover glass with a thickness 0.17 mm.

4 **Mechanical tube length** (**∞**)

∞: The objective can be used with tubes of any length.

Additional label components

Immersion Oil (**Oil**)

The objective can be used with immersion oil.

Phase contrast objective (**Ph**)

The objective can be used in phase contrast microscopy.

Image quality (**i**)

The image quality is improved.

Info

Due to their short working distance, the following objectives have a resilient mount (specimen protection):

- 20x
- 20x/Ph 2
- 40x
- 40x/Ph 2
- 100x dry objective
- 100x/Ph 3 oil
- 100x oil

4 Installation

Perform only the installation work described in this document. All other installation work not described may only be carried out by an authorized ZEISS service representative.

4.1 Safety During Installation

Before installing and starting up the microscope, be sure to carefully read and observe the notes on instrument safety, see *Safety [▸ 9]*.

NOTICE

Pollution of the optics

Dirty optics impair the function of the microscope.

- Do not touch optical surfaces when unpacking the microscope to avoid fingerprints!

4.2 Preparing Installation

The microscope is supplied completely assembled and including accessories that are packed to commercial standards.

Additionally ordered components, such as sliders and the transmitted-light module with illuminating mirror or Reflected-light FL iLED, are delivered in separate packages and must be mounted to the microscope.

Prerequisite

- ✓ The floor on which the support/table stands is subject to the least possible vibration.
- ✓ The surface of the support/table is hard and non-flammable.

Procedure

1. Remove the microscope from the transport case and place it on the work table. The distance of the Microscope to the wall should be at least 9 cm, in order to ensure sufficient air circulation and accessibility of the cabling.
2. Retain the original packaging for longer-period non-use storage of the instrument or for return to the manufacturer.


4.3 Installing Additional Components

4.3.1 Installing or Changing the Tube


Prerequisite ✓ The microscope is unplugged from the mains power.

Procedure

1. Insert the tube to be installed into the stand.
 → The eyepieces must point to the right and its dovetail mount must be slightly inclined, under the two holding elements.

2. Place the tube horizontally onto the stand.
 → The groove on the underside of the tube must be located above the third holding element of the stand.
3. Turn the tube counterclockwise by 90°, align it to the stand.
 → The eyepieces are pointing frontward.

4. Tighten the knurled clamping screw.

5. To detach the tube, loosen the knurled clamping screw.
6. Turn the installed tube clockwise by approximately 90° and detach it on the right side upward.

Info

For space-saving storage of the microscope (e.g. in a cabinet) the tube may also be rotated by 180° to the back.

Info

For reasons of space, the tube may also be clamped by the Allen set screw supplied with the tube.


4.3.2 Installing the Reflected-Light FL iLED onto the Stand

Prerequisite


- ✓ The microscope is unplugged from mains power.
- ✓ Stand without small hand lift.
- ✓ RJ12 cable is available.

Procedure


1. Insert the Reflected-light FL iLED at correct angle and slightly inclined with its dovetail mount into the stand.


2. Position the Reflected-light FL iLED horizontally.


3. Turn the Reflected-light FL iLED clockwise to face backwards with the dovetail mount in the stand.

4. Plug the RJ12 cable into the RJ12 female connector port of the stand.


5. Align the Reflected-light FL iLED to the outer edges of the stand.


6. Tighten the set screw of the stand.

7. Plug the other side of the RJ12 cable into the RJ12 female connector port of the Reflected-light FL iLED.

8. Hide the RJ12 cable at the rear side groove of the Reflected-light FL iLED.
9. Place the tube onto the Reflected-light FL iLED (see *Installing or Changing the Tube [▸ 43]*) and tighten the clamping screw of the Reflected-light FL iLED.

10. If required, place the yellow filter onto the luminous-field diaphragm, see *Assembling the Yellow Filter or the TL Cover Plate (Fixed Köhler Stand) [▸ 47]*.

Info

Special eyecups with light protection (see *Changing the Eyecups [▸ 56]*) are used for fluorescence applications if no darkroom is available. However, they cannot be folded over and are, therefore, not suitable for spectacle wearers. For this reason, spectacle wearers should use the standard eyecups.


Info

For mounting the Reflected-light FL iLED on stands with small hand lift, please contact an authorized ZEISS service representative.


4.3.3 Assembling the Yellow Filter or the TL Cover Plate (Fixed Köhler Stand)

Procedure

1. Unscrew the cover cap from the luminous-field diaphragm and remove it.


2. Place the yellow filter onto the mounting surface of the luminous-field diaphragm or remove it if necessary.
3. Re-attach the cover cap to the luminous-field diaphragm and lock it.
4. If required, place the TL cover plate into the Abbe condenser phase plate slot.


4.3.4 Inserting Color Filters (Full Köhler Stand)

Procedure

1. Move the condenser carrier up as far as it will go by turning the knurled ring for condenser height adjustment.

2. Unscrew the cover cap from the luminous-field diaphragm.


3. Put the desired filter – yellow, green or blue – onto the mounting surface of the luminous-field diaphragm.
4. Screw the cover cap back into place.


4.3.5 Installing the Polarizer Plate (Full Köhler Stand)

Procedure

1. Move the condenser carrier up as far as it will go by turning the knurled ring for condenser height adjustment.

2. Unscrew the cover cap from the luminous-field diaphragm and store it for later use.


3. Put the polarizer plate onto the luminous-field diaphragm.
 → The position stop on the polarizer plate has to match with the position stop on the luminous-field diaphragm.

4.3.6 Installing the Analyzer Plate

Prerequisite
✓ The microscope is unplugged from the mains power.

Procedure
1. Remove the tube, see *Installing or Changing the Tube [▸ 43]*.
2. Place the analyzer plate into the beam path on the top on stand.
 → The position stop on the analyzer plate has to match with the position stop on the stand.


3. Re-attach the tube on the stand.


4.3.7 Installing the Binocular Tube 25°/22 w/int Cam 8 MPx

Prerequisite
✓ The microscope is unplugged from mains power.
✓ Stand without small hand lift
✓ RJ12 cable
✓ If needed, mount the Reflected-light FL iLED onto the stand first, see *Installing the Reflected-Light FL iLED onto the Stand [▸ 44]*. Then, mount the Binocular tube 25°/22 w/int cam 8 MPx onto the Reflected-light FL iLED in the same way as described below.


Procedure
1. Remove the tube, see *Installing or Changing the Tube [▸ 43]*.
2. Unscrew the cover plate.

3. Insert the Binocular tube 25°/22 w/int cam 8 MPx, at the correct angle and slightly inclined with its dovetail mount into the stand.


4. Position the Binocular tube 25°/22 w/int cam 8 MPx horizontally.


5. Turn the Binocular tube 25°/22 w/int cam 8 MPx clockwise to face backwards with the dovetail mount in the stand.

6. Plug the RJ12 cable into the RJ12 female connector port of the stand.
 If the Reflected-light FL iLED is used, plug the RJ12 cable into the upper RJ12 female connector port of the Reflected-light FL iLED.


7. Align the Binocular tube 25°/22 w/int cam 8 MPx to the outer edges of the stand.


8. Tighten the set screw of the stand.

9. Plug the other side of the RJ12 cable into the RJ12 female connector port of the Binocular tube 25°/22 w/int cam 8 MPx.

10. Screw on the cover plate. Make sure the RJ12 cable is invisible.


4.3.8 Installing a Camera to the Trinocular Tube

Prerequisite


- ✓ A trinocular tube (photo tube) is mounted on the microscope.
- ✓ A camera adapter P95-C 2/3" 0.65x or P95-C 1/2" 0.5x

Procedure


1. Loosen the clamping screw and remove the dust cap from the tube.

2. Mount the C-mount camera adapter on the camera.

3. Attach the camera with the adapter to the microscope's C-mount port.

4. Orient the camera to the stand and fix its position by tightening the clamping screw.


Info

With camera/adapter combinations that have not been recommended expressly by ZEISS, it may be quite impossible to obtain an unvignetted image.


4.3.9 Exchanging Objectives

Procedure

1. Turn the focusing drive to move the mechanical stage down as far as it will go.
2. Turn the nosepiece to move the objective to be changed into a lateral position.
3. Unscrew the objective using the rubber strip included in the delivery and remove it downward.

4. Screw the desired objective handtight into the nosepiece as far as it will go.

5. If you intend to insert the objective in a previously unused lens mount, remove the dust cover from the corresponding mount of the nosepiece.

4.3.10 Installing/Removing the Condenser


The mirror serves to illuminate the specimen if no power outlet is available.

Parts and Tools Hex Key


Prerequisite ✓ Full-Köhler stand

Procedure


1. Turn coarse focusing drive to move the stage carrier to the upper stop position.

2. Lower the condenser carrier as far as it will go, using the knurled screw for vertical adjustment.

3. Loosen the clamping screw of the condenser, so that the condenser can be removed frontward.
 → Use an Allen wrench if the clamping screw is an Allen screw.

4. Remove the installed condenser, e.g. Abbe condenser.
5. Insert the condenser to be installed, e.g. turret condenser, with its dovetail ring mount into the condenser carrier.
6. Fix it by tightening the clamping screw.

4.3.11 Installing/Removing the Mirror


The mirror serves to illuminate the specimen if no power outlet is available.

Parts and Tools Hex Key


Prerequisite ✓ Fixed-Köhler stand or Full-Köhler stand

Procedure

1. Remove the cover cap from the luminous-field diaphragm, see *Assembling the Yellow Filter or the TL Cover Plate (Fixed Köhler Stand) [▸ 47]*.
2. Remove the condenser, see *Installing/Removing the Condenser [▸ 53]*.
3. Unscrew the 2 clamping screws from the condenser carrier insert.

4. Push the condenser carrier insert backward against the spring, and inclining it, remove it upward from the condenser carrier.
5. Insert the mirror from the top through the opening of the condenser carrier and place it onto the mount of the luminous-field diaphragm. Take care that the mirror rests horizontally on the mount.

6. Rotate and incline the mirror until the daylight is reflected homogeneously into the light path.


Info

Follow these instructions in reverse order to remove the mirror.

4.3.12 Folding Over the Eyecups

Procedure


1. If required, fold over the rubber eyecups.

4.3.13 Changing the Eyecups

Procedure

1. Remove the existing eyecup from the eyepiece, e.g. the foldover rubber eyecups.

2. Attach the desired eyecups, e.g. the special eyecups with light protection.

4.4 Connecting the Microscope to the Power Supply

4.4.1 Connecting the Microscope to the Mains Power Supply

All microscopes equipped with Fixed-Köhler stands or Full-Köhler stands can be connected directly to the mains power supply.

Prerequisite

✓ The microscope is installed completely with all components to be used.

Procedure

1. Plug the power cable 3 into the mains socket 2 of the microscope.


Fig. 30: Connecting the microscope to the mains power supply

2. Make sure the **DC 5V** switch 1 is in the **OUT** position.
 → The **DC 5V** switch is only available on Fixed-Köhler stands.
3. Plug the other end 5 of the power cable into the mains power supply 4.

4.4.2 Connecting the Microscope to a Power Bank

All microscopes equipped with Fixed-Köhler stands can be connected to a power bank for power supply.

Prerequisite

✓ The microscope is installed completely with all components to be used.

✓ The mobile power bank is charged.

Procedure

1. Plug the USB Type A to A cable 3 into the 5V USB port 2 of the microscope.


Fig. 31: Connecting the microscope to a power bank

2. Switch the **DC 5V** switch 1 to the **IN** position.
3. Plug the other side of the USB Type A to A cable 5 into the appropriate port of the power bank 4.

4.5 Connecting the Binocular Tube 25°/22 w/int Cam 8 MPx to External Devices

The Binocular tube 25°/22 w/int cam 8 MPx can be connected to varios external devices.

The following options are available:

- *Connecting a USB Flash Drive [▸ 57]*
- *Connecting to the PC via USB 3.0 Port [▸ 58]*
- *Connecting to a Display (without PC) [▸ 58]*
- *Integrating the Binocular Tube 25°/22 w/int Cam 8 MPx into a Network [▸ 59]*


4.5.1 Connecting a USB Flash Drive

Prerequisite

✓ The USB flash drive has to be of FAT32 format and have enough free space for storing the data

Procedure

1. Insert the flash drive into one of the USB Type A ports of the Binocular tube 25°/22 w/int cam 8 MPx.

Info

To view the captured images or videos connect the USB flash drive to a PC or laptop.


4.5.2 Connecting to the PC via USB 3.0 Port

Prerequisite

- ✓ USB 3.0 specification defined as maximum: 500 mA / USB port
- ✓ A USB 3.0 (Type B and Type A) cable is required

Procedure

1. Insert the USB 3.0 cable connector into the fast Ethernet port (RJ45).

2. Insert the other USB 3.0 cable connector into the corresponding socket of the PC or laptop.

4.5.3 Connecting to a Display (without PC)


The Binocular tube 25°/22 w/int cam 8 MPx can be connected to a monitor, TV or projector for visualization of the live image data and for operating the On Screen Display menu (OSD) functions.

Prerequisite

- ✓ The monitor/projector supports 1080p and 4K HD formats.
- ✓ To enable the display in HD quality, it is recommended that only HD-capable devices are connected.
- ✓ If connecting to a monitor, use a monitor that supports 4K for live imaging at 4K resolution (e.g. Monitor TFT 32" HP Z32)
- ✓ If a monitor/projector is used, an HDMI/DVI adapter cable or HDMI/DVI adapter is required for the reproduction of image data.

Procedure

1. Connect the monitor/projector to the Binocular tube 25°/22 w/int cam 8 MPx using the HDMI cable.

2. Secure the HDMI plug on the jack of the monitor/projector.
3. Set the display device's aspect ratio to 16:9 or Aspect.

4.5.4 Integrating the Binocular Tube 25°/22 w/int Cam 8 MPx into a Network

The Binocular tube 25°/22 w/int cam 8 MPx is able to communicate and transfer data to network via Ethernet port or Wi-Fi adapter and work with the **ZEISS Labscope** app.

A WLAN router is required to use the built-in camera together with the **ZEISS Labscope** app.

An existing WLAN infrastructure or a separate microscope WLAN setup may be used.

To use an existing WLAN infrastructure, a High-performance 802.11n WLAN with sufficient free bandwidth is required. In the event of an overloaded or slow WLAN the live image of the built-in camera may be delayed or incorrectly displayed on the iPAD.


To create a separate microscope WLAN set up, visit **www.zeiss.com/micro-apps** for further information.

The integrated camera will identify itself automatically to the network on DHCP and will be automatically recognized by **Labscope**, provided the iPAD is in the same network.

4.5.4.1 Connecting the Binocular Tube 25°/22 w/int Cam 8 MPx via Ethernet

Procedure

1. Insert the Ethernet cable into the jack on the reverse side of the Binocular tube 25°/22 w/int cam 8 MPx.


2. Insert the Ethernet cable's opposite connector into the corresponding socket on your WLAN router.
3. Switch on the router.
 → **Labscope** can recognize the network (network name and password are the same as in the WLAN router setting).

4.5.4.2 Connecting the Binocular Tube 25°/22 w/int Cam 8 MPx via Wi-Fi Adapter

Procedure

1. Unscrew the plastic cover plate.

2. Insert the Wi-Fi adapter into the USB type A port.
3. Replace the plastic cover plate and lock it with a screwdriver.
 → **Labscope** can recognize the network (Network name: Primostar 3_last six numbers from MAC address (MAC address label is on top of the Binocular tube 25°/22 w/int cam 8 MPx, for example
 Primostar3_F9A919)
 Password: ZEISS1846

Info

When a compatible USB Wi-Fi adapter is inserted into the Binocular tube 25°/22 w/int cam 8 MPx, the **Wi-Fi Settings** menu will be displayed. It offers advice for connecting the camera to a Wi-Fi device (e.g. iPad or PC).

Info

Consult your network administrator for further guidance.

Further information of all ZEISS Microscopy apps can be found at https://www.zeiss.com/microscopy/int/products/microscope-software/microscopy-apps.html.

4.6 Installing Labscope Software for Using the Binocular Tube 25°/22 w/int Cam 8 MPx

To capture images with the Binocular tube 25°/22 w/int cam 8 MPx, you have the following options:

- Installation of **Labscope** for Windows to PC
- Installation of **Labscope** to iPad

4.6.1 Installing the Labscope Software on PC

Prerequisite ✓ WINDOWS PC is required.

Procedure

1. Go to **www.zeiss.com/labscope** for free download.
2. Install the Labscope App as described.

4.6.2 Installing the Labscope App on iPad

Prerequisite ✓ iPad is required.

Procedure

1. Go to the App Store.
2. Search for **Labscope** by ZEISS.
3. Install the app as decribed.

5 Operation

This chapter describes switching on/off the microscope as well as the first operating steps with the microscope.

Info

For additional information and detailed descriptions, refer to further applicable documents or ask your ZEISS Sales & Service Partner.

Info

Further information on the software and its operation is available in the software's online help.

5.1 Prerequisites for Commissioning and Operation

The following basic prerequisites are necessary for commissioning and operation:

- The Instruction Manual was read prior to commissioning or operation and kept for further use.
- The chapter **Safety** was read and understood.
- The operator is acquainted with the general Windows®-based programs.
- If required: Basic training and safety briefing were successfully completed.


5.2 Switching On the Microscope

5.2.1 Switching On the Microscope with Mains Power Connection


Prerequisite ✓ *The microscope is connected to the mains power supply [▸ 56].*

Procedure

1. At the Fixed-Köhler stand, set the **DC 5V** switch to the **OUT** position.

2. Switch on the **main power ON/OFF** button.


5.2.2 Switching On the Microscope with Power Bank Connection

Prerequisite ✓ *The microscope is connected to the power bank [▸ 56].*


✓ The power bank is charged.

Procedure

1. Set the **DC 5V** to the **IN** position.

2. Switch on the power bank.
3. Switch on the **main power ON/OFF** button.

5.3 Adjusting the Tube

This procedure consists of the following steps:

Prerequisite ✓ The microscope is connected to the mains and is switched on.

Procedure

1. *Adjusting the Position of the Eyepieces [▸ 63]*
2. *Adjusting the Eyepiece with Eyepiece Pointer or Eyepiece Micrometer [▸ 64]*

5.3.1 Adjusting the Position of the Eyepieces

Info


The adjustment of the interpupillary distance is correct when you see only one round image while looking through the two eyepieces.

Procedure

1. Set the interpupillary distance by rotating the eyepiece tubes symmetrically toward or away from one another.


2. Set the viewing height by swivelling the eyepieces a full 180° upwards or downwards.


5.3.2 Adjusting the Eyepiece with Eyepiece Pointer or Eyepiece Micrometer

Procedure


1. Put the eyepiece with pointer or micrometer into the tube and fix it with the set screw.
2. Turn the focusing ring of the eyepiece to focus on the wedge-shaped figure of the eyepiece pointer.

3. Put the specimen onto the mechanical stage.
4. Look at the specimen through the eyepiece with the eyepiece pointer or micrometer.
5. Use the focus drive to bring the microscopic image into focus until both the microscopic image and the eyepiece pointer are sharply defined.

6. Focus the image for the second eye by turning the focusing ring of the second eyepiece.

↳ Having done so, both microscopic images inclusive of the eyepiece pointer or micrometer are focused.

Info

Afterwards, you should focus on the specimen using the focusing drive only.

5.4 Adjusting the Transmitted-Light Brightfield on the Full-Köhler Microscope

Parts and Tools

- High contrast specimen slide with cover glass of 0.17 mm thickness

Info


The knurled knob for the vertical condenser adjustment is always on the side opposite to the stage controls.

Prerequisite

- ✓ The microscope is adapted to the user, see *Adjusting the Tube [▸ 63]*.
- ✓ A 10x objective for use with cover glass of 0.17 mm thickness is mounted.

Procedure


1. Rotate the 10x objective into the light path using the knurled ring of the nosepiece.


2. Place the high-contrast specimen slide on the stage and fix it using the spring lever of the specimen holder.


3. If applicable, move all phase contrast or darkfield elements out of the light path.
4. If the microscope stand is equipped with the turret condenser, move turret to **BF** position.
5. When using the Reflected-light FL iLED, turn the **fluorescence/brightfield** changeover switch to the **brightfield** position (turn the switch fully upward first).

6. Use the rotary knob for illumination intensity to adjust the illumination to a comfortable setting.


7. Use the knurled knob for vertical condenser adjustment to move the condenser to the upper limit stop.


8. Use the control lever of the aperture diaphragm to set the aperture diaphragm to the middle position.


9. Look through one eyepiece and use the focusing drive to bring the specimen into focus.


10. If necessary, readjust the image sharpness for the other eye by turning the adjustable eyepiece.
 → The image of the high-contrast specimen is sharp for both eyes.
11. Close the luminous-field diaphragm until its edges become visible in the field of view.


12. Use the knurled knob for vertical adjustment to slowly move the Abbe condenser downwards until the edge of the luminous-field diaphragm becomes sharp.


13. Use the centering screws to center the image of the luminous-field diaphragm in the field of view.

14. Open the luminous-field diaphragm until its edges just disappear from the field of view.

15. To adjust the contrast, remove one eyepiece from the tube and look through the tube.
16. Use the control lever of the aperture diaphragm to adjust the aperture diaphragm to approximately 2/3 to 4/5 of the diameter of the exit pupil of the objective.
 → The optimum contrast at almost ideal resolution for most applications (best compromise for the human eye) is set.

17. Reinsert the eyepiece into the tube.

Info

The size of the field of view and the objective aperture change with every objective change. Therefore, the adjustment of the luminous-field diaphragm and the aperture diaphragm according to Köhler has to be repeated to obtain optimum microscopy conditions.


5.5 Adjusting the Transmitted-Light Brightfield on the Fixed-Köhler Microscope

The Fixed-Köhler microscope is supplied factory-adjusted. The operation is restricted to a few manipulations.


Prerequisite ✓ The microscope is adapted to the user, see *Adjusting the Tube [▸ 63]*.

Procedure


1. Place the specimen in the specimen holder of the mechanical stage.


2. If the microscope stand is equipped with a dark-field slider, pull this slider to the left up to the lock-in position.
3. When using the Reflected-light FL iLED, turn the **fluorescence/brightfield** changeover switch to the **Brightfield** position (turn the switch fully upward first).


4. Push the slider containing the yellow filter with its filter position into the light path.
5. Adjust to the desired magnification by rotating the corresponding objective into the light path.


6. Set the control lever of the condenser aperture diaphragm to the value of the selected magnification (**10x**, **40x** or **100x**).

7. Focus on the specimen using the focusing drive.

8. Use the rotary knob for illumination intensity to adjust the illumination to a comfortable setting.

If the condenser has been removed (e.g. for installing the illuminating mirror), make sure to reinstall and center it using the two centering screws, see *Installing/Removing the Mirror [▸ 54]*.

5.6 Adjusting the Transmitted-Light Phase Contrast or the Transmitted-Light Darkfield Using a Contrast Slider

Parts and Tools


- 2 x Hex key, 1.5 mm

Prerequisite


- ✓ The phase contrast method on the microscope requires a stand with Full-Köhler equipment.
- ✓ Contrast slider for phase contrast or darkfield.

Procedure


1. Adjust the microscope as you would for transmitted-light brightfield, see *Adjusting the Transmitted-Light Brightfield on the Full-Köhler Microscope [▸ 65]*.
2. Turn the nosepiece to move the phase contrast objective (**Ph 2**) into the light path.


3. Open the luminous-field diaphragm on the stand.

4. Remove the dummy slider from the Abbe condenser.


5. Remove the locking screw from the phase contrast slider (fitting the objective used, **Ph 2**). Push the phase contrast slider from the left into the Abbe condenser until it reaches the lock-in position.


6. Screw the locking screw into the contrast slider.


7. Fully open the aperture diaphragm of the Abbe condenser using its control lever (left limit stop).


8. Adjust the illumination intensity as required.

9. To center the phase stop, remove one eyepiece and replace it with the diopter or telescope.
10. If it is necessary to center the phase stop, insert the two Allen wrenches into the adjusting screws of the phase slider.

11. Center the phase stop from position A to position B by turning the two adjusting screws of the slider using the two Allen wrenches.

12. Afterwards, replace the diopter or telescope with the eyepiece again.

Info

For darkfield application, use the darkfield slider instead of the phase contrast slider.

5.7 Adjusting the Transmitted-Light Phase Contrast or the Transmitted-Light Darkfield Using the Turret Condenser

Prerequisite


- ✓ The phase contrast method on the microscope requires a stand with Full-Köhler equipment..
- ✓ The turret condenser with phase contrast stop or darkfield stop inserted is required.
- ✓ The microscope is adjusted for *transmitted-light brightfield [▸ 65]*.

Procedure


1. Turn the nosepiece to move the phase contrast objective (**Ph 1** or **Ph 2** or **Ph 3**) into the light path.


2. Open the luminous-field diaphragm on the stand.


3. Turn the turret disk of the turret condenser to the phase contrast position (**Ph 1** or **Ph 2** or **Ph 3**) that corresponds to the objective being used.


4. Fully open the aperture diaphragm of the turret condenser using its control lever (left limit stop).

5. Adjust the illumination intensity as required.

6. To center the phase stop, remove one eyepiece and replace it with the diopter or telescope.
7. If necessary, center the phase stop from position A to position B by turning the two knurled screws.

8. Replace the diopter or telescope with the eyepiece.

Info

For darkfield application, use the darkfield position of the turret condenser instead of the phase contrast position.

Info

If the turret condenser is not equipped with a darkfield stop, contact an authorized ZEISS service technician for installation.

5.8 Adjusting the Transmitted-Light Simple Polarization Contrast

Prerequisite

- ✓ The simple polarization contrast on the microscope requires a stand with Full-Köhler equipment.
- ✓ The polarizer plate is inserted into the luminous-field diaphragm in the correct position, see *Installing the Polarizer Plate (Full Köhler Stand) [▸ 48]*.
- ✓ The analyzer plate is inserted into the stand in the correct position, see *Installing the Analyzer Plate [▸ 49]*.

Procedure

1. Adjust the microscope as you would for transmitted-light brightfield, see *Adjusting the Transmitted-Light Brightfield on the Full-Köhler Microscope [▸ 65]*.
2. Place the specimen to be examined under polarized light on the stage and fix it using the spring lever of the specimen holder.


5.9 Adjusting the Reflected-Light Fluorescence

Prerequisite


- ✓ The reflected-light fluorescence method requires the microscope with Reflected-Light FL iLED.

Procedure


1. Adjust the microscope as you would for transmitted-light brightfield, see *Adjusting the Transmitted-Light Brightfield on the Full-Köhler Microscope [▸ 65]*.
2. Turn the nosepiece to move the objective for fluorescence application into the light path (e.g. objective 40x).

3. Set the **fluorescence/brightfield** changeover switch to the**fluorescence** position (first turn it fully upward).


4. Switch on the reflected-light illumination using the rotary knob of the Reflected-light FL iLED and adjust the illumination intensity to a level comfortable for observation.


→ The pilot lamp at the front of the Reflected-light FL iLED is blue. The brightness of the pilot lamp corresponds to the illumination intensity adjusted for reflected light.

5. Focus on the specimen using the focusing drive.

6. To avoid interfering fluorescences (caused by the transmitted-light LED), place the yellow filter onto the luminous-field diaphragm (see *Inserting Color Filters (Full Köhler Stand) [▸ 47]*) and close the aperture diaphragm of the condenser completely.

Info

If the yellow filter is not available, about 90 % of the interfering fluorescence can be eliminated by closing the condenser aperture diaphragm.

It is also possible to insert the TL cover plate into the slot of the Abbe condenser (see *Assembling the Yellow Filter or the TL Cover Plate (Fixed Köhler Stand) [▸ 47]*).

Info

Our special eyecups with light protection can be used for microscopic examinations in rooms which are not darkened. They are, however, not suitable for spectacle wearers and must not be folded over because the required dimensional stability would be lost.


5.10 Capturing Images and Videos Using the Integrated Smart 8 MP Color Camera

5.10.1 Capturing Images and Videos Using the Operator Panel

Prerequisite ✓ The USB flash drive is inserted into the USB A port.

Procedure

1. Position your specimen on the microscope and adjust the microscope to see a focused image through the eyepieces.
2. To take a single image, briefly press the **Snap** button on the camera panel.
 → The image is saved to the USB flash drive in either JPEG or TIFF format.

3. To start video recording, press and hold the **Snap** button on the camera panel for about five seconds.
4. To stop video recording briefly press the **Snap** button on the camera again.
 → The video is saved to the USB flash drive in MP4 format.

5.10.2 Capturing Images and Videos Using the OSD Menu

Prerequisite

- ✓ The USB flash drive is inserted into the USB A port.
- ✓ A mouse and keyboard are plugged into a USB hub, which is connected to the other USB A port.
- ✓ The Binocular tube 25°/22 w/int cam 8 MPx is connected to a monitor via HDMI.

Procedure

1. Position your specimen on the microscope and adjust the microscope to see a focused image on the monitor.
2. To take a single image, click **Snap** in the **OSD** menu.
 → The image is saved to the USB flash drive in either JPEG or TIFF format.

3. To start video recording, click on **Record** in the **OSD** menu.
 → The **Record** button changes to **Stop**.

4. To stop video recording click on **Stop** in the **OSD** menu.
 → The video is saved to the USB flash drive in MP4 format.

5.10.3 Capturing Images with Labscope

Prerequisite

- ✓ Wi-Fi adapter is inserted into the special USB port A; or the binocular tube 25°/22 w/int cam 8 MPx is connected with an Ethernet cable with WLAN. Connection may also be established with a USB 3.0 cable.

Procedure

1. Start **Labscope**.
 → When starting **Labscope** for the first time, each screen starts with an overlaying information screen explaining the functions.
2. Refer to the displayed information for using the software.
3. Disable or re-activate the information screens in the software's **Settings** menu on your iPad or PC.


Info

For further support in using **Labscope**, see www.zeiss.com/microscopy/community. Please check the **Labscope** thread for problem-solving notes.

5.11 Switching Off the Microscope

Procedure

1. After finishing work, switch off the microscope at the **main power ON/OFF** button.

2. If connected, switch off the power bank.
3. Cover the microscope with the dust cover.

6 Care and Maintenance

To ensure the best possible performance of the microscope, maintenance must be performed on a regular basis. Please keep the service logs for your microscope.

To maintain operational safety and reliability of the microscope, we recommend entering into a **ZEISS Protect Service Agreement**.

Info

For additional information and detailed descriptions, refer to further applicable documents or ask your ZEISS Sales & Service Partner.

6.1 Safety During Cleaning and Maintenance

Only conduct preventive measures described here. All tasks of maintenance, service, and cleaning not described must only be performed by an authorized ZEISS service representative.

Any unauthorized intervention or any operation outside the scope of the intended use can lead to injuries and property damage and voids all rights to warranty claims. Only original spare parts from ZEISS may be used.

⚠ DANGER

Electric shock due to live parts

When the microscope is still switched on, coming in contact with live parts can lead to electric shock or burn.

- Switch off the microscope prior to opening or cleaning.
- Disconnect live parts from the power supply.

NOTICE

Functional impairment due to dirt and moisture

Dirt, dust and moisture can impair the microscope functionality and can cause short-circuits.

- Use the dust protection cover if the microscope is not used.
- The ventilation slots must be unobstructed at all times.
- Perform regular maintenance and cleaning according to the instructions in this document and according to the instructions in the applicable documents.
- Make sure that no cleaning liquid or moisture gets inside the microscope.
- In case of damage, the affected parts of the microscope must be taken out of operation.

6.2 Maintenance Schedule

The recommended intervals for maintenance depend on the total uptime of the microscope.

Time Interval	Unit/Component	Activity
If required	6 V/30 W halogen lamp of the halogen module	Replace the 6 V/30 W halogen lamp, see *Exchanging the Halogen Lamp of the Halogen Module* [▸ 82].

Time Interval	Unit/Component	Activity
If required	Firmware update of the Binocular tube 25°/22 w/int cam 8 MPx	*Updating the Firmware of the Binocular Tube 25°/22 w/int Cam 8 MPx [▸ 81].*

Tab. 3: Maintenance Plan

6.3 Maintenance Work

6.3.1 Cleaning an Optical Surface

NOTICE

Damage of optical surfaces due to improper cleaning


- Remove dust from the optical surface slowly and carefully.
- Remove dust on optical surfaces with a natural-hair brush or blow it off with a rubber bellows.
- Avoid touching optical surfaces with fingers.

Parts and Tools

- Clean cloth
- Cotton swab
- Optical cleaning solution (85% n-hexan and 15 vol% isopropyl alcohol (IPA))
- Lint-free cloth

Procedure

1. Moisten a cotton swab or a clean cloth with an optical cleaning solution, if necessary.
2. Wipe optical surfaces in a circular motion towards the edge of the optics with slight pressure.

3. Dry with a lint-free cloth.

6.3.2 Removing Water-soluble Contamination

Parts and Tools

- Clean cloth
- Lint-free cloth

Procedure

1. Moisten a clean cloth with water.
 → A mild detergent may be added to the water (no solvent!).
2. Wipe off the area with the cloth.
3. Dry with a lint-free cloth.

6.3.3 Updating the Firmware of the Binocular Tube 25°/22 w/int Cam 8 MPx

Prerequisite

- ✓ The Binocular tube 25°/22 w/int cam 8 MPx is connected to a PC or laptop.
- ✓ USB flash drive is available (not included in the scope of delivery)

Procedure

1. Download the latest firmware from the website.
2. Save the update file to the root folder of the USB flash drive.
3. Insert the flash drive into the Binocular tube 25°/22 w/int cam 8 MPx's USB port or a connected USB hub.
4. In the **OSD** menu, click **Settings > Operating System Settings > Firmware Update**.
 → The firmware is updated automatically.

6.3.4 Exchanging the Halogen Lamp of the Halogen Module

⚠ DANGER

Electric shock due to live parts

When the microscope is still switched on, coming in contact with live parts can lead to electric shock or burn.

- Switch off the microscope prior to installing components.
- Disconnect live parts from the power supply.

⚠ CAUTION

Burning hazard due to hot halogen lamp

The halogen lamp gets hot during operation.


- Do not touch the halogen lamp before it has cooled down.

Parts and Tools

- Hex Key, 3.0 mm
- Halogen lamp 30W 6V

Procedure

1. Gently press the fastening screw 3 towards the microscope and turn it clockwise by 90°.

 → The halogen module is released.
2. Pull the halogen module 2 out of the stand.
3. Remove the halogen lamp 1 from the lamp holder and insert the new halogen lamp.
 NOTICE Do not touch the new lamp with bare fingers as this will reduce the lifespan of the lamp.
4. Push the halogen module containing the new lamp into the stand.
5. Gently press the fastening screw towards the microscope and turn it counterclockwise by 90°.
 → The halogen module is fixed.

7 Troubleshooting

The following table provides information about solving common problems.

Info

If you cannot solve the problem or if you are unsure about a certain technical difficulty, contact your local ZEISS service representative.

Symptom	Cause	Measure
The field of view is not completely visible.	Nosepiece with objective has not been switched into click-stop position.	Switch nosepiece with objective into click-stop position.
	Condenser has not been set correctly.	Set condenser correctly.
	The aperture diaphragm has not been adjusted correctly.	Adjust aperture diaphragm correctly.
	The luminous-field diaphragm has not been adjusted correctly.	Adjust luminous-field diaphragm correctly.
	The filter has not been inserted correctly into the filter mount.	Insert filter correctly into filter mount.
Low resolving power, poor image contrast	The aperture diaphragm has not been opened to correct size.	Set opening of aperture diaphragm to correct size.
	Condenser not focused correctly.	Focus the condenser.
	Wrong cover glass thickness selected for use of transmitted-light objectives corrected for 0.17 mm cover glass.	Use standard 0.17 mm cover glass.
	Use of no or non-specified immersion oil with immersion objective	Use the supplied immersion oil.
	Air bubbles in immersion oil	Remove the bubbles by applying new oil or moving the objective to and fro.
	Immersion oil on the front lens of a dry objective	Clean the front lens of the dry objective.
	Dirt or dust on the optical surfaces of objectives, eyepieces, condenser, or filters	Clean the respective optical components.
Greater focus differences after changing the objective	The adjustable eyepieces have not been set correctly.	Set the adjustable eyepieces to the corresponding defective vision.

Symptom	Cause	Measure
The 6 V/30 W halogen lamp or the LED source does not light up although the microscope has been switched on.	The power plug has not been plugged into the power outlet.	Connect power plug to power outlet.
	6 V/30 W halogen lamp or LED source is defective.	Replace the defective 6 V/30 W halogen lamp or the LED module.
The 6 V/30 W halogen lamp flickers, its light intensity is unstable, the illumination inhomogeneous	End of average service life of 6 V/30 W halogen lamp has been reached.	Replace the 6 V/30 W halogen lamp.
	Power cable incorrectly installed or broken.	Connect the power cable correctly or replace it.
	The pins of the 6 V/30 W halogen lamp have not been inserted correctly into the lamp holder.	Insert the pins of the 6 V/30 W halogen lamp correctly into the lamp holder.
	The pins of the 6 V/30 W halogen lamp have not been inserted symmetrically into the lamp holder.	Insert the pins of the 6 V/30 W halogen lamp symmetrically into the lamp holder.
Stage comes down by itself, image focus is unstable	Adjusted torque of coarse focusing drive is too low.	Increase the torque of the coarse focusing drive to make motion stiffer.
Transmitted-light illuminator cannot be switched on.	**Fluorescence/brightfield** changeover switch is in **fluorescence** position.	**Fluorescence/brightfield** changeover switch to **brightfield** position.
Reflected-light FL iLED cannot be switched on.	**Fluorescence/brightfield** changeover switch is in **brightfield** position.	**Fluorescence/brightfield** changeover switch to **fluorescence** position.
	Reflected-light FL iLED not sufficiently aligned to microscope stand or displaced, causing interruption of power supply by interlock function.	Align the Reflected-light FL iLED exactly to microscope stand and tighten clamping screw sufficiently.
The X direction range is smaller than expected.	Normal behaviour after long hours of use.	Recover the stage travel range: ▪ Hold the specimen holder's two screws. ▪ Unsing the control knob, move the specimen holder to the left and right until it hits the end stop. The stage travel range is reset.

7.1 Troubleshooting Binocular Tube 25°/22 w/int Cam 8 MPx

Troubleshooting notes for the use of the built-in camera in conjunction with the network/WLAN and Labscope app can be found separately at the ZEISS website www.zeiss.com/labscope.

Symptom	Cause	Measure
LED indicator is off.	The camera is not powered on properly.	Check RJ 12 cable connection to stand.
LED indicator flashes red.	The camera is updating firmware or being reset.	Warning! Do not switch off the main power supply on the stand.
Image/Video cannot be saved to the USB flash drive.	USB flash drive is not in correct format.	Please format the USB flash drive to FAT32 format on a PC.
	USB flash drive is full.	Make sure there is enough free space in the flash drive for storage.
	USB flash drive otherwise not recognized.	Restart the camera by powering off the main power supply on the stand.
Firmware update does not function.	The USB flash drive is not recognized or the firmware is not uploaded properly.	Insert a FAT32-formatted USB flash drive with latest firmware in the root folder and make sure the drive has at least 200 MB free memory space. Be sure to follow the instructions in the firmware update menu.
Date/Time is wrong on the camera.	Date/time is not set correctly.	Date/time can be set in OSD under Operating System setting.
	The buffer battery is empty.	Please contact service to change the battery.
The image is very noisy.	The amplification (gain) is set too high.	Open denoise; reduce the gain.
	The exposure time is set too high.	Open denoise; reduce the exposure time.
	The light intensity is set too low.	Open denoise; increase the light intensity.
The image is too dark or too bright.	Light source intensity too high or too low; or auto exposure time too short	Adjust light source intensity, or change to manual exposure.
The camera forgets the settings (e.g. manual white balance, HDMI resolution) if the power supply is disrupted.	Settings are not stored.	It takes 5 seconds for the settings to be automatically stored.

Symptom	Cause	Measure
Monitor connected via HDMI does not display an image.	The camera is not delivering a signal, or signal is not compatible with the monitor.	Ensure that the camera has been switched on for at least 30 seconds and the LED indicator is blue. Check the plug connections on the camera and monitor. Press factory reset button on the camera, see *Binocular Tube 25°/22 w/int Cam 8 MPx [▸ 28]*.
When switching to 4K, the monitor becomes black and then reverts to 1080p.	The monitor may not support 4K.	Use a 4K monitor such as mode 410350-2405-000 listed on ZEISS micro-shop website.
The image appears distorted on full screen monitor	The image aspect ratio is 16:9. The monitor may be set to a different aspect ratio and causes the distortion.	Set the monitor aspect ratio to 16:9.
The image is blurred on the screen but the sample is in focus through eyepieces	Focus plane of the camera is different from that of the eyepieces	Make sure the sample is focused correctly through eyepieces and calibrate the camera adaptor until image is in focus on the monitor.
No list of network names is shown when trying to connect the camera to a router via Wi-Fi USB adapter.	When the Wi-Fi setting is switching between modes, it takes ~10 s for the new mode to load.	Wait for 10 s and click again to refresh the drop-down list.
Buttons do not react when pressed.	The buttons are temporarily locked, as the camera is being accessed exclusively from the PC or network.	Control the camera via PC software or ZEISS app, or close the PC software or ZEISS app.
The camera shows unexpected behavior.	The camera settings may have been changed unintentionally.	Press factory reset button on the camera, see *Binocular Tube 25°/22 w/int Cam 8 MPx [▸ 28]*.

8 Decommissioning and Disposal

This chapter contains information on the decommissioning and disposal of the microscope and its expansions/components or accessories.

8.1 Decommissioning

If the microscope and its components are not used for an extended period such as several months, they should be shut down completely and secured against unauthorized access.

⚠ DANGER

Electric shock due to live parts

When the microscope is still switched on, coming in contact with live parts can lead to electric shock or burn.

- Switch off the microscope prior to opening or cleaning.
- Disconnect live parts from the power supply.

Procedure

1. Switch off the microscope.
2. Pull the mains plug.

8.2 Transport and Storage

The following regulations must be observed before and during transport:

- The boxes must be secured during transport.
- Avoid rocking the boxes back and forth.
- Note the weight information on the package and on the shipping document.
- Where possible, the original packaging must be used for shipping or transport.

Maximum shock resistance

- Do not drop or bump the boxes during movement or storage. Any acceleration shall be < 10 g.
- Evaluate packaging shock and tilting sensors on delivery and after internal transport.

Allowable temperature

Allowable temperature during on-site storage and transport in packaging:

- Between -40 °C and +70 °C
- Relative humidity less than 93 % at 40 °C

Info

24 hours before installation of the microscope it is required that the boxes be at recommended room temperature to avoid ingress of humidity, which is very harmful to optical paths, and to ensure effective stability of the microscope during installation and testing.

8.3 Disposal

The microscope and its components must not be disposed of as domestic waste or through municipal disposal companies. They must be disposed of in accordance with applicable regulations (WEEE Directive 2012/19/EU). ZEISS has implemented a system for the return and recycling of devices in member states of the European Union that ensures suitable reuse according to the EU Directives mentioned. The customer is responsible for decontamination.

Info

Detailed information on disposal and recycling is available from your ZEISS Sales & Service Partner.

8.4 Decontamination

A decontamination statement must be submitted before returning any used objects to the ZEISS location.

If reliable decontamination cannot be guaranteed, the hazard must be marked according to applicable regulations. In general, a well-visible warning sign must be affixed to the article itself and to the outside of the packaging, together with detailed information on the type of contamination.

9 Technical Data and Conformity

This chapter contains important technical data as well as information on the conformity.

9.1 Performance Data and Specifications

The microscope must only be operated in closed rooms. It is recommended to install the microscope in a dark room where artificial illumination, sunlight or other light sources cannot interfere with image acquisition. The microscope should not be installed near windows with direct sunlight or radiators. Compliance with the installation requirements of the microscope and the availability of the requested supplies is the responsibility of the customer and has to be provided at the time of installation. Due to continuous development, we reserve the right to change specifications without notice.

The microscope must be plugged into a properly installed power socket with protective earth contact using the supplied mains cable. The protective earth connection must not be impaired by the use of extension cables.

Info

Your ZEISS Sales & Service Partner will provide you with the detailed installation requirements.

Weight and Sizes

Main Components: Stand with ...	Length (mm)	Width (mm)	Height (mm)	Weight (kg)
Binocular tube	417	208	395	8.2
Phototube	417	208	395	8.3-9.2*
Binocular tube 25°/22 w/int cam 8 MPx	417	208	395	9.4-10.1*
Tube/phototube turned 180°	306	208	395	
Reflected-light FL iLED	444	208	471	9.6-11.4*

* Depending on the configuration

Location requirements

Installation site	Exclusively inside buildings
Altitude	Max. 2000 m above sea level
Atmospheric pressure	Min. 800 hPa

Air Conditioning and Quality

Temperature range for operation with indicated performance (24 h per day, regardless of whether the microscope is in operation or switched off)	+5 to +40 °C
Relative humidity	< 80 % at 40 °C
Atmospheric pressure	800 to 1060 hPa
Pollution degree	2

Mains connection

Nominal AC voltage	L / N 100 to 240 VAC ± 10 %
Nominal frequency	50 / 60 Hz
Main Power Plug	Local mains plug will be supplied.
Power consumption	max. 100 VA
Protection class	IP20 (IEC 60529)
IEC earth class	Class II of IEC 61140
Overvoltage Category	II

9.2 Performance Data and Specifications of the Optional Components

The customer is responsible for ensuring that the installation conditions for the microscope are met and that the required equipment is already available at the time of installation. Changes are reserved due to continuous technical developments.

Halogen lamp

Voltage	6 V
Power consumption	30 W
Adjustability of light source, continuous	from 1.5 to 6 VDC
Color temperature at 6 V	2800 K
Luminous flux	280 lm
Average service life	1000 h
Luminous area	1.5 x 3 mm

LED illumination

LED	white light
Constant, brightness-independent color temperature	5600 K
Homogeneous field illumination	20 mm (diameter)
Suitable for objectives with magnifications of	4x to 100x
Analogous brightness adjustment from	approx. 15 to 100 %

LED in LED modules

Power consumption	max. 3.5 W
Wavelength	455 / 470 nm

Stand with stage focusing

Coarse focusing drive	42 mm / rev.
Fine focusing drive	0.2 mm / rev.
Total stage lift	20 mm

Nosepiece

Objective change	manual via quadruple objective nosepiece
Objectives	infinity-corrected objective range
Mounting thread	W 0.8

Eyepieces

Tube size	30 mm
Field-of-view number	20 or 22
Magnification	10x
Suitable for spectacle wearers	Br.
Focusable	Foc.

Specimen stage

Type	Rackless
Dimensions (width x depth)	140 x 135 mm
Stage travel (X x Y)	75 x 40 or 75 x 50
Coaxial drive	optional right or left
Vernier scales	readable from the right
Specimen holder	with spring lever, left

Condenser

Abbe condenser 0.9/1.25; Fixed-Köhler	for V_{obj} 4x to 100x
Abbe condenser 0.9/1.25; Full-Köhler	for V_{obj} 4x to 100x
Turret condenser 0.9/1.25; Full-Köhler	for V_{obj} 4x to 100x

Illuminating mirror

Plane surface and spherical surface with f'	75 mm

Binocular Phototube

	Small binocular phototube 25°/20	**Big binocular phototube 25°/22**
Maximum field-of-view number	20	22
Interpupillary distance	adjustable from 48 to 75 mm	adjustable from 48 to 75 mm
Tube angle	25°	25°
Viewing height	380 to 415 mm	380 to 415 mm
Viewing port, tube factor	1x	1x
Photo/video port, tube factor	1x	1x
Photo/video port, mount	60 mm	60 mm
Invariable splitting ratio	50 vis/50 doc %	50 vis/50 doc %

9.3 Performance Data and Specifications of the Binocular Tube 25°/22 w/int Cam 8 MPx

The customer is responsible for ensuring that the installation conditions for the microscope are met and that the required equipment is already available at the time of installation. Changes are reserved due to continuous technical developments.

	Power consumption	12 VDC 1.0 A
Binocular tube 25°/22	Maximum field-of-view number (eyepiece)	22
	Acquired visual field of the camera	14.2 x 8.0 mm 16.3 mm (diagonal)
	Eyepiece (interpupillary) distance	adjustable from 48 to 75 mm
	Viewing angle	25°
	Eyepoint height	from 385 to 420 mm
	Fixed graduation	50 / 50
	Optical adaptation	0.5x
Integrated HD-CMOS camera	Sensor type	Sony CMOS color sensor with back illumination
	Number of pixels	3840 (H) x 2160 (V) = 8.3 megapixels
	Pixel size	1.85 x 1.85 µm
	Sensor size	image diagonal 8.15 mm equivalent to 1/1.9"
	Live image frame rate via HDMI (3840 x 2160)	30 fps (maximal values when using optimal hardware)
	Digitization	8 Bit/pixel
	Exposure time range	61 µs to 1 s
	Interface	USB 3.0 Type B, RJ45 (LAN), HDMI, 3 x USB 3.0 Type A, RJ12 for power supply
	Optical interface	integrated adapter
	Power interface	RJ12
	Cooling	passive cooling by aluminum cooling fins on the top plate

Info

Computer hardware, operating system, and software may decrease the frame rates. All specifications are subject to change without notice.

9.4 Applicable Standards and Regulations

Observe all general and country-specific safety regulations as well as applicable environmental protection laws and regulations.

The microscope is in compliance with the requirements of the following regulations and directives:

2011/65/EU 2015/863/EU	RoHS Directive
2017/746 EU	European Regulation on in-vitro diagnostic medical devices (IVDR)
EN 55011	Industrial, scientific and medical equipment - Radio-frequency disturbance characteristics - Limits and methods of measurement
EN 61010-1 and IEC 61010-2-101	safety requirements for electrical equipment for measurement, control, and laboratory use
EN 61326-1 EN 61326-2-6	electrical equipment for measurement, control and laboratory use – EMC requirements / IVD requirements

According to EU-regulation 2011/65/EU (RoHS) the microscope and its accessories have been classified as instrument category 9 (Monitoring and control instruments including industrial monitoring and control instruments). They also fall under 2012/19/EU (WEEE).

European and International Directives / Standards: For more information on ISO and CSA certificates or CE Declarations of Conformity, contact your ZEISS Sales & Service Partner.

ZEISS works according to a certified Environment Management System according to ISO 14001. The microscope and its components were developed, tested, and produced in accordance with the valid regulations and guidelines for environmental law of the European Union.

10 Accessories and System Expansions

Only the following accessories may be used with the microscope as their safe use has been confirmed by ZEISS. Only original parts from ZEISS may be used. Check in advance whether your microscope can be retrofitted with a system expansion or accessories.

After installation or conversion it must be carefully checked whether the microscope and its system expansions/accessories are in a safe operational state and whether unused ports are closed. For details and safety measures please refer to the associated documents.

Info

Further information on the software and its operation is available in the software's online help.

Info

For additional information and detailed descriptions, refer to further applicable documents or ask your ZEISS Sales & Service Partner.

Name	Description/Info
Objectives	The performance of the microscope objectives affects the image quality of your microscope like no other system component. Whether you work with histological samples, cell samples or entire organisms – the suitability of microscope objectives for your application depends on various factors. More detailed information on available and recommended objectives can be found at https://www.micro-shop.zeiss.com/de/de/shop/objectives or ask your ZEISS Sales & Service Partner.
Sliders	The following sliders are available: ▪ Slider H/Ph2 f/Primostar 3 ▪ Slider H/Ph1/Ph2/Ph3 ▪ Slider DF f/Primostar 3 ▪ Polarizer f/Primostar 3 ▪ Analyzer f/Primostar 3
Color filters	The following color filters are available: ▪ Set of color filters blu/grn/yel d45x1.5
Eyepieces	The following eyepieces and accessories are available: ▪ Eyepiece 10x/20 Br. foc. ▪ Eyepiece 10x/22 Br. foc. ▪ Pair of special eyecups ▪ Pinhole diaphragm D= 30 mm
Condensers	The following condensers are available: ▪ Abbe condenser 0.9/1.25 ▪ Turret condenser with BF/Ph1/Ph2/Ph3/DF
Sample holders	The following sample holders are available: ▪ Specimen holder lh ▪ Specimen holder work f/two slides lh

Name	Description/Info
Light sources	The following light sources are available: ▪ Transmitted-light illumination 3W 5600K ▪ Transmitted-light illumination 6V/30W ▪ Reflected-light FL iLED 455 nm ▪ Reflected-light FL iLED IDC 455 nm ▪ Reflected-light FL iLED 470 nm
Tubes	The following tubes are available: ▪ Binocular tube 25°/20 ▪ Binocular phototube 25°/20 (50.50) ▪ Binocular phototube 25°/22 (50.50)
Cameras	The following cameras and accessories are available: ▪ Axiocam 208 color ▪ WiFi USB adapter for Labscope connection ▪ Binocular tube HD 25°/22 w/int cam 8 MPx demo kit ▪ Camera adapter P95-C 2/3" 0.65x ▪ Camera adapter P95-C 1/2" 0.5x
Transport and storage cases	The following transport and storage cases are available: ▪ Case for transport and storage of Primo Star ▪ Transport+storage case f/Primostar 3

Revision History

Revision	Date of Issue	Introduced Modifications
3	03/2022	Implementation of revision history Adaptation to Regulation (EU) 2017/746 (IVDR)

Glossary

HDMI

High Definition Multimedia Interface

LED

Light-emitting diode

OSD

On Screen Display

TL (transmitted light)

Light used for illuminating a object, where the light is transmitted through the object.

Transmitted-light brightfield

Transmitted-light brightfield microscopy is the most common of all optical microscopic techniques, as it permits high-contrast or stained specimens (e.g. blood smears) to be viewed easily and quickly. Beside the so-called direct bundles of rays, the indirect bundles (i.e. those diffracted and scattered by specimen details) are also of major importance for providing true imaging of the object. The higher the proportion of indirect bundles of rays (aperture), the more realistic the microscopic image according to ABBE. To fully exploit the optical performance of the microscope, particularly that of the objective, the condenser, luminous-field diaphragm and aperture diaphragm should be set based on the rules of the KÖHLER illumination principle.

USB

Universal Serial Bus, an industry standard that defines cables, connectors and communications protocols for connection, communication, and power supply between computers and devices.

User

Person examining a sample under the microscope.

Wi-Fi

Wireless Fidelity

ZEISS

ZEISS is an internationally leading technology enterprise operating in the fields of optics and optoelectronics. Further information about ZEISS can be found at www.zeiss.com.

ZEISS Sales & Service Partner

The Sales & Service Partner is generally in the field for customer support in a regional area and / or a clearly defined customer group.

ZEISS service representative

Specially trained service expert, either ZEISS staff or authorized service partner of ZEISS.

Index

A

Accessories 39, 94
Air Conditioning and Quality 89
Analyzer plate 49

B

Binocular tube 25°/22 w/int cam 8 MPx 28, 30, 49
Brightfield 65, 69

C

Camera 52
Capturing
 Images 77, 78
 Videos 77, 78
Cleaning
 Water-soluble contaminations 81
Components
 Primostar 3 (Fixed-Köhler) 17
 Primostar 3 (Fixed-Köhler) with Binocular tube 25°/22 w/int cam 8 MPx 19
 Primostar 3 (Full-Köhler) 20
 Primostar 3 iLED (Fixed-Köhler) 19
Condenser 54
Connecting the Binocular tube 25°/22 w/int cam 8 MPx 57, 59, 60
Connecting the microscope 56, 57
Connecting to the PC 58
Connections
 Binocular tube 25°/22 w/int cam 8 MPx 29
 Fixed-Köhler stand 23
 Full-Köhler stand 26
 Reflected-light FL iLED 27
Contamination 88
Controls
 Binocular tube 25°/22 w/int cam 8 MPx 28
 Fixed-Köhler stand 21
 Full-Köhler stand 25
 Reflected-light FL iLED 27

D

Darkfield 71, 73
Decontamination 88
Disposal 87

E

Ethernet 60
Eyecup 55
Eyepiece 38
Eyepiece micrometer 64
Eyepiece pointer 64

F

File Options menu 36
Filter 47
Firmware 82
Fixed-Köhler stand 22
Fluorescence 75
Focusing ring of the eyepiece 65
Full-Köhler Stand 26

G

General Safety Information 10

H

Hazard
 Infection hazard 12
Hazards 11
 Prevention 11
Home menu 30

I

Image 77, 78
Image Settings menu 33
Improper use 9
Information labels 13
Installation 42
Installing
 Analyzer plate 49
 Binocular tube 25°/22 w/int cam 8 MPx 49
 Camera 52
 Condenser 54
 Filter 47
 Mirror 54
 Objective 53
 Polarizer plate 48
 Reflected-light FL iLED 44
 Software 61
 Tube 43
Intended Use 9
Interpupillary distance 63

L

Labscope 61, 78

M

Main Components Primostar 3 16
Mains connection 90
Maintenance 80
 interval 80
 schedule 80
Microscope System menu 34
Mirror 54

N

Network 59
Numerical aperture 40

O

Objective 40, 53
Operating System menu 35
Operation
 Prerequisites 62
Optional system expansions 94
 Installation 94
OSD 30, 32, 33, 34, 35, 36, 78

P

Performance data 89
Phase contrast 71, 73
Polarization 75
Polarizer plate 48
Power bank 63
Prerequisites
 Operation 62

R

Reflected-light FL iLED 27, 44
Reflected-light fluorescence 75
Regulations 93
Replacing
 Halogen lamp 82
Requirements
 for Operators 10

S

Safe Operating Condition 10
Safety 9, 80
 Devices 15
 interlocks 15
Settings menu 32
Shutdown 87
Software 61
Spare parts 10
Standards 93
Switching off 79
Switching on 62, 63

T

Tool storage 38
Training 10
Transmitted-light brightfield 65, 69
Transmitted-light darkfield 71, 73
Transmitted-light phase contrast 71, 73
Transmitted-light polarization 75
Troubleshooting 83
Tube 43

U

Update firmware 82

V

Video 77, 78
Viewing height 64

W

Warning
 labels 12
 lights 12
Warning labels 13
Weight and Sizes 89
Wi-Fi adapter 60
Wi-Fi Settings menu 36

Z

ZEISS
 Portal 8
 Service agreements 80
Zoom factor of objective 40

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About Aqua Master Tools

Established in 2018

Aqua Master Tools is the Dutch supplier of user-friendly and reliable pH, EC, and temperature meters. Our meters are used for swimming pools, aquariums, Jacuzzi's, ponds, and liquid plant nutrition.

We understand the importance of convenient, accurate and high-quality meters - and we are on a continuous quest to innovate our products. Aqua Master Tools' wide range of products are customised to our user's needs.

Satisfaction guaranteed

- Easily replaceable electrodes
- Built to last products
- Fully calibrated meters
- Unprecedented quality

Table of Contents

User manual English	6 - 10
Benutzerhandbuch Deutsch	11 - 15
Handleiding Nederlands	16 - 20
Manuel d’utilisation en français	21 - 25
Manuale utente Italiano	26 - 30
Manual del usuario en español	31 - 35

Visit **www.aquamastertools.com** or scan the QR code to watch our **HOW-TO VIDEOS**

Thank you for choosing Aqua Master Tools' Combo pen P110 Pro.

The Combo pen P110 Pro is ideal for measuring pH, EC, and temperature of liquid plant nutrition and water in swimming pools, Jacuzzi's and ponds. The user-friendly buttons can be used to effortlessly switch between the parameters. A rechargeable integrated battery is included - so even when you have no access to power, you will still be able to use the P110 Pro and do not need to purchase batteries.

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
Operation and use

Soak the meter for 12 hours in tap water before first use.

1. Remove the protective cover.
2. Gently clean the electrode with distilled water and dry it off.
3. Press the On/off button to start the meter.
4. Insert the meter upright into the solution up to the edge and stir gently. Wait a few seconds for the correct value to appear on the display.
5. Press the C/F button to switch between Celsius and Fahrenheit.
6. Press the MODE button to switch between pH and EC.
7. Rinse the electrode with water after use and press the On/off button to turn off the meter.

Our electrodes are designed for dry storage. Therefore, storage in liquid is not required for the electrode.

Calibrating the pH meter

Watch the How-to videos at Aqua Master Tools' website.

1. Press the On/off button to start the meter.
2. Press the MODE button to switch to pH on the display.
3. Gently clean the electrode with distilled water and dry it off.
4. Insert the meter in an upright position in pH buffer solution pH 7.01 (make sure this is 25 °C or 77 °F) and stir gently.
5. Wait until the value on the display is stable.
6. Press and hold the °C/°F / CAL button for 5 seconds until the value of the buffer solution is displayed. The meter will automatically recognise the buffer solution and calibrate itself. As soon as the meter switches back to the standard mode, calibration is complete.
7. Gently clean the electrode with distilled water and dry it off.
8. Repeat the process with pH 4.01 buffer solution.
9. It is important to follow the above process order, always start calibrating pH 7.00 and follow with pH 4.01.
10. Carefully clean the meter after this process.
11. Calibration is completed and the meter is ready for usage when all steps have been completed.

Calibrating the EC meter

Watch the How-to videos at Aqua Master Tools' website.

1. Press the On/off button to start the meter.
2. Press the MODE button to switch to EC on the display.
3. Insert the meter in an upright position in calibration solution (µS/cm 1413) EC 1.4 (make sure this is 25 °C or 77 °F) and stir gently.
4. Wait until the value on the display is stable.
5. Press and hold the C/F button for 5 seconds. Then release the button. The value EC 1.4 will now blink on the display. As soon as the meter switches back to the standard mode, calibration is complete.
6. Rinse the meter with clean water.
7. The meter is now ready for use.

When is calibration required?

- If the electrode has been replaced.
- If the meter was last calibrated over a month ago.
- If the meter is used on a regular basis.
- If high precision is needed.

Replacing the electrode

Watch the How-to videos at Aqua Master Tools' website.

1. Properly unscrew the sealing ring in the lower part of the meter.
2. Remove the defected/old electrode by gently pulling it.
3. Carefully place the new electrode with the recess in the correct position so that the connectors properly connect and do not damage the meter or electrode.
4. Check that the rubber seal is properly placed.
5. Tightly tighten the sealing ring.
6. Calibrate the meter as described in this manual.
7. The meter is now ready for use.

Specifications

Range	pH: 0.1 ~ 14.0 pH EC: 0.0 ~ 19.9 µS/cm Temperature: 0°C ~ 50°C (32°F ~ 122°F)
Resolution	pH: 0.1 pH EC: 0.1 µS/cm Temperature: 0,1 °C (1,0 °F)
Accuracy	pH: ± 0.1 pH EC: ± 2% F.S Temp: ±0.5°C Temp. compensation: 0°C - 50°C (32°F - 122°F)
Calibration	pH: 4.0 - 7.0 pH EC: 1.4 EC (factory standard calibrated 1413 µS/cm)
Features	Battery indicator & automatic temp. compensation
Power supply	3 × 1.5 Volt AG-13 button cell battery (included)
Auto. power-off	After 8 minutes
Calibration indicator	With correct calibration
Waterproof	IP67
Dimensions	220 x 40 x 40 mm
Weight	102 grams

Replacing the batteries


If the battery indicator at the top right of the display is activated, or if the display dims, please replace the batteries. To do so, unscrew the top and replace the 3 x 1.5 Volt AG-13 button cell batteries.

Warning or issues

- If your meter is not working properly, watch the troubleshooting videos on the product page of our website aquamastertools.com
- Always make sure that there are no air bubbles present in the solution during calibration. This could affect the measurements.
- If the EC value does not decrease to 0, insert the electrode into ethyl alcohol and shake gently for 5 to 10 seconds.
- Make sure the calibration solution is 25 °C when calibrating the meter. Failure to do so may affect the result of the calibration.
- Never store the electrode in distilled water.
- **DO NOT** immerse the electrode in oil, proteins or suspended solids that leave a layer on the glass bulb.
- **DO NOT** soak in RO (Reverse Osmosis), distilled or de-ionised water. Pure water alters the chemistry of the reference liquid so that it no longer works.
- **DO NOT** place a cold electrode in hot liquid (or vice versa). Sudden changes in temperature may break the glass and permanently damage the pen. Broken glass is not covered by the warranty.
- Make sure not to bump the glass electrode into anything, because this might break the glass bulb on the outside of the internal glass tube. Broken glass is not covered by the warranty.
- Avoid placing the meter near ballasts or any other electrical appliances. Force fields may affect the readings.
- If during calibration the pH/EC value slowly keeps rising or dropping and does not stabilise, please replace the batteries and calibrate the meter. Bad batteries or low power levels may influence the measurement results.

Warranty

This meter is warranted from all defects in material and manufacturing for a period of one year from the date of purchase. The electrode is warranted from all defects in material and manufacturing for a period of six months from the date of purchase. If during this period any parts need reparation or replacement and the damage is not due to incorrect operation by the user, please send the parts to the dealer or to us and the repair will be free of charge. Please refer to the exclusion warning.

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
READ MANUAL

Keep these instructions with the appliance.
You should read this user manual carefully before using the appliance.

Powered by CamScanner

Dear Customer,

Thank you for purchasing this appliance. Please take a few minutes before starting operation of the appliance and read the following operating instructions.

1. Intended use

1.1 This machine is designed for food processing in restaurants, canteens and similar catering establishments.

1.2 The machine is not intended to process for meet which has bone, tough skin, shells or fish scales. Neither use this machine to cutter vegetable or fruit with nut or tough skin. Otherwise, it might lead to damage of the machine, property or personal injury.

1.3 The machine is intended for commercial and indoor use by trained and skilled personnel familiar with these instructions.

The assortment is made up of 4 different sizes:
WB280AF/WB280AV/WB450AF/WB450AV/WB650AF/WB650AV/WB850AF/WB850AV.
The machine body is built in PA66+GF30% and stainless steel for easy cleaning and maintenance. It has safety features compliant with the applicable safety Directives.

2 Safety regulations ⚠

2.1 Safety regulations

2.1.1 This appliance is intended for commercial use only and must not be used for household use.

2.1.2 The appliance must only be used for the purpose for which it was intended and designed. The manufacturer is not liable for any damage caused by incorrect operation and improper use.

2.1.3 Keep the appliance and electrical plug away from water and any other liquids. In the event that the appliance should fall into water, immediately remove plug from the socket and do not use until the appliance has been checked by a certified technician. Failure to follow these instructions could cause a risk to lives.

2.1.4 Never attempt to open the casing of the appliance yourself.

2.1.5 Do not insert any objects in the casing of the appliance.

1.1.6 Do not touch the plug with wet or damp hands.

2.1.7 Danger of electric shock! Do not attempt to repair the appliance yourself. In case of malfunctions, repairs are to be conducted by qualified personnel only.

2.1.8 Never use a damaged appliance! Disconnect the appliance from the electrical outlet and contact the retailer if it is damaged.

2.1.9 Do not immerse the electrical parts of the appliance in water or other liquids. Never hold the appliance under running water.

2.1.10 Regularly check the power plug and cord for any damage. If the power plug or power cord is damaged, it must be replaced by a service agent or similarly qualified persons in order to avoid danger or injury.

2.1.11 Make sure the cord does not come in contact with sharp or hot objects and keep

Powered by CamScanner

it away from open fire. To pull the plug out of the socket, always pull on the plug and not on the cord.

2.1.12 Ensure that the cord (or extension cord) is positioned so that it will not cause a trip hazard.

2.1.13 Always keep an eye on the appliance when in use.

2.1.14 As long as the plug is in the socket the appliance is connected to the power source.

2.1.15 Turn off the appliance before pulling the plug out of the socket.

2.1.16 Never carry the appliance by the cord.

2.1.17 Do not use any extra devices that are not supplied along with the appliance.

2.1.18 Only connect the appliance to an electrical outlet with the voltage and frequency mentioned on the appliance label.

2.1.19 Connect the power plug to an easily accessible electrical outlet so that in case of emergency the appliance can be unplugged immediately. To completely switch off the appliance pull the power plug out of the electrical outlet.

2.1.20 Always turn the appliance off before disconnecting the plug.

2.1.21 Never use accessories other than those recommended by the manufacturer. Failure to do so could pose a safety risk to the user and could damage the appliance. Only use original parts and accessories.

2.1.22 This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge.

2.1.23 This appliance must not be used by children under any circumstances.

2.1.24 Keep the appliance and its cord out of reach of children.

2.1.25 Always disconnect the appliance from the mains if it is left unattended or is not in use, and before assembly, disassembly or cleaning.

2.1.26 Never leave the appliance unattended during use.

2.2 Special Safety Regulations

2.2.1 The appliance is intended for commercial and indoor use by trained and skilled personnel familiar with these instructions.

2.2.2 This machine is designed for food processing in restaurants, canteens and similar catering establishments.

2.2.3 Any other use might lead to damage to the appliance or personal injury.

2.2.4 Do not place the appliance on a heating object (gasoline, electric, charcoal cooker, etc.) Keep the appliance away from any hot surfaces and open flames. Always operate the appliance on a level, stable, clean, heat-resistant and dry surface.

2.2.5 **Danger** of injury! Care is needed when operating or cleaning. Wear protective gloves (not supplied) if necessary.

2.2.6 **Caution!** Securely route the power cord if necessary, in order to prevent unintentional tripping over and falling.


2.2.7 Do not use the appliance before it is properly assembled.

Powered by CamScanner

3 PRODUCT DESCRIPTION


3.1 Parts of description

WB280AF/WB280AV/WB450AF/WB450AV/WB650AF/WB650AV/WB850AF/WB850AV

Powered by CamScanner

3.2 MAIN TECHNICAL PARAMETERS

Model	Voltage Supply	Power	Speed	Net Weight	A (Length)	B (Diameter)	Noise Level
WB280AF	115V~60Hz	280W	20000RPM	1.45KGS	367mm	Ø 88 mm	75~80dB
WB280AV		280W	6000~20000RPM	1.45KGS	367 mm	Ø 88 mm	
WB450AF		450W	16000RPM	2.55KGS	437 mm	Ø 100 mm	80~85dB
WB450AV		450W	4000~16000RPM	2.55KGS	437 mm	Ø 100 mm	
WB650AF		650W	18000RPM	3.10KGS	440 mm	Ø 120 mm	85~90dB
WB650AV		650W	8000~18000RPM	3.10KGS	440 mm	Ø 120 mm	
WB850AF		850W	18000RPM	3.35KGS	440 mm	Ø 120 mm	
WB850AV		850W	8000~18000RPM	3.35KGS	440 mm	Ø 120 mm	

NOTE:
Above machine don't include any following attachments. Below chart you can find available attachments for your choice.

Powered by CamScanner

Model	Net Weight	C Length	D Diameter	Remark
WIK185-T	0.81KGS	185 mm	/	USE FOR WB280AV ONLY
160-T	0.55KGS	240 mm	Ø 28mm	USE FOR WB280AF / WB280AV
WIK250	0.9KGS	250 mm	/	USE FOR WB450AV/WB650AV/WB850AV
250-N	0.93KGS	217 mm	Ø 35 mm	USE FOR WB450AF / WB450AV
300-N	1.05KGS	267 mm	Ø 35 mm	
350-N	1.15KGS	317 mm	Ø 35 mm	
400-N	1.26KGS	367 mm	Ø 35 mm	USE FOR WB450AF / WB450AV / WB650AF / WB650AV / WB850AF / WB850AV
450-N	1.37KGS	417 mm	Ø 35 mm	
500-N	1.60KGS	467 mm	Ø 35 mm	USE FOR WB650AF / WB650AV / WB850AF / WB850AV /
550-N	1.72KGS	517 mm	Ø 35 mm	

4 INTRODUCTION

4.1 Preparation for first use

4.1.1 Check to make sure the appliance is undamaged. In case of any damage, contact your supplier immediately and do NOT use the appliance.

4.1.2 Remove all the packing material and protection film (if applicable). Keep the packaging for future storage.

4.1.3 Clean the appliance, see: **4.6 Cleaning**.

4.1.4 Place the food bowl or other container on a level and steady surface, unless mentioned otherwise. It is forbidden to put the machine above different objects from the working ground used in the food field of height included between 900 – 1100mm from the trampling level.

4.1.5 Position the appliance in such a way that the plug is accessible at all times.

4.1.6 Make sure the power outlet ratings correspond to the specifications printed on the rating plate.

4.1.7 Make sure that the machine is connected to a grounding system if it's needed.

4.1.8 Even if safety systems are installed on the machine, prevent approaching moving parts with hands and do not touch the machine with wet or damp hands.

4.1.9 Warning! When the whisk is used, make sure to start appliance with low speed.

4.1.10 Warning! It is forbidden the use of the whisk in the machines with fixed speed.

Powered by CamScanner

4.1.11 Disconnect the machine's socket from the electric before carrying out any cleaning or maintenance operations.

4.1.12 Do not clean with high-pressure water jet.

4.1.13 It is strongly recommended that the operator contact the after sales center if there is any problem or malfunction occurred. Please do not dismantle the machine or in place of unauthorized parts.

Warning: This machine must be properly immersed and used in water. To avoid overheating or damage, DO NOT allow it to run in air for more than 10 seconds.

4.2 The manufacturer will be not responsible for following cases:

4.2.1 If the machine is tampered with or safety functions are installed by non-authorized staff.

4.2.2 If the components are replaced by unauthorized parts.

4.2.3 If the instruction in this manual is not carefully followed.

4.2.4 If the machine surfaces are treated with unsuitable substance or solvent.

4.3 Safety devices

4.3.1 **Safety Switch**: The safety switch is equipped for IB450AF/IB450AV/IB650AF/IB650AV/IB850AF/IB850AV, to make sure both hands are away from the blade while starting machine.

4.3.2 **Overcurrent protection switch** (Circuit breaker): The switch is designed to protect the electrical circuit from damage caused by excess current as a result of an overload or short circuit. The switch interrupts current flow after a fault is detected. **Note:** To reset the machine, wait a few minutes till the motor cool down, then restart the machine.


4.4 INSTALLATION

4.4.1 Packing and unpacking

4.4.1.1 Remove all packaging materials. Keep the packaging for future storage.

4.4.1.2 After unpacking all accessories, read the manual carefully before installation.

4.4.1.3 Make sure the delivery contents are complete and free of any damage. If you find that parts are missing or show damage, do not use the appliance but contact your dealer.

4.4.2 STARTING AND STOPPING THE BLENDER

4.4.2.1 Check that the on-off button (No.3) is in the released position.

4.4.2.2 Connect the appliance to the power supply.

4.4.2.3 Insert the tube or whisk into the appliance.

4.4.2.4 Press and hold down the safety button (No.1) if

Powered by CamScanner

the appliance has.

4.4.2.5 Press the on-off button (No.3), and the blender will start up. Now, release safety button (No.1).

4.4.2.6 To stop the blender, release the on-off button (No.3), then blender stop.

4.4.3 CONTINUOUS OPERATION:

Follow above steps **4.4.2.1** to **4.4.2.6** above, then:

4.4.3.1 While the on-off button (No.3) is in the depressed position, press the locking button (No.4).

4.4.3.2 Release the on-off button (No.3), and the appliance will function continuously.

To stop the appliance, press the on-off button (No.3) and release it.

4.4.4 VARIABLE-SPEED OPERATION

Follow above steps **4.4.2.1** to **4.4.2.6**, then:

Change the speed of the motor by turning the variable speed button (No.2) towards the maximum or minimum as required. With the self-regulating speed system, once you have selected a speed, it will remain constant, even if the consistency of your mixture changes.

Warning: <u>It is strongly advisable to start a low speed when using the variable-speed blender!</u>

4.4.5 Working position


4.4.5.1 Blender function: For better control, we recommend that you hold the bender by the handle and the bottom of the motor unit. Furthermore, we recommend that you tilt the bender at a slight angle to prevent the bell from touching the bottom of the container. Always make sure that the bell is sufficiently immersed to avoid splashing and that the air vents of the motor unit does not come into contact with any liquid. For optimum efficiency, two thirds of the tube should be immersed in preparation.

4.4.5.2 Whisk function: You can also hold the appliance by the handle with one hand, leaving the other hand free to hold the bowl, if necessary. While you process the preparation, we recommend that you move the whisk around in the bowl, in order to ensure that the mixture is completely homogeneous. We strongly advise you to keep the whisk from touching the sides of the bowl. For Maximum efficiency, at least one-fifth of the whisk length should be immersed. But never immerse the whisk holder into the preparation.

Warning: <u>It is strongly advisable to start a low speed when using the whisk!</u>


Powered by CamScanner

4.5 ASSEMBLING / DISMANTLING: The assembly / disassembly operations of every device of the appliance are performed with the appliance stop and with the plug disconnected from the socket.

Assembly

1. Align and insert the coupling (5) of the tube (3) or whisk (4) in the clutch (2) of the machine body (1). Be sure the concave side of coupling must catch into the convex of the clutch.
2. Screw the nut (6) up to fasten them firmly to the machine body (1).

The disassembly activities perform following in a backgrounds way the assembly phases previously described.

4.6 CLEANING

As a precaution, always unplug your appliance before cleaning it.

4.6.1 After using the immersion blender and whisk in a hot preparation, cool the bottom of the tube in cold water before dismantling the bell from the tube. After completion of preparation, clean immediately to prevent food from sticking to surfaces.

4.6.2 Never put the motor unit and its handle in water, clean them with a damp cloth or a sponge.

4.6.3 Blender tube and bell cleaning: Place the tube and bell into a suitable container with water and switch on for a few seconds to clean the bell and tube.

WARNING: Always dry the blades thoroughly after cleaning to avoid spotting.

4.6.4 Never immerse the whisk holder in water, clean it with a slightly damp cloth or sponge.

4.6.5 It is vital to clean it thoroughly between each job using a detergent or disinfectant, but never use pure bleach. For plastic parts, do not use detergents that are too alkaline.

Powered by CamScanner

Warning 1: The blender tubes and whisks are not suitable for the dishwasher!
Warning 2: Blender tubes and whisk cannot be immersed thoroughly in water!

4.7 Storage

4.7.1 Before reassembling, make sure that all parts are dry. Before storing, remove all residues from the appliance and attachments, disconnect it from the supply and clean the appliance and its accessories thoroughly.

4.7.2 Store the appliance in a dry, protected environment at room temperature and out of the reach of children.

4.7.3 We recommend to put the appliance and attachments onto "wall hanger" for dry and safe storage.

4.8 Troubleshooting

4.8.1 If it is due to a power cut if the appliance had been unplugged, release the on-off button (No.3), check the power supply and then restart the appliance.

4.8.2 If it is due to overheating of the machine, release the on-off button (No.3), and then unplug the appliance, wait a few minutes for the motor to cool down and the thermal safety mechanism to be reset, and then restart the appliance.

4.8.3 If you cannot locate the cause of the problem: release the on-off button, unplug the appliance, check below:


1. the plug.
2. that the blades are free to rotate in the bell, be careful the blade is very sharp.
3. that state of the power cord.
4. the drive shaft is free to rotate. To do this, remove the bell and test the rotation of the end of drive shaft manually.

If a malfunction occurs which is not on our list of possible faults, consult your retailer or return the appliance to the after-sales services.


Powered by CamScanner

5 Circuit diagram


WB280AF


WB280AV

WB450AF, WB650AF, WB850AF


WB450AV, WB650AV, WB850AV

Powered by CamScanner

6 Exploded view and spare parts lists


WB280AF/WB280AV

No.	Part Name
1-1	Bottom Housing 1/2
1-2	Bottom Housing 2/2
1-3	Upper Housing 1/2
1-4	Upper Housing 2/2
2	Motor
3	Air Outlet Spacer
4	Torsion Spring
5	PCB Panel
6	Switch Osculant Staff
7	Switch Sheath
8	Strain Relief
9	Power Cord
10	Pressing Board
11	Screw Plungers
12	Variable Speed Button
13	Silicone Decoration Ring (optional)
14	Air Inlet Spacer
15	Spline Sleeve
16	Aluminium Head
17	Spline
18	Spline Seal Ring
19	Connect Sleeve
20	Fastening Sleeve
21	608 Bearing
22	Output Axis
23	304 S/S Tube
24	Oil Sealing
25	Connect Sleeve Seal Ring
26	Bearing Sleeve
27	Knife Protector
28	Ceramic Ring Sleeve
29	Ceramic Ring
30	Graphite Ring Spring
31	Graphite Ring
32	Convex Ring
33	Smooth Blade

Powered by CamScanner

WB450AF,WB450AV


Spare parts list

No.	Part Name	No.	Part Name	No.	Part Name
1	Power Cord	10	On-off Switch	17	Air Outlet Spacer
2	Strain Relief	11	Switch Osculant Shaft	18	Wire Guard Ring
3	Screw	12	Rubber Interlocking Protector	19	PCB Panel
4	Screw	13	Pressing Board	20	Bottom Housing
5	Screw plunger	14	Safety Button	21	Motor
6	Silicone Decoration Ring	15	Waterproof Cap 1	22	Motor Spacer
7	Upper Housing	16	Variable Button	23	Air Inlet Spacer
8	Waterproof Cap 2	9	Lock Button	24	Spline Sleeve

Powered by CamScanner

WB650AF,WB650AV,WB850AF,WB850AV


Spare parts list

No.	Part Name	No.	Part Name	No.	Part Name
1	Power Cord	10	Safety Button	19	Air Inlet Spacer
2	Screw	11	Waterproof Cap 1	20	Rubber Ring
3	Upper Housing	12	Pressing Board	21	Spline Sleeve
4	Waterproof Cap 2	13	Variable Button	22	Aluminium Head
5	Screw	14	PCB Panel	23	Locking Sleeve
6	Screw Plunger	15	Rubber Interlocking Protector	24	Lock Button
7	On-off Switch	16	Air Outlet Spacer	25	Silicone Decoration Ring
8	Strain Relief	17	Bottom Housing		
9	Switch Osculant Shaft	18	Motor		

Powered by CamScanner

250-N / 300-N / 350-N / 400-N / 450-N / 500-N / 550-N


Spare parts list

No.	Name	No.	Name	No.	Name
1	Smooth blade	6	Bearing	11	Fastening sleeve
2	Oil Seal	7	Blade bell	12	Connect sleeve
3	Bearing fixing seat	8	Seal ring	13	Spline sleeve
4	Seal ring	9	304 S/S tube	4	Seal ring
5	Oil Seal	10	Output axis	5	Oil Seal

Powered by CamScanner

WIK185-T / WIK250

No.	Part Name
1	Silicone rubber ring
2	Spline
3	Φ8 Snap ring
4	Bearing washer
5	Bearing
6	Connect sleeve
7	Gear shaft
8	Screw plunger
9	Housing screw
10	Up seal ring
11	Upper cover
12	Needle bearing
13	Screw
14	Screw
15	Baffle plate
16	Φ16 Snap ring
17	Gear B
18	Cylindrical pin
19	Straight gear
20	Output axis
21	Key
22	Down seal ring
23	Lower cover
24	Oil sealing
25	set screw
26	Cylindrical pin
27	Fixed board
28	Steel ring
29	Steel bar
30	Fastening sleeve
31	Shaft sleeve

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

Z2I-5168-02

SOIL LOAD PENETROMETER
&
GROUND BEARING CAPACITY DETECTOR
FM-204DJ

Installation drawing

I. Product Overview

Universal type penetrometer and filling compactness detector, detector or the foundation bearing capacity is used to detect dry density of backfill soil and uniformity, can save to replace sand ring cutting irrigation irrigation water in the construction site inspect the filled soil bedding face, a full range of quality and reduce the workload of geotechnical tests, shorten the construction period, it is suitable for layering construction units, roller compaction test, It is also suitable for quality inspection and supervision company sampling inspection of any layer, but also suitable for colleges and universities to teach and learn to do geotechnical experiments, it has no radioactive nuclear pollution, low price, small volume, easy to carry on site.

The product has been included in the international standard in 1982, and has been widely used in engineering construction in the United States, Japan, Russia and other countries. It has also been included in the foundation design code, foundation treatment technology code and foundation construction quality acceptance code. At present, it has been widely used in civil and industrial construction, highway, railway, road and bridge, municipal drainage, pipeline line, water conservancy, dam construction, dike, irrigation canal, reservoir construction, airport runway engineering. In the past, it took at least 8 hours to take a set of ring knives for three times of combustion, and the results were far behind the progress of the project, which had no guiding significance for the on-site construction at that time. When the inspection indicators are found to be unqualified or problematic, the site filling has already covered two or three layers of soil. The maximum 15-20 minutes for a group of tests with the penetration method is efficient, fast and productive. Because of its detection depth of 30 cm, but also extended to 1.5 m, so it can timely find soft weakness and hidden trouble at the bottom, to ensure the quality of construction.


Universal type penetro meter is based on the geotechnical experiment: the penetration of soil and the dry density of soil, the relationship between the relationship between the penetration and water content of soil and the relationship between soil water volume and compactness, through a lot of penetration and ring cutting parallel contrast test And establish the regression equation, the rate of penetration resistance between the dry density and the related data sheet. By measuring the penetration resistance and looking up the table, the dry density can be quickly obtained, and then the compaction coefficient can be quickly calculated. Two, main technical parameters

2. main technical parameters


Measuring range: 0kg-100kg (0N-1000N)
Test depth: 0cm-30cm
Material tested: all kinds of backfill soil, clay, sand and mixed soil after rolling
Application environment: highway, railway, reservoir, dike, dam and civil building, etc
Test results: filling resistance, relative density, compaction coefficient and foundation bearing capacity
Power: 2 AA batteries(6 V dc matching)
Weight (excluding battery): instrument 415 g, handle 295 g, probe rod 380 g.
Test accuracy:≤0.25%
Shutdown mode: Close automatically or manually shut down 10 minutes

3. instrument assembly

(I) Battery installation: 2 AA batteries are required for the use of the instrument. Open the back cover, the box is marked with "+" and "-" poles, and assemble the battery according to the polarity marked.

If not used for a long time, you need to take out the battery to avoid damage to the instrument.
(2) mounting the probe rod: assemble the probe rod to the lower part of the instrument. There are three kinds of probe rod in the instrument box, ϕ7mm*240mm, ϕ10*120mm and ϕ10mm*320mm. The two ϕ 10 rods share one probe.

When ϕ 7 probe rod is connected to the main engine, a corresponding joint is required.

Connect the probe rod to the main engine before testing. Users according to different needs, choose different probe rod. Check between the probe and the probe before use, the screw thread buckle between the probe and the host shaft must be tightened, not loose, lest the probe twisted deformation or broken.
(3) Adjust the protection nut: before measuring, screw the two protection nuts at the bottom of the host to the middle, do not touch the protection sleeve. When the instrument is not in use, keep the protective nut close to the protective sheath to protect the force measuring elements

Similarly, when P is any number =725N, plug in ③ to get,

$$\lambda_{\text{Any number}} = \frac{725}{730} = 0.993$$

Conclusion: Firstly, the maximum P value under the maximum Rd state is calculated, and then an arbitrary value is obtained as long as the construction site is selected at any point. By substituting equation ③, λ arbitrary of the corresponding compaction system at this point can be obtained.

Appendix 4-compaction coefficient calculation method

1. According to the known Rd_{MAX} and Rd_{Design}, the pressure coefficient λ design is calculated

①

$$\lambda_{\text{Any number}}\frac{\mathrm{Rd}_{\text{Design}}}{\mathrm{Rd}_{\text{MAX}}}=\frac{1.58}{1.63}=0.97$$

2. Calculate the maximum P in Rdmax state:

②

Equation: $$\frac{P_{\text{MAX}}}{P_{\text{design}}}=\frac{\lambda_{\text{MAX}}}{\lambda_{\text{Design}}}$$

Substitute the specific parameters into:

$$\frac{P_{\text{MAX}}}{709}=\frac{1}{0.97}$$

Equation: $$P_{\text{MAX}}=\frac{1}{0.97}*709=730(\mathrm{N})$$

3.Three, arbitrary calculation of λ :

③


Equation:$$\lambda_{\text{Any number}}\frac{\mathrm{P}_{\text{any number}}}{\mathrm{P}_{\text{MAX}}}$$


When the measured P is arbitrary =700 (N). Substitute in ③ to compute an arbitrary value of λ,

$$\lambda_{\text{Any number}}\frac{700}{730}=0.959$$

Similarly, when P is equal to 715 (N). ① Substitute in ③ and calculate an arbitrary value of λ,

$$\lambda_{\text{Any number}}\frac{715}{730}=0.979$$

measurement **protect**

4. Instrument keys

(1) Instrument switch: on the right side of the screen, used to switch the instrument.

(two) [zero] key: after installing the probe rod, before measuring, the instrument is suspended, because the probe rod is installed, the instrument will have an initial value, press [zero] key to adjust the initial value to zero, equivalent to the meaning of peeling. Zeroing only needs to boot zeroing once.

(3) [follow up] key: press the [follow up] key, the instrument state changes to follow up, the instrument according to the force of the measuring rod, the real-time force value display.

(IV) [Peak] key: After pressing [peak] key, the instrument state changes to the maximum value, and the value displayed by the instrument is the maximum force value after pressing [peak] key. After pressing the [peak] key, the measurement is completed, you need to press the [follow-up] key again, release the record of the maximum value, the next measurement

Press the [peak] key to record the maximum value of this measurement.

(five) [Storage] key: if you feel that the screen is dark and can not see clearly, you can press [storage] key to light the screen backlight, and then press to turn off the backlight.

5. penetration resistance P measurement

(1) Install the probe rod, screw the two protection nuts at the lower part of the main engine to the middle, do not get the protection sleeve.

(two) open the switch on the right of the screen, the screen displays an initial value.

(three) Suspend the instrument and press the [Clear] key to adjust the display value to 0. Reset only need to boot reset once can be.

(4) Drill a hole of suitable depth at the measuring point. For example, if you need to measure penetration resistance at a depth of 30cm, drill a hole at a depth of 20cm. The purpose of drilling is to save the operator the effort to drill the probe to a specified depth.

(5) Insert the probe of the instrument into the depth to be measured. The penetration speed should be uniform, the downward force should be smooth and continuous, the general penetration into 10 cm, need 10-15 seconds, not too fast, do not use impact, do not slow down halfway, to complete the operation in one breath. One person penetration, one person to check the instrument, when nearly penetration to the specified depth, press [peak

Value] key. When penetration reaches the specified depth, the value shown by the instrument is multiplied by 10 to obtain penetration resistance. After recording the data, press the [follow up] key to release the maximum value, so that press the [peak value] key to record the maximum value in the next penetration.

6. Obtaining Rd of relative dry density

(1) Penetration resistance is obtained according to the measurement of penetration resistance P.

(2) According to the soil quality, check the corresponding dry density data table (see Appendix 1- Dry density data for details)

According to the penetration resistance, the relative dry density Rd was obtained. The dry density data table is only of reference value for some special soils. The dry density data table was drawn by a large number of parallel contrast tests of penetration and ring knife method. If you have any questions or concerns, you can draw the dry density data table by comparing with the ring knife method according to the actual soil quality. (See Appendix 2- Calculation method of relative dry density values and Appendix 3- Calibration Method)

7. calculation of compacting coefficient λ

(I) According to the previous step, find the relative dry density dry density Rd.

(2) The relative dry density value divided by the current maximum dry density value of soil, is the required compaction coefficient. (See Appendix 4- Calculation method of compaction coefficient)

8. calculation of foundation bearing capacity

(1) According to the relative dry density value Rd, the relative dry density value Rd was found. (2) Rd multiplied by 100, the value obtained is the foundation bearing capacity, the unit is kPa. If you want to do it in kilograms per square meter, you just multiply it by 100.

9. Storage after measurement

(1) Tighten the protection nut close to the protective sleeve to protect the force measuring original inside the instrument.

(2) The probe rod from the instrument up and down, clean and save.

(3) Remove the battery from the instrument to avoid battery leakage and damage to the instrument.

10. Fault analysis

(I) After the battery is installed and the switch is turned on, the screen of the instrument does not display:

First, make sure the battery is installed correctly. Second, replace the battery and try it. Third, please contact the dealer or manufacturer to solve the problem.

(2) Unstable multiple measurement data:

The data is unstable for many times, which may be caused by insufficient battery power supply. Please try a new battery.

When the relevant data table of the local area is established, and there is no time, no condition to calibrate or do parallel comparison test, this simple judgment method can be used to stipulate: roll four times, then use this instrument for on-site inspection, and record the reading P value of the instrument is higher than or equal to 700N, that is, the dry density is qualified: for example, (750N, 810N) are qualified. If less than 700N is not qualified, must be rolled again - again test. When the water content changes, such as the filler composition changes or changing the material field, a set of ring knife recalibration and its corresponding standard penetration P value should be made. In the above calibration tests, at least three groups of the same soil must be made before curves can be. drawn and regular data tables can be summarized, and regression equations can be established.

Attached table: Times of rolling: four times. Soil classification: cohesive soil, 10cm probe penetration.

Diameter: 10mm. Number of ring knives: 5 groups.

Experimental results loci	one	two	three	four	five	six	seven	eight	average
Ring knife test results Hd	165								
Penetration test	680	695	700	695	715	705	800	390	701
water content WO	18%								

$$Rd_{\text{Any number}} = \frac{720}{709} \times 1.58 = 1.60$$

Measured $\mathbf{P}_{\text{Any number}}$=725N Substitute into equation ② to obtain:

$$Rd_{\text{Any number}} = \frac{725}{709} \times 1.58 = 1.61$$

Conclusion: P design can be calculated as long as Rd design and P design are known. In the case of Rd design and P design, the corresponding dry capacity value can be obtained as long as P is measured at any point on site and substituted into the above formula. Rd is arbitrary, and so on.

Appendix 3- Measurement method

(1) If the conditions permit, the unit shall conduct parallel contrast test by itself with ring knife and penetrometer. The method is to peel off 1/3 of the surface layer and the remaining 1/3 of the bottom layer after rolling three to five times, and test the 1/3 to 2/3 of the middle layer.

(2) Take 3 to 5 ring knives from the position 2/3 deep in the detection layer. When the dry density reaches the qualified requirements, write down the dry density and water content and fill in the attached table as the standard value of ring knives.

(3) At the same time, draw a circle with a radius of about 30cm as the center of the ring knife. Take 7 to 8 points on the circumference to conduct penetration test, and the penetration depth should be consistent with the depth of the ring knife at 1/3 of the cushion* Select between 2/3, eliminate the points with large deviation, and take the average value of the other points as the standard penetration value corresponding to the ring knife.

(4) When the standard penetration resistance value is determined, the penetrometer can be used instead of the ring knife method to enter the construction site to inspect the filling quality in a wide area and in all directions. When the water content is controlled, when the test result is greater than or equal to the standard penetration resistance value, it is considered that the rolling is closed; when the test result is lower than the standard penetration resistance value, it is considered that the dry density is low, and the rolling is continued if it is unqualified. Fill in the results of the parallel tests in the attached table below:

Note: soil classification (cohesive soil, sandy soil, mixed soil, etc.), bedding thickness, rolling times, penetration depth and probe diameter should be indicated in the attached table. Design requirements dry density and maximum dry density. Examples of songhua River embankment construction are listed below:

Dam embankment filling soil composition: just clay. Rolling times: four times. Take a group of ring knives as: 5 soil samples. Experimental results: Rd= 1.65, Wo=18%, the compactness was confirmed to be qualified. At the same time, let the 30cm penetrometer enter the field for the test. Dig out 10cm of the surface and 10cm of the bottom at the same position of the ring cutter, select 10-15 cm of the middle layer for the penetration depth, read the penetration resistance value, and then take seven to eight points of penetration around the ring cutter. The test results of the eight points are:

①P_{20}=680N ②P_{20}=700N ③P_{20}=710N ④P_{20}=695N

11.Question and answer for user operation

(I) Must the depth of the test point taken from 7 to 8 o 'clock of the penometer be consistent?

In the same group of tests on the same measured surface, the test depth must be consistent, for example, when h=10 cm, all measuring points must penetrate 10 cm to read the results.

(2) when is the reading accurate?

When the probe is inserted to a standard depth, the instantaneous peak, the maximum, should be read. The method of reading is: the operator should always keep constant pressure, the recorder should pay attention to observe the insertion depth do not interfere with the operator, when the standard depth boundary is reached, when the soil surface just pressed into the boundary moment, immediately press the button to read the maximum value. If you can't grasp the keys, you can also do not press the keys. Direct observation of the instantaneous peak when the line is pressed. The method is to place your finger on a predetermined standard depth limit, keep your eye on the reading display, and read the instantaneous peak as soon as the soil touches your finger.

(3) What should be paid attention to when detecting the bottom compactness?

It is necessary to calculate the thickness of the interface of each two cushion layers, and avoid detection on the interface, because the high density here is not representative. For example, the first layer is 30cm thick, and the soil is drilled 10cm deep clockwise with a wooden drill with a diameter of 30cm, and then the soil is pulled out counterclockwise. The section from 10cm to 20cm is the best detection layer. When detecting the second layer, attention should be paid to the bottom of the first layer at 20cm-30cm and the top layer at 30-40cm. Also can not test, but also to peel off, with wood drill clockwise to 40 cm and then counterclockwise to pull out the soil, from 40 cm to 50 cm is the second layer of effective detection layer, remember to peel off 1/3 of each cushion layer, leaving 1/3 of the bottom layer, the measured layer in the middle 2/3, and then check a layer, and so on.

12. matters needing attention

(1) To determine the bearing capacity of the foundation at the same depth and at different points on the same working surface, take a point at an interval of 0.5 meters, and take seven to eight measuring points in total, and the penetration depth of each point should be consistent.

(2) to measure the bearing capacity of the different depth, can match with a diameter of 2.5 cm~ 3.0 cm woodworking drill, clockwise, counterclockwise to soil exits, drill a layer of a layer, drilled down another layer, the measured through part of the disturbed soil come back, then a layer of operation in order to detect time, anyhow to grasp a principle, be measured surface and by measuring point must be no disturbance of the earth, Disturb a layer, drill off - - layer, drill off a layer of disturbed soil out and then check down, it is best to keep the same depth of each layer of detection, in order to accurately evaluate the uniformity and change of soil. For example: to select a point in a site, in order to understand the ground endurance of each layer between 0 and 1.5 meters in this hole, take a wooden drill with a diameter of 3cm first, drill into 10cm depth and then pull out the topsoil in reverse. Then the host connection probe is inserted into the hole 10cm, the data is recorded and then drilled down 10cm, the disturbed soil part of the soil out of the test, the instrument is inserted down 10cm. The insertion depth can be calculated by observing the rod scale from the ground.

(4) Remove the probe rod after each use, wipe the rod, butter it, and wrap it in plastic paper to prevent rust. The threaded part of the rear cover of the battery case should also be coated with a little butter to prevent rust.

(5) The interval of a period of time is not used, the battery can be withdrawn, so as not to be affected by moisture leakage.

13. Error analysis

(I) When the water content of dike soil is greatly different, the data of penetration resistance will be affected, such as soil

When the water content increases, the penetration resistance decreases. When the water content decreases, the penetration resistance increases. The best solution to reduce the error is to follow up the test immediately after rolling, because the field air drying and sun exposure after rolling will lose a lot of water every hour, directly affecting the accuracy of the dry density value. In actual operation, there are many factors that affect the compacting uniformity of the penetrometer, such as filling

The composition of the material of cohesive soil, sand soil, mixed soil mixing uniformity), porosity, hole hole, aggregate content, slope, water content change and the degree of compaction, near the center of the operating experience of the author, and will directly influence the resistance value, which requires engineering technical personnel in practice through the penetrometer contrast test with cutting ring, Master the water content identification method of soil, all kinds of soil discrimination method and other means, constantly summarize experience, constantly explore, master rules, overcome shortcomings, perfect the actual application technology of the penetration instrument.

(2) The reading value suddenly increases or decreases in the penetration process: the penetration resistance value will be significantly higher or lower in the case of compaction hard block, edge block hole, root and bark crack during the construction site operation. It can prompt the operator to analyze the cause in time, find out the soft weakness and abnormal point of the bottom layer, and timely find loopholes to eliminate the hidden dangers of the project. The characteristics and uniqueness of the penetrometer just make up for the limitation of local sampling of the ring knife. The penetrometer can pick and fill points randomly, anytime and anywhere, which is very flexible and not limited by time and place.

(3) When the dry capacity is qualified in the ring knife experiment, but the penetration value of each point is low in the perpenter experiment, the reasons should be analyzed: first, whether the sand content of the soil is high or not will directly affect the dry capacity value; Second, the moisture content of the soil is large, and the penetration value is low. At this time, it can be dried by air and then rolled again after half a day, and the penetration value is increased.

Appendix 2- Calculation of relative dry density values

Given that the maximum dry capacity Rd_{MAX}=1.63 and the design requirement Rd_{Design}=1.58, $Rd_{Any\ number}$ value of any measured point is required to be calculated. First, a set of ring knives are taken to measure $Rd_{measured}$= 1.56, and the penetration value $P_{measured}$=700N is measured at the same point.

(1) Equation:

① $$\frac{P_{design}}{P_{measured}} = \frac{R\,d_{Design}}{Rd_{measured}}$$

Substituting the above parameters into formula 1, we can get:

$$\frac{P_{design}}{100} = \frac{1.58}{1.56}$$

Solve the equation: $P_{design} = \frac{1.58}{1.56} \times 700 = 709(N)$

(2) When the measured $P_{Any\ numbe}$=715 (N), $Rd_{Any\ numbe}$ can be calculated according to the following formula.

② Equation: $$\frac{R\,d_{Any\ numbe}}{R\,d_{Design}} = \frac{P_{any\ numbe}}{P_{design}}$$

Substitute the specific parameters into equation ② and get:

$$\frac{Rd_{Any\ numbe}}{1.58} = \frac{715}{709}$$

Solve the equation: $Rd_{Design} = \frac{715}{709} \times 1.58 = 1.59$

When the measured $P_{Any\ numbe}$ =720N, equation ② is also substituted into:

Dry density data sheet　　Northwest yard offers suitable silt fine sand

P	300	330	350	380	400	410	430	440
Rd	1. 35	1. 38	1. 40	1. 46	1. 49	1. 50	1. 51	1. 52
P	450	470	490	500	520	540	550	570
Rd	1. 53	1. 54	1. 55	1. 56	1. 57	1. 58	1. 59	1. 60
P	590	610	630	650	660	670	690	700
Rd	1. 61	1. 62	1. 63	1. 64	1. 65	1. 66	1. 68	1. 70
P	720	730	750	780	800	850	870	900
Rd	1. 72	1. 74	1. 76	1. 78	1. 80	1. 83	1. 85	1. 86

Penetration depth 8cm~20cm probe ϕ10mm　　Water content is 3%~8%

P penetration resistance value in the table (N)　　Rd= relative dry density value


Dry density data sheet　　North China garden provides suitable silty loam

P	140	150	180	200	230	250	280	300
Rd	1. 42	1. 45	1. 47	1. 50	1. 53	1. 55	1. 57	1. 58
P	330	350	380	400	430	450	470	500
Rd	1. 60	1. 61	1. 63	1. 64	1. 65	1. 66	1. 67	1. 68
P	530	550	580	600	630	650	690	740
Rd	1. 69	1. 70	1. 71	1. 72	1. 73	1. 74	1. 75	1. 76
P	780	820	830	860				
Rd	1. 77	1. 78	1. 79	1. 80				

Penetration depth 8cm~20cm probe ϕ10mm　　Water content is 20%~27%

P penetration resistance value in the table (N)　　Rd= relative dry density value

14. Operation buttons

Operation buttons

4-13 (UNIT)

It is a unit conversion key for unit conversion operations.

4-14 (POWER)

It is the power switch, the return key in data browsing mode, and the action key to call up automatic shutdown Settings.


4-15 (PEAK)

It is a switch button between peak holding mode and real-time measurement mode.

4-16 (ZERO)

Instrument zero operation button.

15. Display function description

Screen display

4-9 Battery indicator " "

If the battery voltage is too low, " " appears in the upper left corner of the screen, indicating that the battery voltage is insufficient and the battery needs to be replaced.

4-10 Peak indicator MAX.

When MAX. Is displayed, it indicates the peak holding mode, and the screen displays the peak value. If MAX. Does not appear, it indicates the real-time measurement mode, and the value displayed on the screen changes with the load.

4-11 Units of measurement

Indicates the current measurement unit, including KGF, N, and LBF. Only one of these units is displayed.

Appendix 1- Dry density data table

Dry density data sheet　　Northeast courtyard offers suitable clay soil

P	145	150	170	200	230	250	280	320
Rd	1. 40	1. 42	1. 44	1. 46	1. 48	1. 50	1. 52	1. 54
P	350	380	400	430	450	480	500	540
Rd	1. 55	1. 56	1. 57	1. 58	1. 59	1. 60	1. 61	1. 62
P	560	580	600	630	650	700	750	800
Rd	1. 63	1. 64	1. 65	1. 66	1. 67	1. 68	1. 69	1. 70
P	820	830	840	850	860	880	900	
Rd	1. 74	1. 72	1. 73	1. 74	1. 76	1. 78	1. 80	

Penetration depth 8cm~20cm probe ϕ10mm　　Water content is 18%~26%

P penetration resistance value in the table (N)　　Rd= relative dry density value

Dry density data sheet　　South Central institute provides suitable clay

P	140	150	180	200	230	250	280	300
Rd	1. 42	1. 45	1. 47	1. 50	1. 53	1. 55	1. 57	1. 58
P	330	350	380	400	430	450	470	500
Rd	1. 60	1. 61	1. 63	1. 64	1. 65	1. 66	1. 67	1. 68
P	530	550	580	600	630	650	690	740
Rd	1. 69	1. 70	1. 71	1. 72	1. 73	1. 74	1. 75	1. 76
P	780	820	830	860				
Rd	1. 77	1. 78	1. 79	1. 80				

Penetration depth 8cm~20cm probe ϕ10mm　　Water content is 15%~21%

P penetration resistance value in the table (N)　　Rd= relative dry density value