

BackProductsHeathrow Scientific ProductsHS8881C



Key Documents

COACOO

View All Documentation

Skip To

Properties

Description

Compare Similar Items

Documentation

HS8881C

Share

Vinyl-Coated Lead Ring ("C" shape)

fits beaker size, 500 to 2000 mL, blue

Promotion Available

Sign In to View Organizational & Contract Pricing.

Pack Size	SKU	Availability	Price
1 ea	HS8881C-1EA	Available to ship TODAY from MILWAUKEE	\$43.90

About This Item

NACRES: NB.51

UNSPSC Code: 41121800

HS8881C-1EA

\$43.90

Sign In for Best Price

Available to ship TODAYDetails

1

Add to Cart

Technical Service

Need help? Our team of experienced scientists is here for you.

Let Us Assist

Recommended Products

HS8881A

Vinyl-Coated Lead Ring ("C" shape)

Quick View

HS8881B

Vinyl-Coated Lead Ring ("C" shape)

Quick View

HS8882C

Vinyl-Coated Lead Ring ("O" shape)

Quick View

HS8882A

Vinyl-Coated Lead Ring ("O" shape)

Quick View

Properties

material	blue
feature	fits beaker size 500 to 2000 mL
packaging	pack of 1 ea
manufacturer/tradename	Heathrow Scientific HS8881C

Description

General description

The Lead Ring is perfect for stabilizing vessels and labware or for water bath immersions.

Features and Benefits

- Vinyl-coated lead ring will not mark or scratch benchtops and labware
- O-Shape Style

Compare Similar Items

View Full Comparison

Show Differences

This Item

HS8881C

Vinyl-Coated Lead Ring ("C" shape)

Quick View

manufacturer/tradename

Heathrow Scientific HS8881C

material

blue

HS8881A

Vinyl-Coated Lead Ring ("C" shape)

Quick View

manufacturer/tradename

Heathrow Scientific HS8881A

material

blue

HS8881B

Vinyl-Coated Lead Ring ("C" shape)

Quick View

manufacturer/tradename

Heathrow Scientific HS8881B

material

blue

HS8881D

Vinyl-Coated Lead Ring ("C" shape)

Quick View

manufacturer/tradename

Heathrow Scientific HS8881D

material

blue



This Item packaging pack of 1 ea	packaging pack of 1 ea	packaging pack of 1 ea	packaging pack of 1 ea
feature fits beaker size 500 to 2000 mL	feature fits beaker size 125 to 500 mL	feature fits beaker size 250 to 1000 mL	feature fits beaker size 1000 to 4000 mL

Still not finding the right product?

Explore all of our products under [Vinyl-Coated Lead Ring \("C" shape\)](#)

[View All Products](#)

Documentation

[Certificate Of Analysis](#) [Certificate Of Origin](#)

Choose from one of the most recent versions:

Certificates of Analysis (COA)
Lot/Batch Number
It looks like we've run into a problem, but you can still download Certificates of Analysis from our [Documents](#) section.

If you need assistance, please contact [Customer Support](#)

Already Own This Product?

Find documentation for the products that you have recently purchased in the Document Library.

[Visit the Document Library](#)

Related Categories

[Heathrow Scientific Products](#)
[Labware Partners](#)

Global Trade Item Number

SKU	GTIN
HS8881C-1EA	04065265823050

Recently Viewed Products



CLS430828
Corning® 50 mL centrifuge tubes

Quick View



CLS352070
Corning® Falcon® Centrifuge Tubes

Quick View



55023
Aura Mobile Phase Filter/Degasser...

Quick View



Sigma-Aldrich
8.00133
Barbituric acid

Quick View

Support

[Customer Support](#)
[Contact Us](#)
[FAQ](#)
[Safety Data Sheets \(SDS\)](#)
[Certificates \(COA/COO\)](#)
[Quality & Regulatory](#)
[Calculators & Apps](#)

Orders

[Order Lookup](#)
[Quick Order](#)
[Custom Products](#)
[eCommerce Solutions](#)

Company

[About Us](#)
[Responsibility](#)
[Events](#)
[Press Releases](#)
[Programs](#)
[Careers](#)
[Offices](#)

Social Media



MilliporeSigma

Research. Development. Production.
We are a leading supplier to the global Life Science industry with solutions and services for research, biotechnology development and production, and pharmaceutical drug therapy development and production.



Sigma-Aldrich®
Solutions

BioReliance®
Solutions

Millipore®
Solutions

SAFC®
Solutions

Milli-Q®
Solutions

Supelco®
Solutions