



[Home](#) ▶ [Shop All Products](#) ▶ [pH and Electrochemistry](#) ▶ [Probes and Electrodes](#) ▶ [Combination Ion Selective Electrodes](#) ▶

Orion™ Iodide Ion Selective Elect...



Thermo Scientific™

Orion™ Iodide Ion Selective Electrodes

Get fast, stable and reliable iodide ions readings in aqueous solutions.

Have Questions? [Contact Us](#)

Change view



Catalog Number	Connector Types
9653BNWP	BNC (Waterproof)
9453BN	BNC

Catalog number 9653BNWP



Price (USD)

-

[Request A Quote](#)

Connector Types: BNC (Waterproof)

Product Overview



Ask our AI about this Product

Recommendations

Documents

Get fast and stable readings using Thermo Scientific™ Orion™ Iodide Electrodes with built-in Thermo Scientific™ Sure-Flow™ reference. Combination design decreases the amount of required sample and reduces waste.

Features:

- 9653BNWP is a combination electrode with a Sure-Flow™ junction that allows easy cleaning. When conventional electrode junctions clog, the Sure-Flow junction can be opened and flushed clean simply by pressing the cap.
- Iodide Half-Cell Electrode 9453BN is to be used with the Orion 900200 sure flow double Junction electrode.

Solutions required:

- Standard**
 - 945306 0.1 M iodide standard, 475 mL bottle
- ISA**
 - 940011 ISA for iodide measurements, 475 mL bottle
- Electrode fill solution**
 - 900063 Optimum Results D electrode filling solution, 5 x 60 mL bottles
 - 900002 Inner fill solution for double junction ref electrode, 5 x 60 mL bottles
 - 900003 Outer fill solution for double junction ref electrode, 5 x 60 mL bottles

Items included with:

- 9653BNWP**
 - 900063 Optimum Results D electrode filling solution, 60 mL 1 bottle
 - 948201 Polishing strip. Qty 2.

9453BN

- 948201 Polishing strip. Qty 1.

Warranty period

- 12 months from date of purchase

Specifications

Description	With waterproof BNC connector
Connector Types	BNC (Waterproof)
For Use With (Equipment)	ISE Meter
Type	Combination ISE
Unit Size	Each

Give Feedback

