SMEDENT 上海瑞楼医疗器械有限公司
SHANGHAI SMEDENT MEDICAL INSTRUMENT CO.,LTD

To Whom It May Concern

Ref: Determination of re-using of Dental Diamond Burs

We Shanghai Smedent Medical Instrument Co., Ltd,
with registered address in Room 111, Building E, No. 681, North Huifeng Road, Fengxian, Shanghai, China; certifies that our range of non-sterile **Diamond dental burs** can be subjected to a sterilization process by using vapor sterilizer at 134°C a maximum of three (3) re-process, guaranteeing an adequate cutting efficiency and resistance of the materials.

Shanghai Smedent Medical Instrument Co., Ltd cannot determinate a number of re-uses, It depends on users habit, operating method, polishing target, using environment and other factors. The re-use is normally determined by wear and tear and not through sterilization re-processing. It means that is responsibility of the re-processor routinely inspecting diamond burs before to clean, to sterilize, to use and to re-use, in this case the re-processer is the responsible to determinate the re-use of the device according to parameters inspected and their results.

We are responsible for the quality of said products.

Yours Truly;

Signed: William Fu
For and on behalf of SHANGHAI SMEDENT MEDICAL INSTRUMENT CO.,LTD
Authorized Signature(s)
Date: 2019-07-18

Add : Room 111, Building E, No. 681, North Huifeng Road, Fengxian, Shanghai, China
Tel: 0086 21 34096625 Fax: 0086 21 3409665
Website: http://www.smedentalburs.com