Medical Body Thermometer

KF-HW-001



6 major features

- Precision probe
- One-button
- 50 sets of memory
- Lcd display
- Tri-color backlight
- Multi-purpose




1s infrared thermometer

- ☑ Infrared measurement
- ☑ 1s temperature measurement fast and accurate
- ☑ Temperature measurement at the touch of a button

German imports probe

Smart calibration
Accurate measurement




SMART THREE COLOR
BACKLIGHT REMINDER
红外线体温计 型号：KF-HW-001
SN:P19125520570
体温
36.8
37.8
38.2
Normal
High
Fever



One-button restore
Any key turn on

One-button factory settings can be restored; any key can turn on, the elderly
Children will use, smart keys more peace of mind!



Align the forehead and press the trigger
to measure the temperature

Measuring distancet 1-3cm

MULTIPLE MODE SWITCHING

Multiple uses for one machine, causing harm to the human body from low or high temperatures

SET

Button “SET”

Set other parameters

MEMO

Button "MEMO"

50 sets of memory value query

MODE

Hold “MODE”

Switch body temperature & object temperature



Measure body temperature



Water temperature



Milk temperature




Room temperature

50 GROUPS OF MEMORY VALUES

RECORDING THERMOMETER

Real-time recording of baby's body temperature, large-capacity memory, compare temperature before and after




AUTOMATIC SHUT-DOWN

Automatic shutdown after 15S without operation
Save energy and save power




LCD Display Screen
Long press restore factory settings
Short press to switch measurement mode
Memory value query key
Mute key
Power on/temperature button
Streamlined handle
Battery compartment
红外线体温计 型号：KF-HW-001
SN:P19125520570



PRODUCT INFORMATION			
Brand	Cofoe	Product Name	Cofoe Non-contact IR Digital
Model No	HTD8818D	Battery	DC3V,2*AAA Alkaline Battery
Memory	50 sets	Size	138*95*40mm
Measurement time	≤2 second	Weight	90g
Measurement error	± 0.1- ± 0.3 degrees Celsius		
Measurement Range	34-42.9		
Automatic shutdown	≤18 second		
Measurement Distance	1-3cm		
Different Color	Backlight with Three Different Color Green:34-37.3 Orange:37.4-38.0 Red:38.1-42.9		



Battery installation

Open the thermometer battery cover,

will be two on the 7th

As shown in the figure, the positive and

negative electrodes are placed.

Close the door to measure