DUO-FLOW® SOFT-LINE®

SHORT-TERM HEMODIALYSIS CATHETER



FEATURES & BENEFITS

INTERNAL LUMEN DESIGN
Coaxial

MATERIAL
Polyurethane

TIP DESIGN
Tapered

FRENCH SIZE
7F, 9F, 11.5F

CONFIGURATION
Straight Pre-Curved



DUO-FLOW® SOFT-LINE®

SHORT-TERM HEMODIALYSIS CATHETER

7F BASIC SET		
T73M	7F X 7CM DUO-FLOW® SOFT-LINE®, STRAIGHT	5/BOX
T74M	7F X 10CM DUO-FLOW® SOFT-LINE®, STRAIGHT	5/BOX
9F BASIC SET		
T94M	9F X 12CM DUO-FLOW® SOFT-LINE®, STRAIGHT	5/BOX
T96M	9F X 15CM DUO-FLOW® SOFT-LINE®, STRAIGHT	5/BOX
T98M	9F X 20CM DUO-FLOW® SOFT-LINE®, STRAIGHT	5/BOX
11F BASIC SET		
T114ME	11.5F X 12CM DUO-FLOW® SOFT-LINE®, STRAIGHT	10/BOX
T116ME	11.5F X 15CM DUO-FLOW® SOFT-LINE®, STRAIGHT	10/BOX
T118ME	11.5F X 20CM DUO-FLOW® SOFT-LINE®, STRAIGHT	10/BOX
T119M	11.5F X 24CM DUO-FLOW® SOFT-LINE®, STRAIGHT	5/BOX
T116IJS-2E.	11.5F X 15CM DUO-FLOW® SOFT-LINE®, PRE-CURVED	10/BOX
T118IJS-2E.	11.5F X 20CM DUO-FLOW® SOFT-LINE®, PRE-CURVED	10/BOX

Basic Set Contents: (1) Catheter (1) Guidewire (1) Introducer Needle (1|3*) Dilator (2) End Cap (1) Adhesive Wound Dressing**

*T119M Only
**7F & 9F Sets, T119M Only

IMPORTANT RISK INFORMATION

Indications for Use: The Medcomp® Soft-Line® Catheter is indicated for use in attaining Short-Term vascular access for Hemodialysis and Apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. Alternate insertion sites include subclavian vein or femoral vein as required. The curved Soft- Line® Catheter is intended for internal jugular vein insertion. This catheter is indicated for a duration less than (30) days. For femoral placement, monitor catheter condition closely.

Contraindications: This catheter is intended for Short-Term vascular access only and should not be used for any purpose other than indicated in these instructions.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation

Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438, U.S.A.
P: 215.256.4201 **F**: 215.256.1787
www.**medcomp**.net

ISO13485
Certified Quality System

PN2726BSI Rev. A 5/24