Geuder®
Precision made in Germany

Product description

Product	**Single-use Tubing Sets, sterile**
Article No.	G-19237, G-24203, G-28091, G-28097, G-28098, G-28136, G-28148, G-30114, G-30116, G-30118, G-30146, G-32696, G-46000
Remarks	

GMDN-Code	47625 surgical irrigation, tubing set
UMDNS-Code	15-277 Schlauch, Silikon
PN-Code CA	

Product description
Single-use tubing systems from PVC are used in combination with Geuder ophthalmic surgical systems. Field of application is irrigation and aspiration of aqueous solutions, air injections and silicone oil into the eye.


Product variants

Article No. G-	Description
19237	SINGLE USE TUBING SET with pressure transducer for the intraocular pressure gauge Model Düsseldorf G- 19235, sterile
24203	PHACO-ASPITRON single use tubing set with infusion set, sterile
28091	SINGLE USE HPS I/A TUBING SYSTEM for MEGATRON S3 P / S4 HPS with irrigation set and GEUDER safety-connector 10 pcs. per box, sterile
28097	SINGLE USE TUBING SET for MEGATRON, standard with infusion set, packaging unit 10 pieces per box, sterile
28098	SINGLE USE TUBING SET with infusion set for MEGATRON (external pressure transducer) and MEGATRON S3 P, packaging unit 10 pieces per box, sterile
28136	SINGLE USE TUBING SET with INFUSION SET for MEGATRON (external pressure transducer) and MEGATRON S3 P, bypass with filter to VAA valve, 1.9m length, packaging unit 10 pieces per box, sterile

28148	MEGATRON SINGLE-USE TUBE incl. filter 0.1µ for air injection, length 1.90 m, 5 per box, sterile
30114	MEGATRON S3 V and S3 VIP- SINGLE USE TUBING SET for Venturi cassette, packaging unit 10 pieces per box, sterile
30116	SINGLE-USE HPS I/A TUBING SYSTEM, for MEGATRON S3 V / VIP / S4 HPS with GEUDER safety-connector, 20 pcs. per box, sterile
30118	SINGLE-USE IRRIGATION SET for ophthalmic surgery systems 20 pcs. per box, sterile
30146	SINGLE-USE AIR INJECTION TUBE for MEGATRON S3 / S4 HPS incl. filter 0.1 µ and 2-way-stopcock, length 1.9 m, 6 pcs. per box, sterile
32696	SINGLE-USE PRESSURE TUBE for injection of viscous fluid, Luer-Lock female/male, 10 pcs. per box, sterile
46000	SINGLE-USE ASPIRATION TUBE for vitreous cutters, active backflush handpiece and oil-explantation with GEUDER safety-connector 6 pcs. per box, sterile

Accessories	
Article No. G-	**Description**

Pictures (example)



G- 28136 SINGLE USE TUBING SET with INFUSION SET for MEGATRON (external pressure transducer) and MEGATRON S3 P, bypass with filter to VAA valve, 1.9 m length, packaging unit 10 pieces per box, sterile

Geuder®
Precision made in Germany

Production plants and Marketing authorisation holder (MAH)		
Article No. G-	**Production plant**	**MAH**
	Dimedica GmbH & Co. KG Strasse der Heimkehrer 34 D- 64646 Heppenheim	GEUDER AG Hertzstraße 4 D-69126 Heidelberg

Material in contact with the patient	
Component	**Material**
Tube	PVC

Sterilization	
Sterilization method	Ethylene oxide
Sterilization plant	Sterigenics

Shelf-Life	5 years

Document changes		
Version	**Date**	**Description**
V1.0	December 18 2009	First release
V2.0	August 7 2014	English version

Drafted by	**Reviewed by**	**Released by**
Name: Dr. Anaid Shahbazian, QM/RA	Name: Renate Herbold, QM/RA	Name: Dieter Frauenfeld, head of QM/RA
Date/Signature Anaid Shahbazian 07. AUG. 2014	Date/Signature Renate Herbold 30. SEP. 2014	Date/Signature D. Frauenfeld 13. OKT. 2014